

METROPOLISEDOS_Chennai - Avadi

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
#0002	1p/19q Deletion by FISH Oligodendrogliomas Tissue	Chennai - Avadi	D	14750	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
#0003	3-Methyl Histidine Quantitative Urine	Chennai - Avadi	D	7020	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
#0004	3-Methyl Histidine Quantitative Serum	Chennai - Avadi	R	6365	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
#0006	5-HIAA, Urine 24 hrs by column chormatography	Chennai - Avadi	D	3305	Tue & Fri: 11:30am	Same day : 9:00pm	24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliqout to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
#0007	5q del (MDS) by Fish	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0008	5q deletion/Monosomy 5 by FISH MDS Bone Marrow	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0011	7q deletion/Monosomy 7 by FISH,MDS	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
#0012	7q deletion/Monosomy 7 by FISH MDS Bone Marrow	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0013	Williams Syndrome; Microdeletion of chromosome 7q11.23**	Chennai - Avadi	R	11235	Daily : 11:30am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
#0014	FGFR 1 gene rearrangement t(8p11.2), Blood**	Chennai - Avadi	E	4725	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0015	FGFR 1 gene rearrangement t(8p11.2), Bone Marrow**	Chennai - Avadi	E	4725	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0017	11q (ATM) by FISH	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0018	11q deletion by FISH ATM (CLL) Bone Marrow	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0021	13q del by FISH	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0022	13q deletion by FISH MM, CLL Bone Marrow	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0023	14q32.3 by FISH IGH gene rearrangement (ALL, NHL) Blood**	Chennai - Avadi	D	6490	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0024	14q32.3 by FISH, IGH gene rearrangement (ALL, NHL) Bone Marrow**	Chennai - Avadi	D	6490	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0025	Prader-Willi/Angeman Syndrome (Microdeletion 15q9-13)**	Chennai - Avadi	R	6315	Daily : 11:30am	12th Day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
#0026	17 OH Corticosteroids, Urine by Column Chromatography (#0026)	Chennai - Avadi	D	4130	Tue & Fri: 11:30am	Same day : 9:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
#0027	17-Hydroxy Progesterone Neonatal Screen Dried blood spot	Chennai - Avadi	C	510	Daily : 11:30am	Next day : 5:00pm	accurate results and such history to be conveyed to lab. Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
#0028	17-Hydroxy-Progesterone by Serum	Chennai - Avadi	C	1580	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
#0029	17 Ketosteroids,urine by Column Chromatography (#0029)	Chennai - Avadi	C	3955	Tue & Fri: 11:30am	Same day : 9:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
#0030	Del 17p53 by FISH	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0031	17p deletion by FISH, p53 (MM, CLL) Blood	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0034	20q deletion (MDS) by FISH, Blood FISH probe for chromossome 20q Blood	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0035	Del 20q (MDS) Blood / bone marrow	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0036	DEL 22q (Microdeletion digeorge syndrome) by FISH	Chennai - Avadi	D	5665	Daily : 9:00am	7th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
#1018	5 Nucleotidase activity , Serum**	Chennai - Avadi	R	5350	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
A0001	Abnormal Haemoglobin Studies, Hb Variants, EDTA Blood	Chennai - Avadi	C	1060	Daily : 9:00am	Same day : 7:00pm	3 ml EDTA whole blood
A0003	Absolute Eosinophile Count	Chennai - Avadi	A	205	Daily : 9:00am to 9:00pm	After 4hrs	3 ml EDTA whole blood
A0008	Acetone Ketone by Fully Automated Urine chemistry & Microscopy - Urine	Chennai - Avadi	A	110	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
A0012	AChR Antibodies Acetyl Choline Receptor Antibodies Serum	Chennai - Avadi	C	3505	Daily : 11:30am	Next day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0013	Acid Phosphatase Total Serum	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0014	Acid Phosphatase Total and prostatic fraction Serum	Chennai - Avadi	C	565	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0015	ACTH Stimulation Test - For 17 OH progesterone -Serum	Chennai - Avadi	C	3730	Daily : 11:30am	After 8hrs	3 ml of serum (0 min, 30, 60 min)

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0016	ACTH Stimulation Test - For Aldosterone - Serum	Chennai - Avadi	C	5650	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0 min, 30, 60 min)
A0017	ACTH Stimulation Test - For Androstenedione - Serum	Chennai - Avadi	C	4635	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0 min, 30, 60 min)
A0018	ACTH Stimulation Test - For Cortisol - Serum	Chennai - Avadi	C	1810	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0 min, 30, 60 min)
A0019	ACTH Stimulation Test - For DHEAS-Serum	Chennai - Avadi	C	2825	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0 min, 30, 60 min)
A0020	ACTH Adreno Corticotropic Hormone - Plasma	Chennai - Avadi	C	1750	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of EDTA plasma .
A0021	Actinomyces Culture by Vitek/ MALDITOF - Sputum	Chennai - Avadi	D	1415	Daily : 9:00am to 9:00pm	5th day : 6:00pm	2ml of Sputum in Sterile leak proof container
A0022	Actinomyces Culture by Vitek/ MALDITOF - Tissue	Chennai - Avadi	D	1415	Daily : 9:00am to 9:00pm	5th day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
A0023	Actitest and Fibrotest	Chennai - Avadi	E	15750	Daily : 7:30am	5th day : 6:00pm	5 ml of serum
A0024	Acylcarnitine quantitative Neonatal**	Chennai - Avadi	R	4815	Daily : 11:30am	8th day : 6:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							history and sample to be collected post 72 hrs after transfusion.
A0025	ADA - Adenosine Deaminase by Biochemical Method - Ascitic Fluid	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of Fluid in Sterile leak proof container
A0026	ADA - Adenosine Deaminase by Biochemical Method - CSF	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of CSF in Sterile leak proof tube or container
A0027	ADA (Adenosine Deaminase), Pericardial fluid	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of Fluid in Sterile leak proof container
A0028	ADA - Adenosine Deaminase by Biochemical Method- Peritoneal Fluid	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of Fluid in Sterile leak proof container
A0029	ADA - Adenosine Deaminase by Biochemical Method - Pleural Fluid	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of Fluid in Sterile leak proof container
A0030	ADA - Adenosine Deaminase by Biochemical Method - Pus	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 6hrs	Sample in sterile container
A0031	ADA - Adenosine Deaminase by Biochemical Method - Serum	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0032	Adenovirus DNA Detection by PCR - Respiratory Fluid **	Chennai - Avadi	D	4955	Daily : 11:30am	10th day : 6:00pm	3 ml of Fluid in Sterile leak proof container
A0033	Adenovirus DNA Detection by PCR - Oral swab**	Chennai - Avadi	D	4955	Daily : 11:30am	10th day : 6:00pm	1 Swab with Sufficient Material in Sterile Vacutainer
A0034	Adenovirus DNA Detection by PCR - Eye swab**	Chennai - Avadi	D	4955	Daily : 11:30am	10th day : 6:00pm	1 Swab with Sufficient Material in Sterile Vacutainer
A0035	Adenovirus DNA Detection by PCR - Nasal swab **	Chennai - Avadi	D	4955	Daily : 11:30am	10th day : 6:00pm	1 Swab with Sufficient Material in Sterile Vacutainer
A0036	Adenovirus IgG Antibody by EIA - Serum	Chennai - Avadi	D	1770	Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0037	Adenovirus IgM Antibody by EIA - Serum	Chennai - Avadi	D	1770	Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0038	ADH-Anti Diuretic Hormone Vasopressin Plasma	Chennai - Avadi	D	4840	Mon : 11:30am	Next day : 5:00pm	3 ml of EDTA plasma to be sent in aliquot tube at proper temperature.
A0039	Adiponectin ELISA Serum	Chennai - Avadi	D	2950	Wed : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0041	Adrenaline (Epinephrine) Plasma	Chennai - Avadi	D	3420	Wed : 11:30am	Next day : 5:00pm	Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0042	Adrenaline (Epinephrine] HPLC,Urine24hrs	Chennai - Avadi	D	3420	Wed : 11:30am	Next day : 5:00pm	diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history. 24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
A0052	AFB Rapid culture by Bactec Fx - Blood	Chennai - Avadi	C	1185	Daily : 9:00am	After 6 weeks : 6:00pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
A0053	AFB Rapid culture by MGIT - Body Fluid	Chennai - Avadi	C	1185	Daily : 9:00am	After 6 weeks : 6:00pm	2ml of Body Fluid in Sterile leak proof container
A0054	AFB Rapid culture by Bactec Fx - Bone Marrow	Chennai - Avadi	C	1185	Daily : 9:00am	After 6 weeks : 6:00pm	3ml of Bone Marrow in MycoF/lytic Bactec Bottle or Sodium Heparin
A0055	AFB Rapid culture by MGIT - BAL / Other Respiratory samples	Chennai - Avadi	C	1105	Daily : 9:00am	After 6 weeks : 6:00pm	2ml of BAL / Other Respiratory samples in Sterile leak proof container
A0056	AFB Rapid culture by MGIT - CSF	Chennai - Avadi	C	1185	Daily : 9:00am	After 6 weeks : 6:00pm	2ml of CSF in Sterile leak proof tube or container
A0057	AFB Rapid culture by MGIT - Endometrium	Chennai	C	1185	Daily : 9:00am	After 6	2ml of Endometrium in Sterile leak

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				weeks : 6:00pm	proof container
A0058	AFB Rapid culture by MGIT - Other	Chennai - Avadi	C	1185	Daily : 9:00am	After 6 weeks : 6:00pm	2ml of Specimen in Sterile leak proof container
A0059	AFB Rapid culture by MGIT - Pus, abscess and aspirates	Chennai - Avadi	C	1185	Daily : 9:00am	After 6 weeks : 6:00pm	2ml of Pus / Abscess / Aspirate in Sterile leak proof container
A0060	AFB Rapid Culture by MGIT - Sputum	Chennai - Avadi	C	1185	Daily : 9:00am	After 6 weeks : 6:00pm	2ml of Sputum in Sterile leak proof container
A0061	AFB Rapid culture by MGIT - Tissue/Biopsy	Chennai - Avadi	C	1185	Daily : 9:00am	After 6 weeks : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
A0062	AFB Rapid Culture by MGIT - Urine	Chennai - Avadi	C	1185	Daily : 9:00am	After 6 weeks : 6:00pm	15ml of Spot Urine in Sterile leak proof container
A0063	MOTT Antibiogram for Rapid Growers by Microbroth dilution - Pure Culture	Chennai - Avadi	D	8850	Daily : 9:00am	After 4 weeks : 6:00pm	Isolated colony on culture plates
A0064	MOTT Antibiogram for Slow Growers by Microbroth dilution - Pure Culture	Chennai - Avadi	D	10030	Daily : 9:00am	After 4 weeks : 6:00pm	Isolated colony on culture plates
A0065	AFB (TB) - Antibiogram for MTb complex - Amikacin Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0066	AFB (TB) - Antibiogram for MTb complex - Capreomycin Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0067	AFB (TB) - Antibiogram for MTb complex - Clofazimine Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0068	AFB (TB) - Antibiogram for MTb complex - Ethambutol Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0069	AFB (TB) - Antibiogram for MTb complex - Ethionamide Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0070	AFB (TB) - Antibiogram for MTb complex - Isoniazid Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0071	AFB (TB) - Antibiogram for MTb complex - Kanamycin Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0072	AFB (TB) - Antibiogram for MTb complex - Levofloxacin Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0073	AFB (TB) - Antibiogram for MTb complex - Moxifloxacin Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0074	AFB (TB) - Antibiogram for MTb complex - Ofloxacin Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0075	AFB (TB) - Antibiogram for MTb complex - PAS Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0076	AFB (TB) - Antibiogram for MTb complex - Pyrazinamide Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2940	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0077	AFB (TB) - Antibiogram for MTb complex - Rifabutin Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0078	AFB (TB) - Antibiogram for MTb complex - Rifampicin Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0079	AFB (TB) - Antibiogram for MTb complex - Streptomycin Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0080	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Sputum/ BAL	Chennai - Avadi	C	2205	Daily : 7:30am	Next day : 5:00pm	2ml of Sputum/ BAL/ Other respiratory samples in Sterile leak proof container
A0081	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Body Fluids	Chennai - Avadi	C	2205	Daily : 7:30am	Next day : 5:00pm	2ml of Body Fluid in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0082	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - CSF	Chennai - Avadi	C	2205	Daily : 7:30am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
A0083	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Pus / Abscess / Aspirate	Chennai - Avadi	C	2205	Daily : 7:30am	Next day : 5:00pm	2ml of Pus / Abscess / Aspirate in Sterile leak proof container
A0084	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Other sample type	Chennai - Avadi	C	2205	Daily : 7:30am	Next day : 5:00pm	2ml of Specimen in Sterile leak proof container
A0085	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Tissue / Biopsy /Lymph nodes/Endometrium	Chennai - Avadi	C	2205	Daily : 7:30am	Next day : 5:00pm	Tissue/ Endometrium with Normal Saline/ Biopsy/ Lymph nodes in Sterile leak proof container
A0086	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Urine	Chennai - Avadi	C	2205	Daily : 7:30am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
A0087	AFB Drug susceptibility - 1st and 2nd Line (10 drugs) panel by MGIT - Pure Culture	Chennai - Avadi	D	11210	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0088	AFB Drug susceptibility - 1st Line (4 drugs) SIRE panel by MGIT - Pure Culture	Chennai - Avadi	D	5545	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0089	AFB Drug susceptibility - 2nd Line (4 Drugs) KEPO panel by MGIT - Pure Culture	Chennai - Avadi	D	6845	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0090	AFB Drug susceptibility - Additional drugs 4 Nos. MACC panel by MGIT - Pure Culture	Chennai - Avadi	D	6610	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0091	AFB Drug susceptibility - Comprehensive 13 drugs panel by MGIT - Pure Culture	Chennai - Avadi	D	15340	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0092	AFB Identification & Speciation MTB and MOTT by Hains Line Probe - Pure culture	Chennai - Avadi	D	3480	Tue : 9:00am	Next day : 6:00pm	Pure growth of Mycobacterium spp. on LJ slant/ MGIT Tube/ Bactec alert bottle.
A0093	AFB MDR Screen by Hains Line Probe Assay - Sputum / BAL/ Other respiratory	Chennai - Avadi	D	2125	Tue : 9:00am	Next day : 6:00pm	2ml of Sputum/BAL/other respiratory samples in Sterile leak

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	samples						proof container
A0094	AFB MDR Screen by Hains Line Probe Assay - Extra pulmonary samples	Chennai - Avadi	D	2125	Tue : 9:00am	Next day : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, Synovial fluid)
A0102	AFB-Smear examination by Fluorescent stain - Body fluid	Chennai - Avadi	C	565	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Body Fluid in Sterile leak proof container
A0103	AFB-Smear examination by Fluorescent stain - BAL / Other Respiratory Sample	Chennai - Avadi	C	565	Daily : 9:00am & 3:00pm	After 6hrs	2ml of BAL / Other Respiratory samples in Sterile leak proof container
A0104	AFB-Smear examination by Fluorescent stain - Pus, abscess and aspirates	Chennai - Avadi	C	565	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Pus / Abscess / Aspirate in Sterile leak proof container
A0105_1	AFB-Smear examination by Fluorescent stain - Sputum (1 sample)	Chennai - Avadi	C	510	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Sputum in Sterile leak proof container
A0105_3	AFB-Smear examination by Fluorescent stain - Sputum (3 sample)	Chennai - Avadi	C	1245	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Sputum in Sterile leak proof container
A0106	AFB-Smear examination by Fluorescent stain - Tissue/Biopsy	Chennai - Avadi	C	680	Daily : 9:00am & 3:00pm	After 6hrs	Tissue with Normal Saline / Biopsy / Lymph nodes in Sterile leak proof container
A0107_1	AFB-Smear examination by Fluorescent stain - Urine (1 sample)	Chennai - Avadi	C	510	Daily : 9:00am & 3:00pm	After 6hrs	15ml of Spot Urine in Sterile leak proof container
A0107_3	AFB-Smear examination by Fluorescent stain - Urine (3 sample)	Chennai - Avadi	C	1245	Daily : 9:00am & 3:00pm	After 6hrs	15ml of Spot Urine in Sterile leak proof container
A0108	AFB-Smear examination by ZN Stain - Body fluids	Chennai - Avadi	B	295	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Body Fluid in Sterile leak proof container
A0109	AFB- Smear examination by ZN Stain - BAL / Other Respiratory samples	Chennai - Avadi	B	295	Daily : 9:00am & 3:00pm	After 6hrs	2ml of BAL / Other Respiratory samples in Sterile leak proof container
A0110	AFB-Smear examination by ZN Stain - Pus, abscess and aspirates	Chennai - Avadi	B	295	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Pus / Abscess / Aspirate in Sterile leak proof container
A0111_1	AFB-Smear examination by ZN Stain - Sputum (1 sample)	Chennai - Avadi	B	295	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Sputum in Sterile leak proof container
A0111_3	AFB-Smear examination by ZN Stain - Sputum (3 samples)	Chennai - Avadi	B	755	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Sputum in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0112	AFB- Smear examination by ZN Stain - Tissue/Biopsy	Chennai - Avadi	B	325	Daily : 9:00am & 3:00pm	After 6hrs	Tissue with Normal Saline in Sterile leak proof container
A0113_1	AFB- Smear examination by ZN Stain - Spot Urine (Single Sample)	Chennai - Avadi	B	295	Daily : 9:00am & 3:00pm	After 6hrs	15ml of Spot Urine in Sterile leak proof container
A0113_3	AFB-Smear examination by ZN Stain - Urine (3 samples)	Chennai - Avadi	B	755	Daily : 9:00am & 3:00pm	After 6hrs	15ml of Spot Urine in Sterile leak proof container
A0114	AFB Speciation (MOTT) Additional species by Hains Line Probe Assay - Pure culture	Chennai - Avadi	C	3955	Tue : 9:00am	Next day : 6:00pm	Pure growth of Mycobacterium spp. on LJ slant/ MGIT Tube/ Bactec alert bottle.
A0115	AFB XDR Screen by Hains Line Probe Assay - Sputum / BAL/ Other respiratory samples	Chennai - Avadi	C	5650	Tue : 9:00am	Next day : 6:00pm	2ml of Sputum/BAL/other respiratory samples in Sterile leak proof container
A0116	AFB XDR Screen by Hains Line Probe Assay - Extra pulmonary samples	Chennai - Avadi	C	5650	Tue : 9:00am	Next day : 6:00pm	Extra - Pulmonary Smear positive samples or Culture positive specimen
A0117	AFB-Xpert panel (M.Tb Detection & Rifampicin resistance) by CBNAAT - Pulmonary samples	Chennai - Avadi	D	2400	Daily : 7:30am	Same day : 6:00pm	2ml of Pulmonary specimen in Sterile leak proof container
A0118	AFB-Xpert panel (M.Tb Detection & Rifampicin resistance) by CBNAAT - Extra pulmonary samples	Chennai - Avadi	D	2310	Daily : 7:30am	Same day : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
A0119	AFP-Alpha Feto Protein- CLIA- Serum	Chennai - Avadi	B	910	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0120	Alpha Fetoprotein - Amniotic Fluid (AFP-AF)	Chennai - Avadi	C	1075	Daily : 9:00am to 9:00pm	After 6hrs	2 ml Amniotic Fluid In Sterile container
A0121	AFP-Alpha Feto Protein -CLIA- CSF	Chennai - Avadi	C	1075	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
A0122	Albert stain by smear examination - Throat swab	Chennai - Avadi	B	325	Daily : 9:00am & 3:00pm	After 6hrs	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
A0123	Albumin Qualitative - Urine	Chennai - Avadi	A	110	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0124	Albumin Quantitative Urine 24H	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample.
A0125	Albumin Serum	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0126	Albumin CSF	Chennai - Avadi	C	960	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
A0127	Albumin/Creatinine Ratio, Urine	Chennai - Avadi	A	560	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
A0128	Alcohol Ethanol Urine	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
A0129	Alcohol Ethanol Serum	Chennai - Avadi	C	1075	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0130	Aldehyde Chopra Test for Kala Azar by Chemical Formal Test - Serum**	Chennai - Avadi	R	430	Daily : 11:30am	3rd day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0131	Aldolase enzymatic Serum	Chennai - Avadi	C	1075	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0132	Aldosterone 24H URINE	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							collection, start time to end time and exclude first urine sample and include last voided urine sample.
A0133	Aldosterone ELISA Serum	Chennai - Avadi	C	2205	Daily : 9:00am	Same day after 8 hrs	3 ml Serum (Mention The Position)
A0134	ALK1/EML4 (LUNG CANCER) by FISH	Chennai - Avadi	E	8400	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
A0135	Alkaline Phosphatase Serum	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0136	Alkaline Phosphatase Isoenzymes Electrophoresis - Serum	Chennai - Avadi	D	5900	Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0137	Alkaline Phosphatase with Bone Fraction Ostase Serum	Chennai - Avadi	C	1810	Mon & Thu : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0138	Allergen,Individual-Animal Cat Dander Serum	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0140	Allergen,Individual-Animal Cow DanderSerum	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0141	Allergen,Individual-Animal Dog Dander Serum	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0143	Allergen,Individual-Animal Guinea Pig Epithelium Serum	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0145	Allergen,Individual-Animal Horse Dander Serum	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0152	Allergen,Individual-Drug ACTH	Chennai	C	1580	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi					ml serum for profiles and reflex tests.
A0158	Allergen,Individual-Drug Ciprofloxacinantibiotic	Chennai - Avadi	C	1580	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0173	Allergen, Individual-Drug Sulpha antibiotic	Chennai - Avadi	C	1580	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0175	Allergen,Individual-Dust Dermatophagoides FarinaeHouse Dust Mite	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0176	Allergen,Individual-Dust Dermatophagoides Pteronyssinus House Dust Mite	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0178	Allergen,Individual-Dust House Dust - Hollistier	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0179	Allergen,Individual-DustHouse Dust Greer Labs	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0180	Allergen,Individual-Food Almond Badam (Amygdalus Communis)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0181	Allergen,Individual-Food Alpha - Lactalbumin Milk protein	Chennai - Avadi	C	1580	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0182	Allergen,Individual-Food Apple Safarchand (Malus X Domestica)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0185	Allergen,Individual-Food Baingan	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0186	Allergen,Individual-Food Bakers Yeast (Saccharomyces Cerevisiae)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0187	Allergen,Individual-Food Banana Kela (Musa Spp.)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0189	Allergen,Individual-Food Basil Ocimum basilicum	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0191	Allergen,Individual-Food Beef Cow meat (Bos spp.)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0192	Allergen,Individual-Food Beta - Lactoglobulin- Milk protein	Chennai - Avadi	C	1580	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0193	Allergen,Individual-Food Black Pepper Piper nigrum	Chennai - Avadi	C	1355	Mon & Thu : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0194	Allergen,Individual-Food Blue MusselMytilus Edulis	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0195	Allergen,Individual-Food Brazil Nut Bertholletia excelsa	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0198	Allergen,Individual-Food Cabbage Gobi (Brassica Oleraceavar.Capitata)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0199	Allergen,Individual-Food Cardamon Elettaria cardamomum	Chennai - Avadi	C	1355	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0201	Allergen,Individual-Food Casein(Milk Protein)Dahi	Chennai - Avadi	C	1580	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0202	Allergen,Individual-Food Cashew Nut Kaju (Anacardium Occidentale)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0204	Allergen,Individual-Food Celery Apium Graveolens	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0205	Allergen,Individual-Food Cheese Cheddar	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0207	Allergen,Individual-Food Chick Pea Kabuli Chana (Cicer Arietinus)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tests.
A0208	Allergen,Individual-Food Chicken Meat Murgi	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0209	Allergen,Individual-Food Chilipepper Capsicum frutescens	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0211	Allergen,Individual-Food Clove(Syzygium aromaticum)	Chennai - Avadi	C	1355	Mon & Thu : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0212	Allergen,Individual-Food Cocoa Chocolate (Theobroma Cacao)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0213	Allergen,Individual-Food Coconut Narial (Cocos Nucifera)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0214	Allergen,Individual-Food Cod Fish Machhli (Gadus Morhua)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0215	Allergen,Individual-Food Coffee	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0216	Allergen,Individual-Food Coriander	Chennai - Avadi	C	1355	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0217	Allergen,Individual-Food Corn / Maize Makka (Zea Mays)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0218	Allergen,Individual-Food Cow Milk Doodh	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0219	Allergen,Individual-FoodCrabCancer pagurus	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0224	Allergen,Individual-Food Egg White Anda	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0225	Allergen,Individual-Food Egg Yolk Anda	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0228	Allergen,Individual-Food Garlic Lehsoun (Allium Sativum)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0230	Allergen,Individual-Food Gluten	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0234	Allergen,Individual-Food Green Pea	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0236	Allergen,Individual-Food Hazelnut Corylus avellana	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0238	Allergen,Individual-Food Kiwi Fruit Actinidia deliciosa	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0242	Allergen,Individual-Food Lemon Nimbu (Citrus Limon)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0243	Allergen,Individual-Food Lentil Masoor Dal (Lens Esculenta)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0245	Allergen,Individual-Food Lobster Homarus gammarus	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0249	Allergen,Individual-Food Mango Mangifera Indica	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0252	Allergen,Individual-Food Mushroom champignon / Agaricus hortensis	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0253	Allergen,Individual-Food Mustard Sarson (Brassica / Sinapsis Spp.)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0254	Allergen,Individual-Food Mutton Ovis Spp.	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0256	Allergen,Individual-Food Olive Olea Europaea	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0257	Allergen,Individual-Food Onion Piaz (Allium Cepa)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0258	Allergen,Individual-Food Orangesantra (Citrus sinensis)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0260	Allergen,Individual-Food Papaya Carica Papaya	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0261	Allergen,Individual-Food Paprica Simla Mirchi (Capsicum Annuum)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0263	Allergen,Individual-Food Peach Prunus persica	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0264	Allergen,Individual-Food Peanut Moongfali (Arachis Hypogaea)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0265	Allergen,Individual-Food Pear Pyrus communis	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0266	Allergen,Individual-Food Pecan nut Carya illinoensis	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0267	Allergen,Individual-Food Pineapple Ananas	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0268	Allergen,Individual-Food Pistachio	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0269	Allergen,Individual-Food Plaice Fish (Pleuronectes platessa)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0272	Allergen,Individual-Food Pork Pig meat	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0273	Allergen,Individual-Food Potato Aaloo (Solanum Tuberosum)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0277	Allergen,Individual-Food Rice Chawal (Oryza Sativa)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0278	Allergen,Individual-Food Salmon Rawas fish (Salmo Salar)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0281	Allergen,Individual-Food Shrimp / Prawns	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0282	Allergen,Individual-Food Soyabean Glycine Max	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0283	Allergen,Individual-Food Spinach Palak (Spinachia Oleracea)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0284	Allergen,Individual-Food Strawberry Fragaria vesca	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0289	Allergen,Individual-Food Tomato Tamatar (Lycopersicon Lycopersicum)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0290	Allergen,Individual-Food TunaFish (Pleuronectes platessa)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0292	Allergen,Individual-Food Walnut Juglans californica	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0294	Allergen,Individual-Food Wheat Gehon (Triticum Aestivum)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0295	Allergen,Individual-Food White Bean Phaseolus Vulgaris	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tests.
A0298	Allergen,Individual-Grass pollen Bermuda Grass Durva (Cynodon Grass)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0300	Allergen,Individual-Grass pollen Cultivated Rye Grass Secale Cereale	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0302	Allergen,Individual-Grass pollen Johnson Grass Jowar grass (Sorghum Grass)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0304	Allergen,Individual-Grass pollen Meadow Grass Poa Pratensis	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0308	Allergen,Individual-Grass pollen Timothy Grass Ghass (Phleum Pratense)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0311	Allergen,Individual-Insect Cockroach German Blatella Germanica	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0312	Allergen,Individual-InsectHoney Bee Venom Madhu Makhhi (Apis Mellifera)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0315	Allergen,Individual-Micro organism Alternaria Tenius fungus	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0316	Allergen,Individual-Micro organism Aspergillus Fumigatus fungus	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0317	Allergen,Individual-Micro organism Candida Albicans fungus	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0318	Allergen,Individual-Micro organismCladosporium Herbarum fungus	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0320	Allergen, Individual - Insect Mosquito Machhar (Aedes Communis)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0322	Allergen,Individual-Micro organism Penicillium Notatum fungus (Penicillium chrysogenum)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0338	Allergen,Individual-Tree pollen Acasia Babool (Acasia Longifolia)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0339	Allergen,Individual-Tree pollen Alder Alnus Incana	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0340	Allergen,Individual-Tree pollen Birch Betula Verrucosa	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0342	Allergen,Individual-Tree pollen Eucalyptus Nilgiri (Safeda)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0345	Allergen,Individual-Tree pollen Mesquite Pahadi Keekar (Prosopis Juliflora)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0346	Allergen,Individual-Tree pollen Mulberry Shahtoot (Morus alba)	Chennai - Avadi	C	1355	Daily : 11:30am	5th day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0347	Allergen,Individual-Tree pollen Oak Quercus Alba	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0348	Allergen,Individual-Tree pollen White Pine Pinus strobus	Chennai - Avadi	C	1355	Daily : 11:30am	5th day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0349	Allergen,Individual-Tree pollen Willow Salix Caprea	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0350	Allergen,Individual-Weed pollen Cocklebur Xanthium commune	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0351	Allergen,Individual-Weed pollen Common Pigweed Kaantewali Chauli (Amaranthus Retroflexus)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0352	Allergen,Individual-Weed pollen Common RagweedClose To Parthenium / Congress	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Grass (Ambrosia Elatior)						tests.
A0353	Allergen,Individual-Weed pollen English Plantain Plantago Lanceolata	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0355	Allergen,Individual-Weed pollen Goose Foot Chenopodium Album	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0356	Allergen,Individual-Weed pollen Mugwort Sita Bani (Artemisia Vulgaris / BanoBarna)	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0360	Allergy - Fruit Panel	Chennai - Avadi	D	5430	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0361	Allergy Fx 5 Panel	Chennai - Avadi	D	3775	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0362	Allergy - Milk & Milk Product Panel	Chennai - Avadi	D	4720	Daily : 11:30am	Next day : 6:00pm	5 ml of serum
A0363	Allergy - Non-vegetarian Panel	Chennai - Avadi	D	5900	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0364	Allergy-Nuts Panel	Chennai - Avadi	D	5665	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0367	Allergy - Spices Panel	Chennai - Avadi	D	9675	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0368	Allergy - Vegetable-1 Panel	Chennai - Avadi	D	3540	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0370	Allergy - Vegetarian + Nonvegetarian Panel	Chennai - Avadi	D	11210	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0371	Allergy panel - food Vegetarians Panel	Chennai - Avadi	D	6610	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0372	Allergy - Adult (Comprehensive) Panel	Chennai - Avadi	D	10740	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0373	Allergy - Animal Panel	Chennai - Avadi	D	4600	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0377	Allergy - Drug Panel: (Maxi)	Chennai - Avadi	D	10740	Tue & Fri : 11:30am	After 8hrs	5 ml of serum
A0379	Allergy - Dust Panel	Chennai - Avadi	D	3305	Daily : 11:30am	Next day : 5:00pm	5 ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0382	Allergy - Environmental Panel	Chennai - Avadi	D	14750	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0383	Allergy - GIT Panel	Chennai - Avadi	D	14160	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0384	Allergy-Grass Panel	Chennai - Avadi	D	4130	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0389	Allergy - Mold-2 Panel	Chennai - Avadi	D	3540	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0395	Allergy - Trees pollen Panel	Chennai - Avadi	D	3660	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0396	Allergy - Weeds pollen Panel	Chennai - Avadi	D	3660	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0398	IgETotal antibodySerum	Chennai - Avadi	B	940	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0399	Allergy Screening IgE With Phadiatop (Adult & >5y)	Chennai - Avadi	D	1770	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0399_C	Allergy Screening IgE With PhadiaTop (Child below 5 y)	Chennai - Avadi	D	1600	Daily : 11:30am	After 8hrs	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A0400	Allergy screeningPhadiatop-AdultFor > 5 years	Chennai - Avadi	C	1245	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0401	Phadiatop allergen screening - Infant (For Below 5 yrs)	Chennai - Avadi	C	1245	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0403	Alpha 1 Acid Glycoprotein Serum	Chennai - Avadi	D	2480	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0404	Alpha 1 Antitrypsin Serum	Chennai - Avadi	C	2260	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0406	Alpha 2 Macroglobulin Serum	Chennai - Avadi	D	3305	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0407	Alpha Fucosidase Blood	Chennai - Avadi	R	8025	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0408	Alpha Galactosidase Fabry Blood	Chennai - Avadi	R	6420	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0409	Alpha Glucosidase GSD II / pompe Blood	Chennai - Avadi	R	8240	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0410	Alpha Iduronidase Hurler / MPS I Blood	Chennai - Avadi	R	6955	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0411	Alpha Mannosidase Blood	Chennai - Avadi	R	6420	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0412	Alpha N-Acetyl Galactosaminidase Schindler disease Blood	Chennai - Avadi	R	6420	Daily : 11:30am	After 4 weeks : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0413	Alpha Thalassemia by PCR	Chennai - Avadi	D	4130	Tue & Sat: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Clinical history, referring clinicians contact number, HPLC Report, CBC findings, clinical/family history and TRF mandatory.
A0414	Aluminium in serum	Chennai - Avadi	C	2825	Daily 7:30 AM	Next Day 5:00 PM	HeavyMetalSamplingguidelines.pdf
A0415	Aluminium, Urine 24 Hrs	Chennai - Avadi	C	2825	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0416	Aluminium spot urine	Chennai - Avadi	C	2825	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0417	AMH Mullerian Inhibiting Substance Serum	Chennai - Avadi	B	2050	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0422	Amino Acid Quantitative Analysis by LC MSMS Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
A0423	Amino Acid Quantitative Urine	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
A0424	Amino Levulinic Acid (ALA) 24 H Urine	Chennai - Avadi	D	4840	Thu : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
A0425	Amino Levulinic Acid (ALA) Urine Spot	Chennai - Avadi	D	4840	Thu : 11:30am	Next day : 5:00pm	10 ml of Spot urine (Protect From Light)
A0427	Ammonia Biochemical Plasma	Chennai - Avadi	B	755	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Heparinised Plasma
A0428	Amoebiasis Antibodies-IgG Antibody by IHA - Serum	Chennai - Avadi	C	1470	Mon & Thu : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0429	Amphetamine Quantitative Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of Spot Urine
A0430	Amphiphysin Antibody Neuronal Ag in Paraneoplastic Syndrome Serum	Chennai - Avadi	D	5310	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0431	Amylase Enzymatic Ascitic Fluid	Chennai - Avadi	B	485	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Ascitic Fluid
A0432	Amylase Enzymatic Serum	Chennai - Avadi	B	460	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0433_24	Amylase enzymatic Urine Spot	Chennai - Avadi	B	485	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
A0433_24H	Amylase enzymatic Urine 24H	Chennai - Avadi	B	430	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
A0434	Amyloid A Serum	Chennai - Avadi	D	4305	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0435	ANA Blot 18 Antigen, Serum	Chennai - Avadi	C	3470	Daily:7.30 am	Same day 8:00 PM	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0436	ANA by ELISA Anti Nuclear Antibody Screening Serum	Chennai - Avadi	B	755	Mon, Wed & Fri : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0437	Anti Nuclear Antibody by IFA	Chennai - Avadi	C	1020	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0438	ANA Profile Serum	Chennai - Avadi	D	5500	Tuesday,Thursday and Saturday 7.30	Same day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
					am		tests.
A0440	Anaemia Profile-Haemolytic Type	Chennai - Avadi	C	3505	Daily : 9:00am to 9:00pm	Next day : 5:00pm	4 ml Heparin & EDTA Blood, 4ml Serum
A0441	Anaemia Profile-Maxi	Chennai - Avadi	C	8475	Daily : 9:00am to 9:00pm	Next day : 5:00pm	4 ml Heparin & EDTA Blood, 4ml Serum
A0442	Anaemia Profile-Mini	Chennai - Avadi	C	2035	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of EDTA whole blood, 5 ml of serum
A0443	Anaemia Profile-Nutritional	Chennai - Avadi	C	3955	Daily : 9:00am to 9:00pm	After 6hrs	3ml EDTA blood & Serum
A0444	Anaemia Profile-Pernicious	Chennai - Avadi	C	3335	Mon & Thu : 7:30am	Next day : 5:00pm	3ml EDTA blood & Serum
A0445	Anti neutrophil cytoplasmic antibody (ANCA) Serum	Chennai - Avadi	C	2090	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0446	Androgen Profile	Chennai - Avadi	C	5085	Tue & Fri: 7:30am	Next day : 5:00pm	5 ml of serum
A0447	AndrostenedioneA4Serum	Chennai - Avadi	C	2035	Daily : 9:00am to 9:00pm	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0448	ACE - Angiotensin Converting Enzyme Serum	Chennai - Avadi	C	1185	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0449	ACE-Angiotensin Converting Enzyme, CSF	Chennai - Avadi	D	3305	Daily : 11:30am	4th day : 6:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
A0450	Antenatal Profile-II	Chennai - Avadi	C	1580	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma , 20 ml of spot urine
A0451	Antenatal Profile - I	Chennai - Avadi	C	3275	Daily : 9:00am	3rd day : 6:00pm	2 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma , 20 ml of spot urine
A0453	Anti - A Titre	Chennai - Avadi	D	1770	Daily : 11:30am	After 8hrs	3 ml Serum
A0454	Anti - B Titre	Chennai - Avadi	D	1770	Daily : 11:30am	After 8hrs	3 ml Serum
A0455	Anti DNase B enzymatic Seum	Chennai - Avadi	D	2830	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tests.
A0458	Anti Thrombin III Activity Functional, Citrated plasma	Chennai - Avadi	D	4250	Mon & Thu : 9:00am	Next day : 1:00pm	1 ml of citrated platelet poor plasma in frozen condition.
A0459	Anti Thrombin III Antigen, Citrated plasma	Chennai - Avadi	D	5665	Mon & Thu : 9:00am	Next day : 1:00pm	1 ml of citrated platelet poor plasma in frozen condition.
A0462	Antimony ICPMS Urine Spot	Chennai - Avadi	R	4280	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
A0463_N	Phospholipid Profile - 3	Chennai - Avadi	D	9085	Thu : 9:00am	Same day : 7:00pm	3 ml of Citrate plasma, 3 ml of serum
A0464	APA-IgG Phospholipid Antibody Serum	Chennai - Avadi	C	1015	Daily: 7:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0465	APA-IgM Phospholipid Antibody Serum	Chennai - Avadi	C	1015	Tue, Thu & Sat: 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0466	APCR Activated Protein C Resistance Test, Citrated plasma	Chennai - Avadi	D	4840	Mon & Thu : 9:00am	Next day : 1:00pm	1 ml of citrated platelet poor plasma in frozen condition.
A0467	Apolipoproteins A1 serum	Chennai - Avadi	B	485	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0468	Apolipoproteins B/A1	Chennai - Avadi	C	1245	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0469	Apolipoproteins Bserum	Chennai - Avadi	B	485	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0470	Apolipoproteins Profile, Serum	Chennai - Avadi	C	1920	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0471	APTT Activated Partial Thromboplastin Time, Citrated plasma.	Chennai - Avadi	B	400	Daily : 9:00am to 9:00pm	After 4hrs	1 ml of citrated platelet poor plasma in frozen condition.
A0472	Arginine Quantitative Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
A0473	Arginine Quantitative Urine	Chennai - Avadi	R	6365	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
A0474	Arsenic	Chennai	D	2950	Daily : 11:30am	3rd day :	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0475	Arsenic, Urine 24 Hrs	- Avadi Chennai	D	2950	Daily : 11:30am	6:00pm 3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0476	Arsenic, Urine Spot	- Avadi Chennai	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0477	Arthritis Profile	- Avadi Chennai	C	2710	Daily : 9:00am	Next day : 5:00pm	2 ml of EDTA whole blood, 3 ml of serum, 20 ml of spot urine
A0478	Aryl Sulphatase A metachromatic leucodystrophy Blood	- Avadi Chennai	R	6420	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0479	Aryl Sulphatase B MPS VI Blood	- Avadi Chennai	R	9525	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0480	ASCA-IgA Saccharomyces Cerevisiae Serum	- Avadi Chennai	C	2375	Tue & Fri: 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0481	ASCA-IgG Saccharomyces Cerevisiae Serum	- Avadi Chennai	C	2375	Tue & Fri: 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0482	ASH /NASH / STEATO Profile	- Avadi Chennai	E	15750	Daily : 9:00am	5th day : 6:00pm	5 ml of serum and 2 ml of fluoride plasma . Fasting sample
A0483	ASKA Skeletal (Striated) Muscle Antibody Serum	- Avadi Chennai	D	4010	Mon, Wed & Fri : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0484	ASMA Smooth Muscle Antibody Serum	- Avadi Chennai	C	2145	Daily: 11:30am	Next day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0485	ASO titre Quantitative Serum	- Avadi Chennai	B	540	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0486	Asparagine Quantitative Plasma	- Avadi Chennai	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
A0487	Asparagine Quantitative Urine	- Avadi Chennai	R	6365	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0488	Aspartic Acid Quantitative Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
A0489	Aspartic Acid Quantitative Urine	Chennai - Avadi	R	6365	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
A0490	Aspergillosis IgG antibodies by EIA - Serum	Chennai - Avadi	C	2235	Mon & Thu : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0491	Aspergillosis IgM antibodies by EIA - Serum	Chennai - Avadi	C	2375	Mon & Thu : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0492	Galactomannan, Aspergillus Antigen test by EIA -Serum	Chennai - Avadi	C	5765	Wed & Sat : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0493	Autoimmune Liver Disease Profile -1	Chennai - Avadi	D	6000	Mon, Wed & Fri : 7:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0495	Autoimmune Liver Disease Profile-3	Chennai - Avadi	D	4250	Mon, Wed & Fri : 7:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0496	Autoimmune Liver Disease Profile-4	Chennai - Avadi	D	4840	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0497	Autoimmune Encephalitis Panel	Chennai - Avadi	E	22050	Mon,Wed,Fri 7:30am	Next day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0499	Acid lipase disorder**	Chennai - Avadi	R	7490	Daily : 11:30am	4 weeks	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0500	ACTN3 GENOTYPING_A0500	Chennai - Avadi	R	4280	Thu : 11:30am	10th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
A0509	AFB-Xpert panel by Combo (MTB /RIF	Chennai	D	3480	Daily : 7:30am	After 6	Tissue, Lymph Node, Endometrium,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Detection & AFB Culture) - Extrapulmonary Samples	- Avadi				weeks : 6:00pm	Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
A0510	AFB-Xpert panel by Combo (MTB /RIF Detection & AFB Culture) - Pulmonary Samples	Chennai - Avadi	D	3480	Daily : 7:30am	After 6 weeks : 6:00pm	2ml of Pulmonary specimen in Sterile leak proof container
A0512	ANCA by IFA reflex MPO/PR3 Antibodies Serum	Chennai - Avadi	C	2825	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0513	Apolipoproteins E**	Chennai - Avadi	R	4065	Wed : 11:30am	4th day : 6:00pm	3 ml of serum
A0514	TruHeart Platinum	Chennai - Avadi	W	6000	Mon & Thu : 9:00am	Same day : 5:00pm	5 ml of serum
A0519	Allergen,Individual-Tree pollen Papdi Chibil / Elm_A0519	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0520	Allergen,Individual-Weed pollen SALTWORT (Salsola kali / pestifer)_A0520	Chennai - Avadi	C	1355	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0522	ALK (D5F3) reflex FISH	Chennai - Avadi	D	7080	Daily : 9:00am	7th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
A0525	Allergen,Individual-Animal Pigeon Feathers, Serum_A0525	Chennai - Avadi	C	1355	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0533	ADA - Adenosine Deaminase by Biochemical Method - Synovial fluid	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of Fluid in Sterile leak proof container
A0536	APA-IgA Phospholipid Antibody Serum	Chennai - Avadi	C	1130	Tue, Thu & Sat: 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0537	Atherosclerosis Profile - 1	Chennai - Avadi	D	6370	Mon & Thu : 9:00am	Same day : 5:00pm	5 ml of serum
A0538	Atherosclerosis Profile - 3	Chennai	D	8760	Mon & Thu :	Same day :	5 ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0539	Aldosterone/PRA Ratio	- Avadi Chennai	D	7670	9:00am Tue & Fri: 11:30am	5:00pm Next day : 5:00pm	2 ml of EDTA Plasma and 2 ml Serum
A0560	Galactomannan, Aspergillus Antigen test by EIA - BAL Fluid	- Avadi Chennai	C	5875	Wed & Sat : 11:30am	Next day : 5:00pm	3 ml of Fluid in Sterile leak proof container
A0563	AFB (TB) - Antibioqram for MTb complex - Cycloserine Drug by Conventional LJ - Pure Growth of MTB complex	- Avadi Chennai	C	2090	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0564	Albumin body fluid	- Avadi Chennai	A	820	Daily : 9:00am to 9:00pm	After 6hrs	2 ml body fluid
A0567	ALK/NPM1: t(2;5) translocation Anaplastic Large cell lymphoma FISH, Bone marrow	- Avadi Chennai	D	4720	Daily : 7:30am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
A0568	ALK/NPM1: t(2;5) translocation Anaplastic Large cell lymphoma FISH, BLOOD	- Avadi Chennai	D	5430	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
A0569	ALK/NPM1: t(2;5) translocation Anaplastic Large cell lymphoma FISH, TISSUE	- Avadi Chennai	D	9440	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
A0570	ECP (Eosinophil cationic protein) Serum	- Avadi Chennai	D	1535	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0580	Allergy-GIT Pane I(Veg + Non Veg)	- Avadi Chennai	D	19470	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0584	Allergy screening-Tryptase	- Avadi Chennai	C	2260	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0585	Anaemia Profile - Macrocytic	- Avadi Chennai	C	2035	Daily : 9:00am to 9:00pm	After 8hrs	3ml EDTA blood & Serum
A0586	Anaemia Profile - Microcytic Hypochromic	- Avadi Chennai	C	2035	Daily : 9:00am to 9:00pm	After 8hrs	3ml EDTA blood & Serum
A0587	Anaemia Workup Profile	- Avadi Chennai	C	2260	Daily : 9:00am to 9:00pm	Next day : 5:00pm	3ml EDTA blood & Serum
A0588	Aneuploidy-Sperm 13,18,21,X and Y by	Chennai	D	5665	Daily : 9:00am	7th day :	CytogeneticsandFISH.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0595	FISH Adrenoleukodystrophy (X-linked ALD) sequencing of ABCD1 gene**	- Avadi Chennai - Avadi	R	20330	Daily : 11:30am	6:00pm After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
A0596	Autoimmune encephalitis panel-CSF	Chennai - Avadi	E	19950	Mon,Wed,Fri 7:30am	Next day : 7:00pm	2 ml of CSF in sterile screw capped container.
A0599	APC gene analysis**	Chennai - Avadi	R	17850	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
A0608	AFB (TB) - Antibiogram for MTb complex - Linezolid Drug by MGIT - Pure Growth of MTB complex	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0613	Alpha thalassemia gene analysis (HBA1 & HBA2)	Chennai - Avadi	E	16800	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Clinical history, referring clinicians contact number, HPLC Report, CBC findings, clinical/family history and TRF mandatory.
A0616	Aldosterone/Direct Renin Ratio	Chennai - Avadi	D	8200	Daily : 7:30am	Same day : 5:00pm	2 ml of EDTA Plasma and 2 ml Serum
A0617_CSF	Autoimmune encephalitis panel-2	Chennai - Avadi	E	30450	Tue : 9:00am	Next day : 7:00pm	2 ml of CSF in sterile screw capped container.
A0617_serum	Autoimmune encephalitis panel-2	Chennai	E	34180	Tue : 9:00am	Next day :	3 ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0618	APO E GENOTYPING**	- Avadi Chennai - Avadi	R	26750	Daily : 11:30am	7:00pm 12th Day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
A0623	Alzheimer disease screening profile CSF	Chennai - Avadi	D	14750	Tue : 11:30am	Next day : 5:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
A0624	Arsenic plasma	Chennai - Avadi	D	3245	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0625	Serum Ascites Albumin Gradient (SAAG)	Chennai - Avadi	D	590	Daily : 9:00am to 9:00pm	Next day : 5:00pm	Serum 1 ml, Fluid 1 ml
A0626	Allergy Screening	Chennai - Avadi	D	3305	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0629	Aspergillus broncho pulmonary aspergillosis (ABPA)	Chennai - Avadi	D	3420	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0630	AFB Drug susceptibility - Flexi panel (4 drugs) by MGIT - Pure Culture	Chennai - Avadi	D	7080	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0632	Asthma / Rhinitis Comprehensive	Chennai - Avadi	D	12980	Daily : 11:30am	Next day : 5:00pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A0633	Allergy panel - Asthma / Rhinitis comprehensive with monitoring 4	Chennai - Avadi	D	20120	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0635	Allergy - Adult panel - 2	Chennai - Avadi	D	5780	Daily : 11:30am	Next day : 5:00pm	for 7 days and at Ambient temperature for a day. Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A0636	Allergy - Adult panel - 3	Chennai - Avadi	D	9085	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A0638	Aluminium, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
A0639	Aluminium, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
A0640	Amphetamines confirmation by LCMS, Spot Urine**	Chennai - Avadi	R	9630	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
A0641	Antimony, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
A0642	Arsenic, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							storage and shipping. Plz enquire before sending
A0643	Arsenic, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
A0644	Acute intermittent porphyria (HMBS)-Gene Analysis	Chennai - Avadi	E	21000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
A0645	Alkaptonuria (HGD)-Gene Analysis	Chennai - Avadi	E	19950	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
A0649	Mycobacterial culture identification by MALDI TOF- Pure Culture	Chennai - Avadi	C	2485	Daily : 9:00am	15th day : 6:00pm	Pure growth of Mycobacterium spp. on LJ slant/ MGIT Tube/ Bactec alert bottle.
A0650	SARS-CoV-2 Antibody, qualitative by CLIA - Serum	Chennai - Avadi	C	790	Daily : 9:00am to 9:00pm	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0653	ALPS Screen Panel by Flow cytometry**	Chennai - Avadi	R	9735	Daily : 11:30am	5th day : 6:00pm	sample : 3 ml of EDTA Blood
A0654	SARS-CoV-2 (Spike) Antibody, Quantitative (Coviprotect) by CLIA - Serum	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0656	Aspergillus DNA Detection by Real Time PCR - Body fluids/CSF/EDTA Blood/Pure culture/Sputum/Tissue	Chennai - Avadi	D	7080	Daily : 7:30am	3rd day : 6:00pm	3ml of Body fluids/ CSF/ EDTA Blood/ Pure culture/ Sputum/ Tissue in Sterile leak proof container
A0657	Array Feto	Chennai	E	17640	Daily : 9:00am	15th day :	1)Upto 12 weeks of gestation:

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				6:00pm	50-100mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container 2) In 12-18 weeks of gestation: a) 50-100mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container b) and or fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container 3)IUFD >20 weeks of gestation: fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container. if possible in 20ml normal saline with 2-3 drops of gentamycin and 1-4ml Maternal blood in EDTA tube is required for maternal cell contamination (MCC) studies.
A0658	Amniochrome	Chennai - Avadi	E	19950	Daily : 9:00am	15th day : 6:00pm	Amniotic fluid (20-30ml divided in 2 sterile falcon centrifuge tube) and 1-4ml Maternal blood in EDTA tube is required for maternal cell contamination (MCC) studies. (Completely filled Consent form G with patient and doctor stamp and signature is mandatory as per PCPNDT Act)
A0659	Array Chrome	Chennai - Avadi	E	16800	Daily : 9:00am	15th day : 6:00pm	3ml Peripheral Blood in EDTA Vacutainers
A0660	Anti complement Factor H	Chennai - Avadi	R	5350	Daily : 11:30am	6th day : 6:00pm	3 ml Serum
A0661	ADAMTS 13 Activity	Chennai - Avadi	R	6420	Daily : 11:30am	6th day : 6:00pm	3 ml Citrate Plasma (Frozen mandatory)
A0662	ADAMTS 13 Autoantibody	Chennai - Avadi	E	8400	Daily : 11:30am	5th day : 6:00pm	3 ml Citrate Plasma (Frozen mandatory)
A0680	ALLERGY COMPONENT TESTING (CRD)	Chennai - Avadi	E	21000	Daily : 9:00am	Next day : 6:00pm	Require serum (3ml) sample collected in fasting status

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A1009	Adrenal Gland Antibodies, Serum**	Chennai - Avadi	R	6420	Daily : 11:30am	After 4 weeks : 6:00pm	preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and -20°C for long time. 5 ml Serum
A1901	Albumin, serum	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
A1902	Canine - Alkaline Phosphatase	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
A1903	Canine - Amylase level	Chennai - Avadi	A	355	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
A5583	AMH Premium	Chennai - Avadi	D	3305	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A5592	AMH Supreme	Chennai - Avadi	D	5075	Daily : 9:00am to 9:00pm	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A5702	Comprehensive Allergy Panel	Chennai	D	16000	Daily : 9:00am	Same day :	All tests require serum (3ml)

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				8:00pm	sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A5703	Comprehensive Eczema panel	Chennai - Avadi	D	12980	Daily : 9:00am	Same day : 8:00pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A5704	Common Food panel	Chennai - Avadi	D	7080	Daily : 9:00am	Same day : 8:00pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A8121	Amylase Enzymatic,Pleural Fluid	Chennai - Avadi	B	485	Daily : 9:00am to 9:00pm	After 4hrs	2 ml body fluid
A8285	AFB SUSCEPTIBILITY BEDAQILINE	Chennai - Avadi	C	2710	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8455	AFB Drug susceptibility - 14 Drugs Panel by MGIT - Pure Culture	Chennai - Avadi	D	16520	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8456	AFB Drug susceptibility - First Line 4 Drugs Panel (HREZ) by MGIT - Pure Culture	Chennai - Avadi	D	6845	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8457	AFB Drug susceptibility - 11 Drugs Panel by MGIT - Pure Culture	Chennai - Avadi	D	12330	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8458	AFB Drug susceptibility - Oral Bedaquiline Panel by MGIT - Pure Culture	Chennai - Avadi	D	9440	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8459	AFB Drug susceptibility - Isoniazid Monoresistance Panel by MGIT - Pure	Chennai - Avadi	D	6845	Daily : 9:00am	After 4 weeks :	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Culture					6:00pm	
A8598	Allergen,Individual-Drug CEPHALOSPORIN	Chennai - Avadi	D	1505	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8599	Allergen,Individual-Drug LIDOCAINE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8600	Allergen,Individual-Drug BENZOCAINE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8601	Allergen,Individual-Drug MEPIVACAINE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8602	Allergen,Individual-Drug PRILOCAINE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8603	Allergen,Individual-Drug PARACETAMOL	Chennai - Avadi	D	1505	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8604	Allergen,Individual-Drug PYRAZOLONE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8605	Allergen,Individual-Drug PIROXYCAM	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8606	Allergen,Individual-Drug DICLOFENAC	Chennai - Avadi	D	1505	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8607	Allergen,Individual-Drug IBUPROFEN	Chennai - Avadi	D	1505	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8608	Allergen,Individual-Drug PENICILLOYL G	Chennai - Avadi	D	1505	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8609	Allergen,Individual-Drug PENICILLOYL V	Chennai - Avadi	D	1505	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A8610	Allergen,Individual-Drug DOXYCYCLINE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8611	Allergen,Individual-Drug CLOXACILLIN	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8612	Allergen,Individual-Drug NEOMYCIN	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8613	Allergen,Individual-Drug OXACILLIN	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8614	Allergen,Individual-Drug CHLORAMPHENICOL	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8615	Allergen,Individual-Drug METRONIDAZOLE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8616	Allergen,Individual-Drug AMPICILLIN	Chennai - Avadi	D	1505	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8617	Allergen,Individual-Drug AMOXICILLIN	Chennai - Avadi	D	1505	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8618	Allergen,Individual-Drug TETRACYCLINE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8619	Allergen,Individual-Drug GENTAMYCIN	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8620	Allergen,Individual-Drug ERYTHROMYCIN	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8621	Allergen,Individual-Drug STREPTOMYCIN	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8622	Allergen,Individual-DrugTRIMETHOPRIM	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tests.
A8623	Allergen,Individual-Drug SULFAMETHOXAZOLE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8624	Allergen,Individual-Drug INSULIN HUMAN	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8625	Allergen,Individual-Drug CORTISONE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8626	Allergen,Individual-Drug PREDNISOLONE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8627	Allergen,Individual-Drug THEOPHYLLINE	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8628	Allergen,Individual-Drug AMBROXOL	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8629	Allergen,Individual-Drug PYRIDOXAMINE (VITAMIN B6)	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8630	Allergen,Individual-Drug COBALAMIN (VITAMIN B12)	Chennai - Avadi	D	1500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0007	Bacterial Meningitis Panel (5 antigen) by Latex agglutination - CSF	Chennai - Avadi	D	4955	Daily : 9:00am & 3:00pm	After 6hrs	2ml of CSF in Sterile leak proof tube or container
B0008	Bacterial Meningitis Panel 5 antigen by Latex agglutination - Urine	Chennai - Avadi	D	5665	Daily : 9:00am & 3:00pm	After 6hrs	15ml of Spot Urine in Sterile leak proof container
B0010	Barbiturates Quantitative Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	5 ml of Spot Urine
B0013	Bence Jones Proteins, Urine-Qualitative	Chennai - Avadi	B	1060	Daily : 7:30am	Same day : 5:00pm	10 ml of 24 hrs urine (Preferred) / Spot urine
B0015	Benzodiazapine Diazepam Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	5 ml of Spot Urine
B0017	Beta Galactocerebrosidase Krabbe Blood	Chennai - Avadi	R	8775	Daily : 11:30am	After 4 weeks :	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
						6:00pm	of collection, Avoid collection on weekend & Holidays)
B0018	Beta Galactosidase GM 1 / MPS IV B Blood	Chennai - Avadi	R	5990	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0019	Beta Glucosidase Gaucher Blood	Chennai - Avadi	R	5990	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0020	Beta Hexosaminidase A GM2-1/ taysachs Blood	Chennai - Avadi	R	8025	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0021	Beta Hexosaminidase T GM2-2/ Sandhoff Blood	Chennai - Avadi	R	8025	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0022	Beta Mannosidase Blood	Chennai - Avadi	R	6420	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0023	Beta Thallesemia Prenatal Diagnosis	Chennai - Avadi	D	14160	Daily : 9:00am	7th day : 6:00pm	Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition or Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
B0025	Beta-Thalassemia, Beta Globin (HBB) Gene Sequencing	Chennai - Avadi	D	12980	Wed & Fri: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Clinical history, referring clinicians contact number, HPLC Report, CBC findings, clinical/family history and TRF mandatory.
B0026	Beta-2-Crosslaps CTx-1 Plasma	Chennai - Avadi	C	1920	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
B0027	Beta-2-Glycoprotein 1- IgA serum	Chennai - Avadi	C	1470	Tue & Fri: 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0028	Beta-2-Glycoprotein 1 - IgG Antibody Serum	Chennai - Avadi	C	1185	Tue, Thu & Sat: 9:00am	Same day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0029	Beta-2-Glycoprotein 1- IgM Antibody Serum	Chennai - Avadi	C	1185	Tue, Thu & Sat: 9:00am	Same day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0030	Beta-2-Microglobulin CSF	Chennai - Avadi	C	1810	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of CSF in leakproof container (foilin vu conical tube is recommended to prevent leakage).
B0031	Beta-2-Microglobulin Serum	Chennai - Avadi	C	1640	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0032	Beta-2-Microglobulin Urine Spot	Chennai - Avadi	C	1810	Daily : 9:00am to 9:00pm	After 6hrs	10 ml of Spot urine
B0033	Bicarbonate biochemical Serum	Chennai - Avadi	B	400	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B0034	Bile Acids-Total biochemical Serum	Chennai - Avadi	C	1920	Daily : 9:00am to 9:00pm	After 6hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0035	Bile Salt & Pigments - Urine	Chennai - Avadi	A	110	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
B0036	Bilirubin Direct Serum	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0037	Bilirubin Total Serum	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0038	Bilirubin Total, Direct, Indirect Serum	Chennai - Avadi	A	190	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0039	Biotinidase Neonatal Dried blood spot	Chennai - Avadi	C	780	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B0040	Biotinidase Quantitative Plasma **	Chennai - Avadi	R	5350	Daily : 11:30am	10th day : 6:00pm	4 ml of Heparin plasma
B0044	BK Polyoma DNA Detection by Real Time PCR - Urine	Chennai - Avadi	C	5990	Mon & Thu : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
B0045	Bleeding & Clotting Time (BT,CT)	Chennai - Avadi	B	220	Daily : 9:00am to 9:00pm	After 4hrs	NA
B0046	Blood Group ABO & Rh Typing	Chennai - Avadi	B	180	Daily : 9:00am & 5:00pm	After 6hrs	3 ml of EDTA whole blood and Serum sample required
B0047	Bombay Blood Group Evaluation	Chennai - Avadi	D	1180	Daily : 11:30am	3rd day : 6:00pm	3 ml EDTA whole blood
B0053	BOH (Bad Obstetric History) Profile-Maxi	Chennai - Avadi	C	15255	Daily : 9:00am	15th day : 6:00pm	3 ml EDTA whole blood, 5 ml of serum, 4 ml of heparin whole blood (husband and wife), 3 ml of citrate plasma
B0054	BOH (Bad Obstetric History) Profile-Mini	Chennai - Avadi	C	6100	Daily : 7:30am	Next day : 7:00pm	3 ml EDTA whole blood, 5 ml of serum, 3 ml of citrate plasma
B0055	Bone Marrow Aspiration for morphological examination	Chennai - Avadi	C	960	Daily : 11:30am	Next day : 5:00pm	2 Stained,2 Unstained Fixed BM Smears and 2 Peripheral Smears Required.
B0056	Bone Profile	Chennai - Avadi	C	850	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0057	Bordetella Pertussis IgG antibodies Quantitative by EIA -Serum	Chennai - Avadi	C	2600	Fri : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B0058	BRAF V600 Mutation Analysis, Tissue	Chennai - Avadi	D	6490	Mon, Wed & Fri : 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
B0059	BRAF V600E Mutation Analysis, Hairy Cell Leukemia	Chennai - Avadi	D	6490	Mon & Thu : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B0062	Breast Marker Profile	Chennai - Avadi	D	2185	Daily : 9:00am to 9:00pm	After 8hrs	referring clinicians contact number and TRF mandatory. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0063	Brucella Total antibodies by Tube Agglutination - Serum	Chennai - Avadi	B	970	Daily : 9:00am & 3:00pm	Next day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B0064	Brucella IgG antibody by EIA - Serum	Chennai - Avadi	C	1355	Wed & Sat : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B0065	Brucella IgM antibody by EIA - Serum	Chennai - Avadi	C	1355	Wed & Sat : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B0066	Bullous Pemphigoid 180 Vesiculobullous disorders Serum	Chennai - Avadi	D	4600	Fri : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0067	Bullous Pemphigoid 230 Vesiculobullous disorders Serum	Chennai - Avadi	D	4600	Fri : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0068	BUN Urea Nitrogen Serum	Chennai - Avadi	A	155	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0069	BUN Urea Nitrogen Urine 24H	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
B0070	BUN Urea Nitrogen Pericardial Fluid	Chennai	A	275	Daily : 9:00am to	After 4hrs	2 ml of Pericardial Fluid

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00pm		
B0071	BUN Urea Nitrogen Urine Spot	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
B0074	Bladder Cancer/Urothelial carcinoma Marker Profile (Aneuploidy for chromosomes 3, 7, 17, and loss of the 9p21 **	Chennai - Avadi	R	14015	Daily : 11:30am	10th day : 6:00pm	100 ml Early morning Urine
B0075	BRCA deletion/duplication analysis - EDTA Whole blood	Chennai - Avadi	E	16540	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
B0076	TARGETED MUTATION ANALYSIS (FAMILY MEMBERS)	Chennai - Avadi	E	5250	Daily : 9:00am	7th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
B0077	NextGen BRCA Comprehensive (NGS & MLPA), Blood	Chennai - Avadi	E	30000	Daily : 7:30am	After 4 weeks : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Duly signed BRCA Proforma, BRCA1-2 Consent Form
B0078	Urea Nitrogen Body Fluid	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	2 ml body fluid
B0080	16s rRNA Sequencing Bacterial Identification*	Chennai - Avadi	D	3540	Daily : 9:00am	7th day : 6:00pm	Isolated colony on culture plates. Sample stable for 2 Days at 2 to 8 degree celsius.
B0081	Beta Thalassemia (Prenatal + Parental	Chennai	D	21830	Daily : 9:00am	7th day :	Amniotic Fluid (15-20 ml) in a

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Testing) - 30 Mutaions	- Avadi				6:00pm	sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition or Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
B0083	BCR-ABL Breakpoint Analysis Qualitative by Real Time PCR - Blood	Chennai - Avadi	D	5310	Tue & Fri : 9:00am	3rd day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer
B0084	BCR-ABL Breakpoint Analysis Qualitative by Real Time PCR - Bone Marrow	Chennai - Avadi	D	5310	Tue & Fri : 9:00am	3rd day : 6:00pm	3 ml of Bone Marrow in EDTA Vacutainer
B0090	Beta-Trace protein(BTP) CSF leakage in body fluid	Chennai - Avadi	D	6490	Daily : 11:30am	Next day : 5:00pm	2 ml Nasal Fluid And 3 ml Serum.
B0101	BCL6 Gene rearrangement, FISH, Tissue	Chennai - Avadi	E	8400	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B0104	(1-3)-?-D-Glucan Level by Protease Zymogen Colorimetry by Protease Zymogen Colorimetry - Serum	Chennai - Avadi	D	8510	Daily : 7:30am	Same day : 7:00pm	Collect sample in gel tube directly and send after centrifugation. Kindly do not open tube or separate serum from gel tube to another tube to avoid contamination. Aliquot sample will not accepted. Kindly send minimum 2 ml serum.
B0106	NextGen BRCA Somatic, FFPE tissue	Chennai - Avadi	E	22000	Daily : 7:30am	25th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
B0108	Pneumonia panel by FilmArray (Biofire)- Lower respiratory tract - Sputum/BAL/ Endotracheal secretions	Chennai - Avadi	E	27825	Daily : 9:00am to 9:00pm	After 10hrs	2 ml of BAL/ Sputum/ Endotracheal secretions in Sterile leak proof container
B0109	Barbiturates confirmation by LCMS, Spot Urine**	Chennai - Avadi	R	9630	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
B0111	Barium, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
B0112	Barium, Spot Urine**	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
B0113	Barium, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
B0114	Beryllium, EDTA Whole Blood**	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
B0115	Benzodiazapine confirmation by LCMS,	Chennai	R	9630	Daily : 11:30am	10th day :	10ml of spot urine sample required

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Spot Urine**	- Avadi				6:00pm	to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
B0116	Beryllium, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
B0117	Beryllium, Spot Urine**	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
B0119	Biotinidase deficiency (BTD)-Gene Analysis	Chennai - Avadi	E	19950	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
B0120	Bone Marrow-2 Comprehensive	Chennai - Avadi	D	7670	Daily : 9:00am	5th day : 6:00pm	Bone marrow biopsy (in 10% Neutral buffered formalin) & aspirate (in EDTA) examination with special stain & IHC. Send PBS as well.
B0121	BK Polyoma DNA Detection by Real Time PCR - EDTA Blood	Chennai - Avadi	C	5085	Mon & Thu : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
B0122	Bismuth ICPMS , Blood	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	3 ml EDTA Blood / Metal free EDTA Blood
B0125	Bone Marrow Aspirate Iron stain	Chennai - Avadi	C	510	Daily : 11:30am	Next day : 5:00pm	1 Stained, 2 Unstained air dried Bone marrow smears with good particles and 1 Peripheral smear required._x000D_ _x000D_
B0126	Biofire Blood Panel (BCID)	Chennai	D	19825	Daily : 9:00am to	After 10hrs	Positive blood culture bottle

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B1008	Beta Carotene**	- Avadi Chennai	R	18190	9:00pm Mon : 11:30am	5th day : 6:00pm	3 ml of serum
B1017	Borreliosis-IgG Confirm serology by western blot - Serum **	Chennai - Avadi	R	11450	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B1029	Bacterial Meningitis Panel 5 antigen by Latex agglutination - Serum	Chennai - Avadi	D	5190	Daily : 9:00am & 3:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B1030	Bone formation markers	Chennai - Avadi	D	7080	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B1031	Bone Resorption Markers	Chennai - Avadi	D	4720	Daily : 9:00am	Same day : 5:00pm	3 ml of serum and 20 ml of spot urine
B1032	NextGen BRCA Germline, Blood	Chennai - Avadi	E	16000	Daily : 9:00am	15th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Duly signed BRCA Proforma, BRCA1-2 Consent Form
B1901	Canine - BCP (Biochemistry Panel)	Chennai - Avadi	D	3660	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1902	Bicarbonate, serum by PEPC method	Chennai - Avadi	B	380	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1903	CANINE-Bile Acids-Total	Chennai - Avadi	C	1580	Wed : 9:00am	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1904	Bilirubin Direct, Serum	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1905	Bilirubin Total, Serum	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1906	Canine - Bilirubin Total, Direct, Indirect, Serum	Chennai - Avadi	D	285	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1907	Bone Marrow Aspiration Slides for review	Chennai - Avadi	C	790	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1908	Canine - BUN-Blood Urea Nitrogen	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B8001	Bacteriological Examination of Air by Vitek/ MALDITOF - Blood Agar	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	Blood agar plate , Kindly mention time of plate exposure.
B8046	BK Polyoma DNA Viral load by Real Time PCR - Urine	Chennai - Avadi	C	5085	Mon & Thu : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
B8047	BK Polyoma DNA Viral load by Real Time PCR - EDTA Blood	Chennai - Avadi	C	5085	Mon & Thu : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
B8596	Bone & Joint Infection panel (JI) by FilmArray	Chennai - Avadi	E	24000	Daily : 9:00am to 9:00pm	After 10hrs	3ml of Synovial fluid in sterila container
C0001	C1 Esterase Inhibitor Serum	Chennai - Avadi	D	2950	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0003	C3 - Complement-3 Serum	Chennai - Avadi	C	900	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0004	C4 - Complement-4 Serum	Chennai - Avadi	C	735	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0005	CA-125 Serum	Chennai - Avadi	B	1230	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0006	CA-15.3 Serum	Chennai - Avadi	B	1350	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0007	CA-19.9 Serum	Chennai - Avadi	B	1340	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0009	CA-72.4 Serum	Chennai - Avadi	C	1920	Tue & Fri : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0010	Cadmium Blood	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0011	Cadmium Urine 24H**	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0013	CAH Panel - Congenital Adrenal Hyperplasia	Chennai - Avadi	C	6780	Daily : 9:00am	Next day : 6:00pm	5ml of serum
C0014	Calcitonin Thyrocalcitonin Serum	Chennai - Avadi	C	2375	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0015	Calcium Urine 24H	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0016	Calcium Ionised Serum	Chennai - Avadi	C	565	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0017	Calcium Total Serum	Chennai - Avadi	A	200	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0018	Calcium Profile Serum	Chennai - Avadi	C	620	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0019_24	Calcium/Creatinine Ratio (24 hrs)	Chennai - Avadi	A	465	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0019_SP	Calcium/Creatinine Ratio (Spot sample)	Chennai - Avadi	A	465	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
C0020	Campylobacter antigen detection by Immunochromatography - Stool	Chennai - Avadi	C	1525	Daily : 9:00am & 3:00pm	After 4hrs	5gm of Stool specimen in Sterile leak proof container
C0021	Cancer Detection Profile, Female	Chennai - Avadi	D	9205	Daily : 9:00am to 9:00pm	5th day : 6:00pm	3 ml of EDTA whole blood, 5 ml of serum, 1 thin prep container
C0022	Cancer Detection Profile, Male	Chennai - Avadi	D	8260	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of EDTA whole blood, 5 ml of serum
C0028	Candida Albicans Ab by Indirect Haemagglutination - Serum**	Chennai - Avadi	R	11770	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0029	Cannabis Marijuana, THC Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	5 ml of Spot Urine
C0031	Calprotectin Stool	Chennai - Avadi	C	3100	Daily : 9:00am	Same day : 5:00pm	100 gms of stool in leak proof container. History of medicines, morning sample preferred.
C0032	Carbamazepine Tegretol Serum	Chennai - Avadi	C	960	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0033	Cardiac Injury Profile-Maxi	Chennai - Avadi	D	5310	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of serum
C0034	Cardiac Injury Profile-Mini	Chennai - Avadi	D	1300	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of serum
C0035	Cardiac Risk Profile	Chennai - Avadi	D	4840	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum and 3 ml of citrate plasma
C0036	Cardiac Screen	Chennai - Avadi	C	1245	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma, 10 ml of spot urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0037	Cardiolipin Antibody ACL-IgA antibody Serum	Chennai - Avadi	C	1015	Tue & Fri: 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0038	Cardiolipin Antibody ACL-IgG antibody Serum	Chennai - Avadi	C	905	Daily : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0039	Cardiolipin Antibody ACL-IgM antibody Serum	Chennai - Avadi	C	905	Daily : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0042	Catecholamines in Plasma	Chennai - Avadi	D	8025	Wed : 11:30am	Next day : 7:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
C0043	Catecholamines in urine	Chennai - Avadi	D	5310	Wed : 11:30am	Next day : 7:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0044	CBC Haemogram	Chennai - Avadi	A	380	Daily : 9:00am to 9:00pm	After 4hrs	3 ml EDTA whole blood
C0045	inv (16) CBF Beta Translocation by fish	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
C0046	CBF Beta Inversion (16) by FISH AML-M4E0 Bone Marrow	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
C0047	Cyclic Citrullinated Peptide Antibody (CCP Antibody) Serum	Chennai - Avadi	C	1450	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0048	CD2 by Flow cytometry Blood**	Chennai - Avadi	C	2485	Daily : 7:30am	4th day : 6:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0049	CD 3	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0050	CD 3/CD 4/CD 8	Chennai - Avadi	C	1865	Daily : 9:00am	Same day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0051	CD 3/CD 4/CD 8/CD 45	Chennai - Avadi	C	1865	Daily : 9:00am	Same day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0052	CD4 / CD 8 BY FLOW CYTOMETRY, BODY FLUID**	Chennai - Avadi	R	4280	Daily : 11:30am	5th day : 6:00pm	5ml BAL fluids to be send in Sterile container. Clinical history and clinician number required. (Sample require as soon as possible after collection, Avoid weekend and Holiday)
C0053	CD 4 counts	Chennai - Avadi	C	1130	Daily : 9:00am	Same day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0054	CD 4/CD 8	Chennai - Avadi	C	1865	Daily : 9:00am	Same day : 5:00pm	required. 3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0055	CD 5	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0056	CD 7	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0057	CD 8	Chennai - Avadi	C	1865	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0058	CD 10	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0059	CD13	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							condition(preferably). Clinical history and clinician number required.
C0060	CD14	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0061	CD15	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0062	CD16/56	Chennai - Avadi	C	4405	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0063	CD19	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0064	CD20	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0065	CD22	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0066	CD23	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0067	CD25	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0068	CD33	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0069	CD34 Stem Cells	Chennai - Avadi	C	6215	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0070	CD34	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0071	CD38	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0074	CD45	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0075	CD79A	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0076	CD103	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0077	CEA Carcino Embryonic Antigen Body Fluid	Chennai - Avadi	B	1080	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of Body Fluid
C0078	CEA Carcino Embryonic Antigen Serum	Chennai - Avadi	B	780	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0079	Celiac Disease Profile	Chennai - Avadi	D	5900	Mon & Thu : 9:00am	Same day : 7:00pm	5 ml of serum
C0080	Body fluids Cell Count	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	Body fluid to be send in sterile EDTA container in refrigerated condition. Clinical history and clinician number required. Body

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							fluids (CSF, Pleural, Peritoneal, Pericardial, synovial fluid)
C0081	Centromere Antibody Serum	Chennai - Avadi	C	1750	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0082	Ceruloplasmin (Copper Oxidase)Serum	Chennai - Avadi	C	1245	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0083	Cervical Screening-1 (Conventional PAP & HPV-DNA) - Smear & LBC	Chennai - Avadi	D	3125	Daily : 9:00am	5th day : 6:00pm	Smear & Thin Prep LBC Container
C0084	Cervical Screening-2 (PAP Smear and HPV-DNA) - LBC	Chennai - Avadi	D	3420	Daily : 9:00am	5th day : 6:00pm	10 ml LBC Container
C0086	CH-50 Complement, total Serum	Chennai - Avadi	C	3955	Fri : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0087	Chikungunya RNA Detection by Real Time PCR - EDTA Plasma	Chennai - Avadi	C	4520	Daily : 10:00am	After 24hrs	3ml of EDTA Plasma in Sterile Plain Vacutainer
C0088	Chikungunya IgM antibody by Immunochromatography - Serum	Chennai - Avadi	B	920	Daily : 9:00am & 3:00pm	After 4hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0090	Chitotriosidase storage disorder Plasma	Chennai - Avadi	R	4280	Daily : 11:30am	10th day : 6:00pm	6 ml Heparin Plasma
C0091	Chlamydia Trachomatis DNA Detection by Real Time PCR - Swab	Chennai - Avadi	D	4130	Wed & Sat : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0092	Chlamydia Trachomatis DNA Detection by Real Time PCR - Fluid	Chennai - Avadi	D	4130	Wed & Sat : 9:00am	Next day : 5:00pm	2ml of Body Fluid in Sterile leak proof container
C0093	Chlamydia Trachomatis DNA Detection by Real Time PCR - Urine	Chennai - Avadi	D	4130	Wed & Sat : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
C0094	Chlamydia Trachomatis DNA Detection by Real Time PCR - Vaginal Swab	Chennai - Avadi	D	4130	Wed & Sat : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0096	Chlamydia Pneumoniae IgG Antibody - Serum	Chennai - Avadi	C	2090	Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0097	Chlamydia PneumoniaeIgM Antibody - Serum	Chennai - Avadi	C	2090	Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0098	Chlamydia Trachomatis IgA Antibody by EIA - Serum	Chennai - Avadi	C	1075	Mon, Wed & Fri : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0099	Chlamydia Trachomatis IgG Antibody by EIA - Serum	Chennai - Avadi	C	1075	Mon, Wed & Fri : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0100	Chlamydia Trachomatis IgM Antibody by EIA - Serum	Chennai - Avadi	C	1075	Mon, Wed & Fri : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0101	Chloride (Cl) Serum	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0102	Chlorides Urine 24H	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 6hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0103	Chlorides Urine Spot	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 6hrs	10 ml of Spot urine
C0104	Chlorides Body Fluid	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Body Fluid
C0105	Cholesterol-Total Body Fluid	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Body Fluid
C0106	Cholesterol-Total Serum	Chennai - Avadi	A	180	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0107	Chromium, Blood	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0108	Chromium, Urine 24 Hrs	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0109	ChromiumGraphite Furnace AASUrine Spot	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0110	Chromogranin - A, Plasma	Chennai - Avadi	D	6845	Wed & Sat : 11:30am	Same day : 9:00pm	3 ml of EDTA plasma to be sent in aliquot tube at proper temperature.
C0113	Chromosomal Anaylsis (Fanconis Anemia)**	Chennai - Avadi	R	8025	Daily : 11:30am	20th day : 6:00pm	Clinical history and TRF mandatory. Blood : 5 ml Heparine Whole Blood, Sample sent as soon as possible after collection
C0115	Chylomicron Qualitative Urine	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	After 6hrs	10 ml of Spot urine
C0117	Citrate, Urine 24 H Spectrophotometry	Chennai - Avadi	C	1245	Daily : 11:30am	5th day : 6:00pm	10 mL aliquot of 24-hour Urine. Collect urine with 25 mL 50% concentrated HCl to maintain pH below 3. Record 24-hour urine volume on test request form and urine vial. Ship refrigerated or frozen.
C0118	Citrulline Quantitative Serum	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
C0119	Citrulline Quantitative Urine	Chennai - Avadi	R	6365	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
C0121	C-Kit (GIST) - Exon 9, 11, 13 and 17 by PCR - FFPE Tissue	Chennai - Avadi	D	11210	Mon, Wed & Fri : 9:00am	10th day : 6:00pm	Tumor containing FFPE block (>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours Clinical history with Histopathology report.
C0122	C-Kit - Exon 8 and 17 by PCR - Bone Marrow	Chennai - Avadi	D	6490	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
C0123	C-Kit - Exon 8 and 17 by PCR - Blood	Chennai - Avadi	D	6490	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
C0124	CK-MB Serum	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0128	Clostridium Difficile GDH Reflex Toxin A/B by CLIA - Stool	Chennai - Avadi	C	2485	Daily : 9:00am	Same day : 5:00pm	5gm of Stool specimen in Sterile leak proof container
C0130	Clotting Time (CT)	Chennai - Avadi	A	140	Daily : 9:00am to 9:00pm	After 4hrs	NA
C0131	Clozapine Serum	Chennai - Avadi	D	4130	Wed : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0132	CMV Avidity test by EIA - Serum	Chennai - Avadi	C	1355	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0133	CMV DNA Detection by Real Time PCR - Amniotic Fluid	Chennai - Avadi	C	5085	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
C0134	CMV DNA Detection by Real Time PCR - CSF	Chennai - Avadi	C	5085	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CSF in Sterile leak proof tube or container
C0135	CMV DNA Detection by Real Time PCR - EDTA Plasma	Chennai - Avadi	C	4405	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
C0136	CMV DNA Detection by Real Time PCR - BAL or Saliva	Chennai - Avadi	C	5085	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of BAL/ Saliva in Sterile Leak proof container
C0137	CMV DNA Detection by Real Time PCR - Urine	Chennai - Avadi	C	5085	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	15ml of Spot Urine in Sterile leak proof container
C0138	CMV Combo (CMV Viral load and GCV resistance) by PCR & Sequencing - EDTA Whole Blood**	Chennai - Avadi	R	10700	Tue : 11:30am	12th day : 6:00pm	3 ml of EDTA Whole Blood in EDTA Vacutainer
C0139	CMV IgG & IgM Antibody by CLIA - Serum	Chennai - Avadi	B	1240	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0141	CMV IgG Antibody by CLIA - Serum	Chennai - Avadi	B	650	DAILY : 9:00 AM to 9:00 PM	After 8Hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0142	CMV IgM Antibody by CLIA - Serum	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0143	CMV DNA Viral load by Real Time PCR - EDTA Plasma	Chennai - Avadi	C	7910	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0144	CMV by CISH - FFPE	Chennai - Avadi	D	2950	Daily : 9:00am	4th day : 6:00pm	Tissue In 10% Formalin Or Block / FFPE .
C0145	Cobalt Blood	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0145_CS	Cobalt Serum	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0146	Cobalt Urine 24H**	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0147	Cobalt Urine Spot	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0148	Cocaine Benzoyllecogonine Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	5 ml of Spot Urine
C0149	Cold AgglutininSerum	Chennai - Avadi	B	650	Daily : 11:30am	Next day : 5:00pm	3 ml of Serum
C0150	Cancer marker profile.	Chennai - Avadi	D	5900	Wed & Sat : 9:00am	Next day : 5:00pm	3 ml of Serum and 3 ml of EDTA plasma
C0151	Coombs Test-Direct ,Direct Antiglobulin Test (DAT)-EDTA Blood	Chennai - Avadi	C	565	Daily : 9:00am & 5:00pm	After 6hrs	3 ml EDTA whole blood
C0152	Coombs Test-Indirect, Indirect Antiglobulin Test (IAT)-Serum	Chennai - Avadi	C	680	Daily : 9:00am & 5:00pm	After 6hrs	3 ml Serum sample
C0154	Copper ,24 hr Urine	Chennai - Avadi	D	2300	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0155	Copper, Urine Spot	Chennai - Avadi	D	2240	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0156	Copper Serum	Chennai - Avadi	D	2240	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0157	Cortisol Serum_8.00 AM	Chennai - Avadi	B	640	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0157_4PM	Cortisol Serum_04.00 PM	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0157_R	Cortisol Serum_Random	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0158	Cortisol - Free Urine 24H	Chennai	B	970	Daily : 9:00am to	After 8hrs	24 H urine: 24 hour Urine should be

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00pm		collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0159	Coxiella Burnetti IgG antibody by EIA - Serum	Chennai - Avadi	D	2185	Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0160	Coxiella Burnetti IgM antibody by EIA - Serum	Chennai - Avadi	D	2185	Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0161	Coxsackie IgG Antibody by EIA - Serum	Chennai - Avadi	C	2090	Wed : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0162	Coxsackie IgM Antibody by EIA - Serum	Chennai - Avadi	E	2090	Wed : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0163	C-Peptide, Fasting, Serum	Chennai - Avadi	C	1180	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Serum In Fasting
C0163_PP	C-Peptide Post prandial, Serum	Chennai - Avadi	C	1245	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Serum In Fasting
C0163_R	C-Peptide, Random, Serum	Chennai - Avadi	C	1245	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Serum In Fasting
C0164	CPK Isoenzymes Serum	Chennai - Avadi	D	7670	Fri : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0165	CPK Total Serum	Chennai - Avadi	A	380	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0166	Creatinine Clearance test 24H Urine	Chennai	B	485	Daily : 9:00am to	After 4hrs	3 ml Serum And 24 Hrs Urine . No

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00pm		Preservative
C0167	Creatinine Quantitative Urine 24H	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0168	Creatinine Body Fluid	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Body Fluid
C0169	Creatinine Serum	Chennai - Avadi	A	130	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0170	CRP, C - Reactive Protein Serum	Chennai - Avadi	B	510	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0171	Cryoglobulins Serum	Chennai - Avadi	C	1075	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0172	Cryptococcus Antigen detection by Agglutination - CSF	Chennai - Avadi	C	2260	Daily : 9:00am & 3:00pm	After 6hrs	2ml of CSF in Sterile leak proof tube or container
C0173	Cryptococcus Antigen detection by Agglutination - Serum	Chennai - Avadi	C	2260	Daily : 9:00am & 3:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0179	Cryptosporidium by smear examination - Stool	Chennai - Avadi	C	565	Daily : 9:00am & 3:00pm	After 6hrs	5gm of Stool specimen in Sterile leak proof container
C0180	Cryptosporidium Parvum Antigen Detection by Immunochromatography - Stool	Chennai - Avadi	C	1470	Daily : 9:00am & 3:00pm	After 4hrs	5gm of Stool specimen in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0181	CSF Index (CSF and serum)	Chennai - Avadi	D	2360	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage). 3 ml of Serum.
C0182	Culture & Sensitivity, Aerobic bacteria - Blood single bottle	Chennai - Avadi	C	1000	Daily : 9:00am to 9:00pm	5th day : 6:00pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
C0183	Culture & Sensitivity, Aerobic bacteria - Blood Two bottle set	Chennai - Avadi	C	1695	Daily : 9:00am to 9:00pm	5th day : 6:00pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
C0184	Culture & Sensitivity, Aerobic bacteria - Blood 1 Aerobic 1 Anaerobic	Chennai - Avadi	C	1695	Daily : 9:00am to 9:00pm	5th day : 6:00pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
C0185	Culture & Sensitivity, Aerobic bacteria - Body Fluids	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	2ml of Body Fluid in Sterile leak proof container
C0186	Culture & Sensitivity, Aerobic bacteria - Bone Marrow	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	5th day : 6:00pm	3ml of Bone Marrow in MycoF/lytic Bactec Bottle or Sodium Heparin
C0187	Culture & Sensitivity, Aerobic bacteria - Central Line Catheter Tip	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	Central Line Catheter Tip in sterile container
C0188	Culture & Sensitivity, Aerobic bacteria - Clot	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	5th day : 6:00pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
C0189	Culture & Sensitivity, Aerobic bacteria - CSF	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	2ml of CSF in Sterile leak proof tube or container
C0190	Culture & Sensitivity, Aerobic bacteria - Ear Samples	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	Ear specimens in Sterile leak proof container
C0191	Culture & Sensitivity, Aerobic bacteria - Eye Samples	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	Eye specimens in Sterile leak proof container
C0192	Culture & Sensitivity, Aerobic bacteria - Other samples	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	4th day : 6:00pm	2ml of Specimen in Sterile leak proof container
C0192_C	Culture and Identification, Aerobic bacteria - Other Samples	Chennai - Avadi	C	340	Daily : 9:00am to 9:00pm	4th day : 6:00pm	Pure culture isolate is required
C0193	Culture & Sensitivity, Aerobic bacteria - Pus	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	2ml of PUS in Sterile leak proof container
C0194	Culture & Sensitivity, Aerobic bacteria - Respiratory Samples BAL,ET Secretions	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	2ml of BAL / Other Respiratory samples in Sterile leak proof container
C0195	Culture & Sensitivity, Aerobic bacteria - Semen	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	2ml of Semen in Sterile leak proof container
C0196	Culture & Sensitivity, Aerobic bacteria -	Chennai	C	1015	Daily : 9:00am to	4th day :	2ml of Sputum in Sterile leak proof

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Sputum	- Avadi			9:00pm	6:00pm	container
C0197	Culture & Sensitivity, Aerobic bacteria - Stool	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	5gm of Stool specimen in Sterile leak proof container
C0198	Culture & Sensitivity, Aerobic bacteria - Throat Swab	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0199	Culture & Sensitivity, Aerobic bacteria - Throat Swab For C.Diphtheriae	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0200	Culture & Sensitivity, Aerobic bacteria - Tissue/Biopsy	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
C0201	Culture & Sensitivity, Aerobic bacteria - Urethral Discharge	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	2ml of Urethral Discharge in Sterile leak proof container
C0202	Culture & Sensitivity, Aerobic bacteria - Urine	Chennai - Avadi	C	500	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	15ml of Spot Urine in Sterile leak proof container
C0203	Culture & Sensitivity, Aerobic bacteria - Vaginal Sample Includes Gr B Strepto	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0204	Culture and Identification, Anaerobic bacteria - Blood Single bottle	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	5th day : 6:00pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
C0205	Culture and Identification, Anaerobic bacteria - High Vagina Discharge	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	5th day : 6:00pm	2ml of Vaginal Discharge in Sterile leak proof container
C0206	Culture and Identification, Anaerobic bacteria - Pus	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	5th day : 6:00pm	2ml of PUS in Sterile leak proof container
C0207	Culture and Identification, Anaerobic bacteria (C.Difficile) - Stool	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	5th day : 6:00pm	5gm of Stool specimen in Sterile leak proof container
C0208	Culture and Identification, Anaerobic bacteria - Tissue/Biopsy	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	5th day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
C0209	Cyclospora by smear examination - Stool	Chennai - Avadi	C	565	Daily : 9:00am & 3:00pm	After 6hrs	5gm of Stool specimen in Sterile leak proof container
C0210	Cyclosporine-By LCMS, Blood	Chennai - Avadi	C	2940	Daily:7:30 am	Same day 8:00 PM	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively. Please indicate history of the drug dose taken, time of sample

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
C0212	Cyfra (Ca) 21-1, Serum	Chennai - Avadi	D	3305	Tue & Fri : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0213	CYP2C19 Genotyping	Chennai - Avadi	D	7790	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
C0214	Cystatin C, Serum	Chennai - Avadi	C	1415	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0218	Cystic Fibrosis (CF) Delta F508 Mutation Analysis	Chennai - Avadi	E	4200	Tue & Thu: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
C0219	Cystic Fibrosis IRT (Neonatal screen)Dried blood spot	Chennai - Avadi	C	510	Monday to Friday: 9.00am	Next day 5pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
C0220	Cysticercus (Taenia Solium) IgG Antibody by EIA - CSF	Chennai - Avadi	C	1810	Tue & Fri: 7:30am	Same day : 5:00pm	2 ml of CSF in sterile screw capped container.
C0221	Cysticercus (Taenia Solium) IgG Antibody by EIA - Serum	Chennai - Avadi	C	1695	Tue & Fri: 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0222	Cystine Quantitative Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
C0223	Cystine Quantitative Spot Urine	Chennai - Avadi	R	6365	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
C0224	Cytochemistry Panel	Chennai - Avadi	C	2485	Daily : 9:00am	5th day : 6:00pm	Direct smears & 3 ml of EDTA blood
C0226	Ascitic Fluid by Conventional method	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0227	Bronchoalveolar Lavage (BAL) by LBC	Chennai - Avadi	C	1130	Daily : 9:00am	Next day : 6:00pm	Broncheoalveolar lavage in a plain leak proof container along with brochoscopy findings & relevant clinical/radiological details
C0228	CSF by Cytocentrifuge	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	More than 0.5 ML fluid in plain leak proof container along with relevent clinical details & CT/MRI reports
C0229	Drain Fluid by Conventional method	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0230	Nipple Discharge by Conventional method	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	Prefixed smears along with relevant clinical & radiological details
C0233	Other Fluids by Conventional method	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0234	Pericardial Fluid by Conventional method	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	More than 1 ML fluid in plain leak proof container along with relevant clinical details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0235	Peritoneal washings by Conventional method	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0236	Pleural Fluid by Conventional method	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0238	Sputum by Conventional method	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	Sputum in a plain leak proof container along with relevant clinical/radiological details
C0239	Synovial Fluid by Conventional method	Chennai - Avadi	C	905	Daily : 9:00am	Next day : 6:00pm	0.5 ml Synovial fluid in a plain container
C0241	Urine by Liquid based cytology	Chennai - Avadi	C	1130	Daily : 9:00am	Next day : 6:00pm	More than 1 ML Urine in a leak proof plain container along with relevant clinical & radiological details
C0242	Cytology Reflex to cell block and relevant IHC Body fluid	Chennai - Avadi	C	3390	Daily : 9:00am	5th day : 6:00pm	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0244	Carbohydrate-Deficient Transferrin (CDT) **	Chennai - Avadi	R	4280	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests. International outsource : Category R for referred and I for international. Plz esquire before sending sample for arrangement of proper courier to international lab. Additional sample tubes may be needed for sample adequacy. TAT may vary based on information by concerned lab of delay. Courier charges have to be paid extra based on sample temperature condition and shipment.
C0246	Chromosomal Array CNV+SNP Array Blood	Chennai - Avadi	E	16800	Daily : 9:00am	15th day : 6:00pm	3 ml of EDTA Whole Blood
C0247	Chromosomal Array CNV+SNP Array POC	Chennai - Avadi	E	16800	Daily : 9:00am	15th day : 6:00pm	1)Upto 12 weeks of gestation: 50-100mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container 2) In 12-18 weeks of

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							gestation: a) 50-100mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container b) and or fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container 3)IUFD >20 weeks of gestation: fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container. if possible in 20ml normal saline with 2-3 drops of gentamycin
C0248	Clostridium Difficile DNA Detection by Real Time PCR - Stool	Chennai - Avadi	C	3730	Wed & Sat : 7:30am	Next day : 5:00pm	5gm of Stool specimen in Sterile leak proof container
C0249	Connexin (26 mutations)_C0249	Chennai - Avadi	R	6955	Wed : 11:30am	10th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0250	Copper reflex Ceruloplasmin, Serum	Chennai - Avadi	D	2715	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0251	Calreticulin (CALR) Exon 9 Mutation Analysis, Blood	Chennai - Avadi	D	7080	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
C0252	Calreticulin (CALR) Exon 9 Mutation Analysis, Bone marrow	Chennai - Avadi	D	7080	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0253	Chromium, Serum	Chennai - Avadi	R	2675	Daily : 11:30am	12th Day : 6:00pm	referring clinicians contact number and TRF mandatory. BD RED TOP VACCUTAINER, CENTRIFUGE, SEPARATE SERUM in red top
C0255	Creatinine Urine Spot	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
C0256	CPEO (Chronic Progressive External Ophthalmoplegia) Mitochondrial diseases	Chennai - Avadi	R	7705	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
C0260	NextGen Cancer Hotspot 50 Gene Panel	Chennai - Avadi	E	28000	Daily : 11:30am	After 4 weeks : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
C0264	CD64, Blood by Flow Cytometry	Chennai - Avadi	C	2485	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0265	CEBPA MUTATION ANALYSIS	Chennai - Avadi	E	6300	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood or 3 ml EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
C0267	Cervical Screening-3 (PAP Smear, HPV-	Chennai	D	5780	Daily : 9:00am	7th day :	10 ml LBC Container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0269	DNA and reflex CINtec PLUS) - LBC Cancer marker profile (Brain and Pituitary)	- Avadi Chennai - Avadi	D	3895	Tue & Fri : 9:00am	6:00pm After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0272	Cancer marker profile (Gall Bladder)	Chennai - Avadi	D	1770	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0273	Cancer marker profile Liver-2	Chennai - Avadi	D	6490	Wed & Sat : 9:00am	Next day : 5:00pm	3 ml of Serum and 3 ml of EDTA plasma
C0275	Cancer marker profile (Neuroendocrine tumors)	Chennai - Avadi	D	7315	Wed & Sat : 7:30am	Next day : 5:00pm	3 ml of Serum and 3 ml of EDTA plasma
C0276	Cancer marker profile-Ovarian-Comprehensive	Chennai - Avadi	D	2595	Daily : 9:00am to 9:00pm	After 8hrs	5 ml of serum
C0277	Cancer marker profile - Thyroid-DTC	Chennai - Avadi	D	2595	Daily : 9:00am to 9:00pm	After 8hrs	5 ml of serum
C0279	Cytology CINtec plus	Chennai - Avadi	D	4600	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
C0281	cMYC Gene rearrangement, FISH, BONE MARROW	Chennai - Avadi	D	4720	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
C0282	CMYC Gene rearrangement, FISH, BLOOD	Chennai - Avadi	D	4720	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
C0283	cMYC Gene rearrangement Tissue / FFPE	Chennai - Avadi	D	9440	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
C0284	Coeliac disease profile-1	Chennai - Avadi	D	2950	Mon & Thu : 9:00am	Same day : 7:00pm	5 ml of serum
C0285	Coeliac disease profile-2 Comprehensive	Chennai - Avadi	D	10620	Mon & Thu : 9:00am	3rd day : 6:00pm	5 ml of serum and 4ml EDTA / ACD blood
C0286	Catechol-O-Methyltransferase (COMT)	Chennai	E	6300	Daily : 9:00am	7th day :	3 ml of EDTA whole blood required

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Genotyping	- Avadi				6:00pm	to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0289	CMV DNA Detection by Real Time PCR - Tissue	Chennai - Avadi	C	4520	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	Tissue with Normal Saline in Sterile leak proof container
C0290	Chloride, CSF	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
C0292	CMV DNA Viral load by Real Time PCR - BAL or Saliva	Chennai - Avadi	C	7345	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of BAL/ Saliva in Sterile Leak proof container
C0297	C3 Nephritic Factor	Chennai - Avadi	R	5350	Daily : 11:30am	8th day : 6:00pm	3 ml Serum
C0302	Cryofibrinogen	Chennai - Avadi	D	1060	Daily : 11:30am	4th day : 6:00pm	Collect Blood In Pre-warmed Red Plain Top & Citrate Tube And Pre-warmed Syringe.
C0306	CHIMERISM, PRE-TRANSPLANT (STR ANALYSIS - DONOR & RECIPIENT)*	Chennai - Avadi	E	10395	Daily : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Chimerism Analysis - Information Sheet Clinical details of the recipient Bone Marrow Transplant Date, Donor and Recipient relationship
C0307	CHIMERISM (Post transplant) *	Chennai - Avadi	E	10395	Daily : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Chimerism Analysis Information Sheet Clinical details of the recipient Bone Marrow Transplant Date, Donor and Recipient relationship

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0311	Cardiomyopathy gene panel**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
C0318	Chylomicron Qualitative Body fluid	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	After 8hrs	2 ml body fluid
C0319	Clinical Exome Panel**	Chennai - Avadi	R	22050	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
C0323	Cell free DNA for EGFR*	Chennai - Avadi	E	16800	Daily : 9:00am	7th day : 6:00pm	10-15 ml of blood in streck tube required. Sample stable for 3-7 Days at Ambient Temperature. For streck tubes, enquire before sending the sample. Clinical history with Histopathology report.
C0326	Carba R (Carbapenem Resistance) by CBNAAT - Culure Plate / Anal Swab	Chennai - Avadi	D	10620	Daily : 9:00am & 3:00pm	After 6hrs	Culure Plate / Anal Swab with Sufficient Material in Sterile Falcon Tube
C0332	Cortisol/Creatinine ratio in urine	Chennai - Avadi	B	1190	Daily : 9:00am to 9:00pm	After 6hrs	10 ml of Spot urine
C0337_Adult	Creatinine, advanced (Above18 years) - spot urine	Chennai - Avadi	A	270	Daily : 9:00am to 9:00pm	After 6hrs	10 ml of Spot urine
C0338	Haemogram, advanced	Chennai - Avadi	A	380	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of EDTA whole blood sample
C0340	Chromosomal Array CGH - Prenatal*	Chennai - Avadi	E	16800	Daily : 9:00am	15th day : 6:00pm	Amniotic fluid (20-30ml divided in 2 sterile falcon centrifuge tube) and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0342	Allergy panel - Child panel	Chennai - Avadi	D	14575	Daily : 11:30am	Next day : 5:00pm	1-4ml Maternal blood in EDTA tube is required for maternal cell contamination (MCC) studies. (Completely filled Consent form G with patient and doctor stamp and signature is mandatory as per PCPNDT Act) Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
C0343	Clustered AchR and MuSK and LRP4 Antibodies**	Chennai - Avadi	R	16050	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
C0347	CSF ACPIL**	Chennai - Avadi	R	28890	Daily : 11:30am	After 4 weeks : 6:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage). International outsource : Category R for referred and I for international. Plz esquire before sending sample for arrangement of proper courier to international lab. Additional sample tubes may be needed for sample adequacy. TAT may vary based on information by concerned lab of delay. Courier charges have to be paid extra based on sample temperature condition and shipment.
C0348	Chromosomal Array CGH-Cancer	Chennai - Avadi	E	26250	Daily : 9:00am	After 4 Weeks : 6:00pm	3 ml of EDTA whole blood and 3 ml of Bone marrow
C0349	Cadmium, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
C0350	Cadmium, Water tests by ICPMS**	Chennai	R	1925	Daily : 11:30am	10th day :	2 litter of water in leakproof sterile

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				6:00pm	container. Storage and shipping at 2 to 8 degree celsius.
C0351	Chromium , Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
C0352	Chromium, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
C0353	Cobalt , Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
C0354	Cocaine confirmation by LCMS, Spot Urine**	Chennai - Avadi	R	9630	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
C0356	Copper, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
C0357	Copper, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
C0360	CLL prognostication/risk stratification gene panel**- EDTA Whole Blood	Chennai - Avadi	E	17850	Daily : 11:30am	20th day : 6:00pm	4 ml EDTA Whole blood/bone marrow aspirate
C0361	Carrier ?screening - Gold (500 genes)	Chennai - Avadi	E	43050	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0362	Carrier ?screening - Platinum (2000 genes)	Chennai - Avadi	E	70350	Daily : 7:30am	After 6 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0363	Carrier ?screening - Silver (100 genes)	Chennai - Avadi	E	36855	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0364	Congenital adrenal hyperplasia due to HSD3B2-Gene Analysis	Chennai - Avadi	E	21000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0365	Congenital cataract-Gene Panel	Chennai - Avadi	E	21000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0366	Congenital hepatic fibrosis-Gene Panel	Chennai - Avadi	E	21000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0368	Chromosomal Array CGH-High resolution	Chennai - Avadi	R	26750	Daily : 11:30am	20th day : 6:00pm	3 ml of EDTA Whole Blood
C0369	Chikungunya RNA Detection by Real Time PCR - Serum	Chennai - Avadi	C	5085	Daily : 10:00am	After 24hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0370	CMV DNA Viral load by Real Time PCR - CSF	Chennai - Avadi	C	8360	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CSF in Sterile leak proof tube or container
C0371	CMV DNA Viral load by Real Time PCR - Urine	Chennai - Avadi	C	8360	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	15ml of Spot Urine in Sterile leak proof container
C0372	CT - NG (Chlamydia Trachomatis & Neisseria Gonorrhoeae) Panel - Genital Swab	Chennai - Avadi	D	5900	Wed & Sat : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0373	CT - NG (Chlamydia Trachomatis & Neisseria Gonorrhoeae) Panel - Spot Urine	Chennai - Avadi	D	5900	Wed & Sat : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
C0374	Covid Monitor Initial	Chennai - Avadi	C	6555	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of EDTA whole blood, 3 ml of serum ,3 ml of citrate plasma
C0375	Covid Monitor Maintenance.	Chennai - Avadi	C	4405	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of EDTA whole blood, 3 ml of serum ,3 ml of citrate plasma
C0376	Covid Monitor Recovery.	Chennai - Avadi	C	3165	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of EDTA whole blood, 3 ml of serum ,3 ml of citrate plasma
C0377	Cardiac Screen-1	Chennai - Avadi	C	6780	Mon & Thu : 9:00am	Same day : 5:00pm	5 ml of serum
C0378	Conventional Karyotyping for CLL (B-Cell Stimulation)**	Chennai - Avadi	R	12840	Daily : 11:30am	20th day : 6:00pm	10 ml Heparin Blood
C0379	CoviRisk TBNK by Flow Cytometry - Blood	Chennai	D	5310	Daily : 9:00am to	After 4hrs	3 ml of EDTA Whole Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00pm		
C0381	Covid Associated Biomarker Profile	Chennai - Avadi	C	5085	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of EDTA whole blood, 3 ml of serum ,3 ml of citrate plasma
C0393	Colorectal cancer - Theranostic Panel (Hot Spot) panel	Chennai - Avadi	E	33600	Daily 7.30 am	4 Weeks	Paraffin Block with Histopath report Mandatory
C0394	Cytology - Brush Cytology by conventional method	Chennai - Avadi	C	2825	Daily : 9:00am	3rd day	Unstained slides, Complete clinical & radiology details, clinician contact number
C1001	C1q antibodies test, Serum**	Chennai - Avadi	R	6955	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
C1014	Carnitine free and total plasma	Chennai - Avadi	R	13375	Daily : 11:30am	10th day : 6:00pm	3 ml Heparin Plasma
C1802	Colour Doppler - Carotid	Chennai - Avadi	O	3915			
C1803	Colour Doppler - Limb (Arterial and Venous)	Chennai - Avadi	O	4120			
C1804	Colour Doppler - Limbs - both (Arterial and Venous)	Chennai - Avadi	O	8240			
C1805	Colour Doppler - Obstetric	Chennai - Avadi	O	3915			
C1806	Colour Doppler - Pelvis	Chennai - Avadi	O	2265			
C1807	Colour Doppler - Portal System	Chennai - Avadi	O	3915			
C1808	Colour Doppler - Renal	Chennai - Avadi	O	3915			
C1809	Colour Doppler - Scrotal	Chennai - Avadi	O	3915			
C1818	Cephalogram	Chennai - Avadi	O	515			
C1901	Canine - Calcium	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1902	Canine - CBC Haemogram	Chennai - Avadi	A	250	Daily : 9:00am to 9:00pm	After 24hrs	and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1903	CBF (Complete Body Function)	Chennai - Avadi	C	3615	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1904	Canine - Chlorides	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1905	Canine - Cholesterol-Total	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1906	Canine - Combo Profile	Chennai - Avadi	C	4070	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1907	CANINE-Cortisol	Chennai - Avadi	C	565	Daily : 9:00am to 9:00pm	After 24hrs	and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1908	CANINE-CPK-Creatinine Phospho Kinase	Chennai - Avadi	A	325	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1909	Canine - Creatinine	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1910	CANINE-Culture & Sensitivity, Aerobic bacteria Blood single bottle	Chennai - Avadi	C	905	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1911	CANINE-Culture & Sensitivity, Aerobic bacteria Body Fluids Specify type	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1912	CSF	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1913	Eye Samples	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1914	nasal swab	Chennai - Avadi	#N/A	#N/A	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1915	CANINE-Culture & Sensitivity, Aerobic bacteria Pus	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1916	Semen	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1917	CANINE-Culture & Sensitivity, Aerobic bacteria Sputum	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1918	CANINE-Culture & Sensitivity, Aerobic bacteria Stool	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1919	Throat Swab	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1920	Tissue/Biopsy	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1921	Urethral Discharge	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1922	CANINE-Culture & Sensitivity, Aerobic bacteria Urine	Chennai - Avadi	C	905	Daily : 9:00am to 9:00pm	After 24hrs	and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1923	Vaginal Sample	Chennai - Avadi	C	875	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1924	Second opinion services, Slides	Chennai - Avadi	C	1130	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1925	Ascitic Fluid	Chennai - Avadi	C	905	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1926	CSF	Chennai - Avadi	C	905	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1927	Drain Fluid	Chennai - Avadi	C	905	Daily : 9:00am	After 48hrs	and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1928	Other Fluids	Chennai - Avadi	C	1130	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1929	Other Fluids/Scraping	Chennai - Avadi	C	905	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1930	Pericardial Fluid	Chennai - Avadi	C	905	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1931	Peritoneal washings	Chennai - Avadi	C	905	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1932	Pleural Fluid	Chennai - Avadi	C	905	Daily : 9:00am	After 48hrs	and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1933	Sputum	Chennai - Avadi	C	905	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1934	Synovial Fluid	Chennai - Avadi	C	905	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1936	Metanephrine, Urine spot	Chennai - Avadi	C	3110	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	Spot Urine : 10 ml of Spot sample in Leakproof container. Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
C5054	Cadmium / Creatinine ratio	Chennai - Avadi	D	3245	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C5516	CBC+MP	Chennai - Avadi	B	465	Daily : 9:00am to 9:00pm	After 6hrs	3 ml EDTA whole blood
C5523	CBC, MP, ESR Panel	Chennai	B	515	Daily : 9:00am to	After 4hrs	3 ml EDTA whole blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C5631	COVIPROTECT Plus	- Avadi Chennai - Avadi	C	2710	9:00pm Mon, Wed & Fri : 1:00pm	Next day : 5:00pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
C5643	Comprehensive metabolic panel	Chennai - Avadi	D	5535	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of EDTA whole blood, 2 ml of Flouride plasma, 4 ml of serum, 20 ml of spot urine
C8103	Cytopathology-Expert opinion	Chennai - Avadi	C	2825	Daily : 9:00am	3rd day : 6:00pm	Stained/unstained slides with original cytopathology report & relevant clinical & radiological details
C8369	COVID VARIANT SCREENING	Chennai - Avadi	D	1120	24hr registration	After 24hrs	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
C8551	COMPREHENSIVE OPHTHALMIC GENETIC DISORDER PANEL	Chennai - Avadi	R	21400	Daily 7.30 am	4 Weeks	3 ml of EDTA whole blood
C8606	Critical care antibiotics for drug resistant Gram negative bacteria.	Chennai - Avadi	D	1200	Daily : 9:00am to 9:00pm	4th day : 6:00pm	Pure culture isolate is required
C8623	Celiac Disease Profile - 1	Chennai - Avadi	C	2100	Mon & Thu : 7:30am	Same day : 5:00pm	5 ml of serum
D0001_DCP	DCP- Decarboxy Prothrombin PIVKA II, Serum	Chennai - Avadi	D	3520	Mon, Wed, Fri : 7:30 am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0002	D-Dimer Quant, plasma	Chennai - Avadi	C	1070	Daily : 9:00am to 9:00pm	After 4hrs	1 ml of citrated platelet poor plasma in frozen condition.
D0003	Dengue IgG Antibody by CLIA - Serum	Chennai - Avadi	C	1000	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0004	Dengue IgM Antibody by CLIA - Serum	Chennai - Avadi	C	1000	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0005	Dengue NS1 antigen Detection by CLIA - Serum	Chennai - Avadi	C	1000	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0006	Dengue Profile-1 (Dengue-IgG, IgM, NS1, Platelet count) by EIA - Serum & EDTA blood	Chennai - Avadi	C	2140	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of Serum in Gel Top Vacutainer or Plain Sterile Vacutainer & 3 ml of EDTA Whole

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
D0006_2	Dengue Profile-2 (Dengue-IgG, IgM, NS1, Platelet count & IPF) by EIA - Serum & EDTA blood	Chennai - Avadi	C	2485	Daily : 9:00am to 9:00pm	After 6hrs	Blood in EDTA Vacutainer 3 ml of Serum in Gel Top Vacutainer or Plain Sterile Vacutainer & 3 ml of EDTA Whole Blood in EDTA Vacutainer
D0007	Dengue RNA Detection by Real Time PCR - CSF	Chennai - Avadi	C	4180	Daily : 10:00am	After 24hrs	2ml of CSF in Sterile leak proof tube or container
D0008	Dengue RNA Detection by Real Time PCR - EDTA Plasma	Chennai - Avadi	C	4180	Daily : 10:00am	After 24hrs	3ml of EDTA Plasma in Sterile Plain Vacutainer
D0009	Desmoglein I Antibody Vesiculobullous disorders Serum	Chennai - Avadi	D	2775	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0010	Desmoglein III Antibody Vesiculobullous disorders Serum	Chennai - Avadi	D	2775	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0011	Dexamethasone Suppression test-High dose for ACTH, Plasma	Chennai - Avadi	C	2600	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of EDTA Plasma each
D0012	Dexamethasone Suppression test for Cortisol-High dose, Serum	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	After 8hrs	3 ml Serum each
D0013	Dexamethasone Suppression test-Low dose for ACTH, Serum	Chennai - Avadi	C	2770	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of EDTA Plasma each
D0014	Dexamethasone Suppression test for Cortisol-Low dose, Serum	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	After 8hrs	3 ml Serum each
D0015	Dexamethasone Suppression test Overnight suppression For ACTH, Plasma	Chennai - Avadi	C	2600	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of EDTA Plasma each
D0016	Dexamethasone Suppression test Overnight suppression For Cortisol, Serum	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	After 8hrs	3 ml Serum each
D0017	DHEA- Dehydroepiandrostedione Serum	Chennai - Avadi	C	2825	Tue & Fri: 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0018	DHEAS- Dehydroepiandrostedione Sulphate Serum	Chennai - Avadi	B	1070	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0019	DHT- Dihydrotestosterone Serum	Chennai - Avadi	C	2950	Mon, Wed & Fri : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0020	Diabetes Profile, Mini	Chennai	E	920	Daily : 9:00am to	After 8hrs	3 ml of EDTA whole blood, 3 ml of

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00pm		serum, 2 ml of flouride plasma, 10 ml of spot urine
D0021	Diabetes Profile,Maxi	Chennai - Avadi	E	2890	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma (F, PP), 10 ml of spot urine
D0022	Dialysis Panel	Chennai - Avadi	C	4520	Daily : 9:00am to 9:00pm	After 8hrs	3ml EDTA blood, fluoroide plasma & Serum
D0023	DIC (Coagulation) Profile	Chennai - Avadi	C	2710	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of EDTA whole blood and 1 ml of citrated platelet poor plasma
D0024	Digoxin Lanoxin Serum	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0025	Diphtheria IgG antibodies by EIA - Serum	Chennai - Avadi	D	2830	Fri : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0026	DMD/BMD (18 Exons) for Males Only	Chennai - Avadi	D	6960	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
D0027	DMD/BMD-79 exons Deletion/Duplication Analysis MLPA - Blood	Chennai - Avadi	D	15900	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
D0028	DMD-Test Duchenne Muscular Dystrophy Amniotic Fluid*	Chennai - Avadi	D	17700	Daily : 9:00am	7th day : 6:00pm	Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Duly

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
D0028_c	DMD/ BMD-79 exons Deletion/ Duplication Analysis (Dystrophin gene - Xp21)	Chennai - Avadi	D	20060	Daily : 9:00am	7th day : 6:00pm	Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
D0030	DNA (Double Strand) Antibody By IFA - Reflex end point titre Serum	Chennai - Avadi	C	1525	Tue & Fri: 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0032	DNA (Double Strand) Antibody NcX DNA Serum	Chennai - Avadi	C	1750	Daily : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0033	DNA IsolationBlood	Chennai - Avadi	C	1245	Daily : 9:00am	Next day : 5:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
D0034	DNA Paternity 16 markers *	Chennai - Avadi	E	10500	Daily : 9:00am	10th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Duly signed Paternity Consent form and documents mandatory. Photo Identity proof of father and child 2 Passport size photograph attached to consent form, Birth certificate of child In case of alleged father mother declaration required.
D0036	DNP (Dinitro Phenyl hydralazine) test, urine	Chennai - Avadi	B	325	Daily : 11:30am	After 6hrs	10 ml of Spot urine
D0037	Dopamine ELISA Plasma	Chennai - Avadi	D	2950	Wed : 11:30am	Next day : 7:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
D0038	Dopamine, HPLC, 24 hrs urine	Chennai - Avadi	D	2950	Wed : 11:30am	Next day : 7:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
D0039	Maternal screen-1st Trimester Dual Marker test serum	Chennai - Avadi	C	2600	Daily 07:30 am	same day 09:00 pm	3 ml of Serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
D0039_reflex	Pregascreen Dual (First trimester test-VII, serum Dual Marker test reflex Confirmation Either NIPT or Karyotyping/FISH of CVS)	Chennai - Avadi	D	4130	Daily 07:30 am	same day 09:00 pm	3 ml of serum , If screen is positive reflex confirmation Either NIPT or Karyotyping and FISH of amniotic fluid
D0040	Maternal screen-2nd Trimester- Quadruple Marker test, Serum	Chennai - Avadi	C	3390	Daily 07:30 am	same day 09:00 pm	3 ml of Serum
D0040_reflex	Pregascreen Quadruple (Second trimester test-IV, serum Quadruple marker test reflex confirmation Either NIPT or Karyotyping/ FISH of amniotic fluid).If Positive reflex NIPT	Chennai - Avadi	C	4200	Daily 07:30 am	same day 09:00 pm	3 ml of serum , If screen is positive reflex confirmation Either NIPT or Karyotyping and FISH of amniotic fluid
D0041	Maternal screen (Triple test) 2nd trimester	Chennai - Avadi	C	2825	Daily : 9:00am to 9:00pm	Next day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0042	DPD-Deoxyypyridinoline (Pyrilinks D)Urine	Chennai - Avadi	C	1810	Daily : 9:00am to 9:00pm	After 6hrs	10 ml of Spot Urine (Protect From Light)
D0043	Dihydropyrimidine Dehydrogenase (DPD) Gene Mutation Analysis	Chennai - Avadi	D	8260	Tue & Fri : 9:00am	5th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
D0044	DPT Antibodies by EIA - Serum	Chennai - Avadi	D	7080	Fri : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0045	Drugs Of Abuse Panel-5 Drug Panel, Urine	Chennai - Avadi	C	2375	Daily : 9:00am to 9:00pm	After 8hrs	10 ml of Spot urine
D0046	Drugs Of Abuse Panel-6 Drug Panel, Urine	Chennai - Avadi	C	3050	Daily : 9:00am to 9:00pm	After 8hrs	10 ml of Spot urine
D0047	Drugs Of Abuse Panel - 9 Drug, Urine	Chennai - Avadi	C	4295	Daily : 9:00am to 9:00pm	After 8hrs	10 ml of Spot urine
D0048	DVT-Deep Vein Thrombosis Panel	Chennai - Avadi	D	12980	Mon & Thu : 9:00am	Next day : 1:00pm	1 ml of citrated platelet poor plasma in frozen condition.
D0049	Dysmorphic RBC detection - Spot Urine	Chennai - Avadi	B	380	Daily : 9:00am to 9:00pm	After 6hrs	15ml of Spot Urine in Sterile leak proof container
D0050	D3 Hydroxybutarate (Ranbut), Serum	Chennai - Avadi	C	790	Daily : 9:00am to 9:00pm	After 4hrs	3ml of serum
D0053	Drugs Of Abuse Panel-12 Drug Panel, Urine	Chennai - Avadi	C	6665	Daily : 9:00am to 9:00pm	After 8hrs	10 ml of Spot urine
D0054	DRPLA Gene Analysis	Chennai	E	6300	Mon & Wed:	7th day :	3-5 ml of EDTA whole blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00am	6:00pm	required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
D0059	Diabetes Complications assessment	Chennai - Avadi	E	2890	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of fluoride plasma (F, PP), 10 ml of spot urine
D0060	Diabetes Risk Assessment Profile	Chennai - Avadi	E	2890	Wed : 9:00am	Next day : 5:00pm	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of fluoride plasma
D0061	Diabetes type I profile	Chennai - Avadi	E	11550	Mon & Thu : 9:00am	Next day : 5:00pm	5 ml of serum
D0125	Dystonia gene panel	Chennai - Avadi	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetic history is mandatory.
D0127	DFS-70 Antibody	Chennai - Avadi	C	2375	Wed : 11:30am	Same day : 9:00pm	2 ml serum required
D0128	Deafness panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetic history is mandatory.
D0129	DLBCL (Diffuse large B-cell lymphoma) FISH PANEL Tissue	Chennai - Avadi	E	14175	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
D0130	11 Deoxycorticosterone**	Chennai - Avadi	R	16050	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
D0140	Deletion/ Duplication Analysis	Chennai - Avadi	R	11340	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
D0141	Dengue RNA Detection by Real Time PCR - Serum	Chennai - Avadi	C	4180	Daily : 10:00am	After 24hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D1802	Digital X- Ray- Ankle - AP / Lateral view	Chennai - Avadi	O	565			
D1803	Digital X- Ray- Cervical Spine - Ap / Lateral	Chennai - Avadi	O	650			
D1804	Digital X- Ray- Chest - Decubits	Chennai - Avadi	O	465			
D1805	Digital X- Ray- Chest - Lateral	Chennai - Avadi	O	465			
D1809	Digital X- Ray- Elbow joint - AP view	Chennai - Avadi	O	465			
D1810	Digital X- Ray- Femur bone - AP / Lateral	Chennai - Avadi	O	565			
D1812	Digital X- Ray- Foot - AP / Oblique	Chennai - Avadi	O	550			
D1814	Digital X- Ray- Knee - Both Skyline View	Chennai - Avadi	O	465			
D1816	Digital X- Ray- Knee - joint AP / Lat	Chennai - Avadi	O	550			

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
D1817	Digital X- Ray- Knee Joints (both) - AP/ Lateral view	Chennai - Avadi	0	650			
D1818	Digital X- Ray- KUB	Chennai - Avadi	0	565			
D1819	Digital X- Ray- L S Spine - AP / Lat	Chennai - Avadi	0	650			
D1820	Digital X- Ray- L S spine - AP / Lat / Coned down	Chennai - Avadi	0	670			
D1821	Digital X- Ray- L S Spine - Lateral	Chennai - Avadi	0	465			
D1822	Digital X- Ray- L S spine both - oblique	Chennai - Avadi	0	670			
D1823	Digital X- Ray- Mastoids Schuller s View	Chennai - Avadi	0	465			
D1824	Digital X- Ray- Nasal Bone	Chennai - Avadi	0	465			
D1825	Digital X- Ray- Pelvis with both hip	Chennai - Avadi	0	565			
D1827	Digital X- Ray- PNS Caldwell	Chennai - Avadi	0	565			
D1829	Digital X- Ray- Sacro coccyx Lat	Chennai - Avadi	0	465			
D1830	Digital X- Ray- Sacro iliac joint AP	Chennai - Avadi	0	465			
D1831	Digital X- Ray- Shoulder Joint AP /Lat	Chennai - Avadi	0	565			
D1832	Digital X- Ray- Skull Lat	Chennai - Avadi	0	465			
D1833	Digital X- Ray- Skull AP view	Chennai - Avadi	0	465			
D1834	Digital X- Ray- Tibia / Fibula AP / Lat	Chennai - Avadi	0	565			
D1835	Digital X- Ray- Wrist / Hand AP	Chennai - Avadi	0	465			
D1836	Digital X- Ray- Wrist joint AP / Lat	Chennai - Avadi	0	565			

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
D1854	Digital X- Ray - Right Mastoid	Chennai - Avadi	O	465			
D1901	DIC (Coagulation) Profile	Chennai - Avadi	C	2485	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
D5612	Diabetes Platinum	Chennai - Avadi	E	3570	Daily : 9:00am	Next day : 6:00pm	4 ml Serum ,3 ml EDTA,3 ml Flouride.1Urine container
D5613	Diabetes Gold	Chennai - Avadi	E	1470	Daily : 9:00am	Next day : 6:00pm	4 ml Serum ,3 ml EDTA,3 ml Flouride.1Urine container
D8252	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR and Karyotyping, CVS - Prenatal	Chennai - Avadi	E	15500	Daily : 9:00am	15th day : 6:00pm	â€¢ Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition â€¢ 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius .
D8254	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR and Karyotyping, Amniotic Fluid-Prenatal	Chennai - Avadi	E	15500	Daily : 9:00am	15th day : 6:00pm	â€¢ Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition â€¢ 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius
D8259	Detection of Fetal Aneuploidies (13, 18, 21	Chennai	E	18000	Daily : 9:00am	15th day :	â€¢ Amniotic Fluid (15-20 ml) in a

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	and XY) by QF PCR and Chromosomal Array CGH, Amniotic Fluid - Prenatal	- Avadi				6:00pm	sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius
D8263	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR and Chromosomal Array CGH, POC - Prenatal	Chennai - Avadi	E	18000	Daily : 9:00am	15th day : 6:00pm	POC - 1)Upto 12 weeks of gestation: 50-100 mg placental villi (fetal side) in 20 ml normal saline with 2-3 drops of gentamycin in a sterile container 2) In 12-18 weeks of gestation: a) 50-100 mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container b) and or fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container 3)IUFD >20 weeks of gestation: fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container. if possible in 20ml normal saline with 2-3 drops of gentamycin 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius
D8264	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR, Karyotyping and Chromosomal Array CGH, Amniotic Fluid - Prenatal	Chennai - Avadi	E	19000	Daily : 9:00am	15th day : 6:00pm	Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius
D8266	ALLERGY DRUG PANEL - (COMPREHENSIVE - 35 DRUG)	Chennai - Avadi	E	17325	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D8267	0	Chennai - Avadi	E	11445	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D8268	ALLERGY DRUG PANEL - (ANESTHETIC - 4 DRUG)	Chennai - Avadi	E	4830	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D8269	ALLERGY DRUG PANEL-(PAIN KILLER - 5 DRUG)	Chennai - Avadi	E	5560	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
DL001	Amphetamine, Qualitative	Chennai - Avadi	B	450	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL002	Barbiturates, Qualitative	Chennai - Avadi	B	750	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL003	Benzodiazepine Diazepam	Chennai - Avadi	B	630	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL005	Cannabis Marijuana, THC, Qualitative	Chennai - Avadi	B	630	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL006	Cocaine Benzoylcegonine, Qualitative	Chennai - Avadi	B	630	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL015	Methadone Dolophine, Qualitative	Chennai - Avadi	B	450	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL016	Opiates Morphine, Qualitative, Urine Spot	Chennai - Avadi	B	750	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL017	PCP Phencyclidine Phosphate, Qualitative, Urine Spot_DL017	Chennai - Avadi	B	450	Daily : 9:00am to 9:00pm	After 10 mins	Urine
E0001	E Coli K1 / N. Meningitidis B Antigen detection by Agglutination - CSF	Chennai - Avadi	C	1300	Daily : 9:00am & 3:00pm	After 4hrs	2ml of CSF in Sterile leak proof tube or container
E0002	E Coli K1 / N. Meningitidis B Antigen detection by Agglutination - Serum	Chennai - Avadi	C	1300	Daily : 9:00am & 3:00pm	After 4hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0003	E Coli K1 / N. Meningitidis B Antigen detection by Agglutination - Urine	Chennai - Avadi	C	1300	Daily : 9:00am & 3:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
E0004	E2 Estradiol Serum	Chennai - Avadi	B	640	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0005	E3, unconjugated Estriol Serum	Chennai - Avadi	C	1245	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0006	EBV DNA Detection by Real Time PCR - EDTA Blood	Chennai - Avadi	D	4600	Mon & Thu : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
E0007	EBV DNA Detection by Real Time PCR - CSF	Chennai - Avadi	D	4600	Mon & Thu : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
E0008	EBV DNA Detection by Real Time PCR - Tissue	Chennai - Avadi	D	4600	Mon & Thu : 9:00am	Next day : 5:00pm	Tissue with Normal Saline in Sterile leak proof container
E0009	EBV EA IgG Antibody by CLIA - Serum	Chennai - Avadi	C	2035	Mon : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0010	EBV NA IgG Antibody by CLIA - Serum	Chennai - Avadi	C	2035	Mon : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0011	EBV VCA IgG Antibody by CLIA - Serum	Chennai - Avadi	C	2035	Mon : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0012	EBV NA IgM Antibody by EIA - Serum	Chennai - Avadi	C	2035	Tue & Fri : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0013	EBV VCA IgM Antibody by CLIA - Serum	Chennai - Avadi	C	1810	Daily : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0014	EBV encoded RNA (EBER) by CISH - FFPE	Chennai - Avadi	D	3540	Daily : 9:00am	4th day : 6:00pm	Tissue In 10% Formalin Or Block / FFPE .
E0016	Echinococcus IgG Antibody by EIA - Serum	Chennai - Avadi	C	1750	Tue & Fri: 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0021	Electrolytes Serum	Chennai - Avadi	B	480	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0022	Electrolytes Urine 24H	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 6hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
E0023	Electrolytes Urine Spot	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 4hrs	for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab. 10 ml of Spot urine
E0024	Electrolytes Dialysis Fluid	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Dialysis Fluid
E0027	Endomysial Antibody IgA antibody Serum	Chennai - Avadi	C	2375	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0028	Endomysial Antibody IgG antibody Serum	Chennai - Avadi	C	3390	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0030	Enteroviruses RNA Detection by PCR - CSF**	Chennai - Avadi	R	5350	Daily : 11:30am	5th day : 6:00pm	2 ml of CSF in Sterile leak proof tube or container
E0031	Enteroviruses RNA Detection by PCR - Plasma**	Chennai - Avadi	R	5350	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA Plasma in Sterile Plain Vacutainer
E0032	Eosinophil Detection, URINE	Chennai - Avadi	B	325	Daily : 9:00am to 9:00pm	After 8hrs	5 ml urine sample
E0033	Eosinophil detection (SPUTUM)	Chennai - Avadi	B	325	Daily : 9:00am to 9:00pm	After 8hrs	Sputum In Sterile Container
E0035	Erythropoietin Serum	Chennai - Avadi	C	2260	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0038	ESR Automated Blood	Chennai - Avadi	A	110	Daily : 9:00am to 9:00pm	After 4hrs	3 ml EDTA Whole Blood
E0039	Everolimus, Blood by LCMS method	Chennai - Avadi	D	6100	Daily : 7:30am	Same day : 7:30pm	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
E0042	E1, Estrone, Serum	Chennai - Avadi	R	6420	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
E0043	EGFR MUTATION DETECTION (EXONS 18, 19, 20, 21)	Chennai - Avadi	D	9440	Mon, Wed & Fri : 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>10% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
E0045	EWS gene rearrangement (:Translocation t(11;22) /t(21;22) translocation/ rearrangement **	Chennai - Avadi	R	6955	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
E0051	EBV DNA Viral load by Real Time PCR - EDTA Blood	Chennai - Avadi	D	7080	Mon & Thu : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
E0053	Extended RAS Test (KRAS, BRAF, NRAS, PIK3CA) - FFPE Tissue	Chennai - Avadi	D	17700	Mon, Wed & Fri : 9:00am	10th day : 6:00pm	Tumor containing FFPE block. Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
E0055	Epileptic encephalopathy gene panel**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
E0056	EGFR T790M Mutation (Acquired EGFR TKI Resistance)*	Chennai - Avadi	E	8400	Daily : 9:00am	7th day : 6:00pm	Tumor containing FFPE block (>10% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
E0057	EXPANDED LUNG CANCER PANEL BY FISH	Chennai - Avadi	E	25200	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
E0064	Allergy panel -Common Eczema	Chennai - Avadi	D	8260	Daily : 11:30am	Next day : 5:00pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
E0065	Allergy panel - Eczema Comprehensive Veg	Chennai - Avadi	D	19470	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
E0066	Allergy panel - Eczema Comprehensive Veg + non Veg	Chennai - Avadi	D	24015	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
E0068	ERA-Endometrium Receptivity Assay, Endometrial biopsy**	Chennai - Avadi	R	44940	Daily : 11:30am	15th Day	(2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day. CytogeneticsandFISH.pdf
E0069	EWSR1-FLI1 **	Chennai - Avadi	R	12840	Daily : 11:30am	15th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
E0070	EWSR1(22q12)/NTRK3 (15q25) translocation **	Chennai - Avadi	R	12840	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
E0071	EBV DNA Viral load by Real Time PCR - CSF	Chennai - Avadi	D	7435	Mon & Thu : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
E0072	EBV DNA Viral load by Real Time PCR - Tissue	Chennai - Avadi	D	7435	Mon & Thu : 9:00am	Next day : 5:00pm	Tissue with Normal Saline in Sterile leak proof container
E0075	Endometritis Panel	Chennai - Avadi	D	8495	Wed : 9:00am	Next day : 6:00pm	Genital discharge, Ulcer, Ulcer base, High Vaginal swab, Endocervical swab, Endometrial Pus, Endometrial Fluid, Endometrial Tissue.
E1019	E Histolytica Antigen Detection by Immunochromatography -Stool **	Chennai - Avadi	R	2890	Daily : 11:30am	3rd day : 6:00pm	1gm of Stool specimen in Sterile leak proof container
E1020	Entamoeba histolytica detection by Wet Mount - Stool	Chennai - Avadi	B	325	Daily : 9:00am & 3:00pm	After 6hrs	5gm of Stool specimen in Sterile leak proof container
E1021	Kidney Biopsy by Electron Microscopy	Chennai - Avadi	R	10165	Daily : 7:30am	15 day	Biopsy in 3% Gluteraldehyde along with histopathology report & relevent clinical details
E1902	E2 estradiol, Serum	Chennai	B	565	Daily : 9:00am to	After 24hrs	Please contact lab before sending

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00pm		samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
E1903	ESR Automated, EDTA whole blood	Chennai - Avadi	A	85	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
E8141	Endometrial Cancer panel by NGS	Chennai - Avadi	R	45000	Daily 7.30 am	4 Weeks	Paraffin Block
F0002	Factor V Leiden-Mutant Detection (G1691A Mutation) by Real Time PCR - Blood	Chennai - Avadi	D	5780	Wed & Sat : 9:00am	Next day : 5:00pm	3 ml of EDTA Whole Blood in EDTA Vacutainer
F0003	Factor V Level	Chennai - Avadi	D	6490	4th Mon : 9:00am	Next day : 5:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0004	Factor VII Activity Proconvertin Citrated plasma	Chennai - Avadi	R	3745	Daily : 11:30am	10th day : 6:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0005	Factor VIII Activity Anti Haemophilia factor A, Citrated plasma	Chennai - Avadi	D	2595	Mon : 9:00am	Next day : 1:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0006	Factor IX Activity Anti Haemophilic factor B (Christmas Disease), Citrated plasma	Chennai - Avadi	D	2360	Mon : 9:00am	Same day : 8pm	1 ml of citrated platelet poor plasma in frozen condition.
F0007	Factor X Activity Stuart Power Citrated plasma	Chennai - Avadi	R	4815	Daily : 11:30am	10th day : 6:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0008	Anti Factor X A Assay**	Chennai - Avadi	R	5350	Daily : 11:30am	4th day : 6:00pm	Sample 3 ml Citrate Plasma Frozen with medication history.
F0009	Factor XI Activity Rosenthal Citrated plasma	Chennai - Avadi	R	3745	Daily : 11:30am	10th day : 6:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0010	Factor XII Activity Hageman Citrated plasma	Chennai - Avadi	R	3745	Daily : 11:30am	5th day : 6:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0011	Factor XIII Activity, Citrated plasma	Chennai - Avadi	R	1390	Daily : 11:30am	Next day : 5:00pm	1 ml of citrated platelet poor plasma in frozen condition.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
F0012	Fatty Acids With Very Long Chains C22 To C26 Plasma**	Chennai - Avadi	R	16050	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA Plasma Or Heparin plasma
F0014	FDP DR 70	Chennai - Avadi	R	3850	Daily : 11:30am	7th day : 6:00pm	3 ml of Serum and 3 ml of citrate plasma
F0015	FDP Fibrin degradation products Reflex to end point titre.	Chennai - Avadi	C	1075	Daily : 9:00am, 3:00pm & 9:00pm	After 6hrs	1 ml of citrated platelet poor plasma in frozen condition.
F0016	Febrile by Agglutination - Serum	Chennai - Avadi	B	1405	Daily : 9:00am & 3:00pm	Next day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
F0017	Ferric Chloride, Urine	Chennai - Avadi	C	340	Daily : 11:30am	Next day : 5:00pm	10 ml of Spot urine
F0018	Ferritin Serum	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0019	Fertility Endocrine Panel-Female	Chennai - Avadi	B	1945	Daily : 9:00am to 9:00pm	After 8hrs	5 ml of serum
F0020	Fertility Endocrine Panel-Male	Chennai - Avadi	B	2700	Tue, Thu & Sat: 7:30am	After 8hrs	5 ml of serum
F0021	Fever Profile (Maxi)	Chennai - Avadi	C	2935	Daily: 9am to 9pm	Same Day	1.4 ml of EDTA whole blood 3 ml of serum 10 ml of spot urine
F0022	Fever Profile (Mini)	Chennai - Avadi	C	800	Daily: 9am to 9pm	Report Same Day: 1. Complete Blood Count(CBC Haemogram) 2. ESR - Erythrocyte Sedimentation Rate 3. Malaria Parasite Detection 4. Routine Examination Urine 4 ml of EDTA whole blood 3 ml of serum	4ml of serum 3 ml of EDTA whole blood 20 ml of spot urine Collect blood 2-3 ml in Bactec FX aerobic bottle

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
						10 ml of spot urine F5604 Fever Profile (Super Maxi) Chennai - Avadi C 4745 Daily : 9:00am Report Same Day: 1. Complete Blood Count(CBC Haemogram) 2. Routine Examination Urine 3. CRP - C Reactive Protein 4. Malaria Parasite Detection 5. SGPT 6. Typhi Dot IgM 7. Dengue	
F0023	Fibrinogen Plasma (F0023)	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	After 4hrs	1 ml of citrated platelet poor plasma in frozen condition.
F0028	Fibromax Fibrosis, Activity, STEATO, ASH and NASH of Liver.	Chennai - Avadi	E	13650	Daily : 9:00am to 9:00pm	5th day : 6:00pm	5 ml of serum and 2 ml of fluoride plasma . Fasting sample
F0029	Filaria Antibody detection by Immunochromatography - Serum	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
F0030	Filaria Antigen detection by Immunochromatography - EDTA blood	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer
F0031	Filaria Profile	Chennai - Avadi	B	1240	Daily : 9:00am & 3:00pm	After 6hrs	3 ml EDTA whole blood and serum sample
F0032	FLT-3 (AML) by PCR ITD, D835 - Blood	Chennai - Avadi	D	6490	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
F0033	FLT-3 (AML) by PCR ITD, D835 - Bone Marrow	Chennai - Avadi	D	6490	Mon & Wed: 9:00am	7th day : 6:00pm	Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. 3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
F0034	FMC 7 by Flow cytometry	Chennai - Avadi	D	2715	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
F0039	By Conventional method	Chennai - Avadi	C	1075	Daily : 9:00am	3rd day : 6:00pm	Prefixed smears along with air dry smears, relevant clinical & radiological details mandatory
F0040	By LBC	Chennai - Avadi	C	1470	Daily : 9:00am	3rd day : 6:00pm	More than 1 ML fluid in a leak proof plain container along with relevant clinical & radiological details
F0043	Procedure charge	Chennai - Avadi	C	565	Daily : 9:00am	Same day : 6:00pm	NA
F0046	Folic Acid Serum	Chennai - Avadi	C	1180	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0048	Food intolerance Test	Chennai - Avadi	D	11800	Daily : 11:30am	3rd day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0050	Fragile X (FMR1)CGG Repeat Analysis	Chennai - Avadi	D	8140	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
F0051	Friedreich Ataxia Mutation Analysis	Chennai - Avadi	D	6490	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
F0052	Free Androgen Index FAI	Chennai - Avadi	D	2950	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0053	Free Beta HCG Serum	Chennai - Avadi	C	1415	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0054	Free Kappa Freelite chains Serum	Chennai - Avadi	D	2655	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0055	Free Lambda Freelite chains Serum	Chennai - Avadi	D	2655	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0056	Fructosamine Serum	Chennai - Avadi	C	1015	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0057	Fructose biochemical by Biochemical - Semen	Chennai - Avadi	B	325	Daily : 9:00am & 3:00pm	After 4hrs	2ml of Semen in Sterile leak proof container
F0058	FSH-Follicle Stimulating Hormone ECLIA Serum	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0059	FSH-LH-Prolactin Serum	Chennai - Avadi	C	1450	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0060	FSH-LH-Prolactin-Testosterone	Chennai - Avadi	C	1750	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0061	FSH-LH-Prolactin-TSH	Chennai - Avadi	C	1640	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0062	FSH-LH-Testosterone	Chennai - Avadi	C	1415	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
F0065	Fungal Antibiogram (Yeast) MIC Panel - Pure culture	Chennai - Avadi	C	2825	Daily : 9:00am to 9:00pm	5th day : 6:00pm	tests. Pure culture Growth on Agar slant/ Plate
F0066	Fungal Identification - Pure culture	Chennai - Avadi	C	2260	Daily : 9:00am to 9:00pm	10th day : 6:00pm	Pure culture on SDA plate/ slant
F0067	Fungal stain by KOH & Calcofluor White stain - BAL /E.T Secretion/Respiratory Sample	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	2ml of BAL / Other Respiratory samples in Sterile leak proof container
F0068	Fungal stain by KOH & Calcofluor White stain- Specify type Body fluid	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	3ml of Fluid in Sterile leak proof container
F0069	Fungal stain by KOH & Calcofluor White stain - CSF	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	2ml of CSF in Sterile leak proof tube or container
F0070	Fungal stain by KOH & Calcofluor White stain - Ear Sample	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	Ear specimens should be collected by a medical specialist, an experienced technologist or nurse. Aspirate discharge and collect in a sterile container or collect with a sterile swab. Transport at the earliest. Amies transport medium can be used for transportation.
F0071	Fungal stain by KOH & Calcofluor White stain- Eye Sample	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	Eye specimens in Sterile leak proof container
F0072	Fungal stain by KOH & Calcofluor White stain - Hair	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	Hair sample in Sterile leak proof container
F0073	Fungal stain by KOH & Calcofluor White stain- Nail	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	Nail sample in Sterile leak proof container
F0074	Fungal stain by KOH & Calcofluor White stain- Other samples	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Vaginal swab/ Vaginal discharge in Sterile leak proof container.
F0075	Fungal stain by KOH & Calcofluor White stain- Pus	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	2ml of PUS in Sterile leak proof container
F0076	Fungal stain by KOH & Calcofluor White stain- Skin	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	Skin sample in Sterile leak proof container
F0077	Fungal stain by KOH & Calcofluor White stain- Sputum	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Sputum in Sterile leak proof container
F0078	Fungal stain by KOH & Calcofluor White stain- Tissue	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	Tissue with Normal Saline in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
F0079	Fungal stain by KOH & Calcofluor White stain- Urine	Chennai - Avadi	B	430	Daily : 9:00am & 3:00pm	After 6hrs	15ml of Spot Urine in Sterile leak proof container
F0080	Fungal culture and identification by MALDITOF (Vitek MS) - BAL, E.T Secretion Respiratory Sample	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	2ml of BAL / Other Respiratory samples in Sterile leak proof container
F0081	Fungal culture and identification by MALDITOF (Vitek MS) - Blood	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
F0082	Fungal culture and identification by MALDITOF (Vitek MS) - Body Fluid	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	2ml of Body Fluid in Sterile leak proof container
F0083	Fungal culture and identification by MALDITOF (Vitek MS) - Bone Marrow	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	3ml of Bone Marrow in MycoF/lytic Bactec Bottle or Sodium Heparin
F0084	Fungal culture and identification by MALDITOF (Vitek MS) - CSF	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	2ml of CSF in Sterile leak proof tube or container
F0085	Fungal culture and identification by MALDITOF (Vitek MS) - Ear Sample	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	Ear specimens in Sterile leak proof container
F0086	Fungal culture and identification by MALDITOF (Vitek MS) - Eye Sample	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	Eye specimens in Sterile leak proof container
F0087	Fungal culture and identification by MALDITOF (Vitek MS) - Hair	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	Hair sample in Sterile leak proof container
F0088	Fungal culture and identification by MALDITOF (Vitek MS) - Nail	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	Nail sample in Sterile leak proof container
F0089	Fungal culture and identification by MALDITOF (Vitek MS) - Other Sample	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	2ml of Vaginal swab/ Vaginal discharge in Sterile leak proof container.
F0090	Fungal culture and identification by MALDITOF (Vitek MS) - Pus	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	2ml of PUS in Sterile leak proof container
F0091	Fungal culture and identification by MALDITOF (Vitek MS) - Skin	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	Skin sample in Sterile leak proof container
F0092	Fungal culture and identification by MALDITOF (Vitek MS) - Sputum	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	2ml of Sputum in Sterile leak proof container
F0093	Fungal culture and identification by MALDITOF (Vitek MS) - Stool	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	5gm of Stool specimen in Sterile leak proof container
F0094	Fungal culture and identification by MALDITOF (Vitek MS) - Tissue	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
F0095	Fungal culture and identification by MALDITOF (Vitek MS) - Urine	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	21st day : 6:00pm	15ml of Spot Urine in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
F0099	Fever Profile Detection by Real Time Multiplex PCR - EDTA Blood	Chennai - Avadi	D	5780	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
F0103	FIP1L1-PDGFR A GENE REARRANGEMENT **	Chennai - Avadi	R	11770	Tue : 11:30am	12th day : 6:00pm	4 ML EDTA whole blood / Bone marrow
F0105	Flu Profile-1 Detection by Real Time PCR - Nasopharangeal & Oropharangeal Swab	Chennai - Avadi	D	6960	Daily : 7:00am & 3:00pm	Same Day : 3:00pm & 11:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
F0106	Fibroblast Growth Factor 23, Plasma	Chennai - Avadi	D	11800	Wed 12:00 pm	same day 08:00 pm	3 ml of EDTA plasma to be sent in aliquot tube at proper temperature.
F0107	Flow for Lymphoma/ CLPD / lymphoma MRD**	Chennai - Avadi	R	15515	Daily : 7:30am	5th day : 6:00pm	3 ml of EDTA Bone Marrow
F0109	Flow for B ALL MRD**	Chennai - Avadi	R	15515	Daily : 7:30am	5th day : 6:00pm	3 ml of EDTA Bone Marrow
F0110	Flow for T ALL MRD**	Chennai - Avadi	R	19795	Daily : 7:30am	5th day : 6:00pm	3 ml of EDTA Bone Marrow
F0112	Food Intolerance (Vegetarian panel)	Chennai - Avadi	D	11800	Daily : 11:30am	3rd day : 6:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
F0115	First trimester test-V, Serum Quad Test Free beta HCG, PAPPa, AFP and PLGF	Chennai - Avadi	D	5780	Daily cut off 7:30 AM	same day 9:00 PM	3 ml of serum is required. Mandatory Requirements: Latest Ultrasound (USG) report with a completely filled Maternal screening TRF form including CRL/BPD, nasal bone status, ductus venosus status, and NT value in millimeters. Please note that this test is not valid for triplets and above Preeclampsia: USG findings of mean arterial pressure (MAP) and uterine artery pulsatility index (UAPI), the patients blood pressure (BP), and height should be included as per maternal TRF requirements .
F0116	First trimester test-VI, serum Penta Screen	Chennai	R	6420	Daily : 11:30am	10th day :	3 ml of serum is required.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Free beta HCG, PAPPa, AFP, Inhibin A, PLGF and NT 10 to 13.6 TBA	- Avadi				6:00pm	Mandatory Requirements: Latest Ultrasound (USG) report with a completely filled Maternal screening TRF form including CRL/BPD, nasal bone status, ductus venosus status, and NT value in millimeters. Please note that this test is not valid for triplets and above Preeclampsia: USG findings of mean arterial pressure (MAP) and uterine artery pulsatility index (UAPI), the patients blood pressure (BP), and height should be included as per maternal TRF requirements
F0117	Fragile X mutation analysis+Chromosomal Array profile	Chennai - Avadi	E	22050	Tue & Fri : 9:00am	7th day : 6:00pm	3 ml of Blood In Na-Heparin +3 ml of EDTA Vacutainer
F0118	FOXO1 (13q14) translocation **	Chennai - Avadi	R	17120	Daily : 11:30am	12th Day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
F0119	COVI-INFLUENZA VIRUS PANEL BY REAL TIME PCR (Inf A, Inf B, H1N1(Swine Flu), H3N2, RSV, Covid-19)	Chennai - Avadi	D	6490	Daily : 7:00am & 3:00pm	Same Day : 3:00pm & 11:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
F0120	Mould/Filamentous fungi by microbroth dilution - Agar plate or slant	Chennai - Avadi	C	11750	Daily : 9:00am to 9:00pm	15th day : 6:00pm	Pure culture Growth of Filamentous fungi on Agar slant/ Plate.
F0122	Factor VIII inhibitor assay**	Chennai - Avadi	D	7670	Daily : 11:30am	3rd day : 6:00pm	3 ml Citrate Plasma (Fasting & Frozen mandatory)
F0126	Foundation one CDX	Chennai - Avadi	R	315650	Daily : 7:30am	After 5 weeks : 6:00pm	Paraffin Block (Carcinoma)
F1003	Free Fatty Acid (Non-Esterified)	Chennai - Avadi	D	14750	Daily : 7:30am	10th day : 6:00pm	3 ml Serum
F1901	Second opinion services, Slides	Chennai	C	1130	Daily : 9:00am	After 48hrs	Please contact lab before sending

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi					samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1902	Lymphnode	Chennai - Avadi	C	1075	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1903	CANINE-FNAC-Others, conventional	Chennai - Avadi	C	1075	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1904	CANINE-FNAC-Thyroid, conventional	Chennai - Avadi	C	1075	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1905	Procedure Charges by appointment	Chennai - Avadi	C	565	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1906	Fungal stain	Chennai	C	395	Daily : 9:00am	After 24hrs	Please contact lab before sending

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi					samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F5605	Fatty Liver	Chennai - Avadi	D	6805	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of EDTA whole blood, 2 ml of Flouride plasma, 3 ml of serum
F5606	Fever Profile (Super Mini)	Chennai - Avadi	C	2050	Daily : 9:00am	Same Day	4 ml of EDTA whole blood 3 ml of serum 10 ml of spot urine
F8194	Foundation one Heme**	Chennai - Avadi	R	#N/A	Daily 7.30 am	5 weeks	Paraffin Block (Sarcoma)
F8215	FIT-Fecal Immunochemical Test, quantitative	Chennai - Avadi	D	4200	Daily : 9:00AM	Same day : 7:00 PM	100 gmsof stool is required
G0001	G6PD Neonatal Screen Dried blood spot	Chennai - Avadi	C	510	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0002	G6PD Qualitative Blood	Chennai - Avadi	B	595	Daily : 9:00am to 3:00pm	After 8hrs	3 ml EDTA Whole Blood
G0003	G6PD Quantitative Reflex Qualitative	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	After 8hrs	3 ml EDTA Whole Blood
G0004	GAD Total Antibody Type I Diabetes Serum	Chennai - Avadi	D	6020	Tue, Thu & Sat: 7:30am	Same day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0005	Galactose 6 Sulphate Sulphatase Blood	Chennai - Avadi	R	6955	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
G0006	Galactose Classical Transferase Blood**	Chennai - Avadi	R	2245	Daily : 11:30am	8th day : 6:00pm	6 ml Sodium heparin whole blood required as soon as possible after collection
G0008	Galactose screening,Blood **	Chennai - Avadi	R	855	Daily : 11:30am	7th day : 6:00pm	4 ml of Blood In Sodium Heparin
G0011	Galactose-Total Neonatal Screen Dried blood spot	Chennai - Avadi	C	510	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0012	Ganglioside IgG Antibody,serum	Chennai - Avadi	D	3895	Tue : 11:30am	Next day : 5:00pm	post 72 hrs after transfusion. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0013	Ganglioside IgM Antibody,serum	Chennai - Avadi	D	3895	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0016	Gastrin Serum	Chennai - Avadi	C	1525	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Serum.Fasting Required
G0017	GBM-Antibody Glomerular Basement Membrane serum	Chennai - Avadi	C	1750	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0018	Genetic Hair Loss Test Alopecia Blood	Chennai - Avadi	E	4725	Daily : 7:30am	15th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0019_A	GFR With Creatinine (Adult - Above 18 Years) Serum	Chennai - Avadi	B	320	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0019_C	GFR With Creatinine (Upto 18 Years) Serum	Chennai - Avadi	B	325	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0020	GGTP Gamma GT Serum	Chennai - Avadi	A	260	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0021	Giardia Antigen detection by Immunochromatography - Stool	Chennai - Avadi	B	1240	Daily : 9:00am & 3:00pm	After 4hrs	5gm of Stool specimen in Sterile leak proof container
G0022	Gliadin-IgA Antibodies De-amidated (DGP)Serum	Chennai - Avadi	C	2145	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0023	Gliadin-IgG Antibodies De-amidated (DGP)Serum	Chennai - Avadi	C	2145	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tests.
G0024	Glucagon Stimulation Test - For C-peptide, Serum	Chennai - Avadi	D	4130	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0min, 5, 10, 15 min)
G0025	Glucosamine N Acetyl Transferase Blood	Chennai - Avadi	R	8025	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
G0026	Glucose Urine spot	Chennai - Avadi	A	135	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
G0027	Fasting Glucose, plasma	Chennai - Avadi	A	50	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_100	Glucose Post glucose (100 GMS) Plasma	Chennai - Avadi	A	75	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_120	GLUCOSE 120 MIN	Chennai - Avadi	A	100	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_180	GLUCOSE 180 MIN	Chennai - Avadi	A	100	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_240	GLUCOSE 240 MIN	Chennai - Avadi	A	100	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_300	GLUCOSE 300 MIN	Chennai - Avadi	A	100	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_50	Glucose Post glucose (50 GMS) Plasma	Chennai - Avadi	A	75	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_60	GLUCOSE 60 MIN	Chennai - Avadi	A	100	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_75	Glucose plasma(75 GMS)	Chennai - Avadi	A	75	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_P	Glucose post prandial plasma	Chennai - Avadi	A	50	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0027_R	Random, Glucose	Chennai - Avadi	A	60	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma
G0028	Glucose CSF	Chennai - Avadi	A	135	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Body Fluid
G0028_B	Glucose Body Fluid	Chennai - Avadi	A	135	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Body Fluid
G0028_S	Glucose Synovial Fluid	Chennai	A	135	Daily : 9:00am to	After 4hrs	2 ml of Body Fluid

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0029	Glucose Tolerance Test 2 hours 3 sample	- Avadi Chennai	B	225	Daily : 9:00am to 9:00pm	After 6hrs	0 min, 60 , 120 min -2ml Fluoride plasma
G0030	Glucose Tolerance Test 2.5 hours 5 samples	- Avadi Chennai	B	380	Daily : 9:00am to 9:00pm	After 6hrs	0 min, 30, 60, 120, 150 min - 2ml Fluoride plasma
G0031	Glucose Tolerance Test 5 hours 8 samples	- Avadi Chennai	B	605	Daily : 9:00am to 9:00pm	After 6hrs	0 min, 30, 60, 90, 120, 180, 240, 300 min - 2ml Fluoride plasma
G0032	Glucose Tolerance Test 3 samples for glucose and insulin	- Avadi Chennai	C	2260	Daily : 9:00am to 9:00pm	After 6hrs	0 min, 60 , 120 min -2ml Fluoride plasma and 3ml serum
G0033	Glutamic Acid Quantitative Urine 24H**	- Avadi Chennai	R	6365	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
G0034	Glutamic Acid Quantitative Plasma*	- Avadi Chennai	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
G0035	Glutamine Quantitative Plasma	- Avadi Chennai	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
G0036	Glutamine Quantitative Urine 24H**	- Avadi Chennai	R	6365	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
G0037	Glycine Quantitative CSF	- Avadi Chennai	D	7020	Daily : 11:30am	7th day : 6:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
G0038	Glycine Quantitative Urine 24H**	- Avadi Chennai	R	6365	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
G0039	Glycine Quantitative Plasma*	- Avadi Chennai	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
G0041	Glycophorin-A	- Avadi Chennai	D	2595	Daily : 9:00am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
G0042	GNRH (LHRH) Stimulation Test For FSH and LH	- Avadi Chennai	B	3890	Daily : 9:00am to 9:00pm	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0044	Goodpasture's Syndrome Panel	Chennai - Avadi	D	3540	Fri : 7:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0045	Gram stain For Neisseria sp by Microscopy - Urethral discharge	Chennai - Avadi	B	295	Daily : 9:00am & 3:00pm	After 6hrs	2 ml of Urethral Discharge in Sterile leak proof container
G0046	Gram stain Nugent s score by Microscopy - Vaginal discharge	Chennai - Avadi	B	380	Daily : 9:00am & 3:00pm	After 6hrs	Discharge on sterile swab in Amies transport medium
G0047	Gram stain by Microscopy - Other samples	Chennai - Avadi	B	295	Daily : 9:00am & 3:00pm	After 6hrs	2 ml of Specimen in Sterile leak proof container
G0048	Growth Disorder Panel	Chennai - Avadi	D	6135	Daily : 9:00am to 9:00pm	After 8hrs	5 ml of serum
G0049	Growth Hormone HGH Serum	Chennai - Avadi	C	850	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Serum In Fasting And Resting
G0050	Growth Hormone Stimulation Test- Arginine, Serum	Chennai - Avadi	D	3540	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0051	Growth Hormone Stimulation Test - Clonidine, Serum	Chennai - Avadi	D	3540	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0052	Growth Hormone Stimulation Test- Exercise, Serum	Chennai - Avadi	D	3540	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0053	Growth Hormone Stimulation Test - GHRH, Serum	Chennai - Avadi	D	3895	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0054	Growth Hormone Stimulation Test-Insulin	Chennai - Avadi	D	3540	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0055	Growth hormone stimulation test for L Dopa, Serum	Chennai - Avadi	D	3895	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0056	Growth Hormone Suppression Test- Glucose, Serum	Chennai - Avadi	D	3540	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0058	Anti GAD-65 Paraneoplastic Syndrome Serum	Chennai - Avadi	D	7670	Wed : 1:00pm	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0059	GFR With Cystatin C Serum	Chennai - Avadi	C	1245	Daily : 11:30am	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0059_A	GFR With Cystatin C Serum (Adult - Above 18 Years)	Chennai - Avadi	C	1560	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0059_C	GFR With Cystatin C Serum (Upto 18 Years)	Chennai - Avadi	C	340	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0061	Glycosaminoglycan,Urine	Chennai - Avadi	R	3315	Daily : 11:30am	15th day : 6:00pm	20 ml Spot Urine Sample In Leak Proof Container
G0065	GTT-Gestational ,Two Hour(75 grams glucose)	Chennai - Avadi	B	290	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Fluoride Plasma Each
G0066	GTT-Gestational ,Three Hour(100 grams glucose)	Chennai - Avadi	B	390	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Fluoride Plasma Each
G0073	G6pd Gene Mutation**	Chennai - Avadi	R	4495	Tue : 11:30am	12th day : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability. Amniotic Fluid or Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition.
G0074	GALACTOSEMIA GENE MUTATION**	Chennai - Avadi	R	7385	Tue : 11:30am	12th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0075	Gardenella Vaginalis - Vaginal Discharge/ cervical secretions**	Chennai - Avadi	R	1390	Daily : 9:00am to 9:00pm	7th day : 6:00pm	Collect vaginal discharge on a sterile cotton swab and place in into Amies transport medium with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0078	Councelling Genetics	Chennai - Avadi	E	1260	Daily : 9:00am	Next day : 5:00pm	aseptic precautions. 2. Make vaginal discharge smears on clean glass slide for Gram staining. Prior Appointment
G0079	GIST PDGFRA gene mutation**	Chennai - Avadi	R	10700	Thu & Sat: 11:30am	15th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
G0080	Glucose Tolerance Test 2 hours 2 samples	Chennai - Avadi	B	150	Daily : 9:00am to 9:00pm	After 6hrs	0 min, 120 min -2ml Fluoride plasma
G0082	GLIOMA PANEL (1P/19Q, IDH1/2,MGMT)**	Chennai - Avadi	D	29500	Daily : 7:30am	25th day : 6:00pm	(FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical Diagnosis with Histopathology report, HistologyandIHC.pdf
G0083	Gene analysis Alpha-1 antitrypsin deficiency - SERPINA1**	Chennai - Avadi	R	22470	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0087	Gene analysis Achondroplasia FGFR3 G380R common mutation screening**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0088	Agammaglobulinemia (BTK)-Gene analysis, EDTA Blood	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability. Amniotic Fluid or Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition.
G0089	Gene analysis Alexander disease (GFAP)**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetic history is mandatory.
G0091	Alstrom syndrome (ALMS1) gene analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetic history is mandatory.
G0092	Androgen receptor (AR)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0093	Gene analysis Ataxia-telangiectasia (ATM)**	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0094	Gene analysis ATRX**	Chennai - Avadi	E	19950	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml EDTA Whole Blood
G0095	Brown Vialletto-Van Laere syndrome-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0096	C8orf37 gene-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0097	Canavan disease (ASPA)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0098	Cantu syndrome (ABCC9) gene-Gene Analysis	Chennai - Avadi	R	23540	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0099	Carney complex (PRKAR1A)-Gene Analysis	Chennai - Avadi	R	18190	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0100	Chediak-Higashi syndrome (LYST)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0101	Cohen s syndrome (VPS13B)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0102	Congenital adrenal hyperplasia CYP21A2(21-0H)**	Chennai - Avadi	R	12840	Daily : 11:30am	After 6 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0103	Gene analysis Congenital adrenal hyperplasia CYP21A2 (21-0H) by NGS**	Chennai - Avadi	R	22470	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0104	Gene analysis Cystic fibrosis (CFTR)-Whole gene**	Chennai - Avadi	R	19260	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0105	Gene analysis Cystic megalencephaly (MLC1)**	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0106	Gene analysis DOCK8 (Hyper IgE Syndrome) gene**	Chennai - Avadi	R	17850	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0108	Dravet syndrome (SCN1A) gene panel,NGS**	Chennai - Avadi	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0109	Gene analysis duchenne muscular dystrophy (DMD)**	Chennai - Avadi	R	17000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0110	Fabry disease (GLA)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0111	Factor VII deficiency-Gene analysis, EDTA Blood	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.
G0112	Fanconi bickel syndrome (SLC2A2)-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0113	Gene analysis(Farber lipogranulomatosis - ASAH1 -**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0114	Fraser syndrome genes-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0115	Fucosidosis - Lysosomal storage disorder- Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0116	Galactosialidosis-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0117	GLB1 (GM1 gangliosidosis; mucopolysaccharidosis type IVB) gene analysis**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0118	Gene analysis GLUT1 deficiency (SLC2A1)**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0119	Haemophilia A-Gene analysis, EDTA Blood	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.
G0120	Gene analysis Haemophilia B (F9)**	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.
G0121	Hereditary fructose intolerance-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0122	Hyperargininemia (ARG1)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0123	Krabbe disease (GALC)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0124	Lysosomal acid lipase deficiency (LIPA gene)-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0126	McArdle disease (PYGM)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0127	Menkes disease (ATP7A)-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0128	Mondini defect (SLC26A4)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0129	Gene analysis MYD88 gene L265P mutation**	Chennai - Avadi	R	19260	Daily : 11:30am	12th Day : 6:00pm	3 ml EDTA Whole Blood
G0130	Neurofibromatosis type 1 (NF1) Gene Analysis **	Chennai - Avadi	R	18190	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0131	Gene analysis Neurofibromatosis type 2 (NF2)**	Chennai - Avadi	R	18190	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0132	NOTCH3 (CADASIL) gene analysis -NGS**	Chennai - Avadi	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0133	Pallister Hall syndrome (GLI3)-Gene Analysis	Chennai - Avadi	R	18190	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0134	Phenylketonuria (PAH)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0135	PMP22 Gene analysis**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0136	Polycystic Liver Disease-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0137	Gene analysis Pompe Disease (GAA)**	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0138	PTEN-Gene Analysis	Chennai - Avadi	R	19260	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0139	Pulmonary alveolar microlithiasis (SLC34A2)-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0142	Gene analysis RB1*	Chennai - Avadi	R	18190	Daily : 11:30am	After 4 weeks : 6:00pm	ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0143	Gene analysis RBC membrane disorders and Enzymopathies**	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0144	Sialidosis-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0145	Sjogren-Larsson syndrome (ALDH3A2)-Gene Analysis	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0147	Gene analysis Smith-Lemli-Opitz syndrome (DHCR7)**	Chennai - Avadi	R	18190	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0149	Tay-Sachs disease (HEXA)-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0150	Thrombotic Thrombocytopenic Purpura (TTP) - ADAMTS13 gene analysis**	Chennai - Avadi	R	18900	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0151	TSC1-Gene Analysis	Chennai - Avadi	R	18190	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0152	Gene analysis TSC1 & TSC2**	Chennai - Avadi	R	18190	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0153	TSC2-Gene Analysis	Chennai - Avadi	R	18190	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0154	Tyrosinemia-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0155	Gene analysis Von Hippel-Lindau syndrome (VHL)**	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability.
G0156	Gene analysis Von Willebrand disease** (VWF)	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability.
G0157	WAS (Wiskott Aldrich Syndrome)-Gene analysis, EDTA Blood	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability.
G0158	Wilson disease (ATP7B) gene analysis,NGS**	Chennai - Avadi	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetic history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0159	Xanthinuria (XDH)-Gene Analysis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0162	4H syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0163	Aicardi-Goutieres syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0164	Alagille syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0165	Alport syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0166	Gene Panel Aplastic anaemia**	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0167	Arthrogryposis & congenital myasthenic syndrome Gene Panel**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0168	Bardet-Biedl syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0169	Gene Panel Bartter syndrome**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0170	Benign infantile epilepsy-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0171	Cardiac channelopathy-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0173	Gene Panel Charcot-Marie-Tooth and sensory neuropathies**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0174	Ciliopathy-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0175	Citrullinemia-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0176	Cockayne syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0177	Congenital afibrinogenemia-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0178	Congenital dyserythropoietic anemia gene panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	<p>condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.</p> <p>Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.</p>
G0180	Congenital hypopituitarism-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	<p>3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.</p>
G0181	Gene Panel Congenital neutropenia	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	<p>Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen</p>

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability. Amniotic Fluid or Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition.
G0182	Congenital stationary night blindness-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0183	Gene Panel Connective tissue disorder**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0184	Corneal dystrophy-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0185	Cornelia de Lange Syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0186	Craniosynostosis-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0187	Cutis-laxa-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0188	Diamond blackfan anemia-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0189	Gene Panel Dyskeratosis congenita**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0190	Gene Panel Early-onset juvenile parkinsonism	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0191	Ectodermal dysplasia-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0192	Gene Panel Ehler Danlos syndrome**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0193	Epidermolysis bullosa-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0194	Gene Panel Episodic ataxia**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0195	Familial female mental retardation/epilepsy-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0196	Gene panel Familial hemiplegic migraine**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0197	Gene Panel Fanconi anemia**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	Genetci history is mandatory. Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0198	Fatty acid oxidation disorders-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0199	Glycine encephalopathy-Gene Panel	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0200	Gene Panel Glycogen storage disorder**	Chennai	R	21400	Daily : 11:30am	After 4	3 ml of EDTA whole blood required

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				weeks : 6:00pm	to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0201	Glycosylation (CDG) disorders-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0202	Gene Panel Haemophilia (F8 & F9)Gene Panel, EDTA Blood	Chennai - Avadi	R	19950	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.
G0203	Gene Panel Hemophagocytic lymphohistiocytosis (HLH)**	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0204	Hereditary elliptocytosis-Gene Panel, EDTA Blood	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0205	Hereditary multiple exostoses Panel-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0206	Hereditary pancreatitis gene panel**	Chennai - Avadi	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0207	Hereditary Spherocytosis-Gene Panel, EDTA Blood	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0208	Gene Panel Homocystinuria**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0209	Gene Panel Hypercholesterolemia**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0210	Hyperekplexia-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0211	Hypomyelination syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0212	Icthyosis-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0213	Gene panel Joubert syndrome**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0214	Kallmann syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0215	Leber congenital amaurosis-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0216	Gene Panel Leigh syndrome & mitochondrial encephalopathy**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0217	Gene Panel Leukodystrophy**	Chennai - Avadi	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0218	Lissencephaly-Gene Panel, EDTA Blood	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0219	Maple syrup urine disease-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0220	Meckel Gruber syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0221	Gene Panel Methylmalonic aciduria (MMA)**	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0222	Microcephaly-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0223	Microphthalmia and anophthamia-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0224	Gene Panel Mucopolysaccharidosis**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0225	Myotonia congenita-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0226	Gene Panel Neuronal ceroid lipofuscinosis	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0227	Neuronal migration disorder-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0228	Niemann-Pick disease gene panel**	Chennai - Avadi	R	20330	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0229	Non-immune hydrops-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0230	Gene Panel Noonan syndrome**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0231	Oculocutaneous albinism-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0232	Optic atrophy-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0233	Organic acidemia-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0234	Gene Panel Osteogenesis imperfecta**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0235	Osteopetrosis-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0236	Pachydermoperiostosis & primary hypertrophic osteoarthropathy-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0237	Pantothenate kinase-associated neurodegeneration-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0238	Peroxisomal disorder-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0240	Gene Panel Polycystic kidney disease**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0241	Pontocerebellar hypoplasia-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0242	Gene Panel Primary hyperoxaluria**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0243	Gene Panel (Primary immunodeficiency)**	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability. Amniotic Fluid or Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition.
G0244	Progressive familial intrahepatic cholestasis-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetic history is mandatory.
G0245	Progressive myoclonic epilepsy-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetic history is mandatory.
G0246	Rasopathy-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0247	Retinal Degeneration-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0249	Seckel syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0250	Gene Panel Severe combined immunodeficiency (SCID)**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0251	Sideroblastic anaemia-Gene Panel, EDTA Blood	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition. Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0252	Gene Panel Skeletal dysplasia**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0253	Sotos syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0254	Gene Panel Urea cycle defects**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0255	Usher syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0256	Waardenburg syndrome-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0257	Xeroderma pigmentosum-Gene Panel	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Genetci history is mandatory.
G0323	BTP-Beta Trace Protein GFR*	Chennai - Avadi	D	3895	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0324	Gastrointestinal panel by FilmArray (Biofire) - Stool	Chennai - Avadi	E	38325	Daily : 9:00am to 9:00pm	After 10hrs	5gm of Stool specimen in Sterile leak proof container
G0333	Klippel-Feil syndrome-Gene panel	Chennai - Avadi	E	21000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0334	Nephrotic syndrome-Gene Panel	Chennai - Avadi	R	20000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0335	Genotypic phenotypic correlation / case review by International experts	Chennai - Avadi	R	5350	Daily : 11:30am	Next day : 6:00pm	Prior Appointment
G0336	Glioma Panel-1	Chennai - Avadi	D	21830	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
G0337	Gene panel ALL risk stratification - T-ALL**	Chennai - Avadi	E	17850	Daily : 11:30am	20th day : 6:00pm	4 ml EDTA Whole blood/bone marrow aspirate
G0338	Neurodegeneration with brain iron	Chennai	E	21000	Daily : 7:30am	After 4	3 ml of EDTA whole blood required

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	accumulation-Gene Panel	- Avadi				weeks : 6:00pm	to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0339	FKTN-Gene Sequencing	Chennai - Avadi	E	19950	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0340	Gaucher disease (GBA)-Gene Analysis	Chennai - Avadi	E	19950	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0341	FISH-MDM2 (12q15) GENE AMPLIFICATION**	Chennai - Avadi	E	10185	Daily : 11:30am	15th day : 6:00pm	Formalin fixed paraffin embedded tissue block. Ship at room temperature. Tissue block must contain at least 10% tumor tissue. Histopath report with History mandatory.
G0344	Clinical exome & whole mitochondrial genome sequencing	Chennai - Avadi	E	26000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0345	Giant axonal neuropathy-1 (GAN)-Gene analysis	Chennai - Avadi	E	19950	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0348	Bone marrow failure syndrome gene panel**	Chennai - Avadi	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	4 ml EDTA Whole Blood
G0349	GI mapping biopsies (Upto 5 containers)	Chennai - Avadi	C	5085	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
G1901	CANINE-Gamma GT (GGTP)	Chennai - Avadi	A	240	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
G1903	Canine - Glucose fasting	Chennai - Avadi	A	100	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
G8175	Germline HRR gene panel	Chennai - Avadi	R	24610	Daily 7.30 am	25TH Days	3 ml of EDTA whole blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0001	H1N1 (Swine Flu) RNA Detection by Real Time PCR - Swab	Chennai - Avadi	C	2825	Daily : 7:00am & 3:00pm	Same Day : 3:00pm & 11:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
H0002	Hemochromatosis (HFE) Mutations Analysis (C282Y, H63D)	Chennai - Avadi	D	7670	Tue & Thu: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical history and TRF mandatory.
H0003	Haemoglobin-Free biochemical Plasma	Chennai - Avadi	R	1070	Daily : 11:30am	5th day : 6:00pm	2 ml of Citrated Plasma
H0004	Haemoglobin, free, urine	Chennai - Avadi	C	395	Daily : 11:30am	Same day : 9:00pm	10 ml of Spot urine
H0005	Haemophilia Profile	Chennai - Avadi	D	3895	Mon : 9:00am	Next day : 1:00pm	1 ml of citrated platelet poor plasma in frozen condition.
H0007	Haemophilus Influenza B Antigen Antigen detection by Agglutination - CSF	Chennai - Avadi	D	1415	Daily : 9:00am & 3:00pm	After 4hrs	2ml of CSF in Sterile leak proof tube or container
H0008	Haemophilus Influenza B Antigen Antigen detection by Agglutination - Serum	Chennai - Avadi	D	1415	Daily : 9:00am & 3:00pm	After 4hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0009	Haemophilus Influenza B Antigen Antigen detection by Agglutination - Urine	Chennai - Avadi	D	1415	Daily : 9:00am & 3:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
H0011	Hanging Drop preparation for V Cholerae by Microscopy - Stool	Chennai - Avadi	B	325	Daily : 9:00am & 3:00pm	After 6hrs	5gm of Stool specimen in Sterile leak proof container
H0013	Haptoglobin genotype by PCR**	Chennai - Avadi	R	3960	Wed & Sat : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
H0014	Haptoglobin Serum	Chennai - Avadi	C	2260	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0015	HAV-IgG Antibody by CLIA - Serum	Chennai - Avadi	B	1190	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0016	HAV-IgM Antibody by CLIA - Serum	Chennai	B	1190	Daily : 9:00am to	After 6hrs	3ml of Serum in Gel Vacutainer or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0017	HAV-RNA Detection by PCR - EDTA Blood	- Avadi Chennai - Avadi	D	5545	9:00pm Daily : 7:30am	7th day : 6:00pm	in Plain Sterile Vacutainer 3 ml of Whole Blood in EDTA Vacutainer
H0018	HbA1C- Glycated Haemoglobin, by HPLC method	Chennai - Avadi	B	400	Daily : 9:00am to 9:00pm	After 4hrs	3 ml EDTA whole blood
H0026	Anti HBcAg total by CLIA - Serum	Chennai - Avadi	B	970	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0027	Anti HBcAg-IgM by CLIA - Serum	Chennai - Avadi	B	970	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0028	Anti HBeAg by CLIA - Serum	Chennai - Avadi	B	970	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0029	Anti HbsAg Antibody by CLIA - Serum	Chennai - Avadi	B	970	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0030	HBV DNA Detection by Taqman RT PCR - Serum or Plasma	Chennai - Avadi	C	4465	Daily : 7:30am	Same day : 11:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0031	HBV DNA Detection by PCR - Amniotic fluid**	Chennai - Avadi	D	4660	Daily : 11:30am	7th day : 6:00pm	10 ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
H0032	HBV DNA Detection by PCR - Ascitic fluid**	Chennai - Avadi	D	4660	Daily : 11:30am	7th day : 6:00pm	3 ml of Fluid in Sterile leak proof container
H0033	HBV DNA Detection by PCR - CSF**	Chennai - Avadi	D	4660	Daily : 11:30am	7th day : 6:00pm	3 ml of Body Fluid in Sterile leak proof container
H0034	HBV DNA Detection by PCR - Liver tissue**	Chennai - Avadi	D	4660	Daily : 11:30am	7th day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
H0035	HBV DNA Genotype & Drug Resistance by PCR & Sanger Sequencing - Serum or Plasma	Chennai - Avadi	D	8850	Mon : 9:00am	Wed : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0037	HBeAg by CLIA - Serum	Chennai - Avadi	B	970	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0038	HBsAg-Confirmation by CLIA - Serum	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0039	HBsAg-Quantitative by CLIA - Serum	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0040	HBsAg-Screening by CLIA - Serum	Chennai - Avadi	B	510	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0041	HBV DNA Viral load by Taqman RT PCR - Serum or Plasma	Chennai - Avadi	C	7000	Daily : 7:30am	Same day : 11:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0042	HBV YMDD mutation by PCR & Sanger Sequencing - Serum	Chennai - Avadi	D	5310	Mon : 9:00am	Wed : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0043	HCG Beta Subunit, Serum	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0044	HCG Beta Subunit Urine 24H	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 8hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
H0045	HCG Beta Subunit Urine Spot	Chennai - Avadi	B	755	Daily : 9:00am to 9:00pm	After 8hrs	10 ml of Spot urine
H0045_U	PREGNANCY test by Immunochromatography - Urine	Chennai - Avadi	B	160	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
H0046	HCG Beta Subunit, CSF	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
H0047	HCG Stimulation Test for DHT	Chennai - Avadi	D	4250	Mon, Wed & Fri : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0048	HCG stimulation test for FSH and LH	Chennai - Avadi	B	1730	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0049	HCG Stimulation test for Testosterone	Chennai - Avadi	B	1080	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0050	HCV Genofibro (Fibromax by CLIA , Genotyping and IL28B by PCR & Sequencing - Serum	Chennai - Avadi	C	27120	Daily : 9:00am	10th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0051	HCV RNA Combo 1 (SGPT, Anti HCV by EIA, HCV-RNA Detection & Quantification by RT-PCR, HCV-Genotyping by Sanger Sequencing) - Serum	Chennai - Avadi	C	17515	Daily : 9:00am	4th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0052	HCV Genotyping by PCR & Sanger Sequencing - Serum or Plasma	Chennai - Avadi	D	7905	Tue & Fri : 9:00am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0053	HCV-IgM Antibody by EIA - Serum	Chennai - Avadi	C	1640	Mon & Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0054	HCV RNA Detection by Taqman RT PCR - Serum or Plasma	Chennai - Avadi	C	4405	Daily : 7:30am	Same day : 11:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0055	HCV-total Antibody by CLIA - Serum	Chennai - Avadi	B	1180	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0056	HCV RNA Viral load by Taqman RT PCR - Serum or Plasma	Chennai - Avadi	C	7990	Daily : 7:30am	Same day : 11:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0057	HDL Cholesterol Direct serum	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0058	HDV-IgM Antibody by EIA - Serum	Chennai - Avadi	D	5665	Daily : 7:30am	7th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0059	HDV-RNA Detection by PCR - EDTA Plasma**	Chennai - Avadi	D	6490	Daily : 7:30am	10th day : 6:00pm	3 ml of EDTA Plasma in Sterile Plain Vacutainer
H0060	HDV-total Antibody by EIA - Serum	Chennai - Avadi	C	2145	Mon & Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0061	HE4 with ROMA index	Chennai - Avadi	C	2825	Daily:7:30 am	Same day 8:00 PM	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0062	Heinz BodiesBlood	Chennai - Avadi	B	215	Daily : 11:30am	After 8hrs	3 ml EDTA Whole Blood
H0063	Helicobacter Pylori Antigen detection by CLIA - Stool	Chennai - Avadi	C	1580	Daily : 9:00am & 3:00pm	After 8hrs	5gm of Stool specimen in Sterile leak proof container
H0064	Helicobacter Pylori Culture by Vitek/ MALDITOF - Tissue	Chennai - Avadi	C	1245	Daily : 9:00am to 9:00pm	15th day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
H0065	Helicobacter Pylori IgA antibodies by EIA - Serum	Chennai - Avadi	C	2090	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0066	Helicobacter Pylori IgG antibodies by CLIA - Serum	Chennai - Avadi	C	2090	Daily : 9:00am	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0067	Helicobacter Pylori IgM antibodies by EIA - Serum	Chennai - Avadi	C	2090	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0069	Hemoglobin	Chennai - Avadi	A	80	Daily : 9:00am to 9:00pm	After 4hrs	3 ml EDTA Whole Blood
H0070	Hemosiderin, LBC	Chennai - Avadi	C	1695	Daily : 9:00am	3rd day : 6:00pm	More than 1 ML fluid in a leak proof plain container along with relvent clinical & radiological details
H0071	Heparan Sulphate Sulphamidase Sanfillipo / MPS III A Blood**	Chennai - Avadi	R	10700	Daily : 11:30am	After 4 weeks : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
H0072	HBV Immunity Profile by CLIA - Serum	Chennai - Avadi	C	2260	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0073	HBV Profile by CLIA - Serum	Chennai - Avadi	C	4520	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0074	Hepatitis Marker Profile - Acute by CLIA - Serum	Chennai - Avadi	C	3560	Mon, Wed & Fri : 7:30am	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0075	Hepatitis Profile, Comprehensive by CLIA - Serum	Chennai - Avadi	C	7230	Mon & Thu : 9:00am	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0076	Hepatitis - Atypical Virus Profile by CLIA - Serum	Chennai - Avadi	D	4720	Daily : 7:30am	6th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0077	Hepcidin Serum	Chennai - Avadi	D	2125	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0078	Her-2/Neu(C-Er-B2)(ERBB2)by FISH	Chennai - Avadi	D	11800	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
H0079	HEV IgG Antibody by EIA - Serum	Chennai - Avadi	B	1405	Mon & Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0080	HEV IgM Antibody by EIA - Serum	Chennai - Avadi	B	1405	Daily : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0081	HGV-RNA Detection by PCR - EDTA plasma**	Chennai - Avadi	R	5295	Daily : 11:30am	15th day : 6:00pm	3 ml of EDTA Plasma in Sterile Plain Vacutainer
H0082	Hirsutism Profile-1	Chennai - Avadi	C	3505	Tue & Fri: 7:30am	Same day : 7:00pm	5 ml of serum
H0083	Hirsutism Profile-2	Chennai - Avadi	C	5765	Tue & Fri: 7:30am	Same day : 7:00pm	5 ml of serum
H0084	Hirsutism Profile-3	Chennai - Avadi	C	7910	Tue & Fri: 7:30am	Same day : 7:00pm	5 ml of serum
H0085	Histamine EIA Plasma	Chennai - Avadi	R	6420	Daily : 11:30am	20th day : 6:00pm	3 ml of EDTA plasma
H0086	Histidine Quantitative Serum	Chennai - Avadi	R	6365	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
H0087	Histidine Quantitative Urine 24H**	Chennai - Avadi	R	6365	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
H0088	Histochemistry (special stain with reporting) Alcian Blue Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0089	Histochemistry (special stain with reporting) Bleach Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0090	Histochemistry (special stain with	Chennai	C	1130	Daily : 9:00am	3rd day :	Paraffin blocks (Formalin fixed

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	reporting) Congo Red Tissue/FFPE	- Avadi				6:00pm	paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0091	Histochemistry (special stain with reporting) Fite Farraco Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0092	Histochemistry (special stain with reporting) Fontana Masson Silver Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0093	Histochemistry (special stain with reporting) Giemsa Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0094	Histochemistry (special stain with reporting) Gomori Silver Methanamine (GMS) Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0095	Histochemistry (special stain with reporting) Hale S Colloidal Iron Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0096	Histochemistry (special stain with reporting) Histochemistry Unspecified Special stain Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0097	Histochemistry (special stain with reporting) Masson Trichrome Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0098	Histochemistry (special stain with reporting) Mucicarmine Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0099	Histochemistry (special stain with reporting) Orcein Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0100	Histochemistry (special stain with reporting) PAS with and/or Without Diastase Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0101	Histochemistry (special stain with reporting) PASM Stain Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0102	Histochemistry (special stain with reporting) Prussian Blue Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0103	Histochemistry (special stain with reporting) PTAH Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0104	Histochemistry (special stain with reporting) Reticulin Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0105	Histochemistry (special stain with reporting) Van Gieson Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							detail
H0106	Histochemistry (special stain with reporting) Verhoeff S Vangieson Tissue/ FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0107	Histochemistry (special stain with reporting) Von Kossa Stain Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0108	Histochemistry (special stain with reporting) ZN Tissue/FFPE	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0109	Histochemistry (Special Stain Only) Alcian Blue Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0110	Histochemistry (Special Stain Only) Bleach Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0111	Histochemistry (Special Stain Only) Congo Red Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0112	Histochemistry (Special Stain Only) Fite Farraco Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0113	Histochemistry (Special Stain Only) Fontana Masson Silver Method Tissue/ FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0114	Histochemistry (Special Stain Only) Giemsa Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0115	Histochemistry (Special Stain Only) Gomori Silver Methanamine Stain (GMS) Tissue/ FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0116	Histochemistry (Special Stain Only) Hale S Colloidal Iron Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0117	Histochemistry (Special Stain Only) Masson Trichrome Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0118	Histochemistry (Special Stain Only)	Chennai	C	565	Daily : 9:00am	3rd day :	Single Paraffin block / section on

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Mucicarmine Stain Tissue/FFPE	- Avadi				6:00pm	charged slide
H0120	Histochemistry (Special Stain Only) PAS With and/or Without Diastase Stain Tissue/ FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0121	Histochemistry (Special Stain Only) PASM Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0122	Histochemistry (Special Stain Only) Prussian Blue Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0123	Histochemistry (Special Stain Only) PTAH Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0124	Histochemistry (Special Stain Only) Reticulin Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0125	Histochemistry (Special Stain Only) Unspecified Special Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0126	Histochemistry (Special Stain Only) Van Gieson Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0127	Histochemistry (Special Stain Only) Verhoeff S Vangieson Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0128	Histochemistry (Special Stain Only) Von Kossa Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0129	Histochemistry (Special Stain Only) ZN Stain Tissue/FFPE	Chennai - Avadi	C	565	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
H0130	Histology Frozen Section (Intra Operative) day frozen 8am to 8pm, working days	Chennai - Avadi	D	8850	Daily : 8:00am to 8:00pm	Within 1hr	Fresh tissue to be sent along with relevant clinical & radiological details
H0131	Second opinion	Chennai - Avadi	D	1415	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded)/ Stained unstained slides Please attach original Histopathology report and all relevant clinical and radiological detail
H0132	Histone Antibody Serum	Chennai - Avadi	C	1865	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0133	Histopathological examination of-Amputation of Limb	Chennai - Avadi	C	5085	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0135	Histopathological examination of-Appendix	Chennai - Avadi	C	850	Daily : 9:00am	3rd day : 6:00pm	clinical & radiological details Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0136	Histopathological examination of-Bone Marrow Trepine biopsy	Chennai - Avadi	C	4070	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0137	Histopathological examination of-Bone Specimen	Chennai - Avadi	C	1580	Daily : 9:00am	7th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0138	Histopathological examination of-Brain Tissue	Chennai - Avadi	C	4180	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0139	Histopathological examination of-Breast Lump	Chennai - Avadi	C	2260	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0141	Histopathological examination of- Breast MRM / Breast conservative surgery reflex FISH	Chennai - Avadi	D	7435	Daily : 9:00am	7th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0142	Histopathological examination of-Buccal Mucosal Biopsy	Chennai - Avadi	C	1015	Daily : 9:00am	3rd day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0143	Histopathological examination of- Cervix Biopsy	Chennai - Avadi	C	1015	Daily : 9:00am	3rd day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0144	Histopathological examination of- Cervix Cone Biopsy	Chennai - Avadi	C	4180	Daily : 9:00am	6th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0145	Histopathological examination of-	Chennai	C	5650	Daily : 9:00am	7th day :	Tissue to be sent in adequate 10%

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Colectomy / APR	- Avadi				6:00pm	Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0147	Histopathological examination of-Commando With RND/MND	Chennai - Avadi	C	3390	Daily : 9:00am	7th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0148	Histopathological examination of - Endobronchial biopsy	Chennai - Avadi	C	1580	Daily : 9:00am	3rd day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0149	Histopathological examination of-Endometrium (D&C Material)	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0150	Histopathological examination of- Eye ball Specimen	Chennai - Avadi	C	4070	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0151	Histopathological examination of - Fistulectomy / Sinus tract includes special stains	Chennai - Avadi	C	960	Daily : 9:00am	3rd day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0152	Histopathological examination of- Gall Bladder	Chennai - Avadi	C	850	Daily : 9:00am	3rd day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0153	Histopathological examination of- Gastrectomy	Chennai - Avadi	C	4970	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0154	Histopathological examination of- GI Endoscopic Biopsy / Colonic Biopsy	Chennai - Avadi	C	2035	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0157	Histopathological examination of-Hepatic Resection	Chennai - Avadi	C	4970	Daily : 9:00am	6th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0158	Histopathological examination of-Hirschsprungs Disease	Chennai - Avadi	C	4180	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0159	Histopathological examination of -Kidney Biopsy	Chennai - Avadi	C	1750	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0161	Histopathological examination of-Liver Biopsy	Chennai - Avadi	C	3615	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0162	Histopathological examination of-Lung (Lobectomy)	Chennai - Avadi	C	4970	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0163	Histopathological examination of-Lymph Node	Chennai - Avadi	C	1695	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0164	Histopathological examination of-Mediastinal Mass	Chennai - Avadi	C	5085	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0165	Histopathological examination of-Muscle Biopsy	Chennai - Avadi	C	2035	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0166	Histopathological examination of-Nasal Polyp	Chennai - Avadi	C	905	Daily : 9:00am	3rd day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0167	Histopathological examination of-Nephrectomy (Non Tumorous)	Chennai - Avadi	C	2600	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0168	Histopathological examination of-Oesophagectomy	Chennai - Avadi	C	4180	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							clinical & radiological details
H0169	Histopathological examination of-Orchidectomy	Chennai - Avadi	C	4180	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0170	Histopathological examination of-Ovarian Mass (Large)	Chennai - Avadi	C	4970	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0171	Histopathological examination of-Ovarian Mass (Small To Medium)	Chennai - Avadi	C	2260	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0172	Histopathological examination of-Penis	Chennai - Avadi	C	3110	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0173	Histopathological examination of-Prostate Chips (Tur-P)	Chennai - Avadi	C	1920	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0174	Histopathological examination of-Prostate Trucut Needle Biopsy	Chennai - Avadi	C	3955	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0175	Histopathological examination of-Radical Cystectomy	Chennai - Avadi	C	4180	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0177	Histopathological examination of-Radical Nephrectomy	Chennai - Avadi	C	5085	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0178	Histopathological examination of-Radical Prostatectomy	Chennai - Avadi	C	5085	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0179	Histopathological examination of-	Chennai	C	850	Daily : 9:00am	3rd day :	Tissue to be sent in adequate 10%

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Sebaceous cyst	- Avadi				6:00pm	Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0180	Histopathological examination of-Skin Biopsy	Chennai - Avadi	C	1695	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0182	Histopathological examination of-Small Intestinal Resection	Chennai - Avadi	C	4180	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0183	Histopathological examination of-Soft Tissue Tumor	Chennai - Avadi	C	4970	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0184	Histopathological examination of-Spleen	Chennai - Avadi	C	2485	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0185	Histopathological examination of-Synovial Biopsy	Chennai - Avadi	C	1130	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0186	Histopathological examination of-Thyroidectomy	Chennai - Avadi	C	2035	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0187	Histopathological examination of- Trucut Biopsy (CT/USG Guided) Special stains & Reflex IHC	Chennai - Avadi	C	4070	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0188	Histopathological examination of-TUR Urinary Bladder	Chennai - Avadi	C	2710	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0189	Histopathological examination of-Unspecified specimen-Large	Chennai - Avadi	C	3275	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0190	Histopathological examination of-Unspecified specimen-Large & Complex	Chennai - Avadi	C	5650	Daily : 9:00am	6th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0191	Histopathological examination of-Unspecified specimen-Medium	Chennai - Avadi	C	1695	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0192	Histopathological examination of-Unspecified specimen-Small	Chennai - Avadi	C	1000	Daily : 9:00am	3rd day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0193	Histopathological examination of-Uterus With/Without Adnexae	Chennai - Avadi	C	2035	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0194	Histopathological examination of-Werthiems Hysterectomy	Chennai - Avadi	C	5085	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0195	Histopathological examination of-Whipples Resection	Chennai - Avadi	C	5085	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0198	HIT Antibody Heparin Induced Thrombocytopenia / PF4 Serum**	Chennai - Avadi	R	7500	Daily : 11:30am	4th day : 6:00pm	2ml serum required
H0201	HIV-1 Proviral DNA Detection by RT PCR - EDTA Blood	Chennai - Avadi	C	4070	Daily : 7:30am	Same day : 11:00pm	3 ml of Whole Blood in EDTA Vacutainer
H0202	HIV-1 RNA Detection by Taqman RT PCR - EDTA Plasma	Chennai - Avadi	C	3390	Daily : 7:30am	Same day : 11:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0205	HIV I & II Western Blot Test	Chennai - Avadi	C	2825	Daily : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0206	HIV-1 RNA Combo 1 - EDTA Blood & Plasma	Chennai - Avadi	C	7800	Daily : 7:30am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer & 3ml of EDTA Plasma in Sterile Plain Vacutainer
H0207	HIV1/2 - Human Immunodeficiency Virus - confirmatory serology by immuno-transfer technique - serum	Chennai - Avadi	R	11235	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0208	HIV-2 RNA Detectionby Real Time PCR - EDTA Plasma	Chennai - Avadi	C	4295	Daily : 7:30am	3rd day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0209	HIV-2 RNA Quantitative Viral load by Real Time PCR - EDTA Plasma	Chennai - Avadi	C	6780	Daily : 7:30am	3rd day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0210	HIV-2 RNA Combo by Flow Cytometry & Real Time PCR - EDTA Blood & Plasma	Chennai - Avadi	C	8475	Daily : 7:30am	3rd day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer & 3ml of EDTA Plasma in Sterile Plain Vacutainer
H0211	HIV-DUO by CLIA - Serum	Chennai - Avadi	B	590	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0212	HIV Early Detection Profile-1 - EDTA Plasma & Serum	Chennai - Avadi	C	3955	Daily : 7:30am	3rd day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer & 3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0213	HIV Early Detection Profile-2 - EDTA Blood & Serum	Chennai - Avadi	C	3275	Daily : 7:30am	3rd day : 6:00pm	3 ml of Serum in Gel Top Vacutainer or Plain Sterile Vacutainer & 3 ml of Whole Blood in EDTA Vacutainer
H0214	Donor Specific Antibody (DSA)	Chennai - Avadi	D	10030	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum + 6 ml Heparin/ACD/EDTA whole blood, Donor: 20 ml Heparin/ACD/EDTA Blood
H0215	HLA Typing A Locus	Chennai - Avadi	D	10500	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0216	HLA typing A,B,C,DR,DQ loci (ClassI & Class II)	Chennai - Avadi	D	21240	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0217	HLA Typing A,B,C Locus (ClassI)	Chennai - Avadi	D	14160	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0218	HLA Typing B Locus	Chennai - Avadi	D	10030	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0221	HLA Typing B5 Locus	Chennai - Avadi	D	10030	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0222	HLA Typing B51 Locus	Chennai - Avadi	D	10030	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0224	HLA Typing C Locus	Chennai - Avadi	D	10500	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0225	HLA Typing DQ Locus	Chennai - Avadi	D	6490	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0226	HLA Typing DR Locus	Chennai - Avadi	D	6490	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0227	HLA Typing DR & DQ Loci (ClassII)	Chennai - Avadi	D	12390	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0233	Single Antigen Bead Assay (Class I & Class II)*	Chennai - Avadi	E	42000	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum
H0234	Auto Antibody Detection Test*	Chennai - Avadi	D	1535	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum + 6 ml Heparin/ACD/EDTA Blood
H0235	HLA (Lymphocyte) Cross Match Total (T+B) cells*	Chennai - Avadi	D	2570	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum + 6 ml Heparin/ACD/EDTA Blood, Donor: 6 ml Heparin/ACD/EDTA Blood
H0249	HLA-B*5701 Associated HCP5 rs2395029 Genotyping for Abacavir Hypersensitivity	Chennai - Avadi	D	3540	Tue & Thu: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
H0250	HLA B27 (Human Leucocyte Antigen B27)	Chennai - Avadi	C	2450	Daily : 9:00am	Next day : 5:00pm	3 ml of EDTA whole blood
H0252	HLA-B*27 PCR	Chennai - Avadi	C	3505	Mon, Wed & Fri : 9:00am	Next day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
H0253	HOMA Index Insulin Resistance Test	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	Next day : 5:00pm	3 ml of serum and 2 ml of fluoride plasma. Fasting sample.
H0254	Homocysteine Serum	Chennai - Avadi	C	1300	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0256	Homogentisic Acid (ALKAPTANURIA)	Chennai - Avadi	C	850	Daily : 11:30am	Same day : 9:00pm	10 ml of Spot urine
H0258	HPV-DNA Detectionby PCR and Typing by PCR & Sanger Sequencing - Vaginal swab	Chennai - Avadi	C	2655	Mon & Thu : 9:00am	3rd Day : 6:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
H0259	HPV-DNA Detectionby PCR and Typing PCR & Sanger Sequencing - Cervical scrapping	Chennai - Avadi	C	2655	Mon & Thu : 9:00am	3rd Day : 6:00pm	Thin Prep Container
H0260	HPV-DNA Detectionby PCR and Typing by PCR & Sanger Sequencing - Tissue	Chennai - Avadi	C	3390	Mon & Thurs : 9:00am	3rd day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
H0261	HPV-DNA Detectionby PCR and Typing by PCR & Sanger Sequencing - LBC	Chennai - Avadi	C	2655	Mon & Thu : 9:00am	3rd Day : 6:00pm	Thin Prep Container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0264	HRT Profile	Chennai - Avadi	C	1300	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0265	HsCRP -High Sensitivity CRPSerum	Chennai - Avadi	B	750	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0266	HSV-1 IgG and IgM Antibody by EIA - Serum	Chennai - Avadi	B	1240	Daily : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0268	HSV-1 IgG Antibody by EIA - Serum	Chennai - Avadi	B	650	Daily : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0269	HSV-1 IgM Antibody by EIA - Serum	Chennai - Avadi	B	650	Daily : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0270	HSV 1&2 DNA Detection by Real Time PCR - Blister fluid	Chennai - Avadi	C	4520	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3ml of Fluid in Sterile leak proof container
H0271	HSV 1&2 DNA Detection by Real Time PCR - CSF	Chennai - Avadi	C	4520	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CSF in Sterile leak proof tube or container
H0272	HSV 1&2 DNA Detection by Real Time PCR - Serum	Chennai - Avadi	C	4520	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0273	HSV 1&2 DNA Detection by Real Time PCR - Swab	Chennai - Avadi	C	4520	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
H0274	HSV 1&2 IgG and IgM Antibody by EIA - Serum	Chennai - Avadi	B	1240	Daily : 7:30am	5th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0276	HSV 1+2 IgG Antibody by CLIA - Serum	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0277	HSV 1&2 IgG Antibody by western blot - Serum	Chennai - Avadi	D	2360	Daily: 11:30am	Next day 8pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0278	HSV 1+2 IgM Antibody by CLIA - Serum	Chennai - Avadi	B	650	Daily : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0279	HSV 1&2 IgM Antibody by western blot - Serum	Chennai - Avadi	D	2360	Daily: 11:30am	Next day 8pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0280	HSV-2 IgG and IgM Antibody by EIA - Serum	Chennai - Avadi	B	1240	Daily : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0282	HSV-2 IgG Antibody by EIA - Serum	Chennai - Avadi	B	650	Daily : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0283	HSV-2 IgM Antibody by EIA - Serum	Chennai	B	650	Daily : 7:30am	Same day :	3ml of Serum in Gel Vacutainer or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				5:00pm	in Plain Sterile Vacutainer
H0284	HTLV 1 And 2 Virus Total Antibody by CLIA - Serum	Chennai - Avadi	D	5190	Mon : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0285	Hu AntibodyNeuronal ag in Paraneoplastic syndromesSerum	Chennai - Avadi	D	5310	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0286	Huntington Disease Mutation Analysis	Chennai - Avadi	D	5665	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
H0287	Hydatid Cyst Detection by Microscopy - Fluid	Chennai - Avadi	C	790	Daily : 9:00am & 3:00pm	After 6hrs	3ml of Fluid in Sterile leak proof container
H0288	Hypertension Profile-Routine	Chennai - Avadi	C	1920	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of fluoride plasma, 10 ml of spot urine
H0289	Heart failure profile- maxi	Chennai - Avadi	D	5545	Mon : 7:30am	4th day : 6:00pm	5 ml of serum
H0295	Extra slide set Histopathology cases	Chennai - Avadi	C	285	Daily : 9:00am	4th day : 6:00pm	Stained slide set is issued
H0300	HBV mutation Profile by PCR & Sanger Sequencing - Serum**	Chennai - Avadi	R	10700	Daily : 11:30am	10th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0301	HBV precore mutations by PCR & Sanger Sequencing - Serum**	Chennai - Avadi	R	5350	Wed : 11:30am	10th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0302	HCV Reflex Profile - Serum	Chennai - Avadi	B	1510	Tue & Fri : 9:00am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0304	Histopathological examination of - Fallopian tube	Chennai - Avadi	C	1130	Daily : 9:00am	3rd day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0306	Histopathological examination of-Placenta	Chennai - Avadi	C	2260	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0307	Histopathological examination of-	Chennai	C	1695	Daily : 9:00am	4th day :	Tissue to be sent in adequate 10%

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Testicular biopsy	- Avadi				6:00pm	Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0308	Single Antigen Bead Assay (Class I)*	Chennai - Avadi	E	21000	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum
H0309	Single Antigen Bead Assay (Class II)*	Chennai - Avadi	E	21000	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum
H0310	Homocysteine reflex B12-folate Serum	Chennai - Avadi	C	1920	Daily : 9:00am to 9:00pm	Next day : 1:00pm	4ml of Serum
H0311	HIV-1 RNA Drug Resistance for integrase inhibitor by PCR & Sanger Sequencing - EDTA Plasma	Chennai - Avadi	D	5900	1st & 3rd Mon : 9:00am	15th day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0312	HIV-1 RNA Coreceptor Tropism (Genotype) by PCR & Sanger Sequencing - EDTA Plasma	Chennai - Avadi	D	7080	1st & 3rd Mon : 9:00am	15th day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0313	HBV Basal Core Promoter mutation by PCR & Sanger Sequencing - Serum**	Chennai - Avadi	R	5885	Wed : 11:30am	8th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0314	Homovanillic acid urine	Chennai - Avadi	R	5350	Daily : 11:30am	12th Day : 6:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 30 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
H0316	HIV-1 RNA Drug Resistance-Extended by PCR & Sanger Sequencing - EDTA Plasma	Chennai - Avadi	D	21240	1st & 3rd Mon : 9:00am	15th day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0317	HPFH (Hereditary persistence of fetal haemoglobin) 5 deletion**	Chennai - Avadi	R	12520	Daily : 11:30am	After 4 weeks :	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
						6:00pm	required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability. Amniotic Fluid or Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition.
H0318	HNPCC/ Lynch Syndrome Panel	Chennai - Avadi	R	18000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetic history is mandatory.
H0320	Breast Carcinoma Profile I	Chennai - Avadi	D	10620	Daily : 9:00am	10th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0321	Breast carcinoma Profile II	Chennai - Avadi	D	12980	Daily : 9:00am	10th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0322	Lung Carcinoma (NSCLC) Profile	Chennai - Avadi	D	13570	Daily : 9:00am	After 2 weeks : 6:00pm	(FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical Diagnosis with Histopathology report,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							HistologyandIHC.pdf
H0323	Colorectal Carcinoma Profile	Chennai - Avadi	D	7670	Daily : 9:00am	10th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0324	Head and neck cancer (HNSCC profile)	Chennai - Avadi	D	7670	Daily : 9:00am	5th day : 6:00pm	HistologyandIHC.pdf
H0326	Heavy Metals profile-1, Hair**	Chennai - Avadi	R	5885	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
H0327	Heavy Metals profile-2, Hair**	Chennai - Avadi	R	12305	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
H0328	Heavy Metals profile-3, Hair**	Chennai - Avadi	R	18725	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
H0329	Heavy metals for ayurvedic tablets test	Chennai - Avadi	R	5885	Daily : 11:30am	10th day : 6:00pm	Medicine sample (10 grm), Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clinical symptoms if any. Reported in 7 days
H0357	Histopathological examination of- GI Endoscopic Biopsy with reflex IHC	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0793	H3 Viral marker Profile by CLIA - Serum	Chennai - Avadi	C	1710	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0794	Heavy and Trace Metals profile-1	Chennai - Avadi	D	11800	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
H0795	Heavy and Trace Metals profile-2	Chennai - Avadi	D	7080	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
H0796	CELL BLOCK	Chennai - Avadi	C	4520	Daily : 9:00am	5th day : 6:00pm	HistologyandIHC.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0797	HIV-1 RNA Viral load by Taqman RT PCR - EDTA Plasma	Chennai - Avadi	C	5650	Daily : 7:30am	Same day : 11:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0798	HIV-1 RNA Quantitative Viral load by Taqman RT PCR - CSF	Chennai - Avadi	C	6755	Daily : 7:30am	Same day : 11:00pm	2ml of CSF in Sterile leak proof tube or container
H0799	Panel Reactive Antibody (PRA) Identification by Luminex* (ClassI)	Chennai - Avadi	D	10620	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum
H0800	Panel Reactive Antibody (PRA)Identification by Luminex* (ClassII)	Chennai - Avadi	D	10620	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum
H0801	Panel Reactive Antibody (PRA) Screening by Luminex (ClassI & Class II)	Chennai - Avadi	D	6490	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum
H0808	HEV RNA Detection by Real Time PCR - EDTA Blood	Chennai - Avadi	D	5900	Daily : 9:00am	5th day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer
H0810	Hereditary Cancer panel **	Chennai - Avadi	R	18000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
H0813	HLA typing A,B (ClassI) and DR loci (Class II)	Chennai - Avadi	D	15105	Daily : 9:00am	Same day : 6:00pm	6 ml EDTA/ACD Blood
H0824	Any specimen with reflex IHC	Chennai - Avadi	C	8475	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0825	Histopathological examination of- Nasopharyngeal Biopsy	Chennai - Avadi	D	5310	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0826	Histopathological examination of-Salivary Gland	Chennai - Avadi	C	4070	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0827	Histopathological examination of- Laryngeal Biopsy	Chennai - Avadi	C	2485	Daily : 9:00am	5th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0828	Histopathological examination of-Nasal Polyp with Reflex IHC	Chennai - Avadi	D	2595	Daily : 9:00am	5th day : 6:00pm	proof container along with relevant clinical & radiological details Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0829	Rush Biopsy service	Chennai - Avadi	D	3540	Daily : 9:00am	Next day : 6:00pm	NA
H0841	HPV-DNA Detection by Hybrid Capture reflex Typing by PCR & Sanger Sequencing - LBC or Cervical sampler	Chennai - Avadi	C	3380	Mon & Thu : 9:00am	3rd Day : 6:00pm	Thin Prep Container/ Digene Sampler
H0849	High Resolution HLA Typing (A, B, C, DR, DQ, DP) - Solid Organ Transplant	Chennai - Avadi	E	14980	Mon & Thur : 9.00 am	7 Days	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
H0850	High Resolution HLA Typing (A, B, C, DR, DQ, DP) - BMT/ Stem cell Transplant	Chennai - Avadi	E	14980	Mon & Thur : 9.00 am	7 Days	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
H0856	Sub speciality expert opinion	Chennai - Avadi	D	11565	Daily : 9:00am	7th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container Or paraffin blocks along with relevant clinical & radiological details
H0860	Histology - Frozen Section (Intra Operative) Night frozen 8 pm to 8 am, week off & public holidays	Chennai - Avadi	D	18880	Daily : 8:00pm to 8:00am	Within 1hr	Fresh tissue to be sent along with relevant clinical & radiological details
H0861	HLA-B* 1502 for carbamazepine toxicity**	Chennai - Avadi	D	10030	Daily : 7:30am	20th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
H0862	HLA - T and B Cell Flow Cytometry Cross Match (FCXM)	Chennai - Avadi	D	7100	Daily : 11:30am	3rd day : 6:00pm	Recipient: 5 ml serum, Donor: 15 ml Heparin
H0863	High Resolution HLA Typing B**	Chennai - Avadi	R	6315	Tue & Fri: 11:30am	10th day : 6:00pm	6 ml EDTA Blood
H0864	HCG as Tumour Marker, Serum	Chennai - Avadi	C	905	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0865	High Resolution HLA Typing DQB1**	Chennai - Avadi	R	6315	Tue & Fri: 11:30am	10th day : 6:00pm	6 ml EDTA Blood
H0866	Pre Transplant Workup panel	Chennai - Avadi	D	41300	Daily : 9:00am	Same day : 6:00pm	Recipient : 3 tubes EDTA/ ACD blood & 1 Gel tube, Donor-6 tubes EDTA/ACD Blood
H0867	PRA Identification Class I and Class II*	Chennai - Avadi	D	21240	Daily : 9:00am	Same day : 6:00pm	Recipient: 5 ml serum
H0868	Histopathological examination of- Eyelid Biopsy	Chennai - Avadi	C	1695	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0869	HBV DNA Combo RT-PCR & Sanger Sequencing - Serum or Plasma	Chennai - Avadi	C	15255	Mon : 9:00am	Wed : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0870	HCV RNA Combo 2 by RT-PCR & Sanger Sequencing - Serum or Plasma	Chennai - Avadi	C	10170	Tue & Fri : 9:00am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0871	HIV-1 RNA Combo 2 by PCR & Sanger Sequencing - EDTA Plasma	Chennai - Avadi	D	24780	1st & 3rd Mon : 9:00am	15th day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0872	Anti HMGCR antibodies - serum	Chennai - Avadi	R	10700	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml Serum
H0874	NextGen HRR Gene Panel	Chennai - Avadi	E	25300	Daily : 11:30am	After 4 weeks : 6:00pm	FFPE block with a minimum of 10% (as assessed by the pathologist) of tumor cells should be present

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0876	Hydroxychloroquine (HCQ)	Chennai - Avadi	D	1535	Tue & Fri: 11:30am	Next day : 5:00pm	3 ml EDTA Whole Blood
H0877	Histopathology for fungal infection	Chennai - Avadi	C	1695	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0879	Hemophagocytic lymphohistiocytosis (T-cell Activation HLH Panel)**	Chennai - Avadi	R	10165	Daily : 11:30am	7th day : 6:00pm	Sample required: EDTA Blood 3ml
H1030	HEREDITARY SPASTIC PARAPLEGIA GENE PANEL	Chennai - Avadi	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
H1901	CANINE-Histo-Large specimen (unspecified) More than 5 cm	Chennai - Avadi	C	1695	Daily : 9:00am	8th day : 6:00pm	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
H1902	Second opinion service, Sldies and blocks	Chennai - Avadi	C	1130	Daily : 9:00am	8th day : 6:00pm	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
H1903	CANINE-Histo-Medium specimen (unspecified) 2 to 5 cm	Chennai - Avadi	C	1130	Daily : 9:00am	8th day : 6:00pm	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							and to be interpreted by the clinician.
H5671	Hepatitis B Monitoring profile	Chennai - Avadi	D	11480	Daily : 7:30am	Next day : 5:00pm	3ml of serum, 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H5682	Tru Heart Basic H5682	Chennai - Avadi	W	1399	Daily : 9:00am	Next day : 6:00pm	4 ml Serum
H5692	Hrd track stand-alone test**	Chennai - Avadi	R	139100	Daily : 11:30am	After 4 weeks : 6:00pm	FFPE blocks with 20% tumor content and HP report.
H8335	[Combo] Haemophilia A (F8) gene analysis (along with analysis of the F8-intron 22 inversion**	Chennai - Avadi	R	25680	Daily : 11:30am	After 4 weeks : 6:00pm	4 ml EDTA Whole Blood
H8395	Hypersensitive pneumonitis panel-Extended	Chennai - Avadi	E	11025	Mon, Wed & Fri : 11:30am	Next day : 5:00pm	3ml serum sample required
I0001	IA2 Insulin Autoantibodies serum	Chennai - Avadi	D	2830	Mon & Thu 9 am	Next day 5 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0002	IBD Profile	Chennai - Avadi	C	5875	Fri : 7:30am	Next day : 5:00pm	3 ml of serum, 100 gms of stool in sterile container
I0004	Iduronate 2 Sulphate Sulphatase Hunter/ MPS II Blood	Chennai - Avadi	R	10595	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
I0005	IgA CSF	Chennai - Avadi	C	1470	Mon & Thu : 9:00am	Same day : 5:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
I0006	IgA Serum	Chennai - Avadi	C	790	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0007	IgD by Flow cytometry Blood*	Chennai - Avadi	D	3540	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
I0008	IGF BP3 Serum	Chennai - Avadi	C	3615	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0009	IGF-1 Somatomedin C Serum	Chennai - Avadi	C	3650	Daily : 9:00am to 9:00pm	After 6hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0010	IgG by Flow cytometry Blood**	Chennai - Avadi	R	2945	Daily : 11:30am	4th day : 6:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
I0011	IgG Total antibody CSF	Chennai - Avadi	C	1470	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
I0012	IgG Total antibody Serum	Chennai - Avadi	C	790	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0013	IgG4 Sub Class Serum	Chennai - Avadi	D	7200	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0014	IgM by Flow cytometry Blood**	Chennai - Avadi	R	2945	Daily : 11:30am	4th day : 6:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
I0015	IgM total CSF	Chennai - Avadi	C	1470	Mon & Thu : 9:00am	Same day : 5:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
I0016	IgM total Serum	Chennai - Avadi	C	790	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0017	IHC Panels with reporting- Bone Marrow	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0018	IHC Panels with reporting-Brain Tumor	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0019	IHC Panels with reporting Breast IV Comprehensive diagnosis	Chennai - Avadi	D	7670	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0020	IHC panel-Breast-I, tissue	Chennai - Avadi	C	2825	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0021	IHC Panels with reporting Breast-II ER, PR, c-erb-B2	Chennai - Avadi	C	4180	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0022	IHC Panels with reporting Breast-III ER, PR, c-erb-B2, Ki 67	Chennai - Avadi	C	5085	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0023	IHC Panels with reporting - Cervical Cancer	Chennai - Avadi	C	6555	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0025	IHC Panels with reporting - GastroIntestinal neoplasm	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0026	IHC Panels with reporting - Liver Neoplasm	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0027	IHC Panels with reporting - Lung Neoplasm	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0028	IHC Panels with reporting - Lymph node	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0032	IHC Panels with reporting Ovarian Neoplasm	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0033	IHC Panels with reporting Pancreatic Neoplasm	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0034	IHC Panels with reporting Peritoneal Adenocarcinoma Vs Mesothelioma	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0035	IHC Panels with reporting Pituitary Marker Profile - ACTH, LH, FSH, TSH, PROLACTIN	Chennai - Avadi	C	15255	Daily : 9:00am	12th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0037	IHC Panels with reporting - Round Cell Tumor(Adult/Pediatric)	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0038	IHC Panels with reporting - Carcinoma of unknown primary (CUP)	Chennai - Avadi	C	7910	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0039	IHC Panels with reporting - Salivary Gland Neoplasm	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	all relevant clinical and radiological detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0040	IHC Panels with reporting - Soft Tissue Tumor	Chennai - Avadi	C	9040	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0041	IHC Panels with reporting - Testicular Neoplasm	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0042	IHC Panels with reporting - Thyroid Neoplasm	Chennai - Avadi	C	7345	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0043	IHC Panels with reporting - 34 Beta E12 (HMWCK)	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0046	IHC single marker with reporting Alpha Feto Protein	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0048	IHC single marker with reporting Bcl-2 Oncoprotein	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0049	IHC single marker with reporting Bcl-6 Oncoprotein	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0050	IHC single marker with reporting Beta Catenin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	original Histopathology report and all relevant clinical and radiological detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0051	IHC single marker with reporting Beta HCG	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0053	IHC single marker with reporting C4D	Chennai - Avadi	D	2360	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0054	IHC single marker with reporting CA 125	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0055	IHC single marker with reporting CA-19.9 **	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0056	IHC single marker with reporting Calcitonin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0057	IHC single marker with reporting Calponin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0058	IHC single marker with reporting Carletinin	Chennai	D	2240	Daily : 9:00am	4th day :	Paraffin blocks (Formalin fixed

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				6:00pm	paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0061	IHC multiplex marker with reporting CD 20 + CD 3 Cocktail	Chennai - Avadi	D	3540	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0062	IHC single marker with reporting CD10	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0063	IHC single marker with reporting CD117	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0064	IHC single marker with reporting CD138	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0065	IHC single marker with reporting CD15	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0066	IHC single marker with reporting CD1a	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0067	IHC single marker with reporting CD20	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0068	IHC single marker with reporting CD21	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0070	IHC single marker with reporting CD23	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0071	IHC single marker with reporting CD3	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0072	IHC single marker with reporting CD30	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0073	IHC single marker with reporting CD31	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0074	IHC single marker with reporting CD34	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0075	IHC multiplex marker with reporting CD4 + CD8 Cocktail	Chennai - Avadi	D	3540	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0076	IHC single marker with reporting CD43	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							detail
I0077	IHC single marker with reporting CD45 LCA	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0078	IHC single marker with reporting CD5	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0079	IHC single marker with reporting CD56	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0080	IHC single marker with reporting CD57	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0081	IHC single marker with reporting CD68	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0082	IHC single marker with reporting CD79a	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0083	IHC single marker with reporting CD99 Mic2	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0084	IHC single marker with reporting CDX2	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							all relevant clinical and radiological detail
I0085	IHC single marker with reporting CEA	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0086	IHC-CERB2 (Her 2/neu)-IHC single marker, tissue	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0087	IHC single marker with reporting C-erb-B2 for Gastric carcinoma Her2 Neu	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0088	IHC single marker with reporting Chromogranin A	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0089	IHC single marker with reporting CK10	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0090	IHC single marker with reporting CK19	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0091	IHC single marker with reporting CK20	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0092	IHC single marker with reporting CK5/6	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							original Histopathology report and all relevant clinical and radiological detail
I0093	IHC single marker with reporting CK7	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0094	IHC single marker with reporting CK8	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0095	IHC single marker with reporting p16INK4a	Chennai - Avadi	D	5310	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0096	IHC single marker with reporting Cyclin D1	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0097	IHC single marker with reporting Desmin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0098	IHC multiplex marker with reporting Desmoglein 3 + Napsin A Cocktail	Chennai - Avadi	D	3540	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0099	IHC single marker with reporting DOG-1	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0100	IHC single marker with reporting E	Chennai	D	2240	Daily : 9:00am	4th day :	Paraffin blocks (Formalin fixed

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Cadherin	- Avadi				6:00pm	paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0101	IHC single marker with reporting EGFR	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0102	IHC single marker with reporting EMA Epithelial Membrane Antigen	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0103	IHC-Estrogen Receptor (ER)-single marker, tissue	Chennai - Avadi	D	1770	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0104	IHC single marker with reporting Factor VIII	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0106	IHC single marker with reporting GCDFP15	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0107	IHC single marker with reporting GFAP	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0108	IHC single marker with reporting Glypican-3	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0110	IHC single marker with reporting H.Pylori	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0113	IHC single marker with reporting Hep-Par 1	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0114	IHC single marker with reporting HMB45	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0115	IHC single marker with reporting HPV	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0119	IHC single marker with reporting Inhibin Alpha	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0120	IHC single marker with reporting Insulin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0121	IHC multiplex marker with reporting Kappa + Lambda Cocktail	Chennai - Avadi	D	3540	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0122	IHC single marker with reporting Ki 67	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0126	IHC single marker with reporting Lysoyme	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0127	IHC single marker with reporting Melan A	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0128	IHC single marker with reporting Mesothelin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0129	IHC single marker with reporting MOC 31	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0130	IHC single marker with reporting MPO Myeloperoxidase	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0131	IHC single marker with reporting Mum1 Protein	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0132	IHC single marker with reporting My F4 Myogenin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0133	IHC single marker with reporting Myoglobin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0134	IHC single marker with reporting NSE Neuron Specific Enolase	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	all relevant clinical and radiological detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0135	IHC single marker with reporting Oct.2	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0137	IHC single marker with reporting p53	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0138	IHC single marker with reporting p63	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0139	IHC single marker with reporting Pan CK AE1/AE3	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0140	IHC single marker with reporting Pax 5	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0141	IHC single marker with reporting Pax 8	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0142	IHC multiplex marker with reporting PIN 4 Cocktail Ck5 + Ck14 +P63 + P504S	Chennai - Avadi	D	3540	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0143	IHC single marker with reporting PLAP Placental Alkaline Phosphatase	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	original Histopathology report and all relevant clinical and radiological detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0144	IHC-Progesterone Receptor (PR)single marker,tissue	Chennai - Avadi	D	1770	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0146	IHC single marker with reporting Prostatic Acid Phosphatase	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0147	IHC single marker with reporting PSA Prostate Specific Antigen	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0148	IHC single marker with reporting S100 Protein	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0149	IHC single marker with reporting Smooth Muscle Actin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0150	IHC single marker with reporting Synaptophysin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0151	IHC single marker with reporting TdT	Chennai	D	2240	Daily : 9:00am	4th day :	Paraffin blocks (Formalin fixed

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				6:00pm	paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0152	IHC single marker with reporting Thyroglobulin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0154	IHC single marker with reporting TTF1	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0155	IHC single marker with reporting Vimentin	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0156	IHC single marker with reporting WT1 Wilm S Tumor Protein	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0157	IHC stain only 34 Beta E12	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0158	IHC stain only ALK1	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0159	IHC stain only Alpha Feto Protein	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0160	IHC stain only Bcl-2 Oncoprotein	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0161	IHC stain only Bcl-6 Oncoprotein	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0162	IHC stain only Beta Catenin	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0163	IHC stain only Beta HCG	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0164	IHC stain only CA125	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0165	IHC stain only Calcitonin	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0166	IHC stain only Calponin	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0167	IHC stain only Calretinin	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0170	IHC stain only CD10	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0171	IHC stain only CD117	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0172	IHC stain only CD138	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0173	IHC stain only CD15	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0174	IHC stain only CD1a	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0175	IHC stain only CD20	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0176	IHC stain only CD21	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0178	IHC stain only CD23	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0179	IHC stain only CD3	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0180	IHC stain only CD30	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0181	IHC stain only CD31	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0182	IHC stain only CD34	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0184	IHC stain only CD43	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0185	IHC stain only CD45 LCA	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0186	IHC stain only CD5	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0187	IHC stain only CD79a	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0188	IHC stain only CD99 Mic2	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0189	IHC stain only CDX2	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0190	IHC stain only CEA	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0191	IHC stain only c-erb-B2 Her 2 Neu	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0192	IHC stain only Chromogranin A	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0194	IHC stain only CK 19	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0195	IHC stain only CK 20	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0196	IHC stain only CK 5/6	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0197	IHC stain only CK 7	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0198	IHC stain only CK8	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0200	IHC stain only Cyclin D1	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0201	IHC stain only Desmin	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0202	IHC stain only DOG-1	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0203	IHC stain only E Cadherin	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0204	IHC stain only EGFR	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0205	IHC stain only EMA Epithelial Membrane Antigen	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0206	IHC stain only ER Estrogen Receptor	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0207	IHC stain only Factor VIII	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0208	IHC stain only GCDFP15	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0209	IHC stain only GFAP	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0210	IHC stain only Glypican-3	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0211	IHC stain only H.Pylori	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0212	IHC stain only Hep-Par 1	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0213	IHC stain only HMB45	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0215	IHC stain only Inhibin Alpha	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0216	IHC stain only Insulin	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0217	IHC stain only Ki 67	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0218	IHC stain only Melan A	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0219	IHC stain only Mesothelin	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0220	IHC stain only MOC 31	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0221	IHC stain only MPO Myeloperoxidase	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0222	IHC stain only Multiplex CD 20 + CD3 Cocktail	Chennai - Avadi	D	2240	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0223	IHC stain only Multiplex CD4 + CD8	Chennai - Avadi	D	2240	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0224	IHC stain only Multiplex Desmoglein 3 + Napsin A	Chennai - Avadi	D	2240	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0225	IHC stain only Multiplex Kappa + Lambda	Chennai - Avadi	D	2240	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0226	IHC stain only Multiplex Pin-4 Cocktail Ck5 + Ck14 + P63 + P504S	Chennai - Avadi	D	2240	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0227	IHC stain only Mum1 Protein	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0228	IHC stain only My F4 (Myogenin)	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0230	IHC stain only NSE Neuron Specific Enolase	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0231	IHC stain only OCT.2	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0233	IHC stain only p53	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0234	IHC stain only p63	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0235	IHC stain only Pan CK AE1/AE3	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0236	IHC stain only Pax 5	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0237	IHC stain only Pax 8	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0238	IHC stain only PLAP	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0239	IHC stain only PR Progesterone Receptor	Chennai - Avadi	D	885	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0240	IHC stain only Prostatic Acid Phosphatase	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0241	IHC stain only PSA Prostate Specific Antigen	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0242	IHC stain only S100 Protein	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0243	IHC stain only Single Marker (Any other)	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0244	IHC stain only SMA Smooth Muscle Actin	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0245	IHC stain only Synaptophysin	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0246	IHC stain only TdT	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0247	IHC stain only Thyroglobulin	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0248	IHC stain only TTF1	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0249	IHC stain only Vimentin	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0250	IHC stain only WT1 Wilm S Tumor Protein	Chennai - Avadi	D	1120	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0251	IHC Comprehensive panel for diagnosis	Chennai - Avadi	C	9605	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0258	IL28 B SNP Genotyping by PCR & Sanger Sequencing - EDTA Blood	Chennai - Avadi	D	5310	Daily : 9:00am	10th day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer
I0259	Immune Status Panel	Chennai - Avadi	C	6665	Daily : 9:00am	After 8hrs	5 ml of serum
I0260	Immunofixation Qualitative Serum	Chennai - Avadi	C	5875	Tue, Thu & Sat: 11:30am	Next day : 5:00pm	4 ml of Serum
I0261	Immunofixation Qualitative Urine Profile	Chennai - Avadi	C	6780	Thu : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea ,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0262	Immunofixation- Quantitative Serum	Chennai - Avadi	C	8025	Tue, Thu & Sat: 11:30am	Next day : 5:00pm	coffee , cold drink avoidance for accurate results and such history to be conveyed to lab. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0263	Immunofixation-Quantitative Urine	Chennai - Avadi	C	8925	Mon : 11:30am	Next day : 5:00pm	10 ml 24 Hrs Urine (Preferred) OR Spot urine
I0264	Immunoglobulin Profile IgG, IgM and IgA,CSF	Chennai - Avadi	C	3730	Mon & Thu : 9:00am	Same day : 5:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
I0265	Immunoglobulin Profile IgG, IgM and IgA, Serum	Chennai - Avadi	C	2035	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0266	Impotence Profile	Chennai - Avadi	C	3505	Daily : 9:00am to 9:00pm	After 8hrs	5 ml of serum
I0267	India Ink Preparation by Microscopy - CSF	Chennai - Avadi	B	325	Daily : 9:00am & 3:00pm	After 6hrs	2ml of CSF in Sterile leak proof tube or container
I0268	Infectious Mononucleosis Profile by IFA - Serum	Chennai - Avadi	C	5650	Mon & Thu : 7:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
I0269	Infertility Profile , Female	Chennai - Avadi	C	6780	Daily : 9:00am	5th day : 6:00pm	3 ml of EDTA whole blood, 5 ml of serum, 2 ml of flouride plasma, 10 ml of spot urine
I0270	Infertility Profile, Male	Chennai - Avadi	C	5650	Daily : 9:00am	5th day : 6:00pm	3 ml of EDTA whole blood, 5 ml of serum, 2 ml of flouride plasma, 10 ml of spot urine, Semen sample
I0273	Inhibin A Serum	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0274	Inhibin B Serum	Chennai - Avadi	C	2145	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0275	Insulin (Fasting) Serum	Chennai - Avadi	B	800	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Serum Collected In Fasting.
I0275_PG	Insulin, Post Glucose Serum	Chennai - Avadi	B	865	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of serum collected 2 hours post Glucose Monohydrate 82.5 gm

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0275_PP	Insulin, Post Prandial Serum	Chennai - Avadi	B	865	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of serum collected 2 hours post lunch
I0275_R	Insulin (Random) Serum	Chennai - Avadi	B	865	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of serum
I0276	Insulin Antibody Type I & II diabetes serum	Chennai - Avadi	D	2360	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0277	Insulin Suppression Test - C-peptide Serum	Chennai - Avadi	C	7795	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of Serum (basal, 10,20,30,40,50,60 Minutes)
I0279	Interleukin -6 CLIA Plasma	Chennai - Avadi	C	3335	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of EDTA plasma to be send in aliquot tube at proper temperature.
I0281	Intrinsic Factor Antibody by IFA Serum	Chennai - Avadi	C	2600	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0282	Iron biochemical Serum	Chennai - Avadi	B	380	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0284	Iron Urine 24H**	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
I0285	Iron Serum	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
I0286	Iron Studies, Serum	Chennai - Avadi	B	750	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0287	Islet Cell Antibody ICA-512 serum	Chennai - Avadi	D	2480	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0288	Isoleucine Quantitative Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
I0289	Isoleucine Quantitative Urine 24H**	Chennai - Avadi	R	6365	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
I0290	Isospora by smear examination - Stool	Chennai - Avadi	B	540	Daily : 9:00am & 3:00pm	After 6hrs	5gm of Stool specimen in Sterile leak proof container
I0291	IHC single marker with reporting Unspecified stain	Chennai - Avadi	D	2240	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							all relevant clinical and radiological detail
I0292	Reticulocyte count profile	Chennai - Avadi	B	380	Daily : 9:00am to 9:00pm	After 4hrs	3 ml EDTA Whole Blood
I0293	IHC stain only CD 56_I0293	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0294	IHC stain only CD 57_I0294	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0295	IHC stain only CD68_I0295	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0296	IHC stain only Lysozyme_I0296	Chennai - Avadi	D	1300	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0297	IHC p16INK4a stain only	Chennai - Avadi	D	3420	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0300	IgVH GENE MUTATION**_I0300	Chennai - Avadi	R	9310	Wed : 11:30am	10th day : 6:00pm	3 ml of EDTA whole blood
I0301	BCR ABL1 Mutation Analysis	Chennai - Avadi	E	10290	Daily : 9:00am	8th day : 6:00pm	5-10 ml of EDTA whole blood required. Sample stable for 3 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. Sample should reach to the lab within 48 to 72 hours.
I0305	IGH/CCND1:t(11;14) translocation Mantle Cell Lymphoma FISH, TISSUE	Chennai - Avadi	D	9440	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
I0312	IGH gene rearrangement Translocations involving 14q32 locus	Chennai - Avadi	D	9440	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0315	IDH1 and IDH2 gene mutation - FFPE Tissue	Chennai - Avadi	E	6825	Mon, Wed & Fri : 9:00am	5th day : 6:00pm	FFPE block to be shipped at Ambient temperature. Tumor containing FFPE block (>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
I0316	IGH/BCL2 :t(14;18) translocation Follicular Lymphoma FISH, BONE MARROW	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
I0317	IGH/BCL2 :t(14;18) translocation Follicular Lymphoma FISH, BLOOD	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
I0318	IGH/BCL2 fusion	Chennai - Avadi	D	9440	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
I0319	IHC single marker with reporting p40	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0320	IHC single marker with reporting Uroplakin II	Chennai - Avadi	D	2240	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0321	IHC stain only p40	Chennai - Avadi	D	1120	Daily : 9:00am	4th day : 6:00pm	Single Paraffin block / section on charged slide
I0322	IHC stain only Uroplakin II	Chennai - Avadi	D	1120	Daily : 9:00am	4th day : 6:00pm	Single Paraffin block / section on charged slide
I0323	Immature Platelet fraction (IPF) Count	Chennai - Avadi	B	380	Daily : 9:00am to 9:00pm	After 6hrs	3 ml EDTA whole Blood with direct smear

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				6:00pm	paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0337	PDL1 (SP263) by IHC - Ventana	Chennai - Avadi	D	9440	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0338	IMMUNOHISTOCHEMISTRY C-MYC ANTIBODY	Chennai - Avadi	D	2595	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0339	IMMUNOHISTOCHEMISTRY EBV	Chennai - Avadi	D	2595	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0342	IHC single marker with reporting TFE3	Chennai - Avadi	D	2595	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0347	IHC stain only C-myc antibody	Chennai - Avadi	D	1415	Daily : 9:00am	4th day : 6:00pm	Single Paraffin block / section on charged slide
I0348	IHC stain only EBV	Chennai - Avadi	D	1415	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0350	IHC stain only CD7	Chennai - Avadi	D	2595	Daily : 9:00am	4th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0352	IHC stain only IgG4	Chennai - Avadi	D	2595	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0353	IHC stain only GATA3	Chennai	D	2595	Daily : 9:00am	3rd day :	Paraffin blocks (Formalin fixed

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				6:00pm	paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0354	IHC stain only H caldesmon	Chennai - Avadi	D	2595	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0355	IHC stain only ERG	Chennai - Avadi	D	2595	Daily : 9:00am	3rd day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0363	IHC stain only Androgen Receptor	Chennai - Avadi	D	1770	Daily : 9:00am	3rd day : 6:00pm	Single Paraffin block / section on charged slide
I0366	Hypersensitive pneumonitis panel	Chennai - Avadi	D	8495	Mon, Wed & Fri : 11:30am	Next day : 5:00pm	3ml serum sample required
I0366_1	Penicillium Chrysogenum(Specific IgG)	Chennai - Avadi	C	1695	Mon & Thu : 11:30am	Next day : 6:00pm	Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
I0366_2	Cladosporium herbarum (Specific IgG)	Chennai - Avadi	C	1695	Mon & Thu : 11:30am	Next day : 6:00pm	Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
I0366_3	Aspergillus fumigatus(Specific IgG)	Chennai - Avadi	C	1920	Mon & Thu : 11:30am	Next day : 5:00pm	Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0366_4	Mucor Racemosus(Specific IgG)	Chennai - Avadi	C	1695	Mon & Thu : 11:30am	Next day : 6:00pm	Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
I0366_5	Alternaria alternata(Specific IgG)	Chennai - Avadi	C	1695	Mon & Thu : 11:30am	Next day : 6:00pm	Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
I0368	Tb- Interferon-γ (Gamma) release assay (IGRA)-4 tube	Chennai - Avadi	D	3895	Mon, Wed & Fri : 7:30am	Next day : 5:00pm	5 ml of Lithium Heparin Blood
I0371	Iron, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
I0372	IDH1 and IDH2 gene mutation, AML, Blood	Chennai - Avadi	E	6825	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
I0372_BM	IDH1 and IDH2 gene mutation, AML, Bone Marrow	Chennai - Avadi	E	6825	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
I0373	Tb- Interferon-γ (Gamma) release assay (IGRA)-3 tube	Chennai - Avadi	C	3390	Mon, Wed & Fri : 7:30am	Next day : 5:00pm	5 ml of Lithium Heparin Blood
I0378	Immunofluorescence Analysis in Michels Medium (Kidney / Skin Biopsy)	Chennai - Avadi	D	4130	Daily : 9:00am	5th day : 6:00pm	Biopsy to be sent in Michelâ€™s Medium in a Leak proof container along with relevant clinical &

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							radiological details
I0380	Isavuconazole **	Chennai - Avadi	R	6420	Daily : 7:30am	7th day : 6:00pm	3 ml of EDTA Plasma.Please indicate history of the drug dose taken, time of sample collection and clinical History along with the sample.
I0384	IHC stain only ATRX	Chennai - Avadi	D	4130	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0385	IHC single marker with reporting ATRX	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0386	IHC stain only IDH 1	Chennai - Avadi	D	4130	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0387	IHC single marker with reporting IDH 1	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0388	IHC stain only SOX 10	Chennai - Avadi	D	4130	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0389	IHC single marker with reporting SOX 10	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0391	IHC single marker with reporting HHV 8	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0394	IHC stain only SV 40	Chennai - Avadi	D	4130	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0395	IHC single marker with reporting SV 40	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							detail
I0397	IGLON5 ANTIBODY Serum	Chennai - Avadi	D	14160	Tue & Fri: 11:30am	Next day : 6:00pm	3 ml Serum
I0399	IHC stain only PLA 2 R1	Chennai - Avadi	D	4130	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0400	IHC single marker with reporting PLA 2 R1	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0401	IHC stain only STAT 6	Chennai - Avadi	D	2950	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0402	IHC single marker with reporting STAT 6	Chennai - Avadi	D	3540	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0403	IHC stain only INI 1	Chennai - Avadi	D	4130	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0404	IHC single marker with reporting INI 1	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0405	IHC stain only Sall 4	Chennai - Avadi	D	4130	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0406	IHC single marker with reporting Sall 4	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0407	IHC stain only NKX 3.1	Chennai - Avadi	D	4130	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0408	IHC single marker with reporting NKX 3.1	Chennai - Avadi	D	5310	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0409	IHC stain only CD 4	Chennai - Avadi	D	2360	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0410	IHC single marker with reporting CD 4	Chennai - Avadi	D	3540	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0411	IHC stain only CD 8	Chennai - Avadi	D	2360	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0412	IHC single marker with reporting CD 8	Chennai - Avadi	D	3540	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0413	IHC single marker with reporting CMV	Chennai - Avadi	D	2595	Daily : 9:00am	7th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I1005	IgG Subclass 1(IgG1) serum by Nephelometry**	Chennai - Avadi	R	4280	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum (Age & gender mandatory)
I1006	IgG Subclass 2(IgG2) serum by Nephelometry**	Chennai - Avadi	R	4280	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum (Age & gender mandatory)
I1007	IgG Subclass 3(IgG3) serum by Nephelometry**	Chennai - Avadi	R	4280	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum (Age & gender mandatory)
I1025	IODINE, Serum by ICP-MS method**	Chennai - Avadi	R	5885	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml Serum
I1026	IODINE, Urine by ICP-MS method**	Chennai - Avadi	R	9630	Daily : 11:30am	After 4 weeks : 6:00pm	Individual Metals and Minerals.pdf
I1901	Impression Smears	Chennai - Avadi	#N/A	#N/A	Daily : 9:00am	After 48hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
I1902	Insulin, Serum by CMIA	Chennai - Avadi	B	730	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
I5073	Isoniazid**	Chennai - Avadi	R	7490	Daily : 7:30am	10th day : 6:00pm	2 ml Na-Heparin Plasma
I5649	FLU-XPRT VIRAL PANEL	Chennai - Avadi	D	5310	Daily : 7:00am & 3:00pm	Same day : 3:00pm & 11:00pm	Swab with Sufficient Material in Viral Transport Medium (VTM)
I8151	Immunofixation image	Chennai - Avadi	C	565	Tue, Thu & Sat: 11:30am	Next day : 5:00pm	Pre-requisite: Kindly mention the VID for gel electrophoresis test asked.
I8155	PD-L1 Dako 22C3 IHC	Chennai - Avadi	D	14750	Daily 7.30 am	10th days	Paraffin Block
I8335	ILD Complete	Chennai - Avadi	D	20650	Tue & Fri : 9:00am	Next day : 6:00pm	5 ML OF SERUM
I8336	ILD Screen	Chennai - Avadi	D	8260	Mon & Thu : 9:00am	Next day : 6:00pm	3 ml EDTA WHOLE BLOOD,5 ML OF SERUM,
I8337	ILD Autoimmune Maxi	Chennai - Avadi	D	15930	Tue & Fri : 9:00am	Next day : 6:00pm	5 ML OF SERUM
I8338	ILD Autoimmune Profile Mini	Chennai - Avadi	D	8850	Tue & Fri : 9:00am	Next day : 6:00pm	5 ML OF SERUM
J0001	JAK-2 Mutation (CMPD) Exon 12 Mutation Bone Marrow	Chennai - Avadi	D	7080	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
J0002	JAK-2 Mutation (CMPD) (V617F) by Real	Chennai	D	7465	Tue & Fri : 9:00am	Next day :	3 ml of Bone Marrow in EDTA

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
J0003	Time PCR - Bone Marrow JAK-2 Mutation (CMPD) Exon 12 Mutation Blood	- Avadi Chennai - Avadi	D	7080	Wed & Sat : 9:00am	7:00pm 7th day : 6:00pm	Vacutainer 3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
J0004	JAK-2 Mutation (CMPD) (V617F) by Real Time PCR - Blood	Chennai - Avadi	D	5885	Tue & Fri : 9:00am	Next day : 7:00pm	3 ml of Whole Blood in EDTA Vacutainer
J0005	JC-BK Virus DNA Detection by Real Time PCR - CSF	Chennai - Avadi	D	6490	Mon & Thurs : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
J0006	JC-BK Virus DNA Detection by Real Time PCR - EDTA Plasma	Chennai - Avadi	D	6490	Mon & Thurs : 9:00am	Next day : 5:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
J0007	JC-BK Virus DNA Detection by Real Time PCR - Urine	Chennai - Avadi	D	6490	Mon & Thurs : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
J0008	JEV RNA Detection by Real Time PCR - EDTA Blood	Chennai - Avadi	D	4660	Fri : 7:30am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
J0009	JEV RNA Detection by Real Time PCR - CSF	Chennai - Avadi	D	4660	Fri : 7:30am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
J0010	JEV RNA Detection by Real Time PCR - EDTA Plasma	Chennai - Avadi	D	4660	Fri : 7:30am	Next day : 5:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
J0011	JEV RNA Detection by Real Time PCR - Serum	Chennai - Avadi	D	4660	Fri : 7:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
J0012	Jo-1 Antibody Serum	Chennai - Avadi	C	1640	Wed & Sat : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
J0013	JC DNA Detection by Real Time PCR - CSF	Chennai - Avadi	D	4720	Mon & Thurs : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
J0014	JC DNA Detection by Real Time PCR - EDTA Plasma	Chennai - Avadi	D	4720	Mon & Thurs : 9:00am	Next day : 5:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
J0015	JC DNA Detection by Real Time PCR - Urine	Chennai - Avadi	D	5430	Mon & Thurs : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
J0016	JAK2 V617 Mutation, with Reflex to Jak2 Ex- 12, CALR Ex-9 Mutation and MPL W515, S505 Mutation	Chennai - Avadi	D	15930	Wed & Sat : 9:00am	10th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
J0017	JAK2 V617 Mutation, with Reflex to Jak2 Ex- 12, CALR Ex-9 Mutation and MPL W515, S505 Mutation,Bone Marrow	Chennai - Avadi	D	15930	Wed & Sat : 9:00am	10th day : 6:00pm	Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. 3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
J0020	JC DNA Detection by Real Time PCR - EDTA Blood	Chennai - Avadi	D	5430	Mon & Thurs : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
J0021	JC-BK Virus DNA Detection by Real Time PCR - EDTA Blood	Chennai - Avadi	D	6725	Mon & Thurs : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
K0001	Kappa And Lambda-Freelite, Urine	Chennai - Avadi	D	5310	Daily : 9:00am to 9:00pm	After 8hrs	10 ml of 24 hrs urine (Preferred) / Spot urine
K0002	Kappa And Lambda-Freelite, Serum	Chennai - Avadi	D	4720	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
K0003	Kappa Light Chains by Flow cytometry	Chennai - Avadi	D	2950	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
K0004	Kappa Light Chains-Total Urine 24H	Chennai - Avadi	D	3540	Daily : 9:00am to 9:00pm	After 8hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
K0005	FISH for Aneuploidy in POC - Chromosomes - 13,16,18, 21, 22, X and Y	Chennai - Avadi	D	10620	Daily : 9:00am	6th day : 6:00pm	POC in sterile container with normal saline with few drops of antibiotic. In case of expected delay store at 4 degrees. Never add formalin or freeze. Note: Whole Fetus is not accepted.
K0006	Karyotyping reflex FISH-13,18,21,X,Y AmnioticFluid	Chennai - Avadi	D	13570	Daily : 9:00am	15th day : 6:00pm	20 ml of Amniotic Fluid In sterile centrifuge Falcon 15 ml tubes.
K0007	Karyotyping reflex FISH(13,18,21,X,Y)Chorion Villi	Chennai - Avadi	D	13570	Daily : 9:00am	15th day : 6:00pm	20-30 mg Chorionic Villus in Sterile nutrient medium in 15 ml Falcon tubes(Medium will be provided by Metropolis).
K0010	FISH for Aneuploidy in POC - Chromosomes 13,18,21, X,Y	Chennai - Avadi	D	7670	Daily : 9:00am	6th day : 6:00pm	POC in sterile container with normal saline with few drops of antibiotic. In case of expected delay store at 4 degrees. Never add formalin or freeze.Note: Whole Fetus is not accepted.
K0012	Karyotyping by G-Banding Peripheral Blood	Chennai - Avadi	D	4500	Daily : 9:00am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
K0014	Karyotyping by G-Banding Bone marrow	Chennai - Avadi	D	4720	Daily : 9:00am	10th day : 6:00pm	2-3 ml Bone Marrow/3 ml Blood in Na-heparin vacutainer
K0017	Karyotyping by high resolution banding (HRB) Blood	Chennai - Avadi	D	5310	Daily : 9:00am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
K0018	Kaolin Clotting Time	Chennai - Avadi	B	1350	Daily : 9:00am	After 8hrs	1 ml of citrated platelet poor plasma in frozen condition.
K0019	KRAS-Exon 2 and 3 Mutation Analysis - FFPE Tissue	Chennai - Avadi	D	7080	Mon, Wed & Fri : 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>10% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
K0021	Kennedy disease (CAG repeats) by PCR	Chennai - Avadi	R	6205	Daily : 11:30am	After 4 weeks : 6:00pm	report. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
K0022	Karyotyping from cell pellet /slides	Chennai - Avadi	D	1180	Daily : 9:00am	10th day : 6:00pm	Cell pellet / Slides
K0023	Kala azar	Chennai - Avadi	R	1820	Daily : 11:30am	5th day : 6:00pm	3ml of serum
K0026	Karyotyping G Banding and FISH (Chr 13,18, 21, X and Y),Amniotic Fluid	Chennai - Avadi	D	18290	Daily : 9:00am	15th day : 6:00pm	20 ml of Amniotic Fluid In Sterile Falcon centrifuge tube
K0030	Karyotyping by G-Banding for neonates	Chennai - Avadi	D	5310	Daily : 9:00am	7th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
K0031	Karyotyping by G-Banding for 50 metaphase cells	Chennai - Avadi	D	5310	Daily : 9:00am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
K0033	Karyotyping by G-Banding reflex FISH	Chennai - Avadi	D	11680	Daily : 9:00am	15th day : 6:00pm	3 ml of Prenatal Cord Blood In Na-Heparin Vacutainer
K1901	Canine - Kidney panel	Chennai - Avadi	D	1240	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
K8112	Kidney biopsy for Medical Renal Disease	Chennai - Avadi	D	5900	Daily : 9:00am	7th day	1) Tissue in 10% NB formalin for histopathology 2) Tissue in Michel Medium for Immunofluorescence. Please use this code when you are sending 2 sample vials ; One in Formaldehyde & other in Michels medium.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
K8113	Kidney biopsy for Transplant Pathology	Chennai - Avadi	D	11680	Daily : 9:00am	14th day	1) Tissue in 10% NB formalin for histopathology 2) Tissue in Michel Medium for Immunofluorescence. Please use this code when you are sending 2 sample vials ; One in Formaldehyde & other in Michels medium.
L0001	Lactate CSF	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
L0002	Lactate Biochemical Plasma	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Fluoride Plasma Labile Analyte
L0003	Lactose Tolerance Test 5 Samples for Glucose, Plasma	Chennai - Avadi	C	680	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of flouride plasma (5Samples)
L0004	Lambda Light Chains by Flow cytometry	Chennai - Avadi	D	2950	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
L0005	Lambda Light Chains-Total Urine	Chennai - Avadi	D	3540	Daily : 9:00am to 9:00pm	After 8hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
L0006	Lamotrigine Level Plasma**	Chennai	D	4840	Tue : 11:30am	Next day :	3 ml of EDTA plasma or heparin

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				5:00pm	plasma to be sent in aliquot tube at proper temperature.
L0007	LAP Score	Chennai - Avadi	B	1190	Daily : 11:30am	3rd day : 6:00pm	5 Peripheral Blood Smears and EDTA Blood. Clinical history and clinician no is mandatory
L0008	Liver cytosolic antigen type-1(LC-1)Antibody, Serum	Chennai - Avadi	C	3165	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0009	LD (Leigh`s Disease),Blood **	Chennai - Avadi	R	15515	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
L0010	LD bodies Leishmania Dodovani bodies by Microscopy - Blood	Chennai - Avadi	B	970	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA Bone Marrow in EDTA Vacutainer
L0011	LD bodies Leishmania Dodovani bodies by Microscopy - Bone Marrow	Chennai - Avadi	C	1015	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA Whole Blood in EDTA Vacutainer
L0012	LDH Lactate Dehydrogenase Serum	Chennai - Avadi	A	360	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0013	LDH Lactate Dehydrogenase Ascitic Fluid	Chennai - Avadi	A	415	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Ascitic Fluid
L0014	LDH Lactate Dehydrogenase Peritoneal Fluid	Chennai - Avadi	A	415	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Peritoneal Fluid
L0015	LDH Lactate Dehydrogenase Pleural Fluid	Chennai - Avadi	A	415	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Pleural Fluid
L0016	LDH- Isoenzyme by Elctrophoresis Serum	Chennai - Avadi	D	4720	Wed : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0017	LDL Cholesterol -Direct serum	Chennai - Avadi	A	325	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0019	Lead, Blood	Chennai	C	1980	Daily:7:30 am	Same day	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0020	Lead, Urine 24 Hrs	- Avadi Chennai	C	1980	Daily : 11:30am	8:00 PM 3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
L0021	Lead,Urine Spot	- Avadi Chennai	C	1980	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
L0022	Legionella Culture by Vitek/ MALDITOF - Sputum / BAL**	- Avadi Chennai	R	1175	Daily : 9:00am to 9:00pm	15th day : 6:00pm	2 ml of Sputum in Sterile leak proof container
L0026	Legionella Pneumophila Antigen detection Urine	- Avadi Chennai	C	3165	Daily : 9:00am & 3:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
L0027	Legionella Pneumophila IgG antibodies by EIA - Serum	- Avadi Chennai	C	2825	Mon & Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0028	Legionella Pneumophila IgM antibodies by EIA - Serum	- Avadi Chennai	C	2825	Mon & Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0029	Leishmania Antibody-IgG Antibody for Kala Azar by EIA - Serum	- Avadi Chennai	C	1920	Tue & Fri : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0034	Leptin Human Serum	- Avadi Chennai	D	4600	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0035	Leptospira Detection by smear examination - Blood	- Avadi Chennai	B	595	Daily : 9:00am & 3:00pm	After 6hrs	3 ml of Whole Blood in EDTA Vacutainer
L0036	Leptospira Detection by smear examination - CSF	- Avadi Chennai	B	700	Daily : 9:00am & 3:00pm	After 6hrs	2ml of CSF in Sterile leak proof tube or container
L0037	Leptospira Detection by smear examination - Urine	- Avadi Chennai	B	700	Daily : 9:00am & 3:00pm	After 6hrs	15ml of Spot Urine in Sterile leak proof container
L0038	Leptospira DNA Detection by Real Time PCR - Blood	- Avadi Chennai	D	2360	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
L0039	LeptospiraDNA Detection by Real Time PCR - Urine	- Avadi Chennai	D	2360	Daily : 10:00am	After 24hrs	15ml of Spot Urine in Sterile leak proof container
L0040	Leptospira IgG antibodies by EIA - Serum	- Avadi Chennai	C	1130	Daily : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0041	Leptospira IgM antibodies by EIA - Serum	- Avadi Chennai	C	1130	Daily : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0042	Leptospira Profile	- Avadi Chennai	C	2740	Daily : 11:30am	same day 11:00 pm	3 ml of Whole Blood in EDTA Vacutainer, 3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0043	Leucine Quantitative Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
L0044	Leucine Quantitative Urine 24H**	Chennai - Avadi	R	6365	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
L0046	Leukemia panel by FISH - ALL	Chennai - Avadi	D	14750	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0047	Leukemia panel by FISH - AML	Chennai - Avadi	D	17700	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0051	Leukemia panel by FISH - CLL	Chennai - Avadi	D	16520	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0052	Leukemia panels by FISH - Eosinophilia CBF-Beta, PDGFR-Alpha & beta, FGFR-1 Bone Marrow/Leukemic Blood	Chennai - Avadi	D	14160	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0053	Myelodysplastic Syndrome (MDS) Panel by FISH, Bone marrow	Chennai - Avadi	D	10620	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0054	Comprehensive Multiple Myeloma Panel by FISH	Chennai - Avadi	D	18290	Daily : 9:00am	7th day : 6:00pm	3 ml of Heparin Bone Marrow
L0055	Leukemia -Acute Leukemia Panel, Blood	Chennai - Avadi	D	11210	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA or heparin blood and Direct smears required
L0056	Leukemia - Acute Leukemia Panel, Bone Marrow	Chennai - Avadi	D	11210	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA and heparin bone marrow. Direct smears required
L0057	Leukemia - Chronic Lymphoproliferative Disorder Panel, Blood	Chennai - Avadi	D	14160	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA or heparin blood and Direct smears required
L0058	Leukemia - Chronic Lymphoproliferative Disorder Panel, Bone Marrow	Chennai - Avadi	D	16285	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA or heparin bone marrow and Direct smears required
L0065	Myeloma panel by flowcytometry, Bone marrow**	Chennai - Avadi	D	9440	Daily : 11:30am	5th day : 6:00pm	3 ml of EDTA or heparin bone marrow and Direct smears required
L0066	LH-Luteinizing Hormone ECLIA Serum	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0067	LHON (Leber's hereditary optic neuropathy) by PCR**	Chennai - Avadi	R	17120	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
L0068	Lipase Serum	Chennai - Avadi	B	430	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0068_Fluid	Fluid Lipase	Chennai - Avadi	B	405	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0069	Lipase Urine	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
L0070	Lipid Profile - Maxi	Chennai - Avadi	C	1525	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0071	Lipid Profile-Mini	Chennai - Avadi	B	550	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0071_NF	Lipid profile 2 (Non- fasting)	Chennai - Avadi	B	485	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0072	Lipid ScreenCholesterol and TriglyceridesSerum	Chennai - Avadi	B	430	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0072_NF	Lipid profile 3 (Non- fasting)	Chennai - Avadi	B	430	Daily : 9:00am to 9:00pm	After 6hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0073	Total Lipids, serum	Chennai - Avadi	A	415	Daily : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0076	Lipoprotein Electrophoresis, Serum	Chennai - Avadi	C	2825	Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0077	Lithium Serum	Chennai - Avadi	C	510	Daily : 9:00am & 3:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0079	Liver Function Test-1(Maxi)	Chennai - Avadi	B	1300	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of serum and 3 ml of citrate plasma
L0080	Liver Function Test-2(Mini)	Chennai - Avadi	B	550	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of serum
L0081	Liver Onco Marker Profile	Chennai - Avadi	D	1650	Daily : 9:00am to 9:00pm	After 8hrs	5 ml of serum
L0082	Liver Function Test-3 (LIVER SCREEN)	Chennai - Avadi	B	430	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0083	LKM1 by IFA Liver Kidney Microsomes Reflex to end point titre Serum	Chennai - Avadi	C	2090	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0084	Lp(a)Lipoprotein(a) serum	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0085	Lung Marker Profile	Chennai - Avadi	D	5310	Wed & Sat : 7:30am	Next day : 5:00pm	3 ml of Serum and 3 ml of EDTA plasma
L0086	Lupus Anticoagulants (LAC) profile	Chennai - Avadi	D	2350	Daily : 9:00am & 5:00pm	After 6hrs	1 ml of citrated platelet poor plasma in frozen condition.
L0087	LymeBorrelia Burgdorferi IgG antibodies by EIA - Serum	Chennai - Avadi	C	1920	Sat : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0089	LymeBorrelia Burgdorferi IgM antibodies by EIA - Serum	Chennai - Avadi	C	1920	Sat : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0091	Lysine Quantitative Plasma	Chennai	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi					plasma to be sent in aliquot tube at proper temperature.
L0092	Lysine Quantitative Urine	Chennai - Avadi	R	6365	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
L0097	Leukemia Characterisation Panel for CLPD, Blood	Chennai - Avadi	D	23365	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA and heparin blood. Direct smears required
L0098	Leukemia Characterisation Panel for CLPD , Bone Marrow	Chennai - Avadi	D	23365	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA and heparin bone marrow and Direct smears required
L0099	Leukemia Comprehensive workup panel- Acute Leukemia, Blood	Chennai - Avadi	D	29500	Daily : 11:30am	12th day : 6:00pm	3 ml of EDTA and heparin blood. Direct smears required
L0100	Leukemia Comprehensive workup panel- Acute Leukemia, Bone Marrow	Chennai - Avadi	D	29500	Daily : 11:30am	12th day : 6:00pm	3 ml of EDTA and heparin bone marrow. Direct smears required
L0101	Leukemia CMPD Panel: JAK2 by PCR and BCR-ABL by FISH,Blood	Chennai - Avadi	D	10385	Daily : 9:00am	5th day : 6:00pm	3 ml of EDTA and Na-heparin blood and Direct smears required
L0102	Leukemia CMPD Panel: JAK2 by PCR and BCR-ABL by FISH, Bone Marrow	Chennai - Avadi	D	10385	Daily : 9:00am	5th day : 6:00pm	3 ml of EDTA Bone Marrow
L0103	Lp- PLA2 (PLAC)Cardiac risk marker serum	Chennai - Avadi	D	2950	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0105	Leukemia panel by FISH - ALL, Leukemic Blood	Chennai - Avadi	D	15340	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0106	Leukemia AML panel by FISH, Leukemic Blood	Chennai - Avadi	D	17700	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0107	Leukemia panel by FISH - CLL, Leukemic Blood	Chennai - Avadi	D	16520	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0108	Leukemia panels by FISH - Eosinophilia CBF-Beta, PDGFR-Alpha & beta, FGFR-1 Leukemic Blood	Chennai - Avadi	D	14160	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0109	Myelodysplastic Syndrome (MDS) Panel by FISH	Chennai - Avadi	D	10620	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0111	Lymphocyte subset panel with NK cells	Chennai - Avadi	D	7080	Daily : 11:30am	Next day : 5:00pm	6 ml EDTA whole Blood
L0112	LDH Lactate Dehydrogenase CSF	Chennai - Avadi	A	415	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
L0113	Lyme-Borrelia Burgdorferi IgM by EIA - CSF	Chennai	C	5425	Sat : 11:30am	Next day :	2 ml of CSF in Sterile leak proof

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				5:00pm	tube or container
L0114	Lyme-Borrelia Burgdorferi IgG by EIA - CSF	Chennai - Avadi	C	5425	Sat : 11:30am	Next day : 5:00pm	2 ml of CSF in Sterile leak proof tube or container
L0117	Leukemia translocation Panel (ALL)* - t(1;19), t(4;11), t(9;22), t(12;21)	Chennai - Avadi	D	8850	Daily : 9:00am	7th day : 6:00pm	5-10 ml of EDTA whole blood required. Sample stable for 3 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. Sample should reach to the lab within 48 to 72 hours.
L0118	Leukemia translocation Panel (AML)* - t(8;21), t(15;17), Inv(16)	Chennai - Avadi	D	8850	Daily : 9:00am	7th day : 6:00pm	5-10 ml of EDTA whole blood required. Sample stable for 3 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. Sample should reach to the lab within 48 to 72 hours.
L0120	LDH	Chennai - Avadi	A	390	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
L0123	LDH Body fluid	Chennai - Avadi	A	415	Daily : 9:00am to 9:00pm	After 4hrs	2 ml body fluid
L0133	LUNG CANCER TARGETED THERAPY PANEL (EGFR+ALK+ROS1)	Chennai - Avadi	E	20740	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical Diagnosis with Histopathology report, HistologyandIHC.pdf
L0134	Lung Panel by NGS [EGFR, BRAF, MET (includes Ex-14 skipping), RET, ERBB2, KRAS (SNVs & small INDELS); ALK, ROS1,	Chennai - Avadi	E	32000	Daily : 7:30am	25th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	RET and NTRK1, NTRK2, NTRK3 (Fusions)]						hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
L0135	Lead, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
L0136	Lead, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
L0138	Lithium, Spot Urine**	Chennai - Avadi	C	3110	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
L0139	Leukemia Panel-Ph like B-ALL, Bone Marrow **	Chennai - Avadi	R	20330	Daily : 11:30am	10th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0139_BLOOD	Leukemia Panel-Ph like B-ALL, Blood**	Chennai - Avadi	R	20330	Daily : 9:00am	8th day : 6:00pm	3 ml BLOOD in Na-heparin vacutainer
L0140	Leukemia Panel-T-ALL, Blood**	Chennai - Avadi	R	17120	Daily : 11:30am	8th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0140_BM	Leukemia Panel-T-ALL, Bone marrow**	Chennai - Avadi	R	17120	Daily : 11:30am	10th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0141	Leukemia Panel-AML Mutation, Blood	Chennai - Avadi	D	10620	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
L0141_BM	Leukemia Panel-AML Mutation, Bone marrow	Chennai - Avadi	D	10620	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
L0142	NextGen Oncomine Myeloid 69 Gene Panel	Chennai - Avadi	E	20050	DAILY : 9 : 00 AM	20 Days	6 ml EDTA whole Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0143	Leptin deficiency (LEP) -Gene Analysis	Chennai - Avadi	E	21000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
L0144	Myeloma MRD Flow cytometry Bone marrow **	Chennai - Avadi	R	19795	Daily : 11:30am	5th day : 6:00pm	sample : 3 ml of EDTA Bone Marrow & direct smear
L0145	Flow for AML MRD**	Chennai - Avadi	R	19795	Daily : 11:30am	5th day : 6:00pm	3 ml of EDTA Bone Marrow & General Instruction Local Samples reaching in 24 Hours can be send in Ambient condition and outstation samples to be shipped at Refrigerated condition. Do not keep EDTA or Heparin blood samples in contact with ice packs. Do not freeze. Medical history / Heamato-pathological reports with referring doctors number.
L0146	Lip biopsy for Sjogrens syndrome (Focus score) with WSI	Chennai - Avadi	C	2485	Daily : 9:00am	4th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
L0152	Lung Cancer NGS panel (SNVs, InDels and Fusions)- Liquid Biopsy**	Chennai - Avadi	E	46200	Daily : 11:30am	25th day : 6:00pm	Peripheral Blood in Streck tube (2 tubes of 10 ml each) (Tube suplied upon request)
L1008	Levetiracetam, serum	Chennai - Avadi	D	5900	Wed : 11:30am	Next day : 5:00pm	3 ml of serum(sample collect in, plain vacutainer, not in gel tube)
L1901	Canine - Lipase	Chennai - Avadi	B	565	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L1902	Canine - Liver profile	Chennai - Avadi	D	1005	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
L8015	LH/FSH Ratio Serum	Chennai - Avadi	B	865	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0001	Anti MAG (Myelin associated Glycoprotein) Ab (M0001)	Chennai - Avadi	D	4720	Tue & Fri : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0002	Magnesium Serum	Chennai - Avadi	B	460	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0003	Magnesium, Spot Urine	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
M0004	Magnesium Urine 24H	Chennai - Avadi	B	595	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
M0006	Malaria Parasite detection by Direct Smear examination by Microscopy - EDTA Blood	Chennai - Avadi	B	250	Daily : 9:00am, 1:00pm & 5:00pm	After 4hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0007	Malaria DNA Detection by Real Time PCR - EDTA Blood	Chennai - Avadi	C	3165	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
M0008	Malaria Detection by QBC - EDTA Blood	Chennai - Avadi	B	325	Daily : 9:00am, 1:00pm & 5:00pm	After 4hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer
M0009	Malaria Profile reflex Malaria antigen by Immunochromatography - EDTA Blood	Chennai - Avadi	C	395	Daily : 9:00am, 1:00pm & 5:00pm	After 4hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer
M0010	Malarial Antigen Vivax & Falciparum by Immunochromatography - EDTA Blood	Chennai - Avadi	B	650	Daily : 9:00am, 1:00pm & 5:00pm	After 4hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer
M0011	Manganese Blood	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0012	Manganese Serum	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0013	Manganese Urine	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0014	Manganese Urine 24H**	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0015	Mantoux test Tuberculin test (10TU)	Chennai - Avadi	O	220	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	Mycobacteriology.pdf
M0015_2TU	Mantoux test Tuberculin test (2TU)	Chennai - Avadi	O	165	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	Mycobacteriology.pdf
M0015_5TU	Mantoux test Tuberculin test (5TU)	Chennai - Avadi	O	205	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	Mycobacteriology.pdf
M0019	MDMA (Qualitative) Ecstasy, XTC Urine Spot	Chennai - Avadi	D	765	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
M0020	Measles Virus IgG Antibody by CLIA - Serum	Chennai - Avadi	C	1810	Daily : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0021	Measles Virus IgG Antibody by EIA - CSF	Chennai - Avadi	C	3275	Tue & Fri: 7:30am	Next day : 5:00pm	2 ml of CSF in Sterile leak proof tube or container & 2 mL Serum
M0022	Measles Virus IgM Antibody by CLIA - Serum	Chennai - Avadi	C	1810	Daily : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0023	Medullary Thyroid Carinoma Mutation Detection Sequencing method Blood**	Chennai - Avadi	R	13375	Daily : 11:30am	10th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0024	MELAS-Mutation Detection Blood**	Chennai - Avadi	R	6420	Daily : 11:30am	8th day : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0026	Meningitis (viral) Profile IgG by CLIA - Serum	Chennai - Avadi	D	5430	Wed : 9:00am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0027	Meningitis (viral) Profile IgG & IgM by CLIA - Serum	Chennai - Avadi	D	8850	Wed : 9:00am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0028	Meningitis (viral) Profile IgM by CLIA - Serum	Chennai - Avadi	D	5430	Wed : 9:00am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0029	Menopause	Chennai - Avadi	C	3165	Daily : 9:00am to 9:00pm	After 8hrs	5 ml of serum
M0030	Menopause diagnostic profile	Chennai - Avadi	C	1920	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of serum
M0031	Menopause monitoring profile	Chennai - Avadi	C	2205	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of serum
M0032	Mercury Blood	Chennai - Avadi	D	3540	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0033	Mercury Serum	Chennai - Avadi	D	3540	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0034	Mercury Urine	Chennai - Avadi	D	3540	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0035	Mercury Urine 24H**	Chennai - Avadi	D	3540	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0036	MERRF Mutation Detection Blood**	Chennai - Avadi	R	10700	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0037	Neoxpert by GCMS**	Chennai - Avadi	R	7385	Daily : 11:30am	4th day : 6:00pm	Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory. Spot Urine:10 ml of Spot urine required with clinical history and clinician number.
M0038	Metabolic Screen	Chennai - Avadi	R	3745	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
M0039	Metanephrine- Free Plasma	Chennai - Avadi	D	4720	Tue & Fri : 11:30am	Next day : 5:00pm	Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
M0040	Metanephrine level, HPLC, 24 hrs urine	Chennai - Avadi	D	3500	Mon & Thu : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
M0041	Methadone Dolophine Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
M0042	Methamphetamine Qualitative Urine Spot	Chennai - Avadi	B	680	Daily : 9:00am to 9:00pm	After 10 mins	Urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0042	Methamphetamine Qualitative Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
M0043	Methaqualone Qualitative Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
M0044	Meth-Haemoglobin Blood	Chennai - Avadi	C	1075	Daily : 11:30am	Next day : 5:00pm	3 ml of Heparin Blood. Sample required within 12 to 24 hrs after collection & avoid holiday & weekends.
M0045	Methionine Quantitative Plasma*	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
M0046	Methionine Quantitative Urine SPOT**	Chennai - Avadi	R	6365	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
M0047	Methotrexate Serum	Chennai - Avadi	C	3390	Wed & Sat : 7:30am	Same day : 6:00pm	3 ml of serum(sample collect in, plain vacutainer, not in gel tube)
M0048	Methyl Malonic Acid Qualitative Urine Spot	Chennai - Avadi	R	695	Daily : 11:30am	5th day : 6:00pm	15 ml of Spot Urine
M0049	Methyl Malonic Acid Quantitative Urine SPOT**	Chennai - Avadi	R	10700	Daily : 11:30am	5th day : 6:00pm	15 ml of Spot Urine
M0050	MGMT PCR Temozolomide Resistance Tissue**	Chennai - Avadi	R	9630	Daily : 11:30am	25th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
M0051	MICM By PCR Maternally Inherited Cardiomyopathy Blood**	Chennai - Avadi	R	7490	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0052	Microalbumin Urine 24H	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
M0053	Microalbumin / Creatinine Ratio Urine Spot	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
M0054	Microsporidium Detection by modified ZN Stain - Stool	Chennai - Avadi	B	540	Daily : 9:00am & 3:00pm	After 6hrs	5gm of Stool specimen in Sterile leak proof container
M0055	Mitochondrial (M2) Antibody by IFA Reflex to end point titre Serum	Chennai - Avadi	C	2035	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0059	Mixing Studies (APTT)	Chennai - Avadi	B	2540	Daily : 9:00am & 4:00pm	After 4hrs	1 ml of citrated platelet poor plasma in frozen condition.
M0060	MMC-Mutation Detection Maternal Myopathy With Cardiomyopathy Blood**	Chennai - Avadi	R	7490	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0061	MM-Mutation Detection Mitochondrial Myopathy	Chennai - Avadi	R	9095	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0063	MPO by EIA, p-ANCA Serum	Chennai - Avadi	C	1525	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0064	(cMPO)Myeloid cell marker	Chennai - Avadi	D	2715	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition (preferably). Clinical history and clinician details required.
M0065	MPO (Myeloperoxidase)Stain, Blood	Chennai - Avadi	C	905	Daily : 11:30am	3rd day : 6:00pm	3 Direct Peripheral smear slides, EDTA blood. Clinical history and clinician no is mandatory
M0066	MPO (Myeloperoxidase) Stain, Bone Marrow	Chennai - Avadi	C	960	Daily : 11:30am	3rd day : 6:00pm	3 Direct Bone Marrow Aspirate slides. Clinical history and clinician no is mandatory
M0068	MRSA Screen - Nasal swab	Chennai - Avadi	C	905	Daily : 9:00am to 9:00pm	4th day : 6:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
M0069	MSUD-Maple Syrup Urine Disorder Neonatal Screen Blood	Chennai - Avadi	D	590	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
M0070	MTHFR C677T/A1298C Mutations	Chennai - Avadi	D	7670	Wed & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
M0071	Mucopolysaccharides (MPS) Type-I Blood**	Chennai - Avadi	R	5350	Daily : 11:30am	7th day : 6:00pm	10 ml sodium Heparin OR EDTA Blood
M0072	Mucopolysaccharides (MPS) Type-VI Blood**	Chennai - Avadi	R	2245	Daily : 11:30am	7th day : 6:00pm	10 ml sodium Heparin OR EDTA Blood
M0074	Mucopolysaccharides, urine	Chennai - Avadi	D	590	Daily : 11:30am	Next day : 5:00pm	10 ml of Spot urine
M0075	Multiple Sclerosis Profile, Maxi	Chennai - Avadi	D	10030	Tue & Fri : 9:00am	Next day : 7:00pm	2 ml of CSF in leakproof container (fofin vu conical tube is recommended to prevent leakage). 3 ml of Serum.
M0076	Multiple Sclerosis Profile, Mini	Chennai - Avadi	D	6845	Tue & Fri : 9:00am	Next day : 7:00pm	2 ml of CSF in leakproof container (fofin vu conical tube is recommended to prevent leakage). 3 ml of Serum.
M0078	Mumps virus IgG Antibody by CLIA - Serum	Chennai - Avadi	C	1810	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0079	Mumps virus IgM Antibody by CLIA - Serum	Chennai - Avadi	C	1810	Daily : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0080	MUSK Antibody Myasthenia gravis Serum	Chennai - Avadi	C	7005	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0081	Myasthenia Gravis Profile	Chennai - Avadi	D	11680	Mon : 9:00am	5th day : 6:00pm	5 ml of serum
M0086	Mycoplasma Pneumoniae IgG Antibody by EIA - Serum	Chennai - Avadi	C	2315	Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0087	Mycoplasma Pneumoniae IgM Antibody by	Chennai	C	2315	Mon & Thu :	Next day :	3ml of Serum in Gel Vacutainer or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	EIA - Serum	- Avadi			11:30am	5:00pm	in Plain Sterile Vacutainer
M0091	Myoglobin Serum	Chennai - Avadi	C	2375	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0092	Myoglobin Urine Spot	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	After 8hrs	10 ml of Spot urine
M0095	Myotonic Dystrophy Type 1 (DM1) Mutation Analysis	Chennai - Avadi	E	7350	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
M0099	MRSA screen Methicillin Resistant Staphylococcus aureus - PCR	Chennai - Avadi	C	2260	Daily : 9:00am	Same day : 5:00pm	One swab from Nasal, Groin and Axilla.
M0100	Multiple myeloma workup panel	Chennai - Avadi	D	23540	Daily : 9:00am	7th day : 6:00pm	3 ml of EDTA Bone Marrow, 4ml of heparin Bone marrow, Serum, Direct smears
M0101	MPL W515/S505 Mutation Analysis, Blood	Chennai - Avadi	E	6300	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
M0102	MPL W515/S505 Mutation Analysis, Bone marrow	Chennai - Avadi	E	6300	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
M0103	Microsatellite Instability (MSI), HNPCC/ Lynch Synd	Chennai - Avadi	D	16520	Mon, Wed & Fri : 9:00am	10th day : 6:00pm	â€¢ Clinical history with Histopathology report. FFPE Block - Ambient & EDTA Blood
M0104	MYCN Neuroblastoma (n MYC gene amplification)**	Chennai - Avadi	R	13910	Daily : 11:30am	12th Day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0108	Metachromatic Leukodystrophy gene**	Chennai - Avadi	R	21400	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
M0132	MISMATCH REPAIR (MMR) PROTEIN IMMUNOHISTOCHEMISTRY ONLY	Chennai - Avadi	D	10030	Daily : 9:00am	5th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
M0136	Marfan s syndrome By NGS**	Chennai - Avadi	R	22470	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
M0140	Muscular Dystrophy & Congenital Myopathy gene panel**	Chennai - Avadi	R	17000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
M0141	Methylmalonic acids Quantitative Serum**	Chennai	R	3210	Thu : 11:30am	4th day :	3ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0142	MET Amplification FISH	- Avadi Chennai - Avadi	E	10500	Daily : 9:00am	6:00pm 7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
M0145	Maturity-Onset Diabetes of the Young (MODY) & neonatal diabetes gene panel**	Chennai - Avadi	R	17000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
M0147	MLH1 Promoter Hypermethylation	Chennai - Avadi	E	6825	Mon, Wed & Fri : 9:00am	5th day : 6:00pm	â€¢ Tumor containing FFPE block (>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable.
M0149	MSH2 Sequencing and deletion duplication analysis	Chennai - Avadi	E	18900	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0150	MMR IHC Reflex to BRAF CD600 and MHL1 promoter Hypermethylation	Chennai - Avadi	D	12980	Daily : 9:00am	After 2 weeks : 6:00pm	(FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0151	Metanephrine- Fractionated Plasma	Chennai - Avadi	D	8260	Tue & Fri : 11:30am	Next day : 5:00pm	Clinical Diagnosis with Histopathology report, HistologyandIHC.pdf Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
M0153	Combine screening of Pre-eclampsia & Trisomies	Chennai - Avadi	D	4130	Daily cut off 7:30 AM	same day 9:00 PM	3 ml of serum is required. Mandatory Requirements: Latest Ultrasound (USG) report with a completely filled Maternal screening TRF form including CRL/BPD, nasal bone status, ductus venosus status, and NT value in millimeters. Please note that this test is not valid for triplets and above Preeclampsia: USG findings of mean arterial pressure (MAP) and uterine artery pulsatility index (UAPI), the patients blood pressure (BP), and height should be included as per maternal TRF requirements
M0154	Meningitis/Encephalitis panel by FilmArray (Biofire) - CSF	Chennai - Avadi	E	19425	Daily : 9:00am to 9:00pm	After 10hrs	2ml of CSF in Sterile leak proof tube or container
M0160	MMR IHC Reflex to MSI PCR	Chennai - Avadi	D	9320	Daily : 9:00am	After 2 weeks : 6:00pm	(FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical Diagnosis with Histopathology report, HistologyandIHC.pdf
M0161	Myositis Profile-IgG(16 Antigen)	Chennai - Avadi	D	11680	Daily: 11:30am	Next day 8pm	4 ml of serum
M0163	Magnesium, Hair tests by ICPMS**	Chennai - Avadi	R	5885	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0164	Manganese, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	before sending 1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
M0165	Manganese, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
M0166	Marijuana confirmation by LCMS, Spot Urine**	Chennai - Avadi	R	9630	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
M0167	Mercury, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
M0169	Metals profile-1, EDTA whole blood **	Chennai - Avadi	D	7080	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0170	Metals profile-1, Spot urine	Chennai - Avadi	D	7080	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0172	Metals profile-2, Spot urine**	Chennai - Avadi	R	12305	Daily : 11:30am	10th day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0174	Metals profile-3, Spot urine**	Chennai - Avadi	R	16585	Daily : 11:30am	10th day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0175	Methadone confirmation by LCMS, Spot Urine**	Chennai - Avadi	R	9630	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0176	Methaqualone confirmation by LCMS, Spot Urine**	Chennai - Avadi	R	9630	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
M0178	Minerals profile, Hair**	Chennai - Avadi	R	12305	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
M0179	Minerals profile, Spot Urine**	Chennai - Avadi	R	12305	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
M0181	Molybdenum, Hair tests by ICPMS**	Chennai - Avadi	R	5885	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
M0182	Molybdenum, Spot Urine**	Chennai - Avadi	R	3640	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
M0183	Mucorales (Mucormycosis) DNA Detection by Real Time PCR - Body fluids/CSF/EDTA Blood/Pure culture/Sputum/Tissue	Chennai - Avadi	C	6780	Daily : 7:30am	3rd Day : 6:00pm	3ml of Body fluids/ CSF/ EDTA Blood/ Pure culture/ Sputum/ Tissue in Sterile leak proof container
M0193	MetroHealth Instant Glucose Fasting	Chennai - Avadi	L0	100	Daily : 9:00am to 9:00pm	After 5 mins	Whole Blood
M0194	MetroHealth Instant Glucose Post Prandial	Chennai - Avadi	L0	100	Daily : 9:00am to 9:00pm	After 5 mins	Whole Blood
M0195	MetroHealth Instant Glucose Random	Chennai - Avadi	L0	100	Daily : 9:00am to 9:00pm	After 5 mins	Whole Blood
M0196	MetroHealth Instant Haemoglobin	Chennai - Avadi	L0	100	Daily : 9:00am to 9:00pm	After 5 mins	Whole Blood
M0197	MetroHealth Instant Urine	Chennai - Avadi	L0	150	Daily : 9:00am to 9:00pm	After 5 mins	Urine
M0823	Macro TSH (TSH, Free TSH,% Recovery) Serum	Chennai - Avadi	C	565	Daily : 9:00am to 9:00pm	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tests.
M0824	Anti diabetic drug contents	Chennai - Avadi	R	5885	Daily : 11:30am	10th day : 6:00pm	Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clinical symptoms if any. Reported in 7 days
M0825	Medicine test profile (steroids, heavy metals, anti diabetic, NSAIDS)	Chennai - Avadi	R	21400	Daily : 11:30am	10th day : 6:00pm	Medicine sample (10 grm), Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clinical symptoms if any. Reported in 7 days
M0826	NSAID contents	Chennai - Avadi	R	3210	Daily : 11:30am	10th day : 6:00pm	Medicine sample (10 grm), Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clinical symptoms if any. Reported in 7 days
M0829	Monogenic and syndromic obesity-Gene Panel	Chennai - Avadi	E	21000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0830	MTB Antibiogram (TB NGS) by Next Generation Sequencing - Pure Culture	Chennai - Avadi	R	9095	Daily : 11:30am	15th day : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
M0831	MPA (Mycophenolic Acid)	Chennai - Avadi	C	2545	Mon & Thu : 9:00am	Same day : 5:00pm	require serum (3ml) sample collected in fasting status preferably but not mandatory.
M5085	MPN Comprehensive Panel (Blood)	Chennai - Avadi	D	15930	Tue & Fri : 9:00am	7th day : 6:00pm	3ml whole blood in EDTA Vacutainer & 3ml Heparin blood
M5086	MPN Comprehensive Panel (Bone Marrow)	Chennai	D	15930	Tue & Fri : 9:00am	7th day :	3ml Bone Marrow in EDTA

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M5097	MTB triple combo (For pulmonary samples)	- Avadi Chennai - Avadi	D	7905	Daily : 7:30am	6:00pm After 6 weeks : 6:00pm	Vacutainer & 3ml Bone Marrow 2ml of Pulmonary specimen in Sterile leak proof container
M5098	MTB triple combo (For extra pulmonary samples)	Chennai - Avadi	D	7905	Daily : 7:30am	After 6 weeks : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
M5768	MetroHealth Instant 5 min Screen Pkg	Chennai - Avadi	L0	350	Daily : 9:00am to 9:00pm	After 5 mins	Whole Blood-Finger Prick , 30 to 60 ml of Urine in a leak proof container.
M8027	Myotonic Dystrophy Type 2 (3q 21)	Chennai - Avadi	R	9630	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
M8228	MIC OF COLISTIN BY MICRO BROTH DILUTION	Chennai - Avadi	C	565	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	Pure culture isolate is required with identification upto species level.
M8268	MTB Xpert XDR Rapid (For Pulmonary Samples)	Chennai - Avadi	D	4340	Daily : 7:30am	Same day : 6:00pm	2ml of Pulmonary specimen in Sterile leak proof container
M8269	MTB Xpert XDR Rapid (For Extra pulmonary Samples)	Chennai - Avadi	D	4340	Daily : 7:30am	Same day : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
M8346	MELD package	Chennai - Avadi	R	1285	Daily 9.00am to 9.00pm	Same day at 5.00pm	
N0001	N- Acetyl Glucosamine 6 Sulphate Sulphatase Blood**	Chennai - Avadi	R	9095	Daily : 11:30am	After 4 weeks :	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
						6:00pm	of collection, Avoid collection on weekend & Holidays)
N0002	NARP Neurogenic Ataxia Retinitis Pigmentosa Blood**	Chennai - Avadi	R	6955	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
N0003	Natural Killer cells	Chennai - Avadi	D	4600	Daily : 11:30am	Next day : 5:00pm	6 ml EDTA whole Blood
N0006	Neisseria Gonorrhoeae Culture & Identification by Vitek/ MALDITOF - Urethral discharge	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	4th day : 6:00pm	2 ml of Urethral Discharge in Sterile leak proof container
N0007	Neisseria Gonorrhoeae Culture & Identification Vitek/ MALDITOF - Vaginal discharge	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	4th day : 6:00pm	2 ml of Vaginal Discharge in Sterile leak proof container
N0008	Neisseria Gonorrhoea DNA Detection by Real Time PCR - Genital swab	Chennai - Avadi	C	3615	Wed & Sat : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
N0009	Neisseria Gonorrhoea DNA Detection by Real Time PCR - Urethral discharge	Chennai - Avadi	C	3615	Wed & Sat : 9:00am	Next day : 5:00pm	2ml of Urethral Discharge in Sterile leak proof container
N0010	Neisseria Gonorrhoea DNA Detection by Real Time PCR - Urine	Chennai - Avadi	C	3615	Wed & Sat : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
N0011	Neisseria Meningitidis A,C,Y W135 Antigen detection by Agglutination - CSF	Chennai - Avadi	C	2710	Daily : 9:00am & 3:00pm	After 4hrs	2ml of CSF in Sterile leak proof tube or container
N0012	Neoxpert IAP (Newborn Screening 3 Analytes)	Chennai - Avadi	C	960	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
N0013	Neoxpert 5 (Newborn Screening 5 Analytes)	Chennai - Avadi	C	1750	Daily 9:00am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
N0014	Neoxpert 8 (Newborn Screening 8	Chennai	C	2825	Daily 9:00am	Next day :	Dried blood spot: 1 drop of heel

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Analytes)	- Avadi				5:00pm	prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
N0015	Neoxpert by TMS (Newborn Screening 40+ Analytes)	Chennai - Avadi	C	5535	Daily : 11:30am	5th day : 6:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
N0017	Paraneoplastic syndrome/Neuronal Antibody Profile (IgG)	Chennai - Avadi	E	8400	Daily: 11:30am	Next day 8pm	3 ml of serum
N0018	Nickel EDTA Whole Blood	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
N0019	Nickel Serum	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
N0020	Nickel, Urine Spot	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
N0021	Nicotine Metabolite Quantitative Serum	Chennai - Avadi	C	1470	Daily : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0022	Nicotine Metabolite Quantitative Urine Spot	Chennai - Avadi	C	1470	Daily : 9:00am	Same day : 5:00pm	2 ml of Spot Urine
N0023_NIPS	NEXTGEN NON-INVASIVE PRENATAL SCREENING TEST	Chennai - Avadi	E	16800	Daily : 9:00am	8th day : 6:00pm	Sample collection require in special tube - Streak Tube - 10 ml. We Required filled TRF which include clinician stamp, Ultra Sound report (USG), Dual and Quatraple Marker Report, Photo Identification Proof.
N0027	NMDA Receptor Antibody (NR1), CSF	Chennai - Avadi	D	7080	Mon & Thu : 11:30am	Same day : 9:00pm	2 ml of CSF in sterile screw capped container.
N0028	NMDA Receptor Antibody (NR1) Autoimmune encephalitis Serum	Chennai - Avadi	D	7080	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0029	NMO (Aquaporin 4)Neuromyelitis Optica Antibodies Serum	Chennai - Avadi	C	4860	Tue & Fri: 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0030	NMO (Aquaporin 4)Neuromyelitis Optica Antibodies CSF	Chennai - Avadi	C	4860	Tue & Fri: 11:30am	Same day : 9:00pm	2 ml of CSF in sterile screw capped container.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
N0032	Nocardia Detection by smear examination - Sputum	Chennai - Avadi	B	540	Daily : 9:00am & 3:00pm	After 6hrs	2ml of Sputum in Sterile leak proof container
N0033	Nocardia Detection by smear examination - Tissue	Chennai - Avadi	B	540	Daily : 9:00am & 3:00pm	After 6hrs	Tissue with Normal Saline in Sterile leak proof container
N0034	Nor-Adrenaline (Nor-Epinephrine) Plasma	Chennai - Avadi	D	3420	Wed : 11:30am	Next day : 7:00pm	Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
N0035	Nor-Adrenaline (Nor-Epinephrine) Urine-HPLC	Chennai - Avadi	D	3420	Wed : 11:30am	Next day : 7:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
N0036	Nor-Metanephrine-Free Plasma	Chennai - Avadi	D	4720	Tue & Fri : 11:30am	Next day : 7:00pm	Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
N0037	Nor-Metanephrine, HPLC, urine 24 hrs	Chennai - Avadi	D	2360	Mon & Thu : 11:30am	Next day : 7:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
N0039	NPM1 mutations AML - Blood	Chennai - Avadi	D	4600	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
N0040	NPM1 mutations AML - Bone Marrow	Chennai - Avadi	D	4600	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
N0041	NRAS mutations Exon-2 & 3 FFPE Tissue	Chennai - Avadi	D	7670	Mon, Wed & Fri : 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>10% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
N0042	NSE Neuron Specific Enolase CSF	Chennai - Avadi	D	3775	Tue & Fri: 11:30am	After 8hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
N0043	NSE Neuron Specific Enolase Serum	Chennai - Avadi	D	3540	Tue & Fri: 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0044	NT-Pro BNP Serum	Chennai - Avadi	C	3000	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0045	NTxN-Tellopeptide (Collagen Crosslinked)Urine 24H	Chennai - Avadi	C	3615	Wed : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
N0047	NMDA- CSF & Serum	Chennai - Avadi	D	12980	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests. 2 ml of CSF in sterile screw capped container.
N0048	NMO- CSF & Serum	Chennai - Avadi	D	8850	Tue & Fri: 11:30am	Same day : 9:00pm	2 ml of CSF in sterile screw capped container.
N0056	Paraneoplastic syndrome/Neuronal Antibody Profile-3	Chennai - Avadi	E	13650	Daily: 11:30am	Next day 8pm	3 ml of serum
N0057	Neoxpert Expanded Screening Panel (Newborn Screening 48+ Analytes)	Chennai - Avadi	D	7080	Daily : 11:30am	5th day : 6:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
N0058	Comprehensive Neurology panel**	Chennai - Avadi	R	23540	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
N0059_Serum	NMO with MOG Antibody Profile	Chennai - Avadi	D	7670	Daily : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0060_CSF	NMO with MOG Antibody Profile	Chennai - Avadi	D	7670	Daily : 11:30am	Same day : 9:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
N0061	Nocardia Susceptibility by Microbroth dilution - Pure Culture	Chennai - Avadi	C	3675	Daily : 9:00am	15th day : 6:00pm	Pure growth of Nocardia spp. on Blood/Chocolate Agar plate. Kindly provide Previous identification report of Nocardia.
N0062	Nitrosonaphthol Test Qualitative, Spot Urine	Chennai - Avadi	R	21400	Daily : 11:30am	5th day : 6:00pm	10 ml of Spot sample in Leakproof container. Many analytes like VMA, 5 HIAA, adrenaline may require

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
N0064	Nickel, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
N0065	Nickel, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
N0066	NR4A3 (TEC) 9p22.3 translocation**	Chennai - Avadi	R	12840	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
N0067	NTRK3 and ETV6 translocation**	Chennai - Avadi	R	14980	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
N0068	Neurology comprehensive panel & whole mitochondrial genome sequencing	Chennai - Avadi	E	26250	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
N0069	Neurotransmitter disorders-Gene Panel	Chennai - Avadi	E	21000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
N0070	New born screening panel by NGS	Chennai - Avadi	E	16800	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
N0096_C	Integrated test, standard	Chennai - Avadi	E	5775	Daily : 11:30am	7th day : 6:00pm	3 ml of Serum for all test. Serum to be shipped at 2 to 8 degree Refrigerated. Serum sample to be stored at 2 to 8 degree for 2 days and post 2 days to be send in frozen condition. Latest USG is recommended mentioning CRL or BPD, Single or twin fetus (test not valid for triplets or above), NT value in mm. Maternal screening TRF is mandatory with history of date of birth, weight, Race, Smoking, type 1 diabetes. In IVF pregnancy Egg donor DOB, egg extraction date and embryo transfer date. in Preeclampsia cases history of BP or family histroy additionally required with USG findings of MAP, UAPI values etc as per TRF requirements. Reflex test will be done for positive screen cases and will be conveyed by lab for sample requirements further based on consent with Gynecologist or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							pateint and TAT of report will be based on the test added and registered. For NIPS 10 ml of Whole blood in strek tube to be send in ambinet condition upto 1 day or refrigerated conditio upto 3 days.
N0096_S	Integrated test, sequential	Chennai - Avadi	E	5775	Daily : 11:30am	7th day : 6:00pm	3 ml of Serum for all test. Serum to be shipped at 2 to 8 degree Refrigerated. Serum sample to be stored at 2 to 8 degree for 2 days and post 2 days to be send in frozen condition. Latest USG is recommended mentioning CRL or BPD, Single or twin fetus (test not valid for triplets or above), NT value in mm. Maternal screening TRF is mandatory with history of date of birth, weight, Race, Smoking, type 1 diabetes. In IVF pregnancy Egg donor DOB, egg extraction date and embryo transfer date. in Preeclampsia cases history of BP or family histroy additionally required with USG findings of MAP, UAPI values etc as per TRF requirements. Reflex test will be done for positive screen cases and will be conveyed by lab for sample requirements further based on consent with Gynecologist or pateint and TAT of report will be based on the test added and registered. For NIPS 10 ml of Whole blood in strek tube to be send in ambinet condition upto 1 day or refrigerated conditio upto 3 days.
N5543	NextGen NIPS Advanced	Chennai - Avadi	E	19425	Daily : 9:00am	8th day : 6:00pm	Sample collection require in special tube - Streak Tube - 10 ml. We Required filled TRF which

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							include clinician stamp, Ultra Sound report (USG), Dual and Quatraple Marker Report, Photo Identification Proof.
N8209	NextGen Precision Assay (500 genes +TMB+MSI)**	Chennai - Avadi	R	65000	11.30 am	20 day 6:00pm	EDTA SAMPLE + Paraffin Block
O0001	Obesity Profile	Chennai - Avadi	C	6665	Tue : 7:30am	Next day : 5:00pm	5 ml of serum, 2 ml of flouride plasma , 20 ml of spot urine
O0002	Occult blood by chemical test - Stool	Chennai - Avadi	A	165	Daily : 9:00am to 9:00pm	After 4hrs	5gm of Stool specimen in Sterile leak proof container
O0003	Occult blood by chemical test - Urine	Chennai - Avadi	A	110	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
O0004	Oligoclonal band,CSF	Chennai - Avadi	C	4860	Tue & Fri: 11:30am	Next day : 5:00pm	3 ml of serum & 2 ml of CSF in sterile screw capped container.
O0005	Opiates Morphine Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	5 ml of Spot Urine
O0006	Opportunistic infection panel by Modified ZN smear - Stool	Chennai - Avadi	C	1580	Daily : 9:00am & 3:00pm	After 6hrs	5gm of Stool specimen in Sterile leak proof container
O0008	Ornithine Quantitative Plasma	Chennai - Avadi	R	6365	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
O0009	Ornithine Quantitative Urine 24H**	Chennai - Avadi	R	6365	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
O0010	Osmolality freezing point Serum	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
O0011	Osmolality freezing point Urine	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	After 4hrs	5 ml of Spot Urine
O0012	Osmotic Fragility Test	Chennai - Avadi	C	850	Daily : 9:00am	Same day : 5:00pm	3 ml Sodium heparin sample
O0013	Osteocalcin BGP-Bone Gla Protein Serum	Chennai - Avadi	C	3110	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
O0014	Osteomon Profile	Chennai - Avadi	D	3305	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum and 3 ml of EDTA plasma
O0015	Osteoporosis Profile-I	Chennai - Avadi	C	2090	Daily : 9:00am to 9:00pm	Next day : 5:00pm	2 ml of EDTA Plasma

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
00016	Osteoporosis Profile-II	Chennai - Avadi	C	6440	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum, 2 ml of EDTA plasma and 10 ml of spot urine
00017	Osteoporosis Profile-III	Chennai - Avadi	C	8925	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum, 2 ml of EDTA plasma and 10 ml of spot urine
00018	OT Microbiological Surveillance Panel Bacterial & Fungal by Vitek/ MALDITOF - Per swab	Chennai - Avadi	C	620	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
00020	Ovarian Antibody AOA serum	Chennai - Avadi	C	3730	Mon : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
00021	Ovarian Marker Profile	Chennai - Avadi	D	4600	Daily : 9:00am to 9:00pm	Next day : 5:00pm	5 ml of serum
00022	Ovarian Reserve Profile-Maxi	Chennai - Avadi	C	6665	Mon : 9:00am	Next day : 5:00pm	5 ml of serum
00023	Ovarian Reserve Profile-Mini	Chennai - Avadi	C	3615	Mon & Thu : 9:00am	Same day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
00024	Oxalate, Urine	Chennai - Avadi	C	2090	Daily : 11:30am	Same day : 9:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
00025	OT Microbiological Surveillance Panel Bacterial & Fungal - 5 swabs	Chennai - Avadi	C	2545	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
00026	Oxcarbazepine, Serum**	Chennai	D	5310	Fri : 11:30am	Next day :	3 ml of serum(sample collect in,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
00027	Ova and Parasite Examination,Other Samples	- Avadi Chennai	B	325	Daily : 9:00am & 3:00pm	5:00pm Next day : 6:00pm	plain vacutainer, not in gel tube) 5gm of Stool specimen in Sterile leak proof container
00039	Open Neural tube defect screening, AFP MOM, Second trimester - Serum	- Avadi Chennai	C	1075	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
00040	Opiates confirmation by LCMS, Spot Urine**	- Avadi Chennai	R	9630	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
00041	Oncotrack - Ultima [Liquid biopsy for 56 theranostic genes]**	- Avadi Chennai	R	49400	Daily : 11:30am	After 4 weeks : 6:00pm	Peripheral blood in 2 Streck tubes , (tube supplied upon request)
00042	NextGen Oncomine Focus Assay - 52 Gene Panel	- Avadi Chennai	E	39500	Daily 9:00 am	4 Weeks	Tumor containing FFPE block (>10% of tumor). Clinical history with Histopathology report.
01801	Other diagnostic services- 2 D Echo	- Avadi Chennai	O	2060			
01806	Other diagnostic services- ECG	- Avadi Chennai	O	300			
01810	Other diagnostic services- TMT (Stress test)	- Avadi Chennai	O	1545			
01818	OPG - Orthopantomogram	- Avadi Chennai	O	515			
01901	CANINE-Occult blood Stool	- Avadi Chennai	A	130	Daily: 9:00am, 2:00pm & 8:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
O1902	Occult blood	Chennai - Avadi	A	75	Daily: 9:00am, 2:00pm & 8:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P0001	P1NP Intact N-Terminal Propeptide of Type I Procollagen Serum	Chennai - Avadi	C	1865	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0002	Pacer-26	Chennai - Avadi	C	1920	Daily : 9:00am to 9:00pm	After 8hrs	5 ml of serum and 2 ml of F flouride plasma
P0003	Palmitoyl Protein Thioesterase1	Chennai - Avadi	R	5990	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
P0004	Pancreatic (Acute) Profile	Chennai - Avadi	D	1770	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum and 2 ml of EDTA whole blood
P0005	Pancreatic Marker Profile	Chennai - Avadi	D	4130	Daily : 9:00am to 9:00pm	After 8hrs	5 ml of serum
P0006	Panfungal-DNA detection by PCR & Sanger Sequencing - EDTA Blood	Chennai - Avadi	D	10620	Daily : 7:30am	3rd Day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer
P0007	Panfungal-DNA detection by PCR & Sanger Sequencing - Body fluids	Chennai - Avadi	D	10620	Daily : 7:30am	3rd Day : 6:00pm	2ml of Body Fluid in Sterile leak proof container
P0008	Panfungal-DNA detection by PCR & Sanger Sequencing - CSF	Chennai - Avadi	D	10620	Daily : 7:30am	3rd Day : 6:00pm	2ml of CSF in Sterile leak proof tube or container
P0009	Panfungal-DNA detection by PCR & Sanger Sequencing - Sputum	Chennai - Avadi	D	10620	Daily : 7:30am	3rd Day : 6:00pm	2ml of Sputum in Sterile leak proof container
P0010	Panfungal-DNA detection by PCR & Sanger Sequencing - Tissue	Chennai - Avadi	D	10620	Daily : 7:30am	3rd Day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
P0011	PAP smear by Conventional method	Chennai - Avadi	C	790	Daily : 9:00am	Next day : 6:00pm	Pre fixed smears along with P/S findings
P0012	PAP smear by Liquid based cytology (LBC)	Chennai - Avadi	C	1350	Daily : 9:00am	3rd day : 6:00pm	Cervical scraping in a LBC vial
P0013	PAPPa-Pregnancy Associated Plasma Protein Serum	Chennai - Avadi	D	1535	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P0014	Paracetamol Acetaminophen Serum	Chennai - Avadi	D	1300	Daily : 9:00am to 9:00pm	After 6hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0015	Parainfluenza 1,2 & 3 IgG Antibody by EIA - Serum	Chennai - Avadi	C	1920	Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0016	Parainfluenza 1,2 & 3 IgM Antibody by EIA - Serum	Chennai - Avadi	C	1920	Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0017	Parathyroid Panel	Chennai - Avadi	C	3050	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of EDTA plasma, 3 ml of serum and spot sample urine
P0018	Parietal Cell Antibody GPCA by IFA Serum	Chennai - Avadi	C	2375	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0019	Parvovirus B19 DNA Detection by Real Time PCR - Amniotic fluid	Chennai - Avadi	D	3365	Mon & Thu : 9:00am	Next day : 5:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
P0020	Parvovirus B19 DNA Detection by Real Time PCR - EDTA Blood	Chennai - Avadi	D	3365	Mon & Thu : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
P0021	Parvovirus B19 DNA Detection by Real Time PCR - CSF	Chennai - Avadi	D	3365	Mon & Thu : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
P0022	Parvovirus B19 IgG Antibody by EIA - Serum	Chennai - Avadi	C	2710	Wed : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0023	Parvovirus B19 IgM Antibody by EIA - Serum	Chennai - Avadi	C	2710	Wed & Sat : 1:00pm	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0024	PAS Stain, Blood	Chennai - Avadi	C	790	Mon & Thu : 11:30am	3rd day : 6:00pm	3 Direct Peripheral smear slides, EDTA blood. Clinical history and clinician no is mandatory
P0025	PAS Stain, Bone Marrow	Chennai - Avadi	C	790	Mon & Thu : 11:30am	3rd day : 6:00pm	3 Direct Bone Marrow Aspirate slides. Clinical history and clinician no is mandatory
P0026	Paul bunnel test by Agglutination - Serum	Chennai - Avadi	B	1080	Daily : 9:00am & 3:00pm	3rd day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0028	PCP Phencyclidine Phosphate Urine Spot	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of Spot Urine
P0032	Peripheral smear examination (PBS)	Chennai - Avadi	B	320	Daily : 9:00am & 3:00pm	After 4hrs	3 ml EDTA Whole Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P0033	pH by Indicator Strip - Stool	Chennai - Avadi	A	110	Daily : 9:00am to 9:00pm	After 4hrs	5gm of Stool specimen in Sterile leak proof container
P0034	pH by Indicator Strip - Urine	Chennai - Avadi	A	110	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
P0035	Phenobarbitone Gardinal Serum**	Chennai - Avadi	C	1075	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0036	Phenyl Alanine Neonatal Screen Dried blood spot	Chennai - Avadi	B	485	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
P0037	Phenylalanine Quantitative Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
P0038	Phenylalanine Quantitative Spot Urine	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
P0039	Phenytoin Eptoin, Dilantin Serum	Chennai - Avadi	C	1075	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P0040	Pheochromocytoma Marker Profile-I	Chennai - Avadi	D	8140	Daily : 9:00am	7th day : 6:00pm	20 ml of 24 hrs urine in HCL preservative. Quote Volume.
P0041	Pheochromocytoma Marker Profile-II	Chennai - Avadi	D	19470	Daily : 9:00am	7th day : 6:00pm	4 ml of EDTA plasma and 20 ml of urine (Urine 24H)
P0050	Phospholipid profile IgG	Chennai - Avadi	C	5595	Thu: 1:00pm	Next day : 5:00pm	5 ml of serum
P0051	Phospholipid profile IgM	Chennai - Avadi	C	5595	Thu: 1:00pm	Next day : 5:00pm	5 ml of serum
P0052	Phospholipid Syndrome Profile	Chennai - Avadi	D	4600	Tue, Thu & Sat: 7:30am	Same day : 7:00pm	3 ml of Citrate plasma, 3 ml of serum
P0053	Phosphorus-Inorganic Serum	Chennai - Avadi	A	195	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0054	Phosphorus-Inorganic Urine Spot	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
P0055	Phosphorus-Inorganic Urine 24H	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
P0056	Phosphotidylserine IgG antibodySerum	Chennai - Avadi	C	1245	Thu: 1:00pm	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0057	Phosphotidylserine IgM antibodySerum	Chennai - Avadi	C	1245	Thu: 1:00pm	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P0058	PIK3CA Gene Mutations - FFPE Tissue	Chennai - Avadi	E	5985	Mon, Wed & Fri : 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
P0059	Pituitary Marker Profile	Chennai - Avadi	C	4070	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of EDTA plasma, 3 ml of serum
P0060	PLA2 receptor antibody Primary membranous Glomerulopathy Serum	Chennai - Avadi	D	5310	Tue : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0061	Plate/Slant For Identification & Susceptibility - Automated Plate/Slant	Chennai - Avadi	C	905	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	Pure culture isolate is required
P0062	Platelet Antibodies by IFA Serum	Chennai - Avadi	D	7080	Mon : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0063	Platelet Count ,EDTA Blood	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	3 ml EDTA Whole Blood
P0066	Pneumocystis Carinii Detection by IFA - BAL / Sputum	Chennai - Avadi	C	2825	Daily : 9:00am & 3:00pm	Next day : 6:00pm	2 ml of BAL / Sputum in Sterile leak proof container
P0070	Pneumonia Viral-III Panel by EIA - Serum	Chennai - Avadi	D	7670	Thu : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0071	PNH confirmation test, Blood	Chennai - Avadi	D	7435	Daily : 9:00am	Next day : 5:00pm	3 ml EDTA whole Blood
P0072	PNMA2 (MA2/TA) Antibody Serum	Chennai - Avadi	D	5310	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0073	Synovial Fluid by Polarizing microscopy	Chennai - Avadi	C	1015	Daily : 9:00am	4th day : 6:00pm	0.5 ml Synovial fluid in a plain container
P0074	Pompe Disease Enzyme Assay Blood**	Chennai - Avadi	R	8240	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
P0075	Porphobilinogen (Quantitative), urine 24 Hrs	Chennai - Avadi	C	3955	Thu : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
P0076	Porphyria Profile, 24 hrs Urine	Chennai - Avadi	D	6490	Thu : 11:30am	Next day : 5:00pm	Add 10 ml of 6 M HCL preservative in 24 hrs urine container. Collect 24 hrs urine sample and send 10 ml of 24 Hrs Urine In Dark Container.
P0077	Porphyria Profile, Urine Spot	Chennai - Avadi	D	6490	Thu : 9:00am	Next day : 5:00pm	10 ml of 24 Hours Urine OR spot urine In Dark Container
P0078	Potassium ISE Serum	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0079	Potassium Urine 24H	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea ,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
P0080	Potassium, urine by ISE_P0080	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
P0081	PR3 - c-ANCA Serum	Chennai - Avadi	C	1525	Tue, Thu & Sat: 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0082	PRA- Plasma Renin Activity	Chennai - Avadi	C	4520	Tue & Fri: 11:30am	Next day : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
P0083	Pre Albumin Serum	Chennai - Avadi	D	3540	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0086	Pre-Operative Profile (Maxi)	Chennai - Avadi	C	1865	Daily : 9:00am	Next day : 5:00pm	3 ml of EDTA whole blood, 2 ml of citrate plasma, 3 ml of serum, 20 ml of spot urine
P0087	Pre-Operative Profile (Mini)	Chennai - Avadi	C	1075	Daily : 9:00am	Next day : 5:00pm	3 ml of EDTA whole blood, 3 ml of serum, 20 ml of spot urine
P0088	Procalcitonin Serum	Chennai - Avadi	C	2485	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0089	Progesterone (P4) Serum	Chennai - Avadi	B	640	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0090	Prolactin Serum	Chennai - Avadi	B	540	Daily : 9:00am to 9:00pm	After 6hrs	4 ml of Serum
P0091	Propoxyphene Qualitative Urine Spot	Chennai - Avadi	C	735	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
P0092	Prostate Profile	Chennai - Avadi	D	2240	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0094	Protein C Activity Functional ,Citratd plasma	Chennai - Avadi	C	4405	Mon & Thu : 9:00am	Next day : 1:00pm	1 ml of citrated platelet poor plasma in frozen condition.
P0095	Protein C Antigen Citrated plasma (P0095)	Chennai - Avadi	C	5085	Daily: 7:30am	Same day : 9:00pm	1 ml of citrated platelet poor plasma in frozen condition.
P0096	Serum, Protein Electrophoresis	Chennai	C	780	Daily:9.00 am	Same day	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				8:00 PM	ml serum for profiles and reflex tests.
P0097	Protein Electrophoresis, Urine 24 Hrs	Chennai - Avadi	C	960	Mon & Thu : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
P0098	Protein Electrophoresis reflex immunofixation Serum	Chennai - Avadi	C	2600	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0099	Protein S activity, Plasma	Chennai - Avadi	C	4070	Mon & Thu : 9:00am	Next day : 1:00pm	1 ml of citrated platelet poor plasma in frozen condition.
P0100	Protein S Antigen (Free)Citrated plasma	Chennai - Avadi	C	6215	Mon & Thu : 9:00am	Next day : 1:00pm	1 ml of citrated platelet poor plasma in frozen condition.
P0101	Proteins Serum	Chennai - Avadi	B	215	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0102	Proteins Urine 24H	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
P0103	Proteins Ascitic Fluid	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Ascitic Fluid
P0104	Proteins CSF	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
P0105	Proteins Pedicardial Fluid	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Pericardial Fluid
P0106	Proteins Pleural Fluid	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Pleural Fluid
P0107	Proteins Synovial Fluid	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Synovial Fluid
P0108	Proteins/Creatinine Ratio	Chennai - Avadi	B	560	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
P0109	Prothrombin G20210A (Factor II) Mutation	Chennai - Avadi	D	5900	Wed & Fri: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius .
P0110	Prothrombin Time (PT)	Chennai - Avadi	B	280	Daily : 9:00am to 9:00pm	After 4hrs	Local Labs can send whole Blood within 24 hrs of collection and for Outstations Frozen Plasma required.
P0111	PSA-Free Prostate Specific Antigen Serum	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0112	PSA-Total Prostate Specific Antigen Serum	Chennai - Avadi	B	810	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0113	Pseudo Cholinesterase OPC poisoning Serum	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P0114	PTH-Intact Molecule Parathyroid hormone Serum/Plasma	Chennai - Avadi	B	1610	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of Serum / EDTA Plasma
P0115	PUO Profile	Chennai - Avadi	C	5650	Daily : 9:00am to 9:00pm	After 6hrs	2 Blood culture bottle (Aerobic and anaerobic), 3 ml of serum
P0117	Pyruvate by Biochemical Method	Chennai - Avadi	C	2205	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0119	PSA with Free PSA ratio	Chennai - Avadi	C	1695	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0123	Prolactin, Active & macroprolactin serum	Chennai - Avadi	C	1015	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0124	PSA reflex Free PSA, Serum	Chennai - Avadi	C	1185	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0125	Porphobilinogen (Quantitative), Urine Spot	Chennai - Avadi	C	3955	Thu : 11:30am	Next day : 5:00pm	10 ml of Spot urine
P0126	Protein Body Fluid	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 6hrs	2 ml body fluid
P0127	Potassium Body Fluid	Chennai - Avadi	A	295	Daily : 9:00am to 9:00pm	After 4hrs	2 ml body fluid
P0128	Proteins, Peritoneal Fluid	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	2 ml body fluid
P0135	Neurotropic Virus panel by real time PCR	Chennai - Avadi	D	6845	Daily : 9:00am	3rd day : 6:00pm	2ml of CSF in Sterile leak proof tube or container
P0139	Pancreatic Elastase Test**	Chennai - Avadi	D	7670	Wednesday-9.00AM	same day 5.00 PM	100 gms of stool in leak proof container. History of medicines, morning sample preferred.
P0149	PV panel (Jak-2 V617F Reflex to Jak-2 Exon 12 Mutation)-Blood	Chennai - Avadi	D	9440	Tue & Fri : 9:00am	8th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
P0152	Pre-eclampsia screening Serum	Chennai	D	3775	Mon & Thu :	Next day :	3 ml of Serum for all test. Serum to

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			11:30am	6:00pm	be shipped at 2 to 8 degree Refrigerated. Serum sample to be stored at 2 to 8 degree for 2 days and post 2 days to be send in frozen condition. Latest USG is recommended mentioning CRL or BPD, Single or twin fetus (test not valid for triplets or above), NT value in mm. Maternal screening TRF is mandatory with history of date of birth, weight, Race, Smoking, type 1 diabetes. In IVF pregnancy Egg donor DOB, egg extraction date and embryo transfer date. in Preeclampsia cases history of BP or family histroy additionally required with USG findings of MAP, UAPI values etc as per TRF requirements. Reflex test will be done for positive screen cases and will be conveyed by lab for sample requirements further based on consent with Gynecologist or pateint and TAT of report will be based on the test added and registered. For NIPS 10 ml of Whole blood in strek tube to be send in ambinet condition upto 1 day or refrigerated conditio upto 3 days.
P0154	Posaconazole Serum	Chennai - Avadi	R	6420	Mon & Thu : 11:30am	Next day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0155	Plasma Acylcarnitine profile**	Chennai - Avadi	R	5885	Daily : 11:30am	7th day : 6:00pm	3 ml of separated EDTA plasma required . Sample to be shipped at Frozen condition . Sample stable at refrigerated condition 2 to 8 degree for 1 day. Clinical history, metabolic disorders or inborn errors of metabolsim suspicion or provisional diagnosis of clincian

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							required fasting recommended not mandatory. Amino acid supplementation or medications needs to stoped based on physician recommendations
P0156	PV panel (Jak-2 V617F Reflex to Jak-2 Exon 12 Mutation)-Bone Marrow	Chennai - Avadi	D	9440	Tue & Fri : 9:00am	8th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
P0158	PML/RARA - t(15;17) Qualitative by PCR*	Chennai - Avadi	D	5780	Daily : 9:00am	4th day : 6:00pm	5-10 ml of EDTA whole blood required. Sample stable for 3 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. Sample should reach to the lab within 48 to 72 hours.
P0159	Phencyclidine confirmation by LCMS, Spot Urine**	Chennai - Avadi	R	9630	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
P0160	Platinum, EDTA Whole Blood**	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
P0161	Platinum, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
P0162	Platinum, Spot Urine**	Chennai - Avadi	R	3640	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
P0163	Propoxyphine confirmation by LCMS, Spot	Chennai	R	9630	Daily : 11:30am	10th day :	10ml of spot urine sample required

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Urine**	- Avadi				6:00pm	to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
P0165	POLG-Gene Analysis	Chennai - Avadi	E	17850	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
P0166	Pulmonary surfactant metabolism dysfunction-Gene Panel	Chennai - Avadi	E	21000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
P0167	PGS (Preimplantation genetic screening), Embryo	Chennai - Avadi	R	23540	Daily : 11:30am	15th day : 6:00pm	CytogeneticsandFISH.pdf
P0168	Pre screening for Bleeding Disorders	Chennai - Avadi	C	4745	Daily : 9:00am	Next day : 5:00pm	1 ml of citrated platelet poor plasma in frozen condition.
P0169	Panfungal-DNA detection by PCR & Sanger Sequencing - Pure Culture	Chennai - Avadi	D	10620	Daily : 7:30am	3rd Day : 6:00pm	Pure culture on SDA plate/ slant
P1002	Plasminogen Activator Inhibitor-1 (PAI-1) / SERPINE-1, 4G/5G Genotyping	Chennai - Avadi	R	14445	Wed : 11:30am	10th day : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.
P1003	PCOS Profile - 1	Chennai - Avadi	D	2715	Daily : 7:30am	Same day : 5:00pm	5 ml of serum and 2 ml of F, PP flouride plasma
P1004	PCOS Profile - 2	Chennai - Avadi	D	7790	Mon, Wed & Fri : 7:30am	Same day : 5:00pm	2 ml of Fluoride plasma (F, PP), 5 ml of fasting serum sample
P1005	PCOS profile - 3	Chennai - Avadi	D	11680	Mon, Wed & Fri : 7:30am	Same day : 5:00pm	2 ml of Fluoride plasma (F, PP), 6 ml of fasting serum sample
P1008	Prader-Willi / Angelman Syndrome by MS-PCR	Chennai - Avadi	D	8260	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
P1009	Prolactin (Tripooled Sample) serum	Chennai - Avadi	D	650	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum(0, 20, 40 min)
P1017	Phenol, Urine**	Chennai - Avadi	R	3070	Daily : 11:30am	20th day : 6:00pm	31 ml Urine (collect Urine after end of work shift on last working day of shift)
P1901	CANINE-Peripheral smear examination	Chennai - Avadi	B	325	Daily : 9:00am & 3:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1902	Canine - Phosphorous	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P1903	Platelet Count, EDTA whole blood	Chennai - Avadi	A	205	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1904	Canine - Potassium	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1905	Pre ? Operative Panel	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1906	Progesterone (P4), Serum by CMIA	Chennai - Avadi	B	565	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1907	Protein Electrophoresis, Serum	Chennai - Avadi	C	735	Daily : 7:30am	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P1908	Canine - Proteins, Serum	Chennai - Avadi	C	190	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1909	Prothrombin Time (PT), Citrate plasma	Chennai - Avadi	B	325	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1910	Canine-PTH	Chennai - Avadi	B	1395	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P5629	Prega-Plus	Chennai - Avadi	E	9030	Daily : 9:00am	3rd day : 6:00pm	2 ml of EDTA whole blood in 2 vacutainers, , 5 ml of serum, 2 ml of flouride plasma , 20 ml of spot urine, 10 ml of urine in Sterile Container .(Mid-stream or early morning) . Upon 1st visit sample for all other tests is collected Patient has to be then given 82.5 grams of commercial glucose (75 gram anhydrous Glucose) and Sample has to be collected exactly after two hours.Sample for Pregascreen Dual is collected between 11-13 weeks gestation. Tests for partner to be registered separately under P0630.)

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P5631	Prega-Pro	Chennai - Avadi	E	5880	Daily : 9:00am	Next day : 5:00pm	2 ml of EDTA whole blood in 2 vacutainers, , 5 ml of serum, 2 ml of flouride plasma , 20 ml of spot Urine 1st visit sample for all other tests is collected Patient has to be then given 82.5 grams of commercial glucose (75 gram anhydrous Glucose) and Sample has to be collected exactly after two hours.Sample for Pregascreen Dual is collected between 11-13 weeks gestation...
P8220	Pneumocystis jirovecii (PCP) Detection	Chennai - Avadi	D	7500	Every Tuesday, Thursday and Saturday. 9AM	Same day 11pm	3ml Sputum or BAL
P8229	PLA2 receptor antibody Quantitative , Serum	Chennai - Avadi	D	6255	Tue : 11:30am	Next day : 5:00pm	2 ml Serum
P8277	Liquid Biopsy PIC3CA NGS Mutation	Chennai - Avadi	E	45000	Daily 7.00 am	25 th days	Peripheral Blood in streck tube
P8530	PDL-1 (22C3) Head & Neck squamous cell carcinoma	Chennai - Avadi	D	15020	Daily 9 am	7th day	Paraffin Block
P8531	PDL-1 (22C3) Breast	Chennai - Avadi	D	15020	Daily 9 am	7th day	Paraffin Block
P8532	PDL-1 (22C3) Cervix	Chennai - Avadi	D	15020	Daily 9 am	7th Day	Paraffin Block
P8533	PDL-1 (22C3) Lung	Chennai - Avadi	D	15020	Daily 9 am	7th Day	Paraffin Block
P8534	PDL-1 (22C3) Urothelial	Chennai - Avadi	D	15020	Daily 9 am	7th Day	Paraffin Block
Q8101	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR	Chennai - Avadi	E	5350	Daily : 9:00am	Next day : 6:00pm	â€¢ Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition or â€¢ Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . POC - 1)Upto 12 weeks of gestation: 50-100 mg placental villi (fetal side) in 20 ml normal saline with 2-3 drops of gentamycin in a sterile container 2) In 12-18 weeks of gestation: a) 50-100 mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container b) and or fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container 3)IUFD >20 weeks of gestation: fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container. if possible in 20ml normal saline with 2-3 drops of gentamycin
R0001	RA test Rheumatoid Arthritis Serum	Chennai - Avadi	B	620	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0002	RA test Rheumatoid Arthritis Synovial Fluid	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 4hrs	1ML of Synovial Fluid
R0004	Rabies virus antibodies IgG by EIA	Chennai - Avadi	D	4250	Mon : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0005	RA-IgA Rheumatoid Arthritis Serum	Chennai - Avadi	R	855	Daily : 11:30am	7th day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0007	RBC Folate EDTA Whole Blood	Chennai	C	2145	Daily : 9:00am	Same day	EDTA Blood (2 ml) And Serum (2

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				after 8 hrs	ml)
R0009	Reducing Substances Paediatric - Urine	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
R0010	Reducing Substance by chemical test - Stool	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	5gm of Stool specimen in Sterile leak proof container
R0011	Renal (Kidney) Function Tests-1 RFT Maxi	Chennai - Avadi	C	1500	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of EDTA whole blood, 2 ml of Flouride plasma, 3 ml of serum, 20 ml of spot urine
R0012	Renal (Kidney) Function Tests	Chennai - Avadi	C	800	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum, 20 ml of spot urine
R0012_WithouteGFR	Renal (Kidney) Function Tests-2 RFT Mini	Chennai - Avadi	C	960	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum, 20 ml of spot urine
R0013	Renin Direct Plasma	Chennai - Avadi	C	3505	Daily : 9:00am	Same day : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
R0014	Pneumonia (Atypical) panel-IgG Respiratory Tract Infection Profile-IgG by EIA - Serum	Chennai - Avadi	D	7670	Thu : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0015	Pneumonia (Atypical) panel-IgM Respiratory Tract Infection Profile-IgM by EIA - Serum	Chennai - Avadi	D	7670	Thu : 7:30am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0018	Rett Syndrome MECP2 Sequencing Blood	Chennai - Avadi	E	5775	Daily : 9:00am	10th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
R0020	Rh Antibody, Serum	Chennai - Avadi	B	920	Daily : 9:00am & 5:00pm	After 6hrs	3 ml Serum sample
R0021	Rheumatic Fever Panel	Chennai - Avadi	C	1750	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	3 ml of EDTA blood and serum. And Throat swab
R0022	Rheumatoid Arthritis Panel-1	Chennai - Avadi	C	3505	Wed & Sat : 7:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0023	Rheumatoid Arthritis Panel-2	Chennai - Avadi	C	5535	Wed & Sat : 7:30am	Same day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
R0024	Rheumatoid Arthritis Panel-3	Chennai - Avadi	C	8475	Wed & Sat : 7:30am	Same day : 7:00pm	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0025	Rheumatology Profile	Chennai - Avadi	C	6440	Wed & Sat : 7:30am	Same day : 7:00pm	5 ml of serum
R0026	Ri Antibody, Serum	Chennai - Avadi	D	5310	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0027	RNP-Sm Antibody Serum	Chennai - Avadi	C	1810	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0028	Rota virus Antigen detection - Stool	Chennai - Avadi	C	1470	Daily : 9:00am & 3:00pm	After 4hrs	5gm of Stool specimen in Sterile leak proof container
R0030	Routine examination Comprehensive examination by Microscopy - Stool	Chennai - Avadi	A	270	Daily : 9:00am & 3:00pm	Next day : 6:00pm	5gm of Stool specimen in Sterile leak proof container
R0031	Routine examination - Body fluid	Chennai - Avadi	C	905	Daily : 9:00am & 3:00pm	after 6 hrs	2ml of Body Fluid in Sterile leak proof container
R0032	Routine examination - Ascitic fluid	Chennai - Avadi	C	905	Daily : 9:00am & 3:00pm	after 6 hrs	2ml of Ascitic Fluid in Sterile leak proof container
R0032_2	Routine examination - Ascitic fluid - 2	Chennai - Avadi	C	1120	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Ascitic Fluid in Sterile leak proof container
R0033	Routine examination - CSF	Chennai - Avadi	C	905	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of CSF in Sterile leak proof tube or container
R0033_2	Routine examination - CSF	Chennai - Avadi	C	1415	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of CSF in Sterile leak proof tube or container
R0034	Routine examination (Pericardial Fluid)	Chennai - Avadi	C	905	Daily : 9:00am & 3:00pm	Next day : 10am	2 ml of Pericardial Fluid in Sterile leak proof container
R0035	Routine examination (Peritoneal Fluid)	Chennai - Avadi	C	905	Daily : 9:00am & 3:00pm	Next day : 10am	2 ml of Peritoneal Fluid in Sterile leak proof container
R0036	Routine examination (Pleural Fluid)	Chennai - Avadi	C	905	Daily : 9:00am & 3:00pm	after 6 hrs	2ml of Pleural Fluid in Sterile leak proof container
R0036_2	Routine examination (Pleural Fluid)-2	Chennai - Avadi	C	1570	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Pleural Fluid in Sterile leak proof container
R0037	Routine examination (Prostatic Fluid)	Chennai - Avadi	C	905	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Prostatic Fluid in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
R0038	Semen Analysis by Microscopy - Semen	Chennai - Avadi	B	540	Daily : 9:00am to 7:00pm	After 6hrs	2ml of Semen in Sterile leak proof container
R0039	Routine examination of Sputum	Chennai - Avadi	B	865	Daily : 9:00am & 3:00pm	Next day : 6:00pm	2ml of Sputum in Sterile leak proof container
R0040	Routine Examination of Synovial Fluid	Chennai - Avadi	C	905	Daily : 9:00am & 3:00pm	after 6 hrs	2ml of Synovial Fluid in Sterile leak proof container
R0041	Routine Examination of Urine	Chennai - Avadi	A	150	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
R0042	RSV-IgG Antibody by EIA - Serum	Chennai - Avadi	C	2145	Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0043	RSV-IgM Antibody by EIA - Serum	Chennai - Avadi	C	2145	Thu : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0044	Rubella Virus Avidity test by EIA - Serum	Chennai - Avadi	C	1355	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0045	Rubella Virus IgG and IgM Antibody by CLIA - Serum	Chennai - Avadi	C	1300	Daily : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0047	Rubella Virus IgG Antibody by CLIA - Serum	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0048	Rubella Virus IgM Antibody by CLIA - Serum	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 8hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0049	Rubella Virus RNA Detection by PCR - Amniotic Fluid	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
R0050	Rubella Virus RNA Detection by PCR - EDTA Blood	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3 ml of Whole Blood in EDTA Vacutainer
R0051	Rubella Virus RNA Detection by PCR - CSF	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CSF in Sterile leak proof tube or container
R0052	Rubella Virus RNA Detection by PCR - CVS	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CVS in Sterile leak proof tube or container
R0053	Rubella Virus RNA Detection by PCR - Urine	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	15ml of Spot Urine in Sterile leak proof container
R0055	RVVT-Russell Viper Venom Time	Chennai - Avadi	C	1075	Daily : 9:00am & 3:00pm	After 6hrs	1 ml of citrated platelet poor plasma in frozen condition.
R0056	Retinitis Profile by Real Time PCR - Aqueous/Vitreous Fluid	Chennai - Avadi	D	11800	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3ml of Aqueous/ Vitreous Fluid in Sterile leak proof container
R0059	Retinol binding protein**	Chennai	R	2780	Daily : 11:30am	15th day :	3ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				6:00pm	
R0060	Ricketssia DNA Detection by Real Time PCR - Blood	Chennai - Avadi	D	5900	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
R0063	ROS1 Gene Rearrangement by FISH	Chennai - Avadi	D	13570	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
R0066	Respiratory panel by FilmArray (Biofire) - Nasopharyngeal swab	Chennai - Avadi	E	17640	Daily : 9:00am to 9:00pm	After 10hrs	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
R0067	RET Gene Rearrangement by FISH, Tissue	Chennai - Avadi	E	9450	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
R0069	RNA Polymerase III Antibodies, IgG**	Chennai - Avadi	C	3730	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0077	Allergy panel - Common Rhinitis/ Wheeze	Chennai - Avadi	D	8260	Daily : 11:30am	Next day : 5:00pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
R0080	Renal biopsy profile (LM, IF, EM)	Chennai - Avadi	R	12305	Daily : 9:00am	20 day	1) Tissue in 10% NB formalin 2) Tissue in Michel Medium. Please use this code when you are sending 2 sample vials ; One in Formaldehyde & other in Michels medium. For Electron microscopy please send 3% glutaraldehyde
R0086	Rifampicin**	Chennai	D	7080	Daily : 11:30am	10th day :	2 ml K2-EDTA plasma

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
R1008	Anti Ribosome antibodies, serum**	- Avadi Chennai - Avadi	R	7490	Daily : 11:30am	6:00pm After 4 weeks : 6:00pm	5 ml Serum
R1901	Reticulocyte panel - Automated, EDTA whole blood	Chennai - Avadi	B	380	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
R1902	CANINE-Routine Examination Stool	Chennai - Avadi	A	165	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
R1903	Routine examination reflex to cytology, Body fluids	Chennai - Avadi	C	790	Daily : 9:00am & 3:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
R1904	Canine - Routine Examination Urine	Chennai - Avadi	A	130	Daily : 9:00am to 9:00pm	After 24hrs	10 ml of Spot urine
R5588	Reassure BOH Mini	Chennai - Avadi	D	6960	Daily : 7:00am	3rd day : 6:00pm	5 ml serum, 3 ml citrate plasma
R5589	Reassure BOH Maxi	Chennai - Avadi	D	12980	Daily : 7:00pm	15th day : 6:00pm	5 ml serum, 4 ml Heparin whole Blood (Husband and wife), 3 ml citrate plasma
S0001	Spinocerebellar Ataxia PANEL (SCA1, 2, 3, 6, 7, 12 and 17)	Chennai - Avadi	E	13125	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0002	Spinocerebellar Ataxia (SCA1) DNA TEST	Chennai - Avadi	D	3420	Mon & Wed: 9:00am	7th day : 6:00pm	Clinical/family history, referring clinicians contact number and TRF mandatory. 3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0003	Spinocerebellar Ataxia 2 (SCA2) DNA TEST	Chennai - Avadi	D	3420	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0004	SCA3 (Machado-Joseph Disease) DNA Test	Chennai - Avadi	D	3420	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0005	Spinocerebellar Ataxia 6 (SCA6) DNA TEST	Chennai - Avadi	D	3420	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0006	Spinocerebellar Ataxia TYPE 7 (SCA7) DNA TEST	Chennai - Avadi	D	3420	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0007	Spinocerebellar Ataxia Type 12 (SCA12) DNA Test	Chennai - Avadi	D	3420	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0008	Scl-70 Antibody Serum	Chennai - Avadi	C	1750	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0009	Scleroderma Panel	Chennai - Avadi	C	5650	Wed & Sat : 7:30am	Next day : 5:00pm	5 ml of serum
S0010	Selenium, Blood	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
S0011	Selenium, Serum	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
S0012	Selenium, Urine	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
S0013	Selenium Urine 24H**	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
S0014	Serine Quantitative Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
S0015	Serine Quantitative Urine 24H**	Chennai - Avadi	D	7020	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
S0016	Serotonin5-Hydroxy Tryptamine Serum	Chennai - Avadi	D	4955	Wed : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0017	Serotonin5-Hydroxy Tryptamine,Urine	Chennai - Avadi	D	4955	Wed : 11:30am	Next day : 5:00pm	10 ml of 24 hrs urine (Preferred) / Spot urine
S0018	SGOT AST Serum	Chennai - Avadi	A	200	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0019	SGPT ALT Serum	Chennai - Avadi	A	200	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0020	SHBG Sex hormone binding globulin Serum	Chennai - Avadi	C	2600	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0021	Sickle Cell Anemia Mutation Analysis	Chennai - Avadi	D	3775	Wed & Fri: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Clinical history, referring clinicians contact number, HPLC Report, CBC findings, clinical/family history and TRF mandatory.
S0022	Sickling Test	Chennai - Avadi	B	295	Daily : 9:00am to 9:00pm	after 6 hrs	3 ml EDTA whole blood
S0024	Silver ICPMS Urine	Chennai - Avadi	R	3745	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
S0026	SIROLIMUS, BLOOD BY LC MSMS	Chennai - Avadi	C	5650	Tue & Fri: 11:30am	Next day : 5:00pm	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
S0027	Sjogrens Syndrome Profile	Chennai - Avadi	D	2655	Wed & Sat : 7:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0028	Anti-SLA/LP, IgG	Chennai - Avadi	C	3165	Tues/fri 9am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0029	SLE profile	Chennai - Avadi	C	5085	Wed & Sat : 9:00am	Next day : 5:00pm	3 ml of EDTA whole blood , 3 ml of serum
S0030	Sm Antibody Smith Antibody Serum	Chennai - Avadi	C	1750	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0031	SMA 12+2	Chennai - Avadi	D	1355	Daily : 9:00am to 9:00pm	After 8hrs	2 ml of flouride plasma , 3 ml of serum
S0032	Sodium ISE Serum	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0033	Sodium Urine	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Spot Urine
S0034	Sodium Urine 24H	Chennai - Avadi	A	275	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
S0035	Sodium Body Fluid	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Body Fluid
S0036	Soluble Transferrin Receptor Serum	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0037	Specific gravity by Fully Automated Urine chemistry & Microscopy - Urine	Chennai - Avadi	A	110	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
S0038	Spectral Karyotyping (SKY) (m-FISH) 24 Color Painting Blood**	Chennai - Avadi	D	18880	Daily : 9:00am	15th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
S0039	Spectral Karyotyping (SKY) (m-FISH) 24 Color Painting Bone Marrow**	Chennai - Avadi	R	17120	Daily : 11:30am	15th day : 6:00pm	2-3ml Bone marrow/3ml Blood in Sodium Heparin Vacutainer
S0040	Sperm Antibody-Total ASAB serum	Chennai - Avadi	C	1245	Tue & Fri: 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0041	Sphingomyelinase Niemann Pick Disease Blood**	Chennai - Avadi	R	8025	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
S0044	Spinal Muscular Atrophy Microdeletion By	Chennai	D	3895	Tue & Fri : 9:00am	7th day :	3-5 ml of EDTA whole blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	PCR SMN1 Gene Deletion	- Avadi				6:00pm	required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0045	Sry Deletion FISH	Chennai - Avadi	D	2950	Daily : 9:00am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
S0046	SSA-Antibody Ro Serum	Chennai - Avadi	C	1750	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0047	SSB-Antibody La Serum	Chennai - Avadi	C	1750	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0048	STD Profile -I	Chennai - Avadi	D	3420	Tue, Thu & Sat: 7:30am	5th day : 6:00pm	3 ml of serum
S0049	STD Profile -II	Chennai - Avadi	D	5665	Daily : 8:00pm	5th day : 6:00pm	4 ml Serum
S0050	STD Profile -III (HPV-DNA, Chlamydia Trachomatis DNA, HSV-DNA, HIV-1&2 RNA)	Chennai - Avadi	D	9440	Daily : 9:00am	5th day : 6:00pm	6 ml of EDTA whole blood, 15ml Urine in sterile container. For HPV : Swab, Tissue for male and LBC, Swab, Tissue for female
S0051	Steroid analysis From medicine **	Chennai - Avadi	R	11770	Daily : 11:30am	10th day : 6:00pm	Medicine sample (10 grm)
S0052	Stomach Marker Profile	Chennai - Avadi	D	3895	Tue & Fri : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0053	Calculus (Stone) Analysis by automated FTIR (Gall bladder calculus)	Chennai - Avadi	D	1300	Daily : 11:30am	Next day : 5:00pm	Collect sample in dry plain container, mention stone type and history
S0054	Calculus (Stone) Analysis by automated FTIR (Kidney (Renal) calculus)	Chennai - Avadi	D	1300	Daily : 11:30am	Next day : 5:00pm	Collect sample in dry plain container, mention stone type and history
S0055	Calculus (Stone) Analysis by automated FTIR (Urine calculus)	Chennai - Avadi	D	1300	Daily : 11:30am	Next day : 5:00pm	Collect sample in dry plain container, mention stone type and history
S0056	Stone screening profile	Chennai - Avadi	D	2360	Daily : 9:00am	Next day : 5:00pm	24 hours urine sample with appropriate preservative.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0058	Streptococcus group B Antigen detection by Agglutination - CSF	Chennai - Avadi	C	1355	Daily : 9:00am & 3:00pm	After 4hrs	2ml of CSF in Sterile leak proof tube or container
S0059	Streptococcus group B Antigen detection by Agglutination - Serum	Chennai - Avadi	C	1355	Daily : 9:00am & 3:00pm	After 4hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0060	Streptococcus pneumoniae Antigen detection by Agglutination - CSF	Chennai - Avadi	C	1355	Daily : 9:00am & 3:00pm	After 4hrs	2ml of CSF in Sterile leak proof tube or container
S0061	Streptococcus pneumoniae Antigen detection by Agglutination - Serum	Chennai - Avadi	C	1355	Daily : 9:00am & 3:00pm	After 4hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0062	Streptococcus pneumoniae Antigen detection by Agglutination - Urine	Chennai - Avadi	C	1355	Daily : 9:00am & 3:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
S0063	Sucrose Lysis Test (PNH Screening tests) (S0063)	Chennai - Avadi	B	865	Daily : 9:00am	Next day : 5:00pm	3 ml of EDTA Whole Blood And 3 ml Serum
S0066	Glucose Sugar - Urine	Chennai - Avadi	A	80	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
S0067	Syphilis Antibodies Treponema by EIA - Serum	Chennai - Avadi	C	1245	Tue & Fri: 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0070	Scrub typhus profile Weil felix and scrub typhus-IgM by EIA - Serum	Chennai - Avadi	C	1245	Tue & Fri: 11:30am	Next day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0071	Scrub typhus IgM by EIA - Serum	Chennai - Avadi	C	1015	Tue, Thu & Sat: 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0074	SMA Carrier Detection Genetic Disorders**	Chennai - Avadi	R	9630	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
S0075	SCA17 PCR-Fragment analysis	Chennai - Avadi	D	3420	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0077_Sperm-DFI	SPERM DNA FRAGMENTATION (DFI)	Chennai	D	5300	Daily : 11.30 am	3rd day 1 pm	2-3 ml of freshly collected Semen

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi					sample in a sterile leak proof container. If reaching the lab within 4hrs of collection, transported at ambient temperature.If more than 4 hrs - then transported in frozen condition.3 days abstinence is mandatory.
S0078	Calculus (Stone) Analysis by automated FTIR (other)	Chennai - Avadi	D	1300	Daily : 11:30am	Next day : 5:00pm	Collect sample in dry plain container, mention stone type and history
S0079	SYT gene rearrangement (Translocation: t(X;18) translocation) **	Chennai - Avadi	R	4280	Daily : 11:30am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
S0083	Campylobacter Culture and Identification by Vitek/ MALDITOF - Stool	Chennai - Avadi	C	680	Daily : 9:00am to 9:00pm	5th day : 6:00pm	5gm of Stool specimen in Sterile leak proof container
S0084	SUCCINYLACETONE, BLOOD**	Chennai - Avadi	R	5030	Daily : 11:30am	8th day : 6:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
S0085	Succinylacetone Urine**	Chennai - Avadi	R	4600	Daily : 11:30am	7th day : 6:00pm	Spot Urine:10 ml of Spot urine required with clinical history and clinician number. Baby DOB, sample collection date and relevant clinical history mentioned .
S0096	STRONTIUM ICPMS Blood	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	3 ml EDTA Blood / Metal free EDTA Blood
S0097	Sexually Trnmitted Infections (STI) Urethritis Extended Panel	Chennai - Avadi	D	7670	Wed : 9:00am	Next day : 6:00pm	Urethral discharge, Genital discharge, Blister fluid, Ulcer, Ulcer base, High Vaginal swab, Urethral Swab, Endocervical swab
S0103	Soft Tissue Sarcoma Panel**	Chennai - Avadi	R	42800	Daily 7.30 am	10th days	Paraffin Block
S0547	sFlt-1 (Soluble fms-like tyrosine kinase-1) Serum	Chennai - Avadi	D	4485	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0548	PIGF (Placental Growth Factor) Serum	Chennai - Avadi	D	4485	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0549	sFlt-1/PIGF ratio serum	Chennai - Avadi	D	6845	Daily cut off 7:30 AM	same day 9:00 PM	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0551	Sanger Sequencing / Additional family member screening**	Chennai - Avadi	R	9630	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
S0552	Selenium, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							storage and shipping. Plz enquire before sending
S0553	Selenium, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
S0554	Silver, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
S0555	Silver, Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
S0556	Steroids test by MS method	Chennai - Avadi	R	11770	Daily : 11:30am	10th day : 6:00pm	Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clincial symptoms if any. Reported in 7 days
S0558	SARS- CoV-2 (COVID 19) Detection by Real Time PCR - Nasopharangeal & Oropharangeal Swab	Chennai - Avadi	O	410	24hr registration	After 24hrs	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
S0559	SARS CoV2 Antigen Detection by Immunchromatography - Nasopharangeal Swab	Chennai - Avadi	O	515	Daily : 9:00am	Same day : 5:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
S0560	SARS CoV2 Antigen Reflex to SARS CoV2 real time RT PCR (COVID 19 reflex test)* - Nasopharangeal & Oropharangeal Swab	Chennai - Avadi	O	925	Daily : 9:00am	After 24hrs	For antigen test :at ambient temperature within an hour or if transported at 2 to 8 degree Celsius for 4 to 6 hours.
S0562	SARS- CoV-2 (COVID 19) viral load by Real Time PCR*	Chennai - Avadi	E	2100	Daily : 9:00am	Same day : 5:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
S0565	Soluble CD25 (IL 2 Receptor)	Chennai - Avadi	D	11210	Daily : 11:30am	5th day : 6:00pm	3 ml Serum
S1005	SCC-TA4 ANTIGEN - SERUM **	Chennai - Avadi	R	9415	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
S1006	Schistosomiasis - confirmation by EIA - Serum**	Chennai - Avadi	R	9735	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S1020	Strongyloidiasis serology, Serum by EIA** - Serum	Chennai - Avadi	R	6100	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S1038	Sickle Cell Anemia Mutation Prenatal Diagnosis	Chennai - Avadi	D	7080	Wed & Fri: 9:00am	7th day : 6:00pm	Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition or Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
S1039	Spinal Muscular Atrophy (SMA) Microdeletion	Chennai - Avadi	D	5310	Daily : 9:00am	7th day : 6:00pm	Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Duly

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
S1039_c	Spinal Muscular Atrophy (SMA)-SMN1 Gene deletion	Chennai - Avadi	D	7080	Daily : 9:00am	7th day : 6:00pm	Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G -FORM OF CONSENT (For invasive techniques) Mother 3-5 ml EDTA Blood sample mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius .
S1040	Stone screening profile-1	Chennai - Avadi	D	4130	Daily : 9:00am	Next day : 5:00pm	5 ml of serum and 20 ml of 24 Hr Urine sample. Appropriate preservative needed.
S1041	Systemic sclerosis Profile	Chennai - Avadi	D	4720	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S1901	Canine - SGOT (AST)	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S1902	Canine - SGPT (ALT)	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
S1903	Skin Scraping Examination	Chennai - Avadi	C	565	Daily : 9:00am	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
S1904	Canine - Sodium	Chennai - Avadi	A	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
S1905	CANINE-Calculus (Stone) Analysis by automated FTIR	Chennai - Avadi	C	1130	Daily : 7:30am	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
S8322	SARS-CoV-2 (COVID 19) Detection (Qualitative) by Rapid	Chennai - Avadi	E	2075	Daily : 9:00am	Same day : 5:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
S8399	Streptococcus Group A Antigen Detection	Chennai - Avadi	L4	1200	Daily: 9am to 9pm	After 6 hrs	Throat Swab
S8514	Salivary Cortisol	Chennai - Avadi	E	2940	Tuesday 9 am	next day 5 pm	4-5 ml saliva in sterile lolli sponge tube
T0004	t(1 19)by FISH E2A(ALL)	Chennai	D	4600	Daily : 9:00am	3rd day :	3 ml Blood in Na-heparin

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				6:00pm	vacutainer
T0005	t(1 19)By FISH E2A (ALL) Bone Marrow	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0006	t(11;14) by FISH IGH/CCND1 (ALL, NHL) Bone Marrow**	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0007	t(11;14) by FISH IGH/CCND1 (ALL, NHL) Blood**	Chennai - Avadi	D	5310	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0008	t(11q23) by FISH MLL gene rearrangement (ALL, AML) Blood	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0009	t(11q23) by FISH MLL gene rearrangement (ALL, AML) Bone Marrow	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0010	TEL/AML1(ALL) t(12 21) by FISH	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0011	t(12 21) by FISH TEL/AML1 (ALL) Bone Marrow	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Blood in Na-heparin vacutainer
T0012	PML Ra Ra t(15 17) (AML-M3) by FISH	Chennai - Avadi	D	5780	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0013	t(15 17) by FISH PML Ra Ra (AML-M3) Bone Marrow	Chennai - Avadi	D	5780	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0014	PML Ra Ra (Quanti) by PCR BLOOD**	Chennai - Avadi	D	5780	Daily : 7:30am	7th day : 6:00pm	3 ml of EDTA Whole Blood
T0015	PML Ra Ra (Quanti) by PCR BONE MARROW**	Chennai - Avadi	D	5780	Daily : 7:30am	7th day : 6:00pm	2-3 ml Bone Marrow in EDTA vacutainer
T0016	PDGFR Alpha rearrangement t(4q12) (HES) by FISH	Chennai - Avadi	D	5780	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0017	t(4q12) by FISH PDGFR Alpha Rearrangement (HES) Bone Marrow**	Chennai - Avadi	D	5780	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0018	PDGFR Beta rearrangement t(5q33) (HES) Blood by FISH	Chennai - Avadi	D	5780	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0019	PDGFR Beta rearrangement t(5q33) (HES) Bone marrow by FISH	Chennai - Avadi	D	5780	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0020	ETO/AML1 (AML-M2) t(8;21) by FISH	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0021	t(8 21) by FISH ETO/AML1 (AML-M2) Bone Marrow	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0024	Bcr-Abl (PHILADELPHIA) t(9 22) by FISH,	Chennai	D	4600	Daily : 9:00am	3rd day :	3 ml Blood in Na-heparin

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	BLOOD	- Avadi				6:00pm	vacutainer
T0025	Bcr-Abl (PHILADELPHIA) t(9 22) by FISH, Bone Marrow	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0028	T3-Free Free Tri-Iodothyronine Serum	Chennai - Avadi	B	295	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0029	T3-Total Tri Iodothyronine Serum	Chennai - Avadi	B	215	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0030	T4-Free Free Thyroxine Serum	Chennai - Avadi	B	300	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0031	T4-Total Thyroxine Serum	Chennai - Avadi	B	215	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0032	Tacrolimus By ELCIA Blood	Chennai - Avadi	C	4200	Daily : 1:00pm	After 8hrs	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
T0033	Tacrolimus By LCMS, Blood	Chennai - Avadi	C	4070	Tue, Thu & Sat: 9:00am	Same day : 9:00pm	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
T0036	TCR-Alpha Beta	Chennai - Avadi	D	2715	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							condition(preferably). Clinical history and clinician number required.
T0037	TCR-Gamma Delta	Chennai - Avadi	D	2715	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
T0038	TDT (Nuclear)	Chennai - Avadi	D	2715	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
T0039	Testicular Marker Profile	Chennai - Avadi	D	1415	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0040	Testosterone Bioavailable	Chennai - Avadi	C	2710	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0041	TestosteroneFreeSerum	Chennai - Avadi	C	1610	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0042	Testosterone Total Serum	Chennai - Avadi	B	700	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0043	Testosterone Profile	Chennai - Avadi	C	3955	Tue, Thu & Sat: 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0044	Tetanus Toxoid-IgG Quantitative test by EIA - Serum	Chennai - Avadi	C	4520	Fri : 9:00am	Same day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0045	Thalassemia Studies	Chennai - Avadi	C	1695	Daily : 9:00am	Same day : 7:00pm	3 ml EDTA whole blood and serum sample
T0046	Thallium Blood	Chennai - Avadi	D	3540	Daily : 11:30am	5th day : 6:00pm	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0047	Thallium Serum	Chennai - Avadi	D	3540	Daily : 11:30am	5th day : 6:00pm	HeavyMetalSamplingguidelines.pdf
T0048	Thallium Urine Spot	Chennai - Avadi	D	3540	Daily : 11:30am	5th day : 6:00pm	HeavyMetalSamplingguidelines.pdf
T0049	Theophylline Serum	Chennai - Avadi	C	1120	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0050	Threonine Quantitative Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
T0051	Threonine Quantitative Urine	Chennai - Avadi	R	6365	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
T0052	Thrombin Time (TT)	Chennai - Avadi	B	810	Daily : 9:00am, 3:00pm & 9:00pm	After 6hrs	1 ml of citrated platelet poor plasma in frozen condition.
T0053	Thrombotic Risk DNA Panel (FV, FII, MTHFR)	Chennai - Avadi	D	12980	Wed & Fri: 9:00am	10th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
T0054	Thrombophilia Profile-Maxi	Chennai - Avadi	D	19470	Mon & Thu : 9:00am	Next day : 5:00pm	Serum and 1 ml of citrated platelet poor plasma in frozen condition.
T0055	Thrombophilia Profile-Mini	Chennai - Avadi	D	11210	Mon & Thu : 9:00am	Next day : 5:00pm	Serum and 1 ml of citrated platelet poor plasma in frozen condition.
T0056	Thrombotic Profile	Chennai - Avadi	D	18290	Mon & Thu : 9:00am	Next day : 5:00pm	Serum and 1 ml of citrated platelet poor plasma in frozen condition.
T0057	Thyroglobulin Serum	Chennai - Avadi	C	1750	Daily : 9:00am to 9:00pm	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0058	Thyroglobulin Antibody ATA CSF	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
T0059	Thyroglobulin Antibody ATA Serum	Chennai - Avadi	B	1190	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0060	Thyroid Antibodies-AtAb TPO and ATA,CSF	Chennai - Avadi	C	2375	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
T0061	Thyroid Antibodies-TPO and ATG,Serum	Chennai	C	2040	Daily : 9:00am to	After 6hrs	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00pm		ml serum for profiles and reflex tests.
T0062	Thyroid Comprehensive Profile-1	Chennai - Avadi	B	2810	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0063	Thyroid Comprehensive Profile-2	Chennai - Avadi	B	1460	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0064	Thyroid Marker Profile	Chennai - Avadi	D	6845	Wed & Sat : 9:00am	Next day : 5:00pm	5 ml of serum
T0065	Thyroid panel-1 (T3/T4/TSH)	Chennai - Avadi	B	450	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0066	Thyroid panel - 2(FT3, FT4, TSH)	Chennai - Avadi	B	780	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0067	Thyroid Panel-3 (Reflex)	Chennai - Avadi	B	405	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0068	Thyroid Panel-4	Chennai - Avadi	B	430	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0069	Thyroxine Binding Globulin (TBG) Biochemical Serum	Chennai - Avadi	C	5650	Daily : 11:30am	5th day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0071	TNF Alpha Tumour Necrosis Factor Serum	Chennai - Avadi	C	4405	Wed & Sat : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0072	TORCH PCR (Toxoplasma, CMV, Rubella and HSV 1&2) by PCR - Amniotic Fluid	Chennai - Avadi	C	8475	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
T0073	TORCH PCR (Toxoplasma, CMV, Rubella and HSV 1&2) by PCR - Blood	Chennai - Avadi	C	7345	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3 ml of Whole Blood in EDTA Vacutainer
T0074	TORCH by PCR **	Chennai - Avadi	C	8475	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	Chorionic Villus Sample In Sterile saline container
T0075	TORCH PCR (Toxoplasma, CMV, Rubella and HSV 1&2) by PCR - CSF	Chennai - Avadi	C	8475	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CSF in Sterile leak proof tube or container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0076	TORCH PCR (Toxoplasma, CMV, Rubella and HSV 1&2) by PCR - Urine	Chennai - Avadi	C	8475	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	15ml of Spot Urine in Sterile leak proof container
T0077	TORCH-10	Chennai - Avadi	C	3390	Daily : 9:00am	Same day : 5:00pm	5 ml of serum
T0078	TORCH-13	Chennai - Avadi	C	3955	Mon : 9:00am	Next day : 5:00pm	5 ml of serum
T0079	TORCH-4 IgG	Chennai - Avadi	C	1640	Daily : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0080	Tacrolimus CYP3A5 genotyping**	Chennai - Avadi	R	8560	Tue & Thu: 11:30am	7th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
T0081	TORCH-4 IgM	Chennai - Avadi	C	1640	Daily : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0083	TORCH-5 IgG	Chennai - Avadi	C	1810	Daily : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0084	TORCH-5 IgM	Chennai - Avadi	C	1810	Daily : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0085	TORCH-8	Chennai - Avadi	C	2940	Daily : 9:00am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0086	Toxoplasma Avidity test by EIA- Serum	Chennai - Avadi	C	1355	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0087	Toxoplasma DNA Detection by PCR - Amniotic Fluid	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
T0088	Toxoplasma DNA Detection by PCR - Blood	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3 ml of Whole Blood in EDTA Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0089	Toxoplasma DNA Detection by PCR - Chorionic Villus	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CVS in Sterile leak proof tube or container
T0090	Toxoplasma DNA Detection by PCR - Cord Blood	Chennai - Avadi	C	3390	Mon & Thu : 7:30am	Next day : 6:00pm	3ml of Cord Blood in EDTA Vacutainer
T0091	Toxoplasma DNA Detection by PCR - CSF	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CSF in Sterile leak proof tube or container
T0092	Toxoplasma DNA Detection by PCR - Urine	Chennai - Avadi	C	3390	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	15ml of Spot Urine in Sterile leak proof container
T0094	Toxoplasma-IgG Antibody by CMIA - Serum	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0095	Toxoplasma-IgG & IgM by CLIA - Serum	Chennai - Avadi	B	1240	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0096	Toxoplasma-IgM Antibody by CMIA - Serum	Chennai - Avadi	B	650	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0097	Treponema Pallidum Haemagglutination Assay by Agglutination - Serum	Chennai - Avadi	B	755	Daily : 9:00am	Same Day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0098	TPMT Enzyme ActivityThiopurine Methyl TransferasePlasma	Chennai - Avadi	D	4720	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
T0099	Thiopurine S-Methyltransferase (TPMT) Genotype	Chennai - Avadi	D	5900	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
T0100	TPO (Thyroid Peroxidase) Antibody(AMA) Microsomal Antibody CSF	Chennai - Avadi	C	1355	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of CSF in leakproof container (foin vu conical tube is recommended to prevent leakage).
T0101	TPO (Thyroid Peroxidase) Antibody (AMA) Microsomal Antibody Serum	Chennai - Avadi	B	1180	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0102	TransferrinSerum	Chennai - Avadi	B	1190	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0103	TRH STIMULATION TEST FOR PROLACTIN, SERUM	Chennai - Avadi	C	1525	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum (0, 30, 60 min)
T0104	TRH stimulation test For TSH, Serum	Chennai - Avadi	C	1130	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum (0, 30, 60 min)
T0105	Triglycerides Urine	Chennai	A	325	Daily : 9:00am to	After 4hrs	2 ml of Spot Urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00pm		
T0106	Triglycerides Serum	Chennai - Avadi	A	245	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0107	Triglycerides Body Fluid	Chennai - Avadi	A	325	Daily : 9:00am to 9:00pm	After 4hrs	2 ml of Body Fluid
T0108	Tripeptidyl Peptidase 1 Blood	Chennai - Avadi	R	5990	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
T0109	Trisomy 12 (CLL) by FISH**	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0110	Trisomy 12 by FISH CLL Bone Marrow	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0111	Trisomy 8 by FISH AML, MDS Blood	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0112	Trisomy 8 by FISH AML, MDS Bone Marrow	Chennai - Avadi	D	4600	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0113	Trisomy (Aneuploidy) for Chromosome 13 & 21 By FISH Amniotic Fluid	Chennai - Avadi	D	5190	Daily : 9:00am	3rd day : 6:00pm	20 ml of Amniotic Fluid In Sterile Falcon centrifuge tube
T0114	Trisomy (Aneuploidy) Chromosome 13 & 21 by FISH	Chennai - Avadi	D	5190	Daily : 9:00am	5th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
T0115	Trisomy (Aneuploidy) for Chromosome 13 & 21 By FISH Chorionic Villus	Chennai - Avadi	D	5190	Daily : 9:00am	3rd day : 6:00pm	20-30 mg Chorionic Villus in Sterile nutrient medium in 15 ml Falcon tubes(Medium will be provided by Metropolis).
T0116	Trisomy for Chromosome 13 & 21 By FISH Cord blood	Chennai - Avadi	D	5190	Daily : 9:00am	5th day : 6:00pm	3 ml of Prenatal Cord Blood In Na-Heparin Vacutainer
T0118	Trisomy(Aneuploidy)Chromosome13,18,21, X,Y by FISH	Chennai - Avadi	D	9440	Daily : 9:00am	3rd day : 6:00pm	20 ml of Amniotic Fluid In Sterile Falcon centrifuge tube
T0119	Trisomy for Chromosome 13, 18, 21, X, Y By FISH Blood	Chennai - Avadi	D	9440	Daily : 9:00am	5th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
T0120	Trisomy (Aneuploidy) for Chromosome 13, 18, 21, X, Y By FISH Chorionic Villus	Chennai - Avadi	D	9440	Daily : 9:00am	3rd day : 6:00pm	20-30 mg Chorionic Villus in Sterile nutrient medium in 15 ml Falcon tubes(Medium will be provided by Metropolis).
T0121	Trisomy for Chromosome 13, 18, 21, X, Y	Chennai	D	9440	Daily : 9:00am	5th day :	3 ml of Prenatal Cord Blood In Na-

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	By FISH Cord blood	- Avadi				6:00pm	Heparin Vacutainer
T0123	Trisomy (Aneuploidy) Chromosome 18, X,Y by FISH	Chennai - Avadi	D	5545	Daily : 9:00am	3rd day : 6:00pm	20 ml of Amniotic Fluid In Sterile Falcon centrifuge tube
T0124	Trisomy for Chromosome 18, X & Y By FISH Blood	Chennai - Avadi	D	5545	Daily : 9:00am	5th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
T0125	Trisomy (Aneuploidy) for Chromosome 18, X & Y By FISH Chorionic Villus	Chennai - Avadi	D	5545	Daily : 9:00am	3rd day : 6:00pm	20-30 mg Chorionic Villus in Sterile nutrient medium in 15 ml Falcon tubes(Medium will be provided by Metropolis).
T0126	Trisomy for Chromosome 18, X & Y By FISH Cord blood	Chennai - Avadi	D	5545	Daily : 9:00am	5th day : 6:00pm	3 ml of Prenatal Cord Blood In Na-Heparin Vacutainer
T0128	Troponin-I, Serum	Chennai - Avadi	C	1660	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0129	Troponin-T, Serum	Chennai - Avadi	C	1750	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0130	TSH (Ultrasensitive) ECLIA Serum	Chennai - Avadi	B	350	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0131	TSH Receptor Antibody LATS, TSI - Serum	Chennai - Avadi	C	4600	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0132	TSH-Thyroid Stimulating Hormone Neonatal Screen Dried blood spot	Chennai - Avadi	B	485	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
T0133	tTG Antibody-IgA Tissue Transglutaminase,Celiac disease Serum	Chennai - Avadi	C	1300	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0134	tTG Antibody-IgG Tissue Transglutaminase,Celiac disease Serum	Chennai - Avadi	C	1470	Mon & Thu : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0135	Typhi Dot IgM antibody for Typhoid by Immunochromatography - Serum	Chennai - Avadi	B	550	Daily : 9:00am & 3:00pm	After 4hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0136	Typhoid DNA by Real Time PCR - Blood	Chennai - Avadi	C	6215	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
T0137	Typhoid DNA by Real Time PCR - CSF	Chennai - Avadi	C	6215	Daily : 10:00am	After 24hrs	2ml of CSF in Sterile leak proof tube or container
T0138	Typhoid Profile	Chennai - Avadi	D	2265	Daily : 9:00am to 9:00pm	5th day : 6:00pm	3 ml of EDTA whole blood, 3 ml of serum and 1 blood culture bottle
T0139	Tyrosine Quantitative Plasma	Chennai - Avadi	D	7020	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
T0140	Tyrosine Quantitative Urine**	Chennai - Avadi	R	6365	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
T0141	Tzanck smear (Microscopic examination)	Chennai - Avadi	B	430	Daily : 9:00am	After 8hrs	Air Dried Smears From Vesicular fluid
T0143	TBE-TICKBORNE ENCEPHALITIS Antibodies by EIA - Serum**	Chennai - Avadi	D	8615	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0145	Trypanosoma Cruzi IgG by EIA - Serum**	Chennai - Avadi	R	2245	Daily : 11:30am	20th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0158	t(12;16):FUS rearrangement **	Chennai - Avadi	R	17120	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
T0160	t(X;18):SSX18 translocation **	Chennai - Avadi	R	17120	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
T0162	Tumor marker panel (270 genes), comprehensive, FFPE	Chennai - Avadi	D	47000	Daily : 7:30am	After 4 weeks : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report. 4 ml of EDTA whole blood or collection in special tubes as per test. Sample to be shipped at Ambient temperature for local samples and refrigerated condition for outstation samples. do not freeze or keep EDTA blood samples in contact with ice packs. EDTA samples stable for 2 days at Ambient temperature and 5 days at Refrigerated condition
T0168	NextGen Tumor Mutation Burden	Chennai - Avadi	E	43000	Daily : 7:30am	After 4 weeks : 6:00pm	(FFPE block) Tumor containing tissue fixed in 20% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
T0169	Thallium, hair by ICPMS	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
T0170	Tin, EDTA Whole Blood**	Chennai	R	2945	Daily : 11:30am	10th day :	Individual Metals and Minerals.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0171	Tin, Hair tests by ICPMS**	- Avadi Chennai - Avadi	R	1925	Daily : 11:30am	6:00pm 10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
T0172	Tin, Spot Urine**	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
T0173	TUNGSTEN, EDTA WHOLE BLOOD	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
T0174	Tungsten, Hair tests by ICPMS	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
T0175	TUNGSTEN, SPOT URINE	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
T0182	TB Check (MTB/RIF, AFB Culture, TB NGS) - Pulmonary	Chennai - Avadi	D	11800	Daily : 7:30am	After 6 weeks : 6:00pm	Pulmonary Smear positive samples or Culture positive specimen
T0184	TB Check (MTB/RIF, AFB Culture, TB NGS) - Extra Pulmonary	Chennai - Avadi	D	11800	Daily : 7:30am	After 6 weeks : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
T0192	TB Spectrum (MTB/RIF, AFB Culture, 13 Drugs Panel, TB NGS) - Pulmonary	Chennai - Avadi	D	20650	Daily : 7:30am	After 6 weeks : 6:00pm	Pulmonary Smear positive samples or Culture positive specimen
T0200	TRIO Clinical Exome sequencing - 26MB (80-100x) **	Chennai - Avadi	E	58800	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0201	TB Spectrum (MTB/RIF, AFB Culture, 13 Drugs Panel, TB NGS) - Extra Pulmonary	Chennai - Avadi	D	20650	Daily : 7:30am	After 6 weeks : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
T0502	TBC1D24-Gene Analysis	Chennai - Avadi	E	21000	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
T0504	TRIO-Whole exome sequencing Sure SelectV5-45-50MB (80-100x)	Chennai - Avadi	E	58800	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
T0505	TTN-Gene sequencing	Chennai - Avadi	E	19950	Daily : 7:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
T0506	Trisomy for Chromosome 18, X & Y by FISH, Ovarian follicular fluid **	Chennai - Avadi	D	6255	Daily : 11:30am	5th day : 6:00pm	CytogeneticsandFISH.pdf
T0507	Thrombophilia Super Maxi	Chennai	D	24780	Mon : 9:00am	Wed : 7:00pm	2 ml of citrated platelet poor

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi					plasma in frozen condition ,3 ml serum and 3 ml of EDTA whole blood is required.
T0938	TPMT and NUDT15 gene analysis**	Chennai - Avadi	R	14980	Daily : 11:30am	20th day : 6:00pm	4 ml EDTA Whole blood
T0987	Tricyclic antidepressants Drug(TCA) (Qualitative), Urine	Chennai - Avadi	D	565	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
T0988	Taurine DBS**	Chennai - Avadi	D	8850	Thu : 11:30am	Tue : 5:00pm	ried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. DBS card to be send in Ambient temperature only . Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age, or as close to discharge as possible if discharged prior to 24 hours of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper.The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Plz provide blood transfusion history and sample to be collected post 72 hrs after transfusion. Sample to be shipped at Frozen condition . Sample stable at refrigerated condition 2 to 8 degree for 1 day. Clinical history, metabolic disorders or inborn

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0989	2 step procedure for Screening and diagnosis of Gestational Diabetes Mellitus(ADA)	Chennai - Avadi	E	315	Daily : 9:00am to 9:00pm	After 6hrs	<p>errors of metabolsim suspicion or provisional diagnosis of clinician required fasting recommended not mandatory. Amino acid supplementation or medications needs to stoped based on physician recommendations</p> <p>Glucose Screen in pregnancy. ADA recommends 2 step procedure for Screening and diagnosis of Gestational Diabetes Mellitus (GDM). Step 1: A 50 g glucose challenge test (GCT) is used for screening without regard to the time of last meal or time of the day.Step 2: If 1 hour GCT value is more than 140 mg/dl, 100g Oral Glucose Tolerance Test (OGTT) is recommended and plasma glucose is estimated at 0, 1, 2 and 3 hours.G0027_50(post glucose 1 hr 50 gms) and G0066 (GTT-Gestational 3 hrs 100 gms glucose). If Post Glucose 1hr 50 gms value Above 140 then reflex G0066 GTT-Gestational ,Three Hour(100 grams glucose) taken.</p>
T1032	Tobramycin**	Chennai - Avadi	R	6315	Daily : 11:30am	5th day : 6:00pm	<p>3ml of Serum sample required to be shipped at 2-8C. Do not use gel tubes. Sample is stable for 1day, 7 days & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.</p>
T1170	BCR-ABL IS (International Scale) Transcript fusion Quantification by Real Time PCR - Blood	Chennai - Avadi	D	7080	Tue & Fri : 9:00am	3rd day : 6:00pm	<p>3 ml of Whole Blood in EDTA Vacutainer</p>

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T1171	BCR-ABL IS (International Scale) Transcript fusion Quantification by Real Time PCR - Bone Marrow	Chennai - Avadi	D	7080	Tue & Fri : 9:00am	3rd day : 6:00pm	3 ml of Bone Marrow in EDTA Vacutainer
T1901	Canine - T3 (Total)	Chennai - Avadi	B	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1902	Canine-T4-Free, Serum	Chennai - Avadi	B	270	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1903	Canine-T4 (Total)	Chennai - Avadi	B	185	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1904	Testosterone Total, Serum	Chennai - Avadi	B	565	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1905	CANINE-Thyroid panel-1	Chennai - Avadi	C	510	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							animal. Please correlate clinically and to be interpreted by the clinician.
T1906	CANINE-Thyroid panel - 2	Chennai - Avadi	C	735	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1908	Canine-TSH (Ultrasensitive), Serum	Chennai - Avadi	B	300	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1910	Canine - Triglyceride	Chennai - Avadi	A	220	Daily : 9:00am to 9:00pm	After 24hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1911	Thyroid cancer mutation panel, FFPE Tissue	Chennai - Avadi	E	35000	Daily : 11:30am	After 4 weeks : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
T5041	TruHealth Screen	Chennai - Avadi	W	499			See Individual Test Details
T5585	TruHealth Youth	Chennai - Avadi	W	4000			See Individual Test Details
T5586	TruHealth Master	Chennai - Avadi	W	4500			See Individual Test Details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T5587	TruHealth Senior	Chennai - Avadi	W	5000			See Individual Test Details
T5588	TruHealth Women	Chennai - Avadi	W	5000			See Individual Test Details
T5589	TruHealth Senior Women	Chennai - Avadi	W	5000			See Individual Test Details
T5610	Tru Health Heart Advance	Chennai - Avadi	W	6000	Mon & Thu : 7:30am	Same day : 5:00pm	5 ml of serum
T5611	Tru Health Heart Basic	Chennai - Avadi	W	1500	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of serum,2 ml of flouride plasma
T5612	TruHealth Advance	Chennai - Avadi	W	2800			See Individual Test Details
T5628	Tru Health Heart Pro	Chennai - Avadi	W	3500	Daily : 9:00am to 9:00pm	After 6hrs	5 ml of serum
T5705	T5705 : TruHeart Gold	Chennai - Avadi	W	2499	Daily : 9:00am	Next day : 6:00pm	4 ml Serum
T5709	TB Viram screening	Chennai - Avadi	D	6500	Daily : 9:00am	After 6 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer. 2ml of Fluoride Plasma. 2ml of Pulmonary or Extra Pulmonary specimen in Sterile leak proof container. Also note, kindly refer to individual test for the sample guideline.
T5710	TB Viram susceptibility	Chennai - Avadi	D	22595	Daily : 9:00am	After 4 weeks : 6:00pm	2ml of Sputum/BAL/other respiratory samples in Sterile leak proof container Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
T8006	Tissue processing	Chennai - Avadi	C	565	Daily : 9:00am	Next day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
T8500	Transplant Panel for Viral Infection	Chennai - Avadi	D	15000	Every Tuesday, Thursday and Saturday 9:00 AM	Same day 11:00 PM	2ml EDTA Plasma or Serum
U0001	U1-SnRNP (68 KDa)Serum	Chennai - Avadi	C	1750	Wed & Sat : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
U0002	UGT1A1 Gene Polymorphism (TA Repeat)	Chennai - Avadi	D	7080	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
U0003	Urea, Serum	Chennai - Avadi	A	150	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
U0004	Urea Urine 24H	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
U0005	Urea Breath Test For Helicobacter Pylori	Chennai - Avadi	C	2260	Daily : 9:00am	Same day : 6:00pm	NA
U0006	Urea Clearance Test Serum	Chennai - Avadi	A	545	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of 24 Hrs Urine & 2 ml of serum
U0007	Uric Acid Serum	Chennai - Avadi	A	200	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
U0008	Uric Acid Spot urine	Chennai - Avadi	A	545	Daily : 9:00am to 9:00pm	After 4hrs	10 ml of Spot urine
U0009	Uric Acid Urine 24H	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
U0010	Urobilinogen Urine	Chennai - Avadi	A	110	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
U0014	Urea Urine Spot	Chennai - Avadi	A	300	Daily : 9:00am to 9:00pm	After 4hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
U1801	Ultrasound (USG) by Appointment- Biophysical profile	Chennai - Avadi	O	900			
U1803	Ultrasound (USG) by Appointment- Follicular Study	Chennai - Avadi	O	2060			
U1804	Ultrasound (USG) by Appointment- Full Abdomen with pelvis	Chennai - Avadi	O	2520			
U1805	Ultrasound (USG) by Appointment- K.U.B	Chennai - Avadi	O	1855			
U1806	Ultrasound (USG) by Appointment- Musculoskeletal (Shoulder, Elbow, Ankle, Knee, etc)	Chennai - Avadi	O	925			
U1807	Ultrasound (USG) by Appointment- Pelvis	Chennai - Avadi	O	3915			
U1808	Ultrasound (USG) by Appointment- Pregnancy - 2nd trimester	Chennai - Avadi	O	3915			
U1809	Ultrasound (USG) by Appointment- Pregnancy - 3rd trimester	Chennai - Avadi	O	2265			
U1810	Ultrasound (USG) by Appointment- Pregnancy - early	Chennai - Avadi	O	2265			
U1811	Ultrasound (USG) by Appointment-	Chennai	O	2575			

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Pregnancy - level 3 Scan	- Avadi					
U1812	Ultrasound (USG) by Appointment- Small parts (Thyroid, Scrotum, etc)	Chennai - Avadi	O	1000			
U1813	Ultrasound (USG) by Appointment- Trans rectal with K.U.B	Chennai - Avadi	O	1030			
U1814	Ultrasound (USG) by Appointment- Upper abdomen	Chennai - Avadi	O	1030			
U8003	Urine Sugar Fasting	Chennai - Avadi	A	35	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
U8004	Urine sugar, 2 hrs post prandial	Chennai - Avadi	A	35	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
U8004_PP	Urine sugar, post glucose	Chennai - Avadi	A	35	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
U8005	Urine sugar, random	Chennai - Avadi	A	35	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
V0001	Valproic Acid Sodium Valporate Serum	Chennai - Avadi	C	930	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of Serum (Mention Time of Drug)
V0002	Vancomycin Serum	Chennai - Avadi	R	6420	Daily : 11:30am	4th day : 6:00pm	3 ml of Serum (Mention Time of Drug)
V0003	VDRL RPR by Flocculation - Serum	Chennai - Avadi	B	220	Daily : 9:00am & 3:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
V0004	VDRL by Flocculation - CSF	Chennai - Avadi	C	285	Daily : 9:00am & 3:00pm	After 6hrs	2ml of CSF in Sterile leak proof tube or container
V0005	Voltage Gated Calcium Channel (VGCC) Antibody**	Chennai - Avadi	R	9630	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
V0006	VGKC Antibody, CSF	Chennai - Avadi	D	9440	Mon & Thu : 11:30am	Same day : 9:00pm	2 ml of CSF in sterile screw capped container.
V0007	VGKC Antibody, Serum	Chennai - Avadi	D	9440	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0008	Vitamin A (Retinol), Serum	Chennai - Avadi	D	4720	Mon & Thu : 11:30am	Next day : 5:00pm	2 ml Serum Or EDTA Plasma In aluminium Foil
V0009	Vitamin B1 (Thiamine), EDTA Whole blood by HPLC (V0009)	Chennai - Avadi	D	3895	Wed & Sat : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood covered In Aluminium Foil
V0010	Vitamin B12 Cyanocobalamin Serum	Chennai	B	1070	Daily : 9:00am to	After 6hrs	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			9:00pm		ml serum for profiles and reflex tests.
V0011	Vitamin B12 Active Holo Transcobalamin Serum	Chennai - Avadi	C	1415	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0012	Vitamin B2 Riboflavin EDTA Whole Blood	Chennai - Avadi	D	3540	Wed : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood covered In Aluminium Foil
V0013	Vitamin B6 Pyridoxine EDTA Whole Blood	Chennai - Avadi	D	3500	Sat : 11:30am	Mon : 7:00pm	3 ml of EDTA whole blood covered In Aluminium Foil
V0014	Vitamin C, Serum/Heparin plasma by HPLC	Chennai - Avadi	D	3540	Tue & Fri: 11:30am	Next day : 5:00pm	2 ml Lithium heparin Plasma or serum In Fasting condition covered In Aluminium Foil
V0015	Vitamin D Total-25 Hydroxy Plasma/Serum	Chennai - Avadi	B	1600	Daily : 9:00am to 9:00pm	After 6hrs	2 ml of serum / EDTA Plasma
V0016	Vitamin D plus profile	Chennai - Avadi	A	1635	Daily : 9:00am to 9:00pm	After 6hrs	2 ml Serum / Edta Plasma
V0017	Vitamin D-1,25-Dihydroxy Serum	Chennai - Avadi	C	3615	Daily : 9:00am to 9:00pm	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0018	Vitamin E Tocopherol Serum	Chennai - Avadi	D	4720	Mon & Thu : 11:30am	Next day : 5:00pm	2 ml Serum Or EDTA Plasma In aluminium Foil
V0019	Vitamin profile	Chennai - Avadi	C	3165	Daily : 9:00am to 9:00pm	After 8hrs	3 ml Serum and 3 ml of EDTA plasma
V0020_LC	Vanillyl Mandelic Acid (VMA), LC MSMS, 24hrs urine	Chennai - Avadi	C	3730	Mon & Thu : 1:00pm	Next day : 5:00pm	10 ml 24 Hrs Urine
V0021	VRE Screen Vancomycin Resistant Enterococci by Vitek/ MALDITOF - Perianal swab	Chennai - Avadi	C	1015	Daily : 9:00am to 9:00pm	4th day : 6:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
V0022	Vw-Von Willebrand Antigen, Citrated plasma	Chennai - Avadi	D	10030	Mon : 9:00am	Same day : 8pm	1 ml of citrated platelet poor plasma in frozen condition.
V0023	VZV DNA Detection by Real Time PCR - Amniotic fluid	Chennai - Avadi	C	3390	Mon & Thurs : 9:00am	Next day : 5:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
V0024	VZV DNA Detection by Real Time PCR - Blood	Chennai - Avadi	C	3390	Mon & Thurs : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
V0025	VZV DNA Detection by Real Time PCR - CSF	Chennai - Avadi	C	3390	Mon & Thurs : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
V0026	VZV DNA Detection by Real Time PCR - CVS	Chennai - Avadi	C	3390	Mon & Thurs : 9:00am	Next day : 5:00pm	2ml of CVS in Sterile leak proof tube or container
V0027	VZV DNA Detection by Real Time PCR - Swab	Chennai - Avadi	C	3390	Mon & Thurs : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
V0029	VZV IgG Antibody by CLIA - Serum	Chennai - Avadi	C	1810	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
V0030	VZV IgM Antibody by CLIA - Serum	Chennai - Avadi	C	1810	Daily : 9:00am to 9:00pm	After 6hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
V0031	VDRL reflex Syphilis Total ABS - Serum	Chennai - Avadi	B	380	Mon & Thu : 11:30am	Next day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
V0032	Voriconazole EDTA Plasma	Chennai - Avadi	D	6135	Tue & Fri: 11:30am	Next day : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
V0034	VGKC - CSF & Serum	Chennai - Avadi	E	13125	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.2 ml of CSF in sterile screw capped container.
V0036	Vitamin B Complex profile	Chennai - Avadi	E	6300	Daily : 9:00am	7th day : 6:00pm	3 ml of serum and 3 ml of EDTA whole bloodIn aluminium Foil
V0037	Vitamin K1**	Chennai - Avadi	R	4280	Daily : 11:30am	15th day : 6:00pm	5 ml serum without gel tube. Wrap in aluminium foil. Frozen.
V0038	VANADIUM ICPMS BLOOD	Chennai - Avadi	R	2945	Daily : 11:30am	10th day : 6:00pm	3 ml EDTA Blood / Metal free EDTA Blood
V0144	VITAMIN D GOLD(covers Vitamin D2 , vitamin D3 & D total)	Chennai - Avadi	D	2950	Daily : 11:30am	3rd day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0148	Vasculitis profile, Serum	Chennai - Avadi	D	4720	Tue, Thu & Sat: 7:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0149	Vesicobullous Disorder profile	Chennai - Avadi	E	9450	Mon : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V7005	Vascular Endothelial Growth Factor (VEGF) serum**	Chennai - Avadi	R	6420	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml Serum
V8123	Vitamin B3 (Niacin)	Chennai	D	4250	Daily 7.30 am	7th Days	4 ml Serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
V8124	Vitamin B7 (Biotin)	- Avadi Chennai	D	4250	Daily 7.30 am	7th Days	4 ml Serum
V8125	Vitamin B5	- Avadi Chennai	D	4250	Daily 7.30 am	7th Days	4 ml Serum
W0001	Warfarin Sensitivity CYP2C9*2,CYP2C9*3 &VKORC1 Mutation Blood**	- Avadi Chennai	R	13055	Daily : 11:30am	10th day : 6:00pm	0
W0002	Water analysis Chemical**	- Avadi Chennai	R	1500	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
W0003	Water Culture for coliforms Water	- Avadi Chennai	C	680	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	50ml of Water in Sterile leak proof container
W0006	WBC:Total & Differential Counts, EDTA blood.	- Avadi Chennai	A	195	Daily : 9:00am to 9:00pm	After 4hrs	3 ml EDTA whole Blood
W0007	WBC:Total & Differential Counts, Body Fluid	- Avadi Chennai	A	275	Daily : 9:00am to 9:00pm	After 4hrs	Body fluid to be send in sterile EDTA container in refrigerated condition. Clinical history and clinician number required.Body fulids(CSF, Pleural, Peritoneal, Pericardial, synovial fluid)
W0009	WBC-Total Counts, Body Fluid.	- Avadi Chennai	A	220	Daily : 9:00am to 9:00pm	After 4hrs	Body fluid to be send in sterile EDTA container in refrigerated condition. Clinical history and clinician number required. Body fulids (CSF, Pleural, Peritoneal, Pericardial, synovial fluid)
W0010	Weil Felix Test for Rickettsia by Agglutination - Serum	- Avadi Chennai	B	920	Daily : 9:00am & 3:00pm	Next day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
W0011	Widal Test for Typhoid by Agglutination - Serum	- Avadi Chennai	B	350	Daily : 9:00am & 3:00pm	Next day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
W0013	WilsonS Disease Panel	- Avadi Chennai	D	4955	Mon & Thu : 9:00am	Same day : 5:00pm	3 ml of EDTA whole blood, 2 ml of flouride plasma random, 3 ml of serum, 20 ml aliquot of 24 hrs urine (Volume reqd)
W0014	West nile virus IgG by EIA - Serum	- Avadi Chennai	D	3835	1st & 3rd Fri : 7:30am	4th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
W0015	West nile virus IgM by EIA - Serum	Chennai	D	3835	1st & 3rd Fri :	4th day :	3ml of Serum in Gel Vacutainer or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi			7:30am	6:00pm	in Plain Sterile Vacutainer
W0017	West nile virus DNA Detection by Real Time PCR - EDTA Blood	Chennai - Avadi	D	5900	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
W0018	Whole mitochondrial genome sequencing**	Chennai - Avadi	R	14000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
W0019	WT1 MUTATION DETECTION**	Chennai - Avadi	R	8080	Daily : 11:30am	12th Day : 6:00pm	4 ML EDTA whole blood
W0020	Whole Exome panel**	Chennai - Avadi	R	25000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
W0021	Whole genome sequencing Mean 30x	Chennai - Avadi	R	95550	Daily : 11:30am	After 6 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
W0022	Whole Slide Imaging per slide	Chennai - Avadi	C	565	Daily : 9:00am	Next day : 6:00pm	Well mounted H &E stained Or IHC stained slide, Doctors TRF & E mail ID
W0023	Water chemistry panel - 3**	Chennai	R	4815	Daily : 11:30am	10th day :	2 litter of water in leakproof sterile

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		- Avadi				6:00pm	container. Storage and shipping at 2 to 8 degree celsius.
W0024	Water test profile, Routine	Chennai - Avadi	R	2675	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
W0025	Water test profile, Comprehensive	Chennai - Avadi	R	9630	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
W0026	Water chemistry panel - 1**	Chennai - Avadi	R	640	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
W1802_F	Wellness packages (Body watch) by Appointment- Fitness check-Female Below 45 years	Chennai - Avadi	W	3000			See Individual Test Details
W1802_M	Wellness packages (Body watch) by Appointment- Fitness check-Male Below 45 years	Chennai - Avadi	W	3000			See Individual Test Details
W1803_F	Wellness packages (Body watch) by Appointment- Life style check -Female Below 45 years	Chennai - Avadi	W	7500			See Individual Test Details
W1803_M	Wellness packages (Body watch) by Appointment- Life style check -Male Below 45 years	Chennai - Avadi	W	7500			See Individual Test Details
W1804_F	Wellness packages (Body watch) by Appointment- Vitality check-Female Below 45 years	Chennai - Avadi	W	15000			See Individual Test Details
W1804_M	Wellness packages (Body watch) by Appointment- Vitality check-Male Below 45 years	Chennai - Avadi	W	15000			See Individual Test Details
W1805_F	Wellness packages (Body watch) by Appointment- Fitness check-Female Above 45 years	Chennai - Avadi	W	4500			See Individual Test Details
W1805_M	Wellness packages (Body watch) by Appointment- Fitness check-Male Above 45 years	Chennai - Avadi	W	4500			See Individual Test Details
W1806_F	Wellness packages (Body watch) by Appointment- Life style check-Female Above 45 years	Chennai - Avadi	W	10500			See Individual Test Details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
W1806_M	Wellness packages (Body watch) by Appointment- Life style check-Male Above 45 years	Chennai - Avadi	W	10500			See Individual Test Details
W1807_F	Wellness packages (Body watch) by Appointment- Vitality check- Female Above 45 years	Chennai - Avadi	W	17500			See Individual Test Details
W1807_M	Wellness packages (Body watch) by Appointment- Vitality check-Male Above 45 years	Chennai - Avadi	W	17500			See Individual Test Details
X0001	X And Y (Sex Mismatch) by FISH	Chennai - Avadi	D	4840	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
X0002	X And Y (Sex Mismatch) by FISH Bone Marrow Transplant Bone Marrow	Chennai - Avadi	D	4840	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
X8011	XDR Rapid plus (For extra pulmonary samples)	Chennai - Avadi	D	7200	Daily : 7:30am	Same day : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
X8012	XDR Rapid plus (For pulmonary samples)	Chennai - Avadi	D	7200	Daily : 7:30am	Same day : 6:00pm	2ml of Pulmonary specimen in Sterile leak proof container
Y0001	Y-Chromosome Microdeletion Analysis	Chennai - Avadi	D	8025	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical history and TRF mandatory.
Y0002	Yersinia Enterocolitica Culture by Vitek/ MALDITOF - Stool	Chennai - Avadi	C	1695	Daily : 9:00am to 9:00pm	4th day : 6:00pm	5gm of Stool specimen in Sterile leak proof container
Y0005	Yo Antibody Neuronal antibody for Paraneoplastic syndromes Serum	Chennai - Avadi	D	5310	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
Z0001	ZAP-70	Chennai - Avadi	D	6490	Daily : 11:30am	4th day : 6:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							required.
Z0005	Zinc Urine	Chennai - Avadi	D	3420	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
Z0006	Zinc, Blood	Chennai - Avadi	D	2005	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
Z0007	Zinc, Serum	Chennai - Avadi	D	2300	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
Z0009	Zinc Urine 24H**	Chennai - Avadi	D	2005	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
Z0011	Zika Virus RNA Detection by Real Time PCR - Serum	Chennai - Avadi	D	7670	Daily : 10:00am	After 24hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
Z0012	Zinc , Water tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
Z0013	Zinc, Hair tests by ICPMS**	Chennai - Avadi	R	1925	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
Z0014	Zika Virus RNA Detection by Real Time PCR - Urine	Chennai - Avadi	D	7670	Daily : 10:00am	After 24hrs	15ml of Spot Urine in Sterile leak proof container
Z0015	Zinc Transporter 8 (ZnT8) Antibody	Chennai - Avadi	D	3540	Mon & Thu : 9:00 am	Next day : 5:00pm	3 ml Serum
Â I0390	IHC stain only HHV 8	Chennai - Avadi	D	#N/A	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
Â I0392	IHC stain only MDM 2	Chennai - Avadi	D	#N/A	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
Â I0393	IHC single marker with reporting MDM2	Chennai - Avadi	D	#N/A	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail