

METROPOLISEDOS_Madhya Pradesh

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
#0002	1p/19q Deletion by FISH Oligodendrogliomas Tissue	Madhya Pradesh	D	12305	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
#0003	3-Methyl Histidine Quantitative Urine	Madhya Pradesh	D	6635	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
#0004	3-Methyl Histidine Quantitative Serum	Madhya Pradesh	R	6325	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
#0006	5-HIAA, Urine 24 hrs by column chromatography	Madhya Pradesh	D	3480	Tue & Fri : 11:30am	Same day : 9:00pm	24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
#0007	5q del (MDS) by Fish	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0008	5q deletion/Monosomy 5 by FISH MDS Bone Marrow	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0011	7q deletion/Monosomy 7 by FISH,MDS	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
#0012	7q deletion/Monosomy 7 by FISH MDS Bone Marrow	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0013	Williams Syndrome; Microdeletion of chromosome 7q11.23**	Madhya Pradesh	R	10710	Daily : 11:30am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
#0014	FGFR 1 gene rearrangement t(8p11.2), Blood**	Madhya Pradesh	E	4900	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0015	FGFR 1 gene rearrangement t(8p11.2), Bone Marrow**	Madhya Pradesh	E	4900	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0017	11q (ATM) by FISH	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0018	11q deletion by FISH ATM (CLL) Bone Marrow	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0021	13q del by FISH	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0022	13q deletion by FISH MM, CLL Bone Marrow	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0023	14q32.3 by FISH IGH gene rearrangement (ALL, NHL) Blood**	Madhya Pradesh	D	6090	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0024	14q32.3 by FISH, IGH gene rearrangement (ALL, NHL) Bone Marrow**	Madhya Pradesh	D	5885	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0025	Prader-Willi/Angeman Syndrome (Microdeletion 15q9-13)**	Madhya Pradesh	R	7040	Daily : 11:30am	12th Day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
#0026	17 OH Corticosteroids, Urine by Column Chromatography (#0026)	Madhya Pradesh	D	4385	Tue & Fri : 11:30am	Same day : 9:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
#0027	17-Hydroxy Progesterone Neonatal Screen Dried blood spot	Madhya Pradesh	C	505	Daily : 11:30am	Next day : 5:00pm	accurate results and such history to be conveyed to lab. Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
#0028	17-Hydroxy-Progesterone by Serum	Madhya Pradesh	C	1680	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
#0029	17 Ketosteroids,urine by Column Chromatography (#0029)	Madhya Pradesh	C	4255	Tue & Fri : 11:30am	Same day : 9:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
#0030	Del 17p53 by FISH	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0031	17p deletion by FISH, p53 (MM, CLL) Blood	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0034	20q deletion (MDS) by FISH, Blood FISH probe for chromossome 20q Blood	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
#0035	Del 20q (MDS) Blood / bone marrow	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
#0036	DEL 22q (Microdeletion digeorge syndrome) by FISH	Madhya Pradesh	D	5670	Daily : 9:00am	7th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
#1018	5 Nucleotidase activity , Serum**	Madhya Pradesh	R	6120	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
A0001	Abnormal Haemoglobin Studies, Hb Variants, EDTA Blood	Madhya Pradesh	C	1060	Daily: 9am	Next day 5pm	3 ml EDTA whole blood
A0003	Absolute Eosinophile Count	Madhya Pradesh	A	205	Daily: 9am to 9pm	Absolute Eosinophile Count	3 ml EDTA whole blood
A0008	Acetone Ketone by Fully Automated Urine chemistry & Microscopy - Urine	Madhya Pradesh	A	80	Daily : 9:00am to 9:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
A0012	AChR Antibodies Acetyl Choline Receptor Antibodies Serum	Madhya Pradesh	C	3650	Daily : 11:30am	Next day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0013	Acid Phosphatase Total Serum	Madhya Pradesh	A	325	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0014	Acid Phosphatase Total and prostatic fraction Serum	Madhya Pradesh	C	645	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0015	ACTH Stimulation Test - For 17 OH progesterone -Serum	Madhya Pradesh	C	4060	Daily : 11:30am	After 8hrs	3 ml of serum (0 min, 30, 60 min)
A0016	ACTH Stimulation Test - For Aldosterone -	Madhya	C	6160	Daily: 7.30am	Next day 5pm	3 ml of serum (0 min, 30, 60 min)

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Serum	Pradesh					
A0017	ACTH Stimulation Test - For Androstenedione - Serum	Madhya Pradesh	C	5040	Daily: 7.30am	After 8 hrs	3 ml of serum (0 min, 30, 60 min)
A0018	ACTH Stimulation Test - For Cortisol - Serum	Madhya Pradesh	C	1960	Daily: 9am to 9pm	After 8 hrs	3 ml of serum (0 min, 30, 60 min)
A0019	ACTH Stimulation Test - For DHEAS-Serum	Madhya Pradesh	C	3080	Daily: 9am to 9pm	After 8 hrs	3 ml of serum (0 min, 30, 60 min)
A0020	ACTH Adreno Corticotropic Hormone - Plasma	Madhya Pradesh	C	1790	Daily: 9am to 9pm	After 6 hrs	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
A0021	Actinomyces Culture by Vitek/MALDITOF - Sputum	Madhya Pradesh	D	1390	Daily : 9:00am to 9:00pm	5th day : 5pm	2ml of Sputum in Sterile leak proof container
A0022	Actinomyces Culture by Vitek/MALDITOF - Tissue	Madhya Pradesh	D	1390	Daily : 9:00am to 9:00pm	5th day : 5pm	Tissue with Normal Saline in Sterile leak proof container
A0023	Actitest and Fibrotest	Madhya Pradesh	E	16500	Daily: 7am	5th day	5 ml of serum
A0024	Acylcarnitine quantitative Neonatal**	Madhya Pradesh	R	4590	Daily : 11:30am	8th day : 6:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. DBS card to be send in Ambient temperature only . Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age, or as close to discharge as possible if discharged prior to 24 hours of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper.The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Plz provide blood transfusion histroy and sample to be collected post 72 hrs after transfusion.Clinical history, metabolic disorders or inborn errors of metabolsim suspicion or provisional diagnosis of clincian required fasting recommended not mandatory. Amino acid supplementation or medications needs to stoped based on physician recommendations
A0025	ADA - Adenosine Deaminase by Biochemical Method - Ascitic Fluid	Madhya Pradesh	B	870	Daily : 9:00am to 9:00pm	After 6 hrs	3 ml of Fluid in Sterile leak proof container
A0026	ADA - Adenosine Deaminase by Biochemical Method - CSF	Madhya Pradesh	B	870	Daily : 9:00am to 9:00pm	After 6 hrs	2 ml of CSF in Sterile leak proof tube or container
A0027	ADA (Adenosine Deaminase), Pericardial fluid	Madhya Pradesh	B	870	Daily : 9:00am to 9:00pm	After 6 hrs	3 ml of Fluid in Sterile leak proof container
A0028	ADA - Adenosine Deaminase by Biochemical Method- Peritoneal Fluid	Madhya Pradesh	B	870	Daily : 9:00am to 9:00pm	After 6 hrs	3 ml of Fluid in Sterile leak proof container
A0029	ADA - Adenosine Deaminase by Biochemical Method - PleuralFluid	Madhya Pradesh	B	870	Daily : 9:00am to 9:00pm	After 6 hrs	3 ml of Fluid in Sterile leak proof container
A0030	ADA - Adenosine Deaminase by Biochemical Method - Pus	Madhya Pradesh	B	870	Daily : 9:00am to 9:00pm	After 6 hrs	Sample in sterile container
A0031	ADA - Adenosine Deaminase by Biochemical Method - Serum	Madhya Pradesh	B	870	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0032	Adenovirus DNA Detection by PCR - Respiratory Fluid **	Madhya Pradesh	D	4655	Daily : 11:30am	10th day : 6:00pm	3 ml of Fluid in Sterile leak proof container
A0033	Adenovirus DNA Detection by PCR - Oral swab**	Madhya Pradesh	D	4815	Daily : 11:30am	10th day : 6:00pm	1 Swab with Sufficient Material in Sterile Vacutainer
A0034	Adenovirus DNA Detection by PCR - Eye swab**	Madhya Pradesh	D	4655	Daily : 11:30am	10th day : 6:00pm	1 Swab with Sufficient Material in Sterile Vacutainer
A0035	Adenovirus DNA Detection by PCR - Nasal swab **	Madhya Pradesh	D	4815	Daily : 11:30am	10th day : 6:00pm	1 Swab with Sufficient Material in Sterile Vacutainer
A0036	Adenovirus IgG Antibody by EIA - Serum	Madhya	D	1710	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0037	Adenovirus IgM Antibody by EIA - Serum	Pradesh Madhya Pradesh	D	1710	Thurs : 9:00am	Same day : 5pm	in Plain Sterile Vacutainer 3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0038	ADH-Anti Diuretic Hormone Vasopressin Plasma	Madhya Pradesh	D	4815	Mon : 11:30am	Next day : 5:00pm	3 ml of EDTA plasma to be sent in aliquot tube at proper temperature.
A0039	Adiponectin ELISA Serum	Madhya Pradesh	D	2995	Wed : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0041	Adrenaline (Epinephrine) Plasma	Madhya Pradesh	D	3315	Wed : 11:30am	Next day : 5:00pm	Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
A0042	Adrenaline (Epinephrine] HPLC,Urine24hrs	Madhya Pradesh	D	3315	Wed : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
A0052	AFB Rapid culture by Bactec Fx - Blood	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Upto 6 weeks	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
A0053	AFB Rapid culture by MGIT - Body Fluid	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Upto 6 weeks	2ml of Body Fluid in Sterile leak proof container
A0054	AFB Rapid culture by Bactec Fx - Bone	Madhya Pradesh	C	1345	Daily : 9:00am to	Upto 6 weeks	3ml of Bone Marrow in MycoF/lytic

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0055	Marrow AFB Rapid culture by MGIT - BAL & Other Respiratory samples	Pradesh Madhya Pradesh	C	1345	9:00pm Daily : 9:00am to 9:00pm	Upto 6 weeks	Bactec Bottle or Sodium Heparin 2ml of BAL / Other Respiratory samples in Sterile leak proof container
A0056	AFB Rapid culture by MGIT - CSF	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Upto 6 weeks	2ml of CSF in Sterile leak proof tube or container
A0057	AFB Rapid culture by MGIT - Endometrium	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Upto 6 weeks	2ml of Endometrium in Sterile leak proof container
A0058	AFB Rapid culture by MGIT - Other	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Upto 6 weeks	2ml of Specimen in Sterile leak proof container
A0059	AFB Rapid culture by MGIT - Pus, abscess and aspirates	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Upto 6 weeks	2ml of Pus / Abscess / Aspirate in Sterile leak proof container
A0060	AFB Rapid Culture by MGIT - Sputum	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Upto 6 weeks	2ml of Sputum in Sterile leak proof container
A0061	AFB Rapid culture by MGIT - Tissue/Biopsy	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Upto 6 weeks	Tissue with Normal Saline in Sterile leak proof container
A0062	AFB Rapid Culture by MGIT - Urine	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Upto 6 weeks	15ml of Spot Urine in Sterile leak proof container
A0063	MOTT Antibiogram for Rapid Growers by Microbroth dilution - Pure Culture	Madhya Pradesh	D	8560	Daily : 7:30am	Upto 4 weeks	Isolated colony on culture plates
A0064	MOTT Antibiogram for Slow Growers by Microbroth dilution - Pure Culture	Madhya Pradesh	D	10270	Daily : 7:30am	Upto 4 weeks	Isolated colony on culture plates
A0065	AFB (TB) - Antibiogram for MTb complex - Amikacin Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0066	AFB (TB) - Antibiogram for MTb complex - Capreomycin Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0067	AFB (TB) - Antibiogram for MTb complex - Clofazimine Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0068	AFB (TB) - Antibiogram for MTb complex - Ethambutol Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0069	AFB (TB) - Antibiogram for MTb complex - Ethionamide Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0070	AFB (TB) - Antibiogram for MTb complex - Isoniazid Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0071	AFB (TB) - Antibiogram for MTb complex - Kanamycin Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0072	AFB (TB) - Antibiogram for MTb complex - Levofloxacin Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0073	AFB (TB) - Antibiogram for MTb complex - Moxifloxacin Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0074	AFB (TB) - Antibiogram for MTb complex - Ofloxacin Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0075	AFB (TB) - Antibiogram for MTb complex - PAS Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0076	AFB (TB) - Antibiogram for MTb complex - Pyrazinamide Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2910	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0077	AFB (TB) - Antibiogram for MTb complex - Rifabutin Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0078	AFB (TB) - Antibiogram for MTb complex - Rifampicin Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0079	AFB (TB) - Antibiogram for MTb complex - Streptomycin Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0080	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Other Sample Type	Madhya Pradesh	C	2295	Daily : 7:30am	Next day : 5:00pm	2ml of Sputum/ BAL/ Other respiratory samples in Sterile leak proof container
A0081	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Body Fluids	Madhya Pradesh	C	2295	Daily : 7:30am	Next day : 5:00pm	2ml of Body Fluid in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0082	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - CSF	Madhya Pradesh	C	2295	Daily : 7:30am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
A0083	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Pus / Abscess / Aspirate	Madhya Pradesh	C	2295	Daily : 7:30am	Next day : 5:00pm	2ml of Pus / Abscess / Aspirate in Sterile leak proof container
A0084	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Semen	Madhya Pradesh	C	2295	Daily : 7:30am	Next day : 5:00pm	2ml of Specimen in Sterile leak proof container
A0085	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Tissue / Biopsy / Lymph nodes	Madhya Pradesh	C	2295	Daily : 7:30am	Next day : 5:00pm	Tissue/ Endometrium with Normal Saline/ Biopsy/ Lymph nodes in Sterile leak proof container
A0086	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Urine	Madhya Pradesh	C	2295	Daily : 7:30am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
A0087	AFB Drug susceptibility panel 1st and 2nd Line (10 drugs) by MGIT - Pure Culture	Madhya Pradesh	D	10165	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0088	AFB Drug susceptibility panel 1st Line (4 drugs) SIRE by MGIT - Pure Culture	Madhya Pradesh	D	5030	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0089	AFB Drug susceptibility panel 2nd Line (4 Drugs) KEPO by MGIT - Pure Culture	Madhya Pradesh	D	6205	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0090	AFB Drug susceptibility panel Additional drugs 4 Nos. MACC by MGIT - Pure Culture	Madhya Pradesh	D	5990	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0091	AFB Drug susceptibility panel Comprehensive Panel 13 drugs by MGIT - Pure Culture	Madhya Pradesh	D	13910	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0092	AFB Identification & Speciation MTB and MOTT by Hains Line Probe - Pure culture	Madhya Pradesh	D	3745	Tue : 7:30am	Next day : 5pm	Pure growth of Mycobacterium spp. on LJ slant/ MGIT Tube/ Bactec alert bottle.
A0093	AFB MDR Screen by Hains Line Probe Assay - Sputum, BAL, other respiratory samples	Madhya Pradesh	D	2035	Tue : 7:30am	Next day : 5pm	2ml of Sputum/BAL/other respiratory samples in Sterile leak proof container
A0094	AFB MDR Screen by Hains Line Probe Assay - Extra pulmonary samples	Madhya Pradesh	D	2035	Tue : 7:30am	Next day : 5pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, Synovial fluid)
A0102	AFB-Smear examination by Fluorescent stain - Body fluid	Madhya Pradesh	C	475	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Body Fluid in Sterile leak proof container
A0103	AFB-Smear examination by Fluorescent stain - BAL & Other Respiratory Sample	Madhya Pradesh	C	475	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of BAL / Other Respiratory samples in Sterile leak proof container
A0104	AFB-Smear examination by Fluorescent stain - Pus, abscess and aspirates	Madhya Pradesh	C	475	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Pus / Abscess / Aspirate in Sterile leak proof container
A0105_1	AFB-Smear examination by Fluorescent stain - Sputum (1 sample)	Madhya Pradesh	C	475	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Sputum in Sterile leak proof container
A0105_3	AFB-Smear examination by Fluorescent stain - Sputum (3 sample)	Madhya Pradesh	C	1345	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Sputum in Sterile leak proof container
A0106	AFB-Smear examination by Fluorescent stain - Tissue and Biopsy	Madhya Pradesh	C	475	Daily : 9:00am & 3:00pm	After 6 hrs	Tissue with Normal Saline / Biopsy / Lymph nodes in Sterile leak proof container
A0107_1	AFB-Smear examination by Fluorescent stain - Urine (1 sample)	Madhya Pradesh	C	475	Daily : 9:00am & 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
A0107_3	AFB-Smear examination by Fluorescent stain - Urine (3 sample)	Madhya Pradesh	C	1345	Daily : 9:00am & 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
A0108	AFB-Smear examination by ZN Stain - Body fluids	Madhya Pradesh	B	410	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Body Fluid in Sterile leak proof container
A0109	AFB- Smear examination by ZN Stain - BAL & Other Respiratory samples	Madhya Pradesh	B	410	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of BAL / Other Respiratory samples in Sterile leak proof container
A0110	AFB-Smear examination by ZN Stain - Pus, abscess and aspirates	Madhya Pradesh	B	410	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Pus / Abscess / Aspirate in Sterile leak proof container
A0111_1	AFB-Smear examination by ZN Stain - Sputum (1 sample)	Madhya Pradesh	B	410	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Sputum in Sterile leak proof container
A0111_3	AFB-Smear examination by ZN Stain - Sputum (3 samples)	Madhya Pradesh	B	870	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Sputum in Sterile leak proof container
A0112	AFB- Smear examination by ZN Stain - Tissue / Biopsy	Madhya Pradesh	B	410	Daily : 9:00am & 3:00pm	After 6 hrs	Tissue with Normal Saline in Sterile leak proof container
A0113_1	AFB- Smear examination by ZN Stain - Spot Urine (Single Sample)	Madhya Pradesh	B	410	Daily : 9:00am & 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0113_3	AFB-Smear examination by ZN Stain - Urine (3 samples)	Madhya Pradesh	B	870	Daily : 9:00am & 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
A0114	AFB Speciation (MOTT) Additional species by Hains Line Probe Assay - Pure culture	Madhya Pradesh	C	4815	Tue : 7:30am	Next day : 5pm	Pure growth of Mycobacterium spp. on LJ slant/ MGIT Tube/ Bactec alert bottle.
A0115	AFB XDR Screen by Hains Line Probe Assay - Sputum / BAL/ Other respiratory samples	Madhya Pradesh	C	6330	Tue : 7:30am	Next day : 5pm	2ml of Sputum/BAL/other respiratory samples in Sterile leak proof container
A0116	AFB XDR Screen by Hains Line Probe Assay - Extra pulmonary samples	Madhya Pradesh	C	5770	Tue : 7:30am	Next day : 5pm	Extra - Pulmonary Smear positive samples or Culture positive specimen
A0117	AFB-Xpert panel (M.Tb Detection & Rifampicin resistance) by Genexpert - Pulmonary samples	Madhya Pradesh	D	2200	Daily : 7:30am	Same day : 6:00pm	2ml of Pulmonary specimen in Sterile leak proof container
A0118	AFB-Xpert panel (M.Tb Detection & Rifampicin resistance) by Genexpert - Extra pulmonary samples	Madhya Pradesh	D	2400	Daily : 7:30am	Same day : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
A0119	AFP-Alpha Feto Protein- CLIA- Serum	Madhya Pradesh	B	970	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0120	Alpha Fetoprotein - Amniotic Fluid (AFP-AF)	Madhya Pradesh	C	1230	Daily: 9am to 9pm	After 6 hrs	2 ml Amniotic Fluid In Sterile container
A0121	AFP-Alpha Feto Protein -CLIA- CSF	Madhya Pradesh	C	1120	Daily: 9am to 9pm	After 8 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
A0122	Albert stain by smear examination - Throat swab	Madhya Pradesh	B	380	Daily : 9:00am to 3:00pm	After 6 hrs	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
A0123	Albumin Qualitative - Urine	Madhya Pradesh	A	80	Daily : 9:00am to 9:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
A0124	Albumin Quantitative Urine 24H	Madhya Pradesh	B	630	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
A0125	Albumin Serum	Madhya Pradesh	A	195	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0126	Albumin CSF	Madhya Pradesh	C	1345	Daily: 9am to 9pm	After 8 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
A0127	Albumin/Creatinine Ratio, Urine	Madhya Pradesh	A	590	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
A0128	Alcohol Ethanol Urine	Madhya Pradesh	C	1455	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
A0129	Alcohol Ethanol Serum	Madhya Pradesh	C	1120	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0130	Aldehyde Chopra Test for Kala Azar by Chemical Formal Test - Serum**	Madhya Pradesh	R	460	Daily : 11:30am	3rd day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0131	Aldolase enzymatic Serum	Madhya Pradesh	C	1110	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0132	Aldosterone 24H URINE	Madhya Pradesh	C	2465	Daily : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
A0133	Aldosterone ELISA Serum	Madhya Pradesh	C	2185	Daily : 9:00am	Same day after 8 hrs	3 ml Serum (Mention The Position)
A0134	ALK1/EML4 (LUNG CANCER) by FISH	Madhya Pradesh	E	8500	Daily 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
A0135	Alkaline Phosphatase Serum	Madhya Pradesh	A	195	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0136	Alkaline Phosphatase Isoenzymes Electrophoresis - Serum	Madhya Pradesh	D	5350	Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0137	Alkaline Phosphatase with Bone Fraction Ostase Serum	Madhya Pradesh	C	2015	Mon & Thu : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0138	Allergen,Individual-AnimalCat DanderSerum	Madhya Pradesh	C	1455	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0140	Allergen,Individual-AnimalCow DanderSerum	Madhya Pradesh	C	1455	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0141	Allergen,Individual-AnimalDog DanderSerum	Madhya Pradesh	C	1455	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0143	Allergen,Individual-Animal Guinea Pig Epithelium Serum	Madhya Pradesh	C	1455	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0145	Allergen,Individual-AnimalHorse Dander Serum	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0152	Allergen,Individual-DrugACTH	Madhya Pradesh	C	1850	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0158	Allergen,Individual-DrugCiprofloxacinantibiotic	Madhya Pradesh	C	1680	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0173	Allergen, Individual-Drug Sulpha antibiotic	Madhya Pradesh	C	1680	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0175	Allergen,Individual-DustDermatophagoides FarinaeHouse Dust Mite	Madhya Pradesh	C	1455	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0176	Allergen,Individual-DustDermatophagoides PteronyssinusHouse Dust Mite	Madhya Pradesh	C	1455	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0178	Allergen,Individual-Dust House Dust - Hollistier	Madhya Pradesh	C	1455	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0179	Allergen,Individual-DustHouse Dust Greer Labs	Madhya Pradesh	C	1455	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0180	Allergen,Individual-FoodAlmondBadam (Amygdalus Communis)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0181	Allergen,Individual-FoodAlpha - LactalbuminMilkMilk protein	Madhya Pradesh	C	1570	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0182	Allergen,Individual-FoodAppleSafarchand (Malus X Domestica)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0185	Allergen,Individual-FoodBaingan	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0186	Allergen,Individual-FoodBakers YeastYeast (Saccharomyces Cerevisiae)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0187	Allergen,Individual-FoodBananaKela (Musa Spp.)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0189	Allergen,Individual-Food Basil Ocimum basilicum	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0191	Allergen,Individual-FoodBeefCow meat (Bos spp.)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0192	Allergen,Individual-FoodBeta - Lactoglobulin- MilkMilk protein	Madhya Pradesh	C	1570	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0193	Allergen,Individual-FoodBlack PepperPiper nigrum	Madhya Pradesh	C	1570	Mon & Thu : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0194	Allergen,Individual-FoodBlue MusselMytilus Edulis	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0195	Allergen,Individual-FoodBrazil NutBertholletia excelsa	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0198	Allergen,Individual-FoodCabbageGobi (Brassica Oleraceavar.Capitata)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0199	Allergen,Individual-FoodCardamonElettaria cardamomum	Madhya Pradesh	C	1570	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0201	Allergen,Individual-FoodCasein(Milk Protein)Dahi	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0202	Allergen,Individual-FoodCashew NutKaju (Anacardium Occidentale)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0204	Allergen,Individual-Food Celery Apium Graveolens	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0205	Allergen,Individual-FoodCheese Cheddar	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0207	Allergen,Individual-FoodChick PeaKabuli Chana (Cicer Arietinus)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0208	Allergen,Individual-FoodChicken MeatMurgi	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0209	Allergen,Individual-FoodChilipepperCapsicum frutescens	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0211	Allergen,Individual-FoodClove(Syzygium aromaticum)	Madhya Pradesh	C	1570	Mon & Thu : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0212	Allergen,Individual-FoodCocoaChocolate (Theobroma Cacao)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0213	Allergen,Individual-FoodCoconutNarial (Cocos Nucifera)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0214	Allergen,Individual-FoodCod FishMachhli (Gadus Morhua)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0215	Allergen,Individual-FoodCoffee	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0216	Allergen,Individual-FoodCoriander	Madhya Pradesh	C	1570	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0217	Allergen,Individual-FoodCorn / MaizeMakka (Zea Mays)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0218	Allergen,Individual-FoodCow MilkDoodh	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0219	Allergen,Individual-FoodCrabCancer pagurus	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0224	Allergen,Individual-FoodEgg WhiteAnda	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0225	Allergen,Individual-FoodEgg YolkAnda	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0228	Allergen,Individual-FoodGarlicLehsoon (Allium Sativum)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0230	Allergen,Individual-FoodGluten	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0234	Allergen,Individual-FoodGreen Pea	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0236	Allergen,Individual-FoodHazelnutCorylus avellana	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0238	Allergen,Individual-FoodKiwi FruitActinidia deliciosa	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0242	Allergen,Individual-FoodLemonNimbu (Citrus Limon)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0243	Allergen,Individual-FoodLentilMasoor Dal (Lens Esculenta)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0245	Allergen,Individual-FoodLobsterHomarus gammarus	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0249	Allergen,Individual-FoodMangoMangifera Indica	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0252	Allergen,Individual-FoodMushroomchampignon / Agaricus hortensis	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0253	Allergen,Individual-FoodMustardSarson (Brassica / Sinapsis Spp.)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0254	Allergen,Individual-FoodMuttonOvis Spp.	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0256	Allergen,Individual-FoodOliveOlea Europaea	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0257	Allergen,Individual-FoodOnionPiaz (Allium Cepa)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0258	Allergen,Individual-FoodOrangesantra (Citrus sinensis)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0260	Allergen,Individual-FoodPapayaCarica Papaya	Madhya Pradesh	C	1680	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0261	Allergen,Individual-FoodPapricaSimla Mirchi (Capsicum Annuum)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0263	Allergen,Individual-FoodPeachPrunus persica	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0264	Allergen,Individual-FoodPeanutMoongfali (Arachis Hypogaea)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0265	Allergen,Individual-FoodPearPyrus communis	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0266	Allergen,Individual-Food Pecan nut Carya illinoensis	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0267	Allergen,Individual-FoodPineappleAnanas	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0268	Allergen,Individual-FoodPistachio	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0269	Allergen,Individual-FoodPlaiceFish (Pleuronectes platessa)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0272	Allergen,Individual-FoodPorkPig meat	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0273	Allergen,Individual-FoodPotatoAaloo (Solanum Tuberosum)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0277	Allergen,Individual-FoodRiceChawal (Oryza Sativa)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0278	Allergen,Individual-FoodSalmonRawas fish (Salmo Salar)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0281	Allergen,Individual-FoodShrimp / Prawns	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0282	Allergen,Individual-FoodSoyabeanGlycine Max	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0283	Allergen,Individual-FoodSpinachPalak (Spinachia Oleracea)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0284	Allergen,Individual-FoodStrawberryFragaria vesca	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0289	Allergen,Individual-FoodTomatoTamatar (Lycopersicon Lycopresicum)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0290	Allergen,Individual-FoodTunaFish (Pleuronectes platessa)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0292	Allergen,Individual-FoodWalnutJuglans californica	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0294	Allergen,Individual-FoodWheatGehon (Triticum Aestivum)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0295	Allergen,Individual-FoodWhite BeanPhaseolus Vulgaris	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0298	Allergen,Individual-Grass pollenBermuda GrassDurva (Cynodon Grass)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0300	Allergen,Individual-Grass pollen Cultivated RyeGrass Secale Cereale	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0302	Allergen,Individual-Grass pollenJohnson GrassJowar grass (Sorghum Grass)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0304	Allergen,Individual-Grass pollen Meadow Grass Poa Pratensis	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0308	Allergen,Individual-Grass pollenTimothy GrassGhass (Phleum Pratense)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0311	Allergen,Individual-InsectCockroach GermanBlatella Germanica	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0312	Allergen,Individual-InsectHoney Bee VenomMadhu Makhhi (Apis Mellifera)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0315	Allergen,Individual-Micro organismAlternaria Teniusfungus	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0316	Allergen,Individual-Micro organismAspergillus Fumigatusfungus	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0317	Allergen,Individual-Micro organismCandida Albicansfungus	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0318	Allergen,Individual-Micro organismCladosporium Herbarumfungus	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0320	Allergen, Individual - Insect Mosquito Machhar (Aedes Communis)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0322	Allergen,Individual-Micro organismPenicillium Notatumfungus (Penicillium chrysogenum)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0338	Allergen,Individual-Tree pollenAcasiaBabool (Acasia Longifolia)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0339	Allergen,Individual-Tree pollenAlderAlnus Incana	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0340	Allergen,Individual-Tree pollenBirchBetula Verrucosa	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0342	Allergen,Individual-Tree pollenEucalyptusNilgiri (Safeda)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0345	Allergen,Individual-Tree pollenMesquitePahadi Keekar (Prosopis Juliflora)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0346	Allergen,Individual-Tree pollen Mulberry Shahtoot (Morus alba)	Madhya Pradesh	C	1570	Daily : 11:30am	5th day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0347	Allergen,Individual-Tree pollen Oak Quercus Alba	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0348	Allergen,Individual-Tree pollenWhite PinePinus strobus	Madhya Pradesh	C	1570	Daily : 11:30am	5th day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0349	Allergen,Individual-Tree pollen Willow Salix Caprea	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0350	Allergen,Individual-Weed pollenCockleburXanthium commune	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0351	Allergen,Individual-Weed pollen Common Pigweed Kaantewali Chauli (Amaranthus Retroflexus)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0352	Allergen,Individual-Weed pollenCommon RagweedClose To Parthenium / Congress Grass (Ambrosia Elatior)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0353	Allergen,Individual-Weed pollenEnglish PlantainPlantago Lanceolata	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0355	Allergen,Individual-Weed pollenGoose FootChenopodium Album	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0356	Allergen,Individual-Weed pollenMugwortSita Bani (Artemisia Vulgaris / BanoBarna)	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0360	Allergy - Fruit Panel	Madhya Pradesh	D	6740	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0361	Allergy Fx 5 Panel	Madhya Pradesh	D	5245	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0362	Allergy - Milk & Milk Product Panel	Madhya Pradesh	D	4495	Daily : 11:30am	Next day : 6:00pm	5 ml of serum
A0363	Allergy - Non-vegetarian Panel	Madhya Pradesh	D	7650	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0364	Allergy-Nuts Panel	Madhya Pradesh	D	6675	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0367	Allergy - Spices Panel	Madhya Pradesh	D	5990	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0368	Allergy - Vegetable-1 Panel	Madhya Pradesh	D	6740	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0370	Allergy - Vegetarian + Nonvegetarian Panel	Madhya Pradesh	D	17655	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0371	Allergy panel-foodVegetarians Panel	Madhya Pradesh	D	8345	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0372	Allergy - Adult (Comprehensive) Panel	Madhya Pradesh	D	19155	Daily : 11:30am	Next day : 5:00pm	5 ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
A0373	Allergy - Animal Panel	Madhya Pradesh	D	4495	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0377	Allergy - Drug Panel: (Maxi)	Madhya Pradesh	D	7000	Tue & Fri : 11:30am	After 8hrs	5 ml of serum
A0379	Allergy - Dust Panel	Madhya Pradesh	D	4280	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0382	Allergy - Environmental Panel	Madhya Pradesh	D	16690	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0383	Allergy - GIT Panel	Madhya Pradesh	D	19420	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0384	Allergy-Grass Panel	Madhya Pradesh	D	3745	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0389	Allergy - Mold-2 Panel	Madhya Pradesh	D	4815	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0395	Allergy - Trees pollen Panel	Madhya Pradesh	D	4065	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0396	Allergy - Weeds pollen Panel	Madhya Pradesh	D	5085	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0398	IgETotal antibodySerum	Madhya Pradesh	B	970	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0399	Allergy Screening IgE With Phadiatop (Adult & >5y)	Madhya Pradesh	D	1750	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0399_C	Allergy Screening IgE With PhadiaTop (Child below 5 y)	Madhya Pradesh	D	1600	Daily : 11:30am	After 8hrs	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A0400	Allergy screeningPhadiatop-AdultFor > 5 years	Madhya Pradesh	C	1465	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0401	Phadiatop allergen screening - Infant (For Below 5 yrs)	Madhya Pradesh	C	1455	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0403	Alpha 1 Acid Glycoprotein Serum	Madhya Pradesh	D	2515	Mon and Thu: 7.30am	Same day 5pm	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0404	Alpha 1 Antitrypsin Serum	Madhya Pradesh	C	2465	Mon and Thu: 7.30am	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0405	alpha 1 anti trypsin (Z, X & M mutations) by PCR**	Madhya Pradesh	R	7140	Daily: 7:30am	1 month	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
A0406	Alpha 2 Macroglobulin Serum	Madhya Pradesh	D	3315	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0407	Alpha Fucosidase Blood	Madhya Pradesh	R	7650	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0408	Alpha Galactosidase Fabry Blood	Madhya Pradesh	R	6120	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0409	Alpha Glucosidase GSD II / pompe Blood	Madhya Pradesh	R	7855	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0410	Alpha Iduronidase Hurler / MPS I Blood	Madhya Pradesh	R	6630	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0411	Alpha Mannosidase Blood	Madhya Pradesh	R	6120	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0412	Alpha N-Acetyl Galactosaminidase Schindler disease Blood	Madhya Pradesh	R	6120	Daily : 11:30am	After 4 weeks : 6:00pm	weekend & Holidays) sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0413	Alpha Thalassemia by PCR	Madhya Pradesh	D	5885	Tue & Sat: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Clinical history, referring clinicians contact number, HPLC Report, CBC findings, clinical/family history and TRF mandatory.
A0414	Aluminium in serum	Madhya Pradesh	C	3025	Daily: 7:30am	Next day 5:00pm	HeavyMetalSamplingguidelines.pdf
A0415	Aluminium, Urine 24 Hrs	Madhya Pradesh	C	3025	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0416	Aluminium spot urine	Madhya Pradesh	C	3025	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0417	AMH Mullerian Inhibiting Substance Serum	Madhya Pradesh	B	2100	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0422	Amino Acid Quantitative Analysis by LC MSMS Plasma	Madhya Pradesh	D	7760	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
A0423	Amino Acid Quantitative Urine	Madhya Pradesh	D	7630	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
A0424	Amino Levulinic Acid (ALA) 24 H Urine	Madhya Pradesh	D	5030	Thu : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
A0425	Amino Levulinic Acid (ALA) Urine Spot	Madhya Pradesh	D	4495	Thu : 11:30am	Next day : 5:00pm	10 ml of Spot urine (Protect From Light)
A0427	Ammonia Biochemical Plasma	Madhya Pradesh	B	655	Daily: 9am to 9pm	After 3 hrs	2 ml of Heparinised Plasma
A0428	Amoebiasis Antibodies-IgG Antibody by IHA - Serum	Madhya Pradesh	C	1510	Mon & Thu : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0429	Amphetamine Quantitative Urine Spot	Madhya Pradesh	B	740	Daily: 9am to 9pm	After 6 hrs	5 ml of Spot Urine
A0430	Amphiphysin Antibody Neuronal Ag in Paraneoplastic Syndrome Serum	Madhya Pradesh	D	5030	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0431	Amylase Enzymatic Ascitic Fluid	Madhya Pradesh	B	545	Daily: 9am to 9pm	After 6 hrs	2 ml of Ascitic Fluid
A0432	Amylase Enzymatic Serum	Madhya Pradesh	B	450	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0433_24	Amylase enzymatic Urine Spot	Madhya Pradesh	B	570	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
A0433_24H	Amylase enzymatic Urine 24H	Madhya Pradesh	B	570	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0434	Amyloid A Serum	Madhya Pradesh	D	4335	Tue : 11:30am	Next day : 5:00pm	coffee , cold drink avoidance for accurate results and such history to be conveyed to lab. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0435	ANA Blot 18 Antigen, Serum	Madhya Pradesh	C	3700	Daily:7.30 am	Same day 8:00 PM	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0436	ANA by ELISA Anti Nuclear Antibody Screening Serum	Madhya Pradesh	B	820	M,W,F: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0437	Anti Nuclear Antibody by IFA	Madhya Pradesh	C	1060	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0438	ANA Profile Serum	Madhya Pradesh	D	5500	Tue,Thur and Sat 7.30am	Same day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0440	Anaemia Profile-Haemolytic Type	Madhya Pradesh	C	3585	Daily: 9am to 9pm	Next day 5pm	4 ml Heparin & EDTA Blood, 4ml Serum
A0441	Anaemia Profile-Maxi	Madhya Pradesh	C	8680	Daily: 9am to 9pm	Next day 5pm	4 ml Heparin & EDTA Blood, 4ml Serum
A0442	Anaemia Profile-Mini	Madhya Pradesh	C	2070	Daily: 9am to 9pm	After 6 hrs	3 ml of EDTA whole blood, 5 ml of serum
A0443	Anaemia Profile-Nutritional	Madhya Pradesh	C	4200	Daily: 9am to 9pm	After 6 hrs	3ml EDTA blood & Serum
A0444	Anaemia Profile-Pernicious	Madhya Pradesh	C	3470	Mon/Thurs - 7:30 am	Next day 5:00 pm	3ml EDTA blood & Serum
A0445	Anti neutrophil cytoplasmic antibody (ANCA) Serum	Madhya Pradesh	C	2220	Daily: 7.30am	Next day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0446	Androgen Profile	Madhya Pradesh	C	5040	Tue and Fri: 7am	Next day 5pm	5 ml of serum
A0447	AndrostenedioneA4Serum	Madhya Pradesh	C	2115	Daily: 9am to 9pm	Next day 5:00 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0448	ACE - Angiotensin Converting Enzyme	Madhya	C	1230	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Serum	Pradesh					ml serum for profiles and reflex tests.
A0449	ACE-Angiotensin Converting Enzyme, CSF	Madhya Pradesh	D	3425	Daily : 11:30am	4th day : 6:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
A0450	Antenatal Profile-II	Madhya Pradesh	C	1175	Daily: 9am to 9pm	After 8 hrs	2 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma , 20 ml of spot urine
A0451	Antenatal Profile - I	Madhya Pradesh	C	3000	Daily: 8pm	3rd day	2 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma , 20 ml of spot urine
A0453	Anti - A Titre	Madhya Pradesh	D	1605	Daily : 11:30am	After 8hrs	3 ml Serum
A0454	Anti - B Titre	Madhya Pradesh	D	1605	Daily : 11:30am	After 8hrs	3 ml Serum
A0455	Anti DNase B enzymatic Seum	Madhya Pradesh	D	2225	Daily: 7.30am	3rd day	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0458	Anti Thrombin III Activity Functional, Citrated plasma	Madhya Pradesh	D	4065	Mon and Thurs: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
A0459	Anti Thrombin III Antigen, Citrated plasma	Madhya Pradesh	D	5135	Mon and Thurs: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
A0462	Antimony ICPMS Urine Spot	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
A0463_N	Phospholipid Profile - 3	Madhya Pradesh	D	8455	Thursday 9.0 am	Same day 7.0 pm	3 ml of Citrate plasma, 3 ml of serum
A0464	APA-IgG Phospholipid Antibody Serum	Madhya Pradesh	C	1120	Daily: 7:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0465	APA-IgM Phospholipid Antibody Serum	Madhya Pradesh	C	1120	Tue,Thu, S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0466	APCR Activated Protein C Resistance Test, Citrated plasma	Madhya Pradesh	D	4495	Mon and Thurs: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
A0467	Apolipoproteins A1serum	Madhya Pradesh	B	680	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0468	Apolipoproteins B/A1	Madhya Pradesh	C	1230	Daily: 9am to 9pm	After 4 hour	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0469	Apolipoproteins Bserum	Madhya Pradesh	B	680	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0470	Apolipoproteins Profile, Serum	Madhya Pradesh	C	2015	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0471	APTT Activated Partial Thromboplastin Time, Citrated plasma.	Madhya Pradesh	B	500	Daily: 9am to 9pm	After 6 hrs	1 ml of citrated platelet poor plasma in frozen condition.
A0472	Arginine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
A0473	Arginine Quantitative Urine	Madhya Pradesh	R	6325	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
A0474	Arsenic	Madhya Pradesh	D	3210	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0475	Arsenic, Urine 24 Hrs	Madhya Pradesh	D	3210	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0476	Arsenic, Urine Spot	Madhya Pradesh	D	3210	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0477	Arthritis Profile	Madhya Pradesh	C	2855	Daily: 9am	Next day 5pm	2 ml of EDTA whole blood, 3 ml of serum, 20 ml of spot urine
A0478	Aryl Sulphatase A metachromatic leucodystrophy Blood	Madhya Pradesh	R	6120	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0479	Aryl Sulphatase B MPS VI Blood	Madhya Pradesh	R	9080	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0480	ASCA-IgA Saccharomyces Cerevisiae Serum	Madhya Pradesh	C	2350	Tue,Fri: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0481	ASCA-IgG Saccharomyces Cerevisiae Serum	Madhya Pradesh	C	2350	Tue,Fri: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0482	ASH /NASH / STEATO Profile	Madhya Pradesh	E	16000	Daily: 9am	5th day	tests. 5 ml of serum and 2 ml of fluoride plasma . Fasting sample
A0483	ASKA Skeletal (Striated) Muscle Antibody Serum	Madhya Pradesh	D	3640	Mon, Wed & Fri : 11:30am	Same day: 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0484	ASMA Smooth Muscle Antibody Serum	Madhya Pradesh	C	2130	Daily: 11:30am	Next day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0485	ASO titre Quantitative Serum	Madhya Pradesh	B	545	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0486	Asparagine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
A0487	Asparagine Quantitative Urine	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
A0488	Aspartic Acid Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
A0489	Aspartic Acid Quantitative Urine	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
A0490	Aspergillosis IgG antibodies by EIA - Serum	Madhya Pradesh	C	2350	Mon & Thu : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0491	Aspergillosis IgM antibodies by EIA - Serum	Madhya Pradesh	C	2350	Mon & Thu : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0492	Galactomannan, Aspergillus Antigen test by EIA -Serum	Madhya Pradesh	C	6160	Wed & Sat : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0493	Autoimmune Liver Disease Profile -1	Madhya Pradesh	D	5000	Mon/Wed - 7:30 am	Next day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0495	Autoimmune Liver Disease Profile-3	Madhya Pradesh	D	4065	Mon/Wed - 7:30 am	Next day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0496	Autoimmune Liver Disease Profile-4	Madhya Pradesh	D	5780	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0497	Autoimmune Encephalitis Panel	Madhya Pradesh	E	23000	Mon, Wed, Fri : 7:30am	Next Day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0499	Acid lipase disorder**	Madhya Pradesh	R	7140	Daily : 11:30am	4 weeks	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
A0500	ACTN3 GENOTYPING_A0500	Madhya Pradesh	R	4080	Thu : 11:30am	10th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
A0509	AFB-Xpert panel by Genexpert-Combo - Extrapulmonary Samples	Madhya Pradesh	D	3210	Daily : 7:30am	After 6 weeks : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
A0510	AFB-Xpert panel by Genexpert-Combo - Pulmonary Samples	Madhya Pradesh	D	3210	Daily : 7:30am	After 6 weeks : 6:00pm	2ml of Pulmonary specimen in Sterile leak proof container
A0512	ANCA by IFA reflex MPO/PR3 Antibodies Serum	Madhya Pradesh	C	2800	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0513	Apolipoproteins E**	Madhya Pradesh	R	4180	Wed : 11:30am	4th day : 6:00pm	3 ml of serum
A0514	TruHeart Platinum	Madhya Pradesh	W	5800	Mon and Thu: 9am	Same day 5pm	5 ml of serum
A0519	Allergen,Individual-Tree pollen Papdi Chibil / Elm_A0519	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0520	Allergen,Individual-Weed pollen SALTWORT (Salsola kali / pestifer)_A0520	Madhya Pradesh	C	1570	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0522	ALK (D5F3) reflex FISH	Madhya Pradesh	D	6420	Daily : 9:00am	7th day : 6:00pm	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
A0525	Allergen,Individual-Animal Pigeon Feathers, Serum_A0525	Madhya Pradesh	C	1455	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0533	ADA - Adenosine Deaminase by Biochemical Method - Synovial fluid	Madhya Pradesh	B	870	Daily : 9:00am to 9:00pm	After 6 hrs	3 ml of Fluid in Sterile leak proof container
A0536	APA-IgA Phospholipid Antibody Serum	Madhya Pradesh	C	1120	Tue,Thu, S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0537	Atherosclerosis Profile - 1	Madhya Pradesh	D	5295	Mon and Thu: 9am	Same day 5pm	5 ml of serum
A0538	Atherosclerosis Profile - 3	Madhya Pradesh	D	7435	Mon and Thu: 9am	Same day 5pm	5 ml of serum
A0539	Aldosterone/PRA Ratio	Madhya Pradesh	D	7300	Tue & Fri: 11:30am	Next day : 5:00pm	2 ml of EDTA Plasma and 2 ml Serum
A0560	Galactomannan, Aspergillus Antigen test by EIA - BAL Fluid	Madhya Pradesh	C	6160	Wed & Sat : 11:30am	Next day : 5:00pm	3 ml of Fluid in Sterile leak proof container
A0563	AFB (TB) - Antibiogram for MTb complex - Cycloserine Drug by Conventional LJ - Pure Growth of MTB complex	Madhya Pradesh	C	2690	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0564	Albumin body fluid	Madhya Pradesh	A	220	Daily: 9am to 9pm	After 6 hrs	2 ml body fluid
A0567	ALK/NPM1: t(2;5) translocation Anaplastic Large cell lymphoma FISH, Bone marrow	Madhya Pradesh	D	4920	Daily: 7:30 am	3rd day	2-3 ml Bone Marrow in Na-heparin vacutainer
A0568	ALK/NPM1: t(2;5) translocation Anaplastic Large cell lymphoma FISH, BLOOD	Madhya Pradesh	D	4920	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
A0569	ALK/NPM1: t(2;5) translocation Anaplastic Large cell lymphoma FISH, TISSUE	Madhya Pradesh	D	8560	Daily 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0570	ECP (Eosinophil cationic protein) Serum	Madhya Pradesh	D	1605	Mon,Thu: 9am	Same day: 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0580	Allergy-GIT Pane I(Veg + Non Veg)	Madhya Pradesh	D	26750	Daily : 11:30am	Next day : 5:00pm	5 ml of serum
A0584	Allergy screening-Tryptase	Madhya Pradesh	C	2240	Mon,Thu: 9am	Same day: 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0585	Anaemia Profile - Macrocytic	Madhya Pradesh	C	2130	Daily: 9am to 9pm	After 8 hrs	3ml EDTA blood & Serum
A0586	Anaemia Profile - Microcytic Hypochromic	Madhya Pradesh	C	2130	Daily: 9am to 9pm	After 8 hrs	3ml EDTA blood & Serum
A0587	Anaemia Workup Profile	Madhya Pradesh	C	2350	Daily: 9am to 9pm	Next day 5pm	3ml EDTA blood & Serum
A0588	Aneuploidy-Sperm 13,18,21,X and Y by FISH	Madhya Pradesh	D	5135	Daily : 9:00am	7th day : 6:00pm	CytogeneticsandFISH.pdf
A0595	Adrenoleukodystrophy (X-linked ALD) sequencing of ABCD1 gene**	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
A0596	Autoimmune encephalitis panel-CSF	Madhya Pradesh	E	19900	Mon, Wed, Fri : 7:30am	Next Day : 7:00pm	2 ml of CSF in sterile screw capped container.
A0599	APC gene analysis**	Madhya Pradesh	R	17850	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0608	AFB (TB) - Antibiogram for MTb complex - Linezolid Drug by MGIT - Pure Growth of MTB complex	Madhya Pradesh	C	2910	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0613	Alpha thalassemia gene analysis (HBA1 & HBA2)	Madhya Pradesh	E	16000	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Clinical history, referring clinicians contact number, HPLC Report, CBC findings, clinical/family history and TRF mandatory.
A0616	Aldosterone/Direct Renin Ratio	Madhya Pradesh	D	6420	Daily: 7.30am	Same day 5pm	2 ml of EDTA Plasma and 2 ml Serum
A0617_CSF	Autoimmune encephalitis panel-2	Madhya Pradesh	E	31000	Tue: 9am	Next Day 7pm	2 ml of CSF in sterile screw capped container.
A0617_serum	Autoimmune encephalitis panel-2	Madhya Pradesh	E	31000	Tue: 9am	Next Day 7pm	3 ml of serum
A0618	APO E GENOTYPING**	Madhya Pradesh	R	25500	Daily : 11:30am	12th Day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
A0623	Alzheimer disease screening profile CSF	Madhya Pradesh	D	13375	Tue : 11:30am	Next day : 5:00pm	2 ml of CSF in leakproof container. (folin vu conical tube is recommended to prevent leakage).
A0624	Arsenic plasma	Madhya Pradesh	D	3210	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
A0625	Serum Ascites Albumin Gradient (SAAG)	Madhya Pradesh	D	535	Daily: 9am to 9pm	After 6 hrs	Serum 1 ml, Fluid 1 ml
A0626	Allergy Screening	Madhya Pradesh	D	5990	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A0629	Aspergillus broncho pulmonary	Madhya	D	2900	Mon & Thu :	Next day : 5:00pm	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	aspergillosis (ABPA)	Pradesh			11:30am		ml serum for profiles and reflex tests.
A0630	AFB Drug susceptibility Flexi panel (4 drugs) by MGIT - Pure Culture	Madhya Pradesh	D	6420	Daily : 7:30am	Upto 4 weeks	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A0632	Asthma / Rhinitis Comprehensive	Madhya Pradesh	D	12305	Daily : 11:30am	Next day : 5:00pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A0633	Allergy panel - Asthma / Rhinitis comprehensive with monitoring 4	Madhya Pradesh	D	18245	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A0635	Allergy - Adult panel - 2	Madhya Pradesh	D	5245	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A0636	Allergy - Adult panel - 3	Madhya Pradesh	D	8240	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A0638	Aluminium, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
A0639	Aluminium, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A0640	Amphetamines confirmation by LCMS, Spot Urine**	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
A0641	Antimony, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
A0642	Arsenic, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
A0643	Arsenic, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
A0644	Acute intermittent porphyria (HMBS)-Gene Analysis	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
A0645	Alkaptonuria (HGD)-Gene Analysis	Madhya Pradesh	E	19000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							history is mandatory.
A0647	TRUNAT MTB/RIF Resistance by TAQMAN	Madhya Pradesh	C	2295	Daily	Same Day (12 hr TAT)	Mycobacteriology.pdf
A0648	AFB-DNA (TB-PCR) Detection by Real Time PCR , Reflex to Rifampicin resistance by Ultra CBNAAT (If - MTB Positive) - Sputum / BAL	Madhya Pradesh	C	2295	Daily : 7:30am	Next day : 5pm	Mycobacteriology.pdf
A0649	Mycobacterial culture identification by MALDI TOF- Pure Culture	Madhya Pradesh	C	2465	Tue : 7:30am	Next day : 5pm	Pure growth of Mycobacterium spp. on LJ slant/ MGIT Tube/ Bactec alert bottle.
A0650	SARS-CoV-2 Antibody, qualitative by CLIA - Serum	Madhya Pradesh	C	785	Daily : 9:00am to 9:00pm	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0653	ALPS Screen Panel by Flow cytometry**	Madhya Pradesh	R	9280	Daily : 11:30am	5th day : 6:00pm	sample : 3 ml of EDTA Blood
A0654	SARS-CoV-2 (Spike) Antibody, Quantitative (Coviprotect) by CLIA - Serum	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	Report after 24 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
A0656	Aspergillus DNA Detection by Real Time PCR - Body fluids/CSF/EDTA Blood/Pure culture/Sputum/Tissue	Madhya Pradesh	D	8345	Daily : 7:30am	3rd day : 6:00pm	3ml of Body fluids/ CSF/ EDTA Blood/ Pure culture/ Sputum/ Tissue in Sterile leak proof container
A0657	Array Feto	Madhya Pradesh	E	16000	Daily : 9:00am	15th day : 6:00pm	1)Upto 12 weeks of gestation: 50-100mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container 2) In 12-18 weeks of gestation: a) 50-100mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container b) and or fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container 3)IUFD >20 weeks of gestation: fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container. if possible in 20ml normal saline with 2-3 drops of gentamycin and 1-4ml Maternal

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							blood in EDTA tube is required for maternal cell contamination (MCC) studies.
A0658	Amniochrome	Madhya Pradesh	E	19000	Daily : 9:00am	15th day : 6:00pm	Amniotic fluid (20-30ml divided in 2 sterile falcon centrifuge tube) and 1-4ml Maternal blood in EDTA tube is required for maternal cell contamination (MCC) studies. (Completely filled Consent form G with patient and doctor stamp and signature is mandatory as per PCPNDT Act)
A0659	Array Chrome	Madhya Pradesh	E	16000	Daily : 9:00am	15th day : 6:00pm	3ml Peripheral Blood in EDTA Vacutainers
A0660	Anti complement Factor H	Madhya Pradesh	R	5100	Daily : 11:30am	6th day : 6:00pm	3 ml Serum
A0661	ADAMTS 13 Activity	Madhya Pradesh	R	7900	Daily : 11:30am	6th day : 6:00pm	3 ml Citrate Plasma (Frozen mandatory)
A0662	ADAMTS 13 Autoantibody	Madhya Pradesh	E	8000	Daily : 11:30am	5th day : 6:00pm	3 ml Citrate Plasma (Frozen mandatory)
A0680	ALLERGY COMPONENT TESTING (CRD)	Madhya Pradesh	E	20000	Daily:- 9.00 am	Next day 5 pm	Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and -20Â°C for long time.
A1009	Adrenal Gland Antibodies, Serum**	Madhya Pradesh	R	6120	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
A1901	Albumin, serum	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
A1902	Canine - Alkaline Phosphatase	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A1903	Canine - Amylase level	Madhya Pradesh	A	520	Daily: 9am to 9pm	Report after 24 hrs	for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
A5583	AMH Premium	Madhya Pradesh	D	3150	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A5592	AMH Supreme	Madhya Pradesh	D	4600	Daily: 9am to 9pm	Next day at 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A5702	Comprehensive Allergy Panel	Madhya Pradesh	D	15000	Daily: 9am	Same day 8pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A5703	Comprehensive Eczema panel	Madhya Pradesh	D	11770	Daily: 9am	Same day 8pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
A5704	Common Food panel	Madhya Pradesh	D	6420	Daily: 9am	Same day 8pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A8121	Amylase Enzymatic,Pleural Fluid	Madhya Pradesh	B	570	Daily: 9am to 9pm	After 6 hrs	temperature for a day. 2 ml body fluid
A8285	AFB SUSCEPTIBILITY BEDAQILINE	Madhya Pradesh	C	2690	Daily : 9:00am	After 4 weeks : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8455	AFB DRUG SUSCEPTIBILITY PANEL 14 DRUGS	Madhya Pradesh	D	14980	Daily: 7:30am	Upto 4 weeks 5 PM	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8456	AFB DRUG SUSCEPTIBILITY FIRST LINE PANEL (HREZ), 4 DRUG BY MGIT	Madhya Pradesh	D	6205	Daily: 7:30am	Upto 4 weeks 5 PM	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8457	AFB DRUG SUSCEPTIBILITY PANEL 11 DRUGS	Madhya Pradesh	D	11180	Daily: 7:30am	Upto 4 weeks 5 PM	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8458	AFB DRUG SUSCEPTIBILITY ORAL BEDAQUILINE PANEL, 8 DRUGS BY MGIT	Madhya Pradesh	D	8560	Daily: 7:30am	Upto 4 weeks 5 PM	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8459	AFB DRUG SUSCEPTIBILITY ISONIAZID MONORESITANCE PANEL (LFxREZ), 4 DRUGS BY MGIT	Madhya Pradesh	D	6205	Daily: 7:30am	Upto 4 weeks 5 PM	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
A8598	Allergen,Individual-Drug CEPHALOSPORIN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8599	Allergen,Individual-Drug LIDOCAINE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8600	Allergen,Individual-Drug BENZOCAINE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8601	Allergen,Individual-Drug MEPIVACAINE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8602	Allergen,Individual-Drug PRILOCAINE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8603	Allergen,Individual-Drug PARACETAMOL	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8604	Allergen,Individual-Drug PYRAZOLONE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A8605	Allergen,Individual-Drug PIROXYCAM	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8606	Allergen,Individual-Drug DICLOFENAC	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8607	Allergen,Individual-Drug IBUPROFEN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8608	Allergen,Individual-Drug PENICILLOYL G	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8609	Allergen,Individual-Drug PENICILLOYL V	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8610	Allergen,Individual-Drug DOXYCYCLINE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8611	Allergen,Individual-Drug CLOXACILLIN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8612	Allergen,Individual-Drug NEOMYCIN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8613	Allergen,Individual-Drug OXACILLIN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8614	Allergen,Individual-Drug CHLORAMPHENICOL	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8615	Allergen,Individual-Drug METRONIDAZOLE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8616	Allergen,Individual-Drug AMPICILLIN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8617	Allergen,Individual-Drug AMOXICILLIN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A8618	Allergen,Individual-Drug TETRACYCLINE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8619	Allergen,Individual-Drug GENTAMYCIN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8620	Allergen,Individual-Drug ERYTHROMYCIN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8621	Allergen,Individual-Drug STREPTOMYCIN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8622	Allergen,Individual-DrugTRIMETHOPRIM	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8623	Allergen,Individual-Drug SULFAMETHOXAZOLE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8624	Allergen,Individual-Drug INSULIN HUMAN	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8625	Allergen,Individual-Drug CORTISONE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8626	Allergen,Individual-Drug PREDNISOLONE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8627	Allergen,Individual-Drug THEOPHYLLINE	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8628	Allergen,Individual-Drug AMBROXOL	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
A8629	Allergen,Individual-Drug PYRIDOXAMINE (VITAMIN B6)	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
A8630	Allergen,Individual-DrugCOBALAMIN (VITAMIN B12)	Madhya Pradesh	D	1365	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0007	Bacterial Meningitis Panel 5 antigen (Ecoli K1, H.influenza, Streptococcus pneumoniae, Neisseria meningitidis, Streptococcus agalactiae) by Latex agglutination - CSF	Madhya Pradesh	D	4710	Daily : 9:00am to 3:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
B0008	Bacterial Meningitis Panel 5 Full Panel by Latex agglutination - Serum	Madhya Pradesh	D	5990	Daily : 9:00am to 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
B0010	Barbiturates Quantitative Urine Spot	Madhya Pradesh	B	740	Daily: 9am to 9pm	After 6 hrs	5 ml of Spot Urine
B0013	Bence Jones Proteins, Urine-Qualitative	Madhya Pradesh	B	1070	Daily: 7.30am	Same day 5pm	10 ml of 24 hrs urine (Preferred) / Spot urine
B0015	Benzodiazapine Diazepam Urine Spot	Madhya Pradesh	B	740	Daily: 9am to 9pm	After 6 hrs	5 ml of Spot Urine
B0017	Beta Galactocerebrosidase Krabbe Blood	Madhya Pradesh	R	8365	Daily : 11:30am	After 4 weeks : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0018	Beta Galactosidase GM 1 / MPS IV B Blood	Madhya Pradesh	R	5710	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0019	Beta Glucosidase Gaucher Blood	Madhya Pradesh	R	5710	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0020	Beta Hexosaminidase A GM2-1/ taysachs Blood	Madhya Pradesh	R	7650	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0021	Beta Hexosaminidase T GM2-2/ Sandhoff Blood	Madhya Pradesh	R	7650	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
B0022	Beta Mannosidase Blood	Madhya Pradesh	R	6120	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B0023	Beta Thallesemia Prenatal Diagnosis	Madhya Pradesh	D	17655	Daily : 9:00am	7th day : 6:00pm	of collection, Avoid collection on weekend & Holidays) Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition or Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
B0025	Beta-Thalassemia, Beta Globin (HBB) Gene Sequencing	Madhya Pradesh	D	9415	Wed & Fri: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Clinical history, referring clinicians contact number, HPLC Report, CBC findings, clinical/family history and TRF mandatory.
B0026	Beta-2-Crosslaps CTx-1 Plasma	Madhya Pradesh	C	2015	Mon and Thu: 7.30am	Same day 5pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B0027	Beta-2-Glycoprotein 1- IgA serum	Madhya Pradesh	C	1455	Tue,Fri: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0028	Beta-2-Glycoprotein 1 - IgG Antibody Serum	Madhya Pradesh	C	1230	Tuesday, Thursday and Saturday 7:30 am	Same day 7pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0029	Beta-2-Glycoprotein 1- IgM Antibody Serum	Madhya Pradesh	C	1230	Tuesday, Thursday and Saturday 7:30 am	Same day 7pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0030	Beta-2-Microglobulin CSF	Madhya Pradesh	C	1905	Daily: 9am to 9pm	After 8 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
B0031	Beta-2-Microglobulin Serum	Madhya Pradesh	C	1680	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0032	Beta-2-Microglobulin Urine Spot	Madhya Pradesh	C	1905	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
B0033	Bicarbonate biochemical Serum	Madhya Pradesh	B	400	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0034	Bile Acids-Total biochemical Serum	Madhya Pradesh	C	1905	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0035	Bile Salt & Pigments - Urine	Madhya Pradesh	A	100	Daily : 9:00am to 9:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
B0036	Bilirubin Direct Serum	Madhya Pradesh	A	165	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0037	Bilirubin Total Serum	Madhya Pradesh	A	165	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0038	Bilirubin Total, Direct, Indirect Serum	Madhya Pradesh	A	200	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0039	Biotinidase Neonatal Dried blood spot	Madhya Pradesh	C	775	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. DBS card to be send in Ambient

Test Code	Test Name	Location	Category	Price(Rs.)	Test Schedule	Reported on	Sample Guideline
							<p>temperature only . Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age, or as close to discharge as possible if discharged prior to 24 hours of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Plz provide blood transfusion history and sample to be collected post 72 hrs after transfusion. Spot Urine: 10 ml of Spot urine required with clinical history and clinician number. Amino acid Plasma: 3 ml of separated EDTA plasma required . Sample to be shipped at Frozen condition . Sample stable at refrigerated condition 2 to 8 degree for 1 day. Clinical history, metabolic disorders or inborn errors of metabolism suspicion or provisional diagnosis of clinician required fasting recommended not mandatory. Amino acid supplementation or medications needs to be stopped based on physician recommendations</p>

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B0040	Biotinidase Quantitative Plasma **	Madhya Pradesh	R	5305	Daily : 11:30am	10th day : 6:00pm	4 ml of Heparin plasma
B0044	BK Polyoma DNA Detection by Real Time PCR - Urine	Madhya Pradesh	C	5275	Mon, Thursday: 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
B0045	Bleeding & Clotting Time (BT,CT)	Madhya Pradesh	B	150	Daily: 9am to 9pm	After 6 hrs	NA
B0046	Blood Group ABO & Rh Typing	Madhya Pradesh	B	120	Daily: 9am and 5pm	After 6 hrs	3 ml of EDTA whole blood and Serum sample required
B0047	Bombay Blood Group Evaluation	Madhya Pradesh	D	1070	Daily : 11:30am	3rd day : 6:00pm	3 ml EDTA whole blood
B0053	BOH (Bad Obstetric History) Profile-Maxi	Madhya Pradesh	C	16240	Daily : 9:00am	15th day : 6:00pm	3 ml EDTA whole blood, 5 ml of serum, 4 ml of heparin whole blood (husband and wife), 3 ml of citrate plasma
B0054	BOH (Bad Obstetric History) Profile-Mini	Madhya Pradesh	C	6270	Daily: 7.30am	Next day 7pm	3 ml EDTA whole blood, 5 ml of serum, 3 ml of citrate plasma
B0055	Bone Marrow Aspiration for morphological examination	Madhya Pradesh	C	1455	Daily : 11:30am	Next day : 5:00pm	2 Stained,2 Unstained Fixed BM Smears and 2 Peripheral Smears Required.
B0056	Bone Profile	Madhya Pradesh	C	1230	Daily 9:00 am	After 7 days	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0057	Bordetella Pertussis IgG antibodies Quantitative by EIA -Serum	Madhya Pradesh	C	2800	Fri : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B0058	BRAF V600 Mutation Analysis, Tissue	Madhya Pradesh	D	6420	Mon, Wed & Fri: 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
B0059	BRAF V600E Mutation Analysis, Hairy Cell Leukemia	Madhya Pradesh	D	5885	Mon & Thu: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							referring clinicians contact number and TRF mandatory.
B0062	Breast Marker Profile	Madhya Pradesh	D	2055	Daily: 8pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0063	Brucella Total antibodies by Tube Agglutination - Serum	Madhya Pradesh	B	1090	Daily : 7:30pm	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B0064	Brucella IgG antibody by EIA - Serum	Madhya Pradesh	C	1455	Wed & Sat : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B0065	Brucella IgM antibody by EIA - Serum	Madhya Pradesh	C	1455	Wed & Sat : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B0066	Bullous Pemphigoid 180 Vesiculobullous disorders Serum	Madhya Pradesh	D	4710	Fri : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0067	Bullous Pemphigoid 230 Vesiculobullous disorders Serum	Madhya Pradesh	D	4710	Fri : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0068	BUN Urea Nitrogen Serum	Madhya Pradesh	A	175	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B0069	BUN Urea Nitrogen Urine 24H	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
B0070	BUN Urea NitrogenPericardial Fluid	Madhya	A	275	Daily: 9am to 9pm	After 6 hrs	2 ml of Pericardial Fluid

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B0071	BUN Urea Nitrogen Urine Spot	Pradesh Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
B0074	Bladder Cancer/Urothelial carcinoma Marker Profile (Aneuploidy for chromosomes 3, 7, 17, and loss of the 9p21 **)	Madhya Pradesh	R	14280	Daily : 11:30am	10th day : 6:00pm	100 ml Early morning Urine
B0075	BRCA deletion/duplication analysis - EDTA Whole blood	Madhya Pradesh	E	16540	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
B0076	TARGETED MUTATION ANALYSIS (FAMILY MEMBERS)	Madhya Pradesh	E	5000	Daily : 9:00am	7th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
B0077	NextGen BRCA Comprehensive (NGS & MLPA), Blood	Madhya Pradesh	E	30000	Daily: 9 am	4 weeks	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Duly signed BRCA Proforma, BRCA1-2 Consent Form
B0078	Urea Nitrogen Body Fluid	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	2 ml body fluid
B0080	16s rRNA Sequencing Bacterial Identification*	Madhya Pradesh	D	3210	Daily : 9:00am	7th day : 6:00pm	Isolated colony on culture plates. Sample stable for 2 Days at 2 to 8 degree celsius.
B0081	Beta Thalassemia (Prenatal + Parental	Madhya	D	19795	Daily : 9:00am	7th day : 6:00pm	Amniotic Fluid (15-20 ml) in a

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Testing) - 30 Mutaions	Pradesh					sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition or Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
B0083	BCR-ABL Breakpoint Analysis Qualitative by Real Time PCR - Blood	Madhya Pradesh	D	4815	Tue & Fri : 9:00am	3rd Day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer
B0084	BCR-ABL Breakpoint Analysis Qualitative by Real Time PCR - Bone Marrow	Madhya Pradesh	D	4815	Tue & Fri : 9:00am	3rd Day : 6:00pm	3 ml of Bone Marrow in EDTA Vacutainer
B0090	Beta-Trace protein(BTP) CSF leakage in body fluid	Madhya Pradesh	D	5885	Daily : 11:30am	Next day : 5:00pm	2 ml Nasal Fluid And 3 ml Serum.
B0101	BCL6 Gene rearrangement, FISH, Tissue	Madhya Pradesh	E	8000	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B0104	(1-3)-?-D-Glucan Level by Protease Zymogen Colorimetry by Protease Zymogen Colorimetry - Serum	Madhya Pradesh	D	7715	Tue, Thurs & Sat : 9:00am	Same day : 7pm	Collect sample in gel tube directly and send after centrifugation. Kindly do not open tube or separate serum from gel tube to another tube to avoid contamination. Aliquot sample will not accepted. Kindly send minimum 2 ml serum.
B0106	NextGen BRCA Somatic, FFPE tissue	Madhya Pradesh	E	22000		25 Days	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
B0108	Pneumonia panel by FilmArray (Biofire)- Lower respiratory tract - Sputum/BAL/ Endotracheal secretions	Madhya Pradesh	E	26500	Daily : 9:00am to 9:00pm	After 10 hrs	2 ml of BAL/ Sputum/ Endotracheal secretions in Sterile leak proof container
B0109	Barbiturates confirmation by LCMS, Spot Urine**	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
B0111	Barium, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
B0112	Barium, Spot Urine**	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
B0113	Barium, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
B0114	Beryllium, EDTA Whole Blood**	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
B0115	Benzodiazapine confirmation by LCMS,	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Spot Urine**	Pradesh					to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
B0116	Beryllium, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
B0117	Beryllium, Spot Urine**	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
B0119	Biotinidase deficiency (BTD)-Gene Analysis	Madhya Pradesh	E	19000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
B0120	Bone Marrow-2 Comprehensive	Madhya Pradesh	D	6955		5th Day	Bone marrow biopsy (in 10% Neutral buffered formalin) & aspirate (in EDTA) examination with special stain & IHC. Send PBS as well.
B0121	BK Polyoma DNA Detection by Real Time PCR - EDTA Blood	Madhya Pradesh	C	5275	Mon & Thurs : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
B0122	Bismuth ICPMS , Blood	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	3 ml EDTA Blood / Metal free EDTA Blood
B0125	Bone Marrow Aspirate Iron stain	Madhya Pradesh	C	505	Daily : 11:30am	Next day : 5:00pm	1 Stained, 2 Unstained air dried Bone marrow smears with good particles and 1 Peripheral smear required._x000D_ _x000D_
B0126	Biofire Blood Panel (BCID)	Madhya	D	17975	Daily : 9:00am to	After 10 hrs	Positive blood culture bottle

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B1008	Beta Carotene**	Pradesh Madhya Pradesh	R	2960	9:00pm Mon : 11:30am	5th day : 6:00pm	3 ml of serum
B1017	Borreliosis-IgG Confirm serology by western blot - Serum **	Madhya Pradesh	R	13055	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B1029	Bacterial Meningitis Panel 5 antigen (Ecoli K1, H.influenza, Streptococcus pneumoniae, Neisseria meningitidis, Streptococcus agalactiae) by Latex agglutination - Serum	Madhya Pradesh	D	4710	Daily : 9:00am to 3:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
B1030	Bone formation markers	Madhya Pradesh	D	6690	Mon & Thur 9 am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
B1031	Bone Resorption Markers	Madhya Pradesh	D	4550	Daily: 9 am	After 5 days	3 ml of serum and 20 ml of spot urine
B1032	NextGen BRCA Germline, Blood	Madhya Pradesh	E	16000	Daily : 9:00am	15th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Duly signed BRCA Proforma, BRCA1-2 Consent Form
B1901	Canine - BCP (Biochemistry Panel)	Madhya Pradesh	D	3315	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1902	Bicarbonate, serum by PEPC method	Madhya Pradesh	B	655	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1903	CANINE-Bile Acids-Total	Madhya Pradesh	C	1905	Wed: 9am	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
B1904	Bilirubin Direct, Serum	Madhya Pradesh	A	195	Daily: 9am to 9pm	Report after 24 hrs	for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1905	Bilirubin Total, Serum	Madhya Pradesh	A	195	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1906	Canine - Bilirubin Total, Direct, Indirect, Serum	Madhya Pradesh	D	300	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1907	Bone Marrow Aspiration Slides for review	Madhya Pradesh	C	1455	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B1908	Canine - BUN-Blood Urea Nitrogen	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
B8001	Bacteriological Examination of Air by Vitek/ MALDITOF - Blood Agar	Madhya Pradesh	B	545	Daily : 9:00am to 9:00pm	3rd day : 5pm	Blood agar plate , Kindly mention time of plate exposure.
B8046	BK Polyoma DNA Viral load by Real Time PCR - Urine	Madhya Pradesh	C	6160	Mon, Thursday: 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
B8047	BK Polyoma DNA Viral load by Real Time PCR - EDTA Blood	Madhya Pradesh	C	6160	Mon, Thursday: 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
B8596	Bone & Joint Infection panel (JI) by FilmArray	Madhya Pradesh	E	24000	Daily : 9:00am to 9:00pm	After 10hrs	3ml of Synovial fluid in sterila container
C0001	C1 Esterase Inhibitor Serum	Madhya Pradesh	D	3530	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0003	C3 - Complement-3 Serum	Madhya Pradesh	C	900	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0004	C4 - Complement-4 Serum	Madhya Pradesh	C	785	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0005	CA-125 Serum	Madhya Pradesh	B	1290	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0006	CA-15.3 Serum	Madhya Pradesh	B	1415	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0007	CA-19.9 Serum	Madhya Pradesh	B	1400	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0009	CA-72.4 Serum	Madhya Pradesh	C	2015	Tue & Fri : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0010	Cadmium Blood	Madhya Pradesh	D	3640	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0011	Cadmium Urine 24H**	Madhya Pradesh	D	3415	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0013	CAH Panel - Congenital Adrenal Hyperplasia	Madhya Pradesh	C	8625	Daily: 9am	Next day 5pm	5ml of serum
C0014	Calcitonin Thyrocalcitonin Serum	Madhya Pradesh	C	2350	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0015	Calcium Urine 24H	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0016	Calcium Ionised Serum	Madhya Pradesh	C	650	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0017	Calcium Total Serum	Madhya Pradesh	A	160	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0018	Calcium Profile Serum	Madhya Pradesh	C	615	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0019_24	Calcium/Creatinine Ratio (24 hrs)	Madhya Pradesh	A	435	Daily: 9am to 9pm	Calcium/Creatinine Ratio (24 hrs)	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
C0019_SP	Calcium/Creatinine Ratio (Spot sample)	Madhya Pradesh	A	435	Daily: 9am to 9pm	Calcium/Creatinine Ratio (Spot sample)	10 ml of Spot urine
C0020	Campylobacter antigen detection by Immunochromatography - Stool	Madhya Pradesh	C	1625	Daily : 9:00am to 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
C0021	Cancer Detection Profile, Female	Madhya Pradesh	D	9900	Daily: 9am to 9pm	5th day	3 ml of EDTA whole blood, 5 ml of serum, 1 thin prep container
C0022	Cancer Detection Profile, Male	Madhya Pradesh	D	8025	Daily: 9am to 9pm	After 8 hrs	3 ml of EDTA whole blood, 5 ml of serum
C0028	Candida Albicans Ab by Indirect Haemagglutination - Serum**	Madhya Pradesh	R	11220	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0029	Cannabis Marijuana, THC Urine Spot	Madhya Pradesh	B	765	Daily: 9am to 9pm	After 6 hrs	5 ml of Spot Urine
C0030	Capecitabine (induced toxicity) DPD,TS,MTHFR Blood**	Madhya Pradesh	R	7140	Daily: 7.30am	5th day	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
C0031	Calprotectin Stool	Madhya Pradesh	C	3100	Daily: 7.30am	Same day 5pm	100 gms of stool in leak proof container. History of medicines, morning sample preferred.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0032	Carbamazepine Tegretol Serum	Madhya Pradesh	C	950	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0033	Cardiac Injury Profile-Maxi	Madhya Pradesh	D	5030	Daily: 9am to 9pm	After 6 hrs	5 ml of serum
C0034	Cardiac Injury Profile-Mini	Madhya Pradesh	D	1285	Daily: 9am to 9pm	After 6 hrs	5 ml of serum
C0035	Cardiac Risk Profile	Madhya Pradesh	D	4600	Daily: 9am to 9pm	After 6 hrs	3 ml of serum and 3 ml of citrate plasma
C0036	Cardiac Screen	Madhya Pradesh	C	1570	Daily: 9am to 9pm	After 6 hrs	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma, 10 ml of spot urine
C0037	Cardiolipin Antibody ACL-IgA antibody Serum	Madhya Pradesh	C	1065	Tue & Fri : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0038	Cardiolipin Antibody ACL-IgG antibody Serum	Madhya Pradesh	C	985	Daily: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0039	Cardiolipin Antibody ACL-IgM antibody Serum	Madhya Pradesh	C	985	Daily: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0042	Catecholamines in Plasma	Madhya Pradesh	D	5885	Wed : 11:30am	Next day : 7:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
C0043	Catecholamines in urine	Madhya Pradesh	D	4815	Wed : 11:30am	Next day : 7:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea ,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0044	CBC Haemogram	Madhya Pradesh	A	280	Daily: 9am to 9pm	After 6 hrs	3 ml EDTA whole blood
C0045	inv (16) CBF Beta Translocation by fish	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
C0046	CBF Beta Inversion (16) by FISH AML-M4E0 Bone Marrow	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
C0047	Cyclic Citrullinated Peptide Antibody (CCP Antibody) Serum	Madhya Pradesh	C	1550	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0048	CD2 by Flow cytometry Blood**	Madhya Pradesh	C	2575	Daily: 9am	Next day 5pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0049	CD 3	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0050	CD 3/CD 4/CD 8	Madhya Pradesh	C	2130	Daily: 9am	Same day 5pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0051	CD 3/CD 4/CD 8/CD 45	Madhya Pradesh	C	2130	Daily: 9am	Same day 5pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0052	CD4 / CD 8 BY FLOW CYTOMETRY, BODY FLUID**	Madhya Pradesh	R	4285	Daily : 11:30am	5th day : 6:00pm	required. 5ml BAL fluids to be send in Sterile container. Clinical history and clinician number required. (Sample require as soon as possible after collection, Avoid weekend and Holiday)
C0053	CD 4 counts	Madhya Pradesh	C	1290	Daily: 9am	Same day 5pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0054	CD 4/CD 8	Madhya Pradesh	C	2130	Daily: 9am	Same day 5pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0055	CD 5	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0056	CD 7	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0057	CD 8	Madhya Pradesh	C	2130	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0058	CD 10	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	history and clinician number required. 3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0059	CD13	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0060	CD14	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0061	CD15	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0062	CD16/56	Madhya Pradesh	C	5040	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0063	CD19	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0064	CD20	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0065	CD22	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0066	CD23	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0067	CD25	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0068	CD33	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0069	CD34 Stem Cells	Madhya	C	6160	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0070	CD34	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0071	CD38	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0074	CD45	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0075	CD79A	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0076	CD103	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0077	CEA Carcino Embryonic Antigen Body Fluid	Madhya Pradesh	B	1090	Daily: 9am to 9pm	After 8 hrs	required. 2 ml of Body Fluid
C0078	CEA Carcino Embryonic Antigen Serum	Madhya Pradesh	B	840	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0079	Celiac Disease Profile	Madhya Pradesh	D	5565	Mon and Thu: 9am	Same day 7pm	5 ml of serum
C0080	Body fluids Cell Count	Madhya Pradesh	A	305	Daily: 9am to 9pm	After 6 hrs	Body fluid to be send in sterile EDTA container in refrigerated condition. Clinical history and clinician number required. Body fulids (CSF, Pleural, Peritoneal, Pericardial, synovial fluid)
C0081	Centromere Antibody Serum	Madhya Pradesh	C	1790	W,S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0082	Ceruloplasmin (Copper Oxidase)Serum	Madhya Pradesh	C	1310	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0083	Cervical Screening-1 (Conventional PAP & HPV-DNA) - Smear & LBC	Madhya Pradesh	D	3105	Daily : 9:00am	5th day : 6:00pm	Smear & Thin Prep LBC Container
C0084	Cervical Screening-2 (PAP Smear and HPV-DNA) - LBC	Madhya Pradesh	D	3500	Daily : 9:00am	5th day : 6:00pm	10 ml LBC Container
C0086	CH-50 Complement, total Serum	Madhya Pradesh	C	5265	Fri : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0087	Chikungunya RNA Detection by Real Time PCR - EDTA Plasma	Madhya Pradesh	C	5040	Daily : 10:00am	After 24hrs	3ml of EDTA Plasma in Sterile Plain Vacutainer
C0088	Chikungunya IgM antibody by Immunochromatography - Serum	Madhya Pradesh	B	1035	Daily : 9:00am	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0090	Chitotriosidase storage disorder Plasma	Madhya Pradesh	R	4080	Daily : 11:30am	10th day : 6:00pm	6 ml Heparin Plasma
C0091	Chlamydia Trachomatis DNA Detection by Real Time PCR - Eye Swab	Madhya Pradesh	D	4300	Wed & Sat : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0092	Chlamydia Trachomatis DNA Detection by	Madhya Pradesh	D	4300	Wed & Sat : 9:00am	Next day : 5:00pm	2ml of Body Fluid in Sterile leak

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Real Time PCR - Fluid	Pradesh					proof container
C0093	Chlamydia Trachomatis DNA Detection by Real Time PCR - Urine	Madhya Pradesh	D	4300	Wed & Sat : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
C0094	Chlamydia Trachomatis DNA Detection by Real Time PCR - Vaginal Swab	Madhya Pradesh	D	4300	Wed & Sat : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0096	Chlamydia Pneumoniae IgG Antibody - Serum	Madhya Pradesh	C	2185	Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0097	Chlamydia Pneumoniae IgM Antibody - Serum	Madhya Pradesh	C	2185	Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0098	Chlamydia Trachomatis IgA Antibody by EIA - Serum	Madhya Pradesh	C	1140	Mon, Wed & Fri : 7:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0099	Chlamydia Trachomatis IgG Antibody by EIA - Serum	Madhya Pradesh	C	1140	Mon, Wed & Fri : 7:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0100	Chlamydia Trachomatis IgM Antibody by EIA - Serum	Madhya Pradesh	C	1140	Mon, Wed & Fri : 7:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0101	Chloride (Cl) Serum	Madhya Pradesh	A	175	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0102	Chlorides Urine 24H	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
C0103	Chlorides Urine Spot	Madhya	A	275	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
C0104	Chlorides Body Fluid	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	2 ml of Body Fluid
C0105	Cholesterol-Total Body Fluid	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	2 ml of Body Fluid
C0106	Cholesterol-Total Serum	Madhya Pradesh	A	170	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0107	Chromium, Blood	Madhya Pradesh	D	3380	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0108	Chromium, Urine 24 Hrs	Madhya Pradesh	D	3530	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0109	ChromiumGraphite Furnace AASUrine Spot	Madhya Pradesh	D	3640	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0110	Chromogranin - A, Plasma	Madhya Pradesh	D	6900	Wed & Sat : 11:30am	Same day : 9:00pm	3 ml of EDTA plasma to be sent in aliquot tube at proper temperature.
C0113	Chromosomal Anaylsis (Fanconis Anemia)**	Madhya Pradesh	R	8935	Daily : 11:30am	20th day : 6:00pm	Clinical history and TRF mandatory. Blood : 5 ml Heparine Whole Blood, Sample sent as soon as possible after collection
C0115	Chylomicron Qualitative Urine	Madhya Pradesh	B	545	Daily: 7.30am	Next day 5pm	10 ml of Spot urine
C0117	Citrate, Urine 24 H Spectrophotometry	Madhya Pradesh	C	1290	Daily : 11:30am	5th day : 6:00pm	10 mL aliquot of 24-hour Urine. Collect urine with 25 mL 50% concentrated HCl to maintain pH below 3. Record 24-hour urine volume on test request form and urine vial. Ship refrigerated or frozen.
C0118	Citrulline Quantitative Serum	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
C0119	Citrulline Quantitative Urine	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
C0121	C-Kit (GIST) - Exon 9, 11, 13 and 17 by PCR - FFPE Tissue	Madhya Pradesh	D	11770	Mon, Wed & Fri: 9:00am	10th day : 6:00pm	Tumor containing FFPE block (>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							formalin with duration of 6 to 72 hours Clinical history with Histopathology report.
C0122	C-Kit - Exon 8 and 17 by PCR - Bone Marrow	Madhya Pradesh	D	5885	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.
C0123	C-Kit - Exon 8 and 17 by PCR - Blood	Madhya Pradesh	D	5885	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.
C0124	CK-MB Serum	Madhya Pradesh	B	710	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0128	Clostridium Difficile GDH Reflex Toxin A/B by CLIA - Stool	Madhya Pradesh	C	2465	Daily : 9:00am	Same day : 5pm	Category R: Referred (Please enquire before sending sample), Category I: International (Please enquire before sending sample)
C0130	Clotting Time (CT)	Madhya Pradesh	A	140	Daily: 9am to 9pm	After 6 hrs	NA
C0131	Clozapine Serum	Madhya Pradesh	D	4175	Wed : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0132	CMV Avidity test by EIA - Serum	Madhya Pradesh	C	1455	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0133	CMV DNA Detection by Real Time PCR - Amniotic Fluid	Madhya Pradesh	C	5375	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
C0134	CMV DNA Detection by Real Time PCR - CSF	Madhya Pradesh	C	5375	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	2ml of CSF in Sterile leak proof tube or container
C0135	CMV DNA Detection by Real Time PCR - EDTA Plasma	Madhya Pradesh	C	4590	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0136	CMV DNA Detection by Real Time PCR - BAL or Saliva	Madhya Pradesh	C	5375	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	2ml of BAL/ Saliva in Sterile Leak proof container
C0137	CMV DNA Detection by Real Time PCR - Urine	Madhya Pradesh	C	5375	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	15ml of Spot Urine in Sterile leak proof container
C0138	CMV Combo (CMV Viral load and GCV resistance) by PCR & Sequencing - EDTA Whole Blood**	Madhya Pradesh	R	11220	Tue : 11:30am	12th day : 6:00pm	3 ml of EDTA Whole Blood in EDTA Vacutainer
C0139	CMV IgG & IgM Antibody by CLIA - Serum	Madhya Pradesh	B	1310	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0141	CMV IgG Antibody by CLIA - Serum	Madhya Pradesh	B	720	DAILY : 9:00 AM to 9:00 PM	After 8Hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0142	CMV IgM Antibody by CLIA - Serum	Madhya Pradesh	B	720	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0143	CMV DNA Viral load by Real Time PCR - EDTA Plasma	Madhya Pradesh	C	8250	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
C0144	CMV by CISH - FFPE	Madhya Pradesh	D	2945	Daily : 7:30pm	Same day : 5pm	Tissue In 10% Formalin Or Block / FFPE .
C0145	Cobalt Blood	Madhya Pradesh	D	3530	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0145_CS	Cobalt Serum	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0146	Cobalt Urine 24H**	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0147	Cobalt Urine Spot	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0148	Cocaine Benzoylcegonine Urine Spot	Madhya Pradesh	B	740	Daily: 9am to 9pm	After 6 hrs	5 ml of Spot Urine
C0149	Cold AgglutininSerum	Madhya Pradesh	B	655	Daily : 11:30am	Next day : 5:00pm	3 ml of Serum
C0150	Cancer marker profile.	Madhya Pradesh	D	6400	Wed and Sat: 9am	next day 5.0 pm	3 ml of Serum and 3 ml of EDTA plasma
C0151	Coombs Test-Direct ,Direct Antiglobulin Test (DAT)-EDTA Blood	Madhya Pradesh	C	450	Daily: 9am and 5pm	After 6 hrs	3 ml EDTA whole blood
C0152	Coombs Test-Indirect, Indirect Antiglobulin Test (IAT)-Serum	Madhya Pradesh	C	450	Daily: 9am and 5pm	After 6 hrs	3 ml Serum sample
C0154	Copper ,24 hr Urine	Madhya Pradesh	D	2300	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
C0155	Copper, Urine Spot	Madhya Pradesh	D	2245	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0156	Copper Serum	Madhya Pradesh	D	2245	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0157	Cortisol Serum_8.00 AM	Madhya Pradesh	B	670	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0157_4PM	Cortisol Serum_04.00 PM	Madhya Pradesh	B	680	Daily: 9am to 9pm	Cortisol CLIA Serum_04.00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0157_R	Cortisol Serum_Random	Madhya Pradesh	B	680	Daily: 9am to 9pm	Cortisol CLIA Serum_Random	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0158	Cortisol - Free Urine 24H	Madhya Pradesh	B	980	Daily: 9am to 9pm	After 8 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0159	Coxiella Burnetti IgG antibody by EIA - Serum	Madhya Pradesh	D	2085	Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0160	Coxiella Burnetti IgM antibody by EIA - Serum	Madhya Pradesh	D	2085	Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0161	Coxsackie IgG Antibody by EIA - Serum	Madhya Pradesh	C	2800	Wed : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0162	Coxsackie IgM Antibody by EIA - Serum	Madhya Pradesh	E	2800	Wed : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0163	C-Peptide, Fasting, Serum	Madhya Pradesh	C	1210	Daily: 9am to 9pm	After 6 hrs	2 ml of Serum In Fasting
C0163_PP	C-Peptide Post prandial, Serum	Madhya Pradesh	C	1265	Daily: 9am to 9pm	C-PeptidePost prandial	2 ml of Serum In Fasting
C0163_R	C-Peptide, Random, Serum	Madhya Pradesh	C	1265	Daily: 9am to 9pm	C-Peptide, Random	2 ml of Serum In Fasting
C0164	CPK Isoenzymes Serum	Madhya Pradesh	D	8990	Fri : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0165	CPK Total Serum	Madhya Pradesh	A	380	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0166	Creatinine Clearance test 24H Urine	Madhya Pradesh	B	525	Daily: 9am to 9pm	After 6 hrs	3 ml Serum And 24 Hrs Urine . No Preservative
C0167	Creatinine Quantitative Urine 24H	Madhya Pradesh	A	315	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
C0168	Creatinine Body Fluid	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	2 ml of Body Fluid
C0169	Creatinine Serum	Madhya Pradesh	A	170	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0170	CRP, C - Reactive Protein Serum	Madhya Pradesh	B	250	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0171	Cryoglobulins Serum	Madhya Pradesh	C	1230	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0172	Cryptococcus Antigen detection by Agglutination - CSF	Madhya Pradesh	C	2520	Daily : 9:00am,1:00pm & 6:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
C0173	Cryptococcus Antigen detection by Agglutination - Serum	Madhya Pradesh	C	2520	Daily : 9:00am,1:00pm & 6:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0179	Cryptosporidium by smear examination - Stool	Madhya Pradesh	C	615	Daily : 9:00am to 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
C0180	Cryptosporidium Parvum Antigen Detection by Immunochromatography - Stool	Madhya Pradesh	C	1570	Daily : 9:00am	Same day : 5pm	5gm of Stool specimen in Sterile leak proof container
C0181	CSF Index (CSF and serum)	Madhya Pradesh	D	2780	Daily: 9am to 9pm	After 8 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage). 3 ml of Serum.
C0182	Culture & Sensitivity, Aerobic bacteria - Blood single bottle	Madhya Pradesh	C	1000	Daily : 9:00am to 9:00pm	After 5 days : 5pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
C0183	Culture & Sensitivity, Aerobic bacteria - Blood Two bottle set	Madhya Pradesh	C	3135	Daily : 9:00am to 9:00pm	After 5 days : 5pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
C0184	Culture & Sensitivity, Aerobic bacteria - Blood 1 Aerobic 1 Anaerobic	Madhya Pradesh	C	1790	Daily : 9:00am to 9:00pm	After 5 days : 5pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
C0185	Culture & Sensitivity, Aerobic bacteria - Body Fluids	Madhya Pradesh	C	1010	Daily : 9:00am to 9:00pm	Upto 72 hours	2ml of Body Fluid in Sterile leak proof container
C0186	Culture & Sensitivity, Aerobic bacteria - Bone Marrow	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	After 5 days : 5pm	3ml of Bone Marrow in MycoF/lytic Bactec Bottle or Sodium Heparin
C0187	Culture & Sensitivity, Aerobic bacteria - Central Line Catheter Tip	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	5th day : 5pm	Central Line Catheter Tip in sterile container
C0188	Culture & Sensitivity, Aerobic bacteria - Clot	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	After 5 days : 5pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
C0189	Culture & Sensitivity, Aerobic bacteria - CSF	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	2ml of CSF in Sterile leak proof tube or container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0190	Culture & Sensitivity, Aerobic bacteria - Ear Samples	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	Ear specimens in Sterile leak proof container
C0191	Culture & Sensitivity, Aerobic bacteria - Eye Samples	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	Eye specimens in Sterile leak proof container
C0192	Culture & Sensitivity, Aerobic bacteria - Other samples	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	2ml of Specimen in Sterile leak proof container
C0192_C	Culture and Identification, Aerobic bacteria - Other Samples	Madhya Pradesh	C	1400	Daily : 9:00am to 9:00pm	3rd day : 5pm	Pure culture isolate is required
C0193	Culture & Sensitivity, Aerobic bacteria - Pus	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	2ml of PUS in Sterile leak proof container
C0194	Culture & Sensitivity, Aerobic bacteria - Respiratory Samples BAL,ET Secretions	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	2ml of BAL / Other Respiratory samplesin Sterile leak proof container
C0195	Culture & Sensitivity, Aerobic bacteria - Semen	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	2ml of Semen in Sterile leak proof container
C0196	Culture & Sensitivity, Aerobic bacteria - Sputum	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	2ml of Sputum in Sterile leak proof container
C0197	Culture & Sensitivity, Aerobic bacteria - Stool	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	5gm of Stool specimen in Sterile leak proof container
C0198	Culture & Sensitivity, Aerobic bacteria - Throat Swab	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0199	Culture & Sensitivity, Aerobic bacteria - Throat Swab For C.Diphtheriae	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0200	Culture & Sensitivity, Aerobic bacteria - Tissue/Biopsy	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	Tissue with Normal Saline in Sterile leak proof container
C0201	Culture & Sensitivity, Aerobic bacteria - Urethral Discharge	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	2ml of Urethral Discharge in Sterile leak proof container
C0202	Culture & Sensitivity, Aerobic bacteria - Urine	Madhya Pradesh	C	860	Daily : 9:00am to 9:00pm	Upto 72 hours	15ml of Spot Urine in Sterile leak proof container
C0203	Culture & Sensitivity, Aerobic bacteria - Vaginal Sample Includes Gr B Strepto	Madhya Pradesh	C	1175	Daily : 9:00am to 9:00pm	Upto 72 hours	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0204	Culture and Identification, Anaerobic bacteria - Blood Single bottle	Madhya Pradesh	C	1570	Daily : 9:00am to 9:00pm	After 5 days : 5pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0205	Culture and Identification, Anaerobic bacteria - High Vagina Discharge	Madhya Pradesh	C	1570	Daily : 9:00am to 9:00pm	5th day : 5pm	2ml of Vaginal Discharge in Sterile leak proof container
C0206	Culture and Identification, Anaerobic bacteria - Pus	Madhya Pradesh	C	1570	Daily : 9:00am to 9:00pm	5th day : 5pm	2ml of PUS in Sterile leak proof container
C0207	Culture and Identification, Anaerobic bacteria (C.Difficile) - Stool	Madhya Pradesh	C	1570	Daily : 9:00am to 9:00pm	5th day : 5pm	5gm of Stool specimen in Sterile leak proof container
C0208	Culture and Identification, Anaerobic bacteria - Tissue/Biopsy	Madhya Pradesh	C	1570	Daily : 9:00am to 9:00pm	5th day : 5pm	Tissue with Normal Saline in Sterile leak proof container
C0209	Cyclospora by smear examination - Stool	Madhya Pradesh	C	615	Daily : 9:00am to 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
C0210	Cyclosporine-By LCMS, Blood	Madhya Pradesh	C	3135	Daily : 7:30am	Same day : 8:00pm	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
C0212	Cyfra (Ca) 21-1, Serum	Madhya Pradesh	D	3155	Tue & Fri : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0213	CYP2C19 Genotyping	Madhya Pradesh	D	8345	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius .
C0214	Cystatin C, Serum	Madhya Pradesh	C	1455	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0218	Cystic Fibrosis (CF) Delta F508 Mutation Analysis	Madhya Pradesh	E	4300	Tue & Thu: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
C0219	Cystic Fibrosis IRT (Neonatal screen)Dried	Madhya	C	505	Monday to Friday:	Next day 5pm	Dried blood spot: 1 drop of heel

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	blood spot	Pradesh			9.00am		prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
C0220	Cysticercus (Taenia Solium) IgG Antibody by EIA - CSF	Madhya Pradesh	C	1905	Tue & Fri : 7:30am	Same day : 5pm	2 ml of CSF in sterile screw capped container.
C0221	Cysticercus (Taenia Solium) IgG Antibody by EIA - Serum	Madhya Pradesh	C	1765	Tue & Fri : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0222	Cystine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
C0223	Cystine Quantitative Spot Urine	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
C0224	Cytochemistry Panel	Madhya Pradesh	C	2465	Daily: 9am	5th day	Direct smears & 3 ml of EDTA blood
C0226	Ascitic Fluid by Conventional method	Madhya Pradesh	C	1010		2nd day	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0227	Bronchoalveolar Lavage (BAL) by LBC	Madhya	C	1230		2nd day	Broncheoalveolar lavage in a plain

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					leak proof container along with brochscopy findings & relevant clinical/radiological details
C0228	CSF by Cyto centrifuge	Madhya Pradesh	C	1010		2nd day	More than 0.5 ML fluid in plain leak proof container along with relevant clinical details & CT/MRI reports
C0229	Drain Fluid by Conventional method	Madhya Pradesh	C	1010		2nd day	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0230	Nipple Discharge by Conventional method	Madhya Pradesh	C	1010		3rd day	Prefixed smears along with relevant clinical & radiological details
C0233	Other Fluids by Conventional method	Madhya Pradesh	C	1010		2nd day	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0234	Pericardial Fluid by Conventional method	Madhya Pradesh	C	1010		2nd day	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0235	Peritoneal washings by Conventional method	Madhya Pradesh	C	1010		2nd day	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0236	Pleural Fluid by Conventional method	Madhya Pradesh	C	1010		2nd day	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0238	Sputum by Conventional method	Madhya Pradesh	C	1010		2nd day	Sputum in a plain leak proof container along with relevant clinical/radiological details
C0239	Synovial Fluid by Conventional method	Madhya Pradesh	C	1010		2nd day	0.5 ml Synovial fluid in a plain container
C0241	Urine by Liquid based cytology	Madhya Pradesh	C	1230		2nd day	More than 1 ML Urine in a leak proof plain container along with relevant clinical & radiological details
C0242	Cytology Reflex to cell block and relevant IHC Body fluid	Madhya Pradesh	C	4480		5th day	More than 1 ML fluid in plain leak proof container along with relevant clinical details
C0244	Carbohydrate-Deficient Transferrin (CDT) **	Madhya Pradesh	R	4080	Daily: 9am	Same day 5 PM	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tests. International outsource : Category R for referred and I for international. Plz esquire before sending sample for arrangement of proper courier to international lab. Additional sample tubes may be needed for sample adequacy. TAT may vary based on information by concerned lab of delay. Courier charges have to be paid extra based on sample temperature condition and shipment.
C0246	Chromosomal Array CNV+SNP Array Blood	Madhya Pradesh	E	16000	Daily : 9:00am	15th day : 6:00pm	3 ml of EDTA Whole Blood
C0247	Chromosomal Array CNV+SNP Array POC	Madhya Pradesh	E	16800	Daily: 9am	15th Day	1)Upto 12 weeks of gestation: 50-100mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container 2) In 12-18 weeks of gestation: a) 50-100mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container b) and or fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container 3)IUFD >20 weeks of gestation: fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container. if possible in 20ml normal saline with 2-3 drops of gentamycin
C0248	Clostridium Difficile DNA Detection by Real Time PCR - Stool	Madhya Pradesh	C	3695	Wed & Sat : 7:30am	Next day : 5:00pm	5gm of Stool specimen in Sterile leak proof container
C0249	Connexin (26 mutations)_C0249	Madhya Pradesh	R	6630	Wed : 11:30am	10th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0250	Copper reflex Ceruloplasmin, Serum	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C0251	Calreticulin (CALR) Exon 9 Mutation Analysis, Blood	Madhya Pradesh	D	6740	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
C0252	Calreticulin (CALR) Exon 9 Mutation Analysis, Bone marrow	Madhya Pradesh	D	6740	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
C0253	Chromium, Serum	Madhya Pradesh	R	3570	Daily : 11:30am	12th Day : 6:00pm	BD RED TOP VACCUTAINER, CENTRIFUGE, SEPARATE SERUM in red top
C0255	Creatinine Urine Spot	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
C0256	CPEO (Chronic Progressive External Ophthalmoplegia) Mitochondrial diseases	Madhya Pradesh	R	7345	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
C0260	NextGen Cancer Hotspot 50 Gene Panel	Madhya Pradesh	E	28560	Daily : 11:30am	After 4 weeks : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0264	CD64, Blood by Flow Cytometry	Madhya Pradesh	C	2575	Daily : 11:30am	Next day : 5:00pm	hours and paraffin embedded. Clinical history with Histopath report. 3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
C0265	CEBPA MUTATION ANALYSIS	Madhya Pradesh	E	6000	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood or 3 ml EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
C0267	Cervical Screening-3 (PAP Smear, HPV-DNA and reflex CINtec PLUS) - LBC	Madhya Pradesh	D	5885	Daily : 9:00am	7th day : 6:00pm	10 ml LBC Container
C0269	Cancer marker profile (Brain and Pituitary)	Madhya Pradesh	D	3690	Tue and Fri: 9am	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0272	Cancer marker profile (Gall Bladder)	Madhya Pradesh	D	1710	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
C0273	Cancer marker profile Liver-2	Madhya Pradesh	D	6155	Wed and Sat: 9am	next day 5.0 pm	3 ml of Serum and 3 ml of EDTA plasma
C0275	Cancer marker profile (Neuroendocrine tumors)	Madhya Pradesh	D	6635	Wed/Sat- 7:30 am	Next day 5:00 pm	3 ml of Serum and 3 ml of EDTA plasma
C0276	Cancer marker profile-Ovarian-Comprehensive	Madhya Pradesh	D	2410	Daily: 9am to 9pm	After 8 hrs	5 ml of serum
C0277	Cancer marker profile - Thyroid-DTC	Madhya Pradesh	D	2355	Daily: 9am to 9pm	After 8 hrs	5 ml of serum
C0279	Cytology CINtec plus	Madhya Pradesh	D	4175		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0281	cMYC Gene rearrangement, FISH, BONE MARROW	Madhya Pradesh	D	4710	Daily : 9:00am	3rd day : 6:00pm	detail 2-3 ml Bone Marrow in Na-heparin vacutainer
C0282	CMYC Gene rearrangement, FISH, BLOOD	Madhya Pradesh	D	4710	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
C0283	cMYC Gene rearrangement Tissue / FFPE	Madhya Pradesh	D	6420	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
C0284	Coeliac disease profile-1	Madhya Pradesh	D	2780	Mon and Thu: 9am	Same day 7pm	5 ml of serum
C0285	Coeliac disease profile-2 Comprehensive	Madhya Pradesh	D	10595	Mon and Thu: 9am	3rd day	5 ml of serum and 4ml EDTA / ACD blood
C0286	Catechol-O-Methyltransferase (COMT) Genotyping	Madhya Pradesh	E	6000	Daily : 9:00am	7th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0289	CMV DNA Detection by Real Time PCR - Tissue	Madhya Pradesh	C	5375	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	Tissue with Normal Saline in Sterile leak proof container
C0290	Chloride, CSF	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
C0292	CMV DNA Viral load by Real Time PCR - BAL or Saliva	Madhya Pradesh	C	8290	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	2ml of BAL/ Saliva in Sterile Leak proof container
C0297	C3 Nephritic Factor	Madhya Pradesh	R	5100	Daily : 11:30am	8th day : 6:00pm	3 ml Serum
C0302	Cryofibrinogen	Madhya Pradesh	D	965	Daily : 11:30am	4th day : 6:00pm	Collect Blood In Pre-warmed Red Plain Top & Citrate Tube And Pre-

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0306	CHIMERISM, PRE-TRANSPLANT (STR ANALYSIS - DONOR & RECIPIENT)*	Madhya Pradesh	E	9900	Daily : 9:00am	7th day : 6:00pm	warmed Syringe. 3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Chimerism Analysis - Information Sheet Clinical details of the recipient Bone Marrow Transplant Date, Donor and Recipient relationship
C0307	CHIMERISM (Post transplant) *	Madhya Pradesh	E	9900	Daily : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Chimerism Analysis Information Sheet Clinical details of the recipient Bone Marrow Transplant Date, Donor and Recipient relationship
C0311	Cardiomyopathy gene panel**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
C0318	Chylomicron Qualitative Body fluid	Madhya Pradesh	B	545	Daily: 9am to 9pm	After 8 hrs	2 ml body fluid
C0319	Clinical Exome Panel**	Madhya Pradesh	R	22050	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0323	Cell free DNA for EGFR*	Madhya Pradesh	E	16000	Daily : 9:00am	7th day : 6:00pm	Genetci history is mandatory. 10-15 ml of blood in streck tube required. Sample stable for 3-7 Days at Ambient Temperature. For streck tubes, enquire before sending the sample. Clinical history with Histopathology report.
C0326	Carba R (Carbapenem Resistance) by GeneXpert - Culure Plate / Anal Swab	Madhya Pradesh	D	10380	Daily : 9:00am & 3:00pm	After 6hrs	Culure Plate / Anal Swab with Sufficient Material in Sterile Falcon Tube
C0332	Cortisol/Creatinine ratio in urine	Madhya Pradesh	B	1200	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
C0337_Adult	Creatinine, advanced (Above18 years) - spot urine	Madhya Pradesh	A	285	Daily: 9am to 9pm	Creatinine, advanced(Above18 years)	10 ml of Spot urine
C0338	Haemogram, advanced	Madhya Pradesh	A	350	Daily: 9am to 9pm	After 6 hrs	3 ml of EDTA whole blood sample
C0340	Chromosomal Array CGH - Prenatal*	Madhya Pradesh	E	16000	Daily : 9:00am	15th day : 6:00pm	Amniotic fluid (20-30ml divided in 2 sterile falcon centrifuge tube) and 1-4ml Maternal blood in EDTA tube is required for maternal cell contamination (MCC) studies. (Completely filled Consent form G with patient and doctor stamp and signature is mandetory as per PCPNDT Act)
C0342	Allergy panel - Child panel	Madhya Pradesh	D	13215	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
C0343	Clustered AchR and MuSK and LRP4 Antibodies**	Madhya Pradesh	R	15300	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
C0347	CSF ACPII**	Madhya Pradesh	R	27540	Daily : 11:30am	After 4 weeks : 6:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage). International outsource : Category

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							R for referred and I for international. Plz esquire before sending sample for arrangement of proper courier to international lab. Additional sample tubes may be needed for sample adequacy. TAT may vary based on information by concerned lab of delay. Courier charges have to be paid extra based on sample temperature condition and shipment.
C0348	Chromosomal Array CGH-Cancer	Madhya Pradesh	E	25000	Daily : 9:00am	After 4 Weeks : 6:00pm	3 ml of EDTA whole blood and 3 ml of Bone marrow
C0349	Cadmium, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
C0350	Cadmium, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
C0351	Chromium , Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
C0352	Chromium, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
C0353	Cobalt , Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
C0354	Cocaine confirmation by LCMS, Spot Urine**	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C0356	Copper, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	e.g C0 or C2) and clinical History along with the sample. 1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
C0357	Copper, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
C0359	COCKROACH, GERMAN BLATELLA GERMANICA (SPECIFIC IGG4)	Madhya Pradesh	D	1605	Mon,Thu: 9am	Next day: 5pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
C0360	CLL prognostication/risk stratification gene panel**- EDTA Whole Blood	Madhya Pradesh	E	17000	Daily : 11:30am	20th day : 6:00pm	4 ml EDTA Whole blood/bone marrow aspirate
C0361	Carrier ?screening - Gold (500 genes)	Madhya Pradesh	E	41000	Daily: 7:30am	5 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0362	Carrier ?screening - Platinum (2000 genes)	Madhya Pradesh	E	67000	Daily: 7:30am	6 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0363	Carrier ?screening - Silver (100 genes)	Madhya	E	35100	Daily: 7:30am	5 weeks	3 ml of EDTA whole blood required

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0364	Congenital adrenal hyperplasia due to HSD3B2-Gene Analysis	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0365	Congenital cataract-Gene Panel	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0366	Congenital hepatic fibrosis-Gene Panel	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
C0368	Chromosomal Array CGH-High resolution	Madhya	R	25500	Daily : 11:30am	20th day : 6:00pm	3 ml of EDTA Whole Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
C0369	Chikungunya RNA Detection by Real Time PCR - Serum	Madhya Pradesh	C	5040	Daily : 10:00am	24 Hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
C0370	CMV DNA Viral load by Real Time PCR - CSF	Madhya Pradesh	C	8290	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CSF in Sterile leak proof tube or container
C0371	CMV DNA Viral load by Real Time PCR - Urine	Madhya Pradesh	C	8290	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	15ml of Spot Urine in Sterile leak proof container
C0372	CT - NG (Chlamydia Trachomatis & Neisseria Gonorrhoeae) Panel - Genital Swab	Madhya Pradesh	D	5780	Wed & Sat : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
C0373	CT - NG (Chlamydia Trachomatis & Neisseria Gonorrhoeae) Panel - Urine spot	Madhya Pradesh	D	5780	Wed & Sat : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
C0374	Covid Monitor Initial	Madhya Pradesh	C	6495	Daily: 9am	Same day 5pm	2 ml of EDTA whole blood, 3 ml of serum ,3 ml of citrate plasma
C0375	Covid Monitor Maintenance.	Madhya Pradesh	C	4370	Daily: 9am	Same day 5pm	2 ml of EDTA whole blood, 3 ml of serum ,3 ml of citrate plasma
C0376	Covid Monitor Recovery.	Madhya Pradesh	C	3135	Daily: 9am to 9 pm	Report after 24 hours	2 ml of EDTA whole blood, 3 ml of serum ,3 ml of citrate plasma
C0377	Cardiac Screen-1	Madhya Pradesh	C	6720	Mon and Thu: 9am	Same day 5pm	5 ml of serum
C0378	Conventional Karyotyping for CLL (B-Cell Stimulation)**	Madhya Pradesh	R	12240	Daily : 11:30am	20th day : 6:00pm	10 ml Heparin Blood
C0379	CoviRisk TBNK by Flow Cytometry - Blood	Madhya Pradesh	D	4815	Daily : 9:00am to 9:00pm	After 6 hrs	3 ml of EDTA Whole Blood
C0381	Covid Associated Biomarker Profile	Madhya Pradesh	C	5040	9 am to 9 pm	Report after 24 hrs	2 ml of EDTA whole blood, 3 ml of serum ,3 ml of citrate plasma
C0393	Colorectal cancer - Theranostic Panel (Hot Spot) panel	Madhya Pradesh	E	32000	Daily 7.30 am	4 Weeks	Paraffin Block with Histopath report Mandatory
C0394	Cytology - Brush Cytology by conventional method	Madhya Pradesh	C	2800	Daily : 9:00am	3rd day	Unstained slides, Complete clinical & radiology details, clinician contact number
C1001	C1q antibodies test, Serum	Madhya Pradesh	R	6630	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
C1014	Carnitine free and total plasma**	Madhya Pradesh	R	12750	Daily : 11:30am	10th day : 6:00pm	3 ml Heparin Plasma
C1802	Colour Doppler - Carotid	Madhya	O	2500			

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
C1803	Colour Doppler - Limb (Arterial and Venous)	Madhya Pradesh	O	4000			
C1804	Colour Doppler - Limbs - both (Arterial and Venous)	Madhya Pradesh	O	8000			
C1805	Colour Doppler - Obstetric	Madhya Pradesh	O	2000			
C1806	Colour Doppler - Pelvis	Madhya Pradesh	O	2200			
C1807	Colour Doppler - Portal System	Madhya Pradesh	O	2200			
C1808	Colour Doppler - Renal	Madhya Pradesh	O	2200			
C1809	Colour Doppler - Scrotal	Madhya Pradesh	O	2200			
C1818	Cephalogram	Madhya Pradesh	O	500			
C1901	Canine - Calcium	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1902	Canine - CBC Haemogram	Madhya Pradesh	A	340	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1903	CBF (Complete Body Function)	Madhya Pradesh	C	3585	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							animal. Please correlate clinically and to be interpreted by the clinician.
C1904	Canine - Chlorides	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1905	Canine - Cholesterol-Total	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1906	Canine - Combo Profile	Madhya Pradesh	C	4030	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1907	CANINE-Cortisol	Madhya Pradesh	C	700	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1908	CANINE-CPK-Creatinine Phospho Kinase	Madhya Pradesh	A	380	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1909	Canine - Creatinine	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1910	CANINE-Culture & Sensitivity, Aerobic bacteria Blood single bottle	Madhya Pradesh	C	895	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1911	CANINE-Culture & Sensitivity, Aerobic bacteria Body Fluids Specify type	Madhya Pradesh	C	870	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1912	CSF	Madhya Pradesh	C	1175	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1913	Eye Samples	Madhya Pradesh	C	1175	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1914	nasal swab	Madhya Pradesh	C	1050	Daily: 9am to 9pm	Report after 24 hrs	animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1915	CANINE-Culture & Sensitivity, Aerobic bacteria Pus	Madhya Pradesh	C	1175	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1916	Semen	Madhya Pradesh	C	1175	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1917	CANINE-Culture & Sensitivity, Aerobic bacteria Sputum	Madhya Pradesh	C	870	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1918	CANINE-Culture & Sensitivity, Aerobic bacteria Stool	Madhya Pradesh	C	870	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1919	Throat Swab	Madhya Pradesh	C	1175	Daily: 9am to 9pm	Report after 24 hrs	animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1920	Tissue/Biopsy	Madhya Pradesh	C	1175	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1921	Urethral Discharge	Madhya Pradesh	C	1175	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1922	CANINE-Culture & Sensitivity, Aerobic bacteria Urine	Madhya Pradesh	C	895	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1923	Vaginal Sample	Madhya Pradesh	C	1175	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1924	Second opinion services, Slides	Madhya Pradesh	C	1120	Daily: 9am	Report after 48 hrs	animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1925	Ascitic Fluid	Madhya Pradesh	C	1010	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1926	CSF	Madhya Pradesh	C	1010	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1927	Drain Fluid	Madhya Pradesh	C	1010	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1928	Other Fluids	Madhya Pradesh	C	1120	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C1929	Other Fluids/Scraping	Madhya Pradesh	C	1010	Daily: 9am	Report after 48 hrs	animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1930	Pericardial Fluid	Madhya Pradesh	C	1010	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1931	Peritoneal washings	Madhya Pradesh	C	1010	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1932	Pleural Fluid	Madhya Pradesh	C	1010	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1933	Sputum	Madhya Pradesh	C	1010	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							animal. Please correlate clinically and to be interpreted by the clinician.
C1934	Synovial Fluid	Madhya Pradesh	C	1010	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
C1936	Metanephrine, Urine spot	Madhya Pradesh	C	3080	Daily: 9am to 9pm	Report after 5 days	Spot Urine : 10 ml of Spot sample in Leakproof container. Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
C5054	Cadmium / Creatinine ratio	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
C5516	CBC+MP	Madhya Pradesh	B	435	Daily: 9am to 9pm	After 6 hrs	3 ml EDTA whole blood
C5523	CBC, MP, ESR Panel	Madhya Pradesh	B	520	Daily: 9am to 9pm	After 6 hrs	3 ml EDTA whole blood
C5631	COVIPROTECT Plus	Madhya Pradesh	C	2690	Mon, Wed, Fri 1:00 PM	Next Day 5:00 PM	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
C5643	Comprehensive metabolic panel	Madhya Pradesh	D	5020	Daily: 9am to 9pm	After 6hrs	3 ml of EDTA whole blood, 2 ml of Flouride plasma, 4 ml of serum, 20 ml of spot urine
C8103	Cytopathology-Expert opinion	Madhya Pradesh	C	2800		7th day	Stained/unstained slides with original cytopathology report & relevant clinical & radiological details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
C8369	COVID VARIANT SCREENING	Madhya Pradesh	D	805	24hr registration	After 24hrs	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
C8551	COMPREHENSIVE OPHTHALMIC GENETIC DISORDER PANEL	Madhya Pradesh	R	20400	Daily 7.30 am	4 Weeks	3 ml of EDTA whole blood
C8606	Critical care antibiotics for drug resistant Gram negative bacteria.	Madhya Pradesh	D	1200	Daily : 9:00am to 9:00pm	4th day : 6:00pm	Pure culture isolate is required
C8623	Celiac Disease Profile - 1	Madhya Pradesh	C	2100	Mon & Thu : 7:30am	Same day : 5:00pm	5 ml of serum
D0001_DCP	DCP- Decarboxy Prothrombin PIVKA II, Serum	Madhya Pradesh	D	3550	Mon, Wed, Fri : 7:30am	Same day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0002	D-Dimer Quant, plasma	Madhya Pradesh	C	1000	Daily: 9am and 3pm and 9pm	After 6 hrs	1 ml of citrated platelet poor plasma in frozen condition.
D0003	Dengue IgG Antibody by CLIA - Serum	Madhya Pradesh	C	1010	Daily : 6:30pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0004	Dengue IgM Antibody by CLIA - Serum	Madhya Pradesh	C	670	Daily : 6:30pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0005	Dengue NS1 antigen Detection by EIA - Serum	Madhya Pradesh	C	630	Daily : 6:30pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0006	Dengue Profile-1 (Dengue-IgG, IgM, NS1, Platelet count) by EIA - Serum & EDTA blood	Madhya Pradesh	C	2240	Daily : 6:30pm	After 8 hrs	3 ml of Serum in Gel Top Vacutainer or Plain Sterile Vacutainer & 3 ml of EDTA Whole Blood in EDTA Vacutainer
D0006_2	Dengue Profile-2 (Dengue-IgG, IgM, NS1, Platelet count & IPF) by EIA - Serum & EDTA blood	Madhya Pradesh	C	2465	Daily : 6:30pm	After 8 hrs	3 ml of Serum in Gel Top Vacutainer or Plain Sterile Vacutainer & 3 ml of EDTA Whole Blood in EDTA Vacutainer
D0007	Dengue RNA Detection by Real Time PCR - CSF	Madhya Pradesh	C	4145	Daily : 10:00am	After 24hrs	2ml of CSF in Sterile leak proof tube or container
D0008	Dengue RNA Detection by Real Time PCR - EDTA Plasma	Madhya Pradesh	C	4145	Daily : 10:00am	After 24hrs	3ml of EDTA Plasma in Sterile Plain Vacutainer
D0009	Desmoglein I Antibody Vesiculobullous disorders Serum	Madhya Pradesh	D	3000	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0010	Desmoglein III Antibody Vesiculobullous disorders Serum	Madhya Pradesh	D	3000	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
D0011	Dexamethasone Suppression test-High dose for ACTH, Plasma	Madhya Pradesh	C	2745	Daily: 9am to 9pm	Same day 5pm	3 ml of EDTA Plasma each
D0012	Dexamethasone Suppression test for Cortisol-High dose, Serum	Madhya Pradesh	C	1455	Daily: 9am to 9pm	Same day 5pm	3 ml Serum each
D0013	Dexamethasone Suppression test-Low dose for ACTH, Serum	Madhya Pradesh	C	2745	Daily: 9am to 9pm	Same day 5pm	3 ml of EDTA Plasma each
D0014	Dexamethasone Suppression test for Cortisol-Low dose, Serum	Madhya Pradesh	C	1455	Daily: 9am to 9pm	Same day 5pm	3 ml Serum each
D0015	Dexamethasone Suppression test Overnight suppression For ACTH, Plasma	Madhya Pradesh	C	2745	Daily: 9am to 9pm	Same day 5pm	3 ml of EDTA Plasma each
D0016	Dexamethasone Suppression test Overnight suppression For Cortisol, Serum	Madhya Pradesh	C	1455	Daily: 9am to 9pm	Same day 5pm	3 ml Serum each
D0017	DHEA- Dehydroepiandrostedione Serum	Madhya Pradesh	C	3080	Tue,Fri: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0018	DHEAS- Dehydroepiandrostedione Sulphate Serum	Madhya Pradesh	B	1150	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0019	DHT- Dihydrotestosterone Serum	Madhya Pradesh	C	2950	M,W,F: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0020	Diabetes Profile, Mini	Madhya Pradesh	E	925	Daily: 8pm	After 8 hrs	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma, 10 ml of spot urine
D0021	Diabetes Profile,Maxi	Madhya Pradesh	E	3000	Wed: 7.30am	Same day 5pm	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma (F, PP), 10 ml of spot urine
D0022	Dialysis Panel	Madhya Pradesh	C	4480	Daily: 9am to 9pm	After 8 hrs	3ml EDTA blood, fluoroide plasma & Serum
D0023	DIC (Coagulation) Profile	Madhya Pradesh	C	2910	Daily: 9am to 9pm	After 8 hrs	3 ml of EDTA whole blood and 1 ml of citrated platelet poor plasma
D0024	Digoxin Lanoxin Serum	Madhya Pradesh	C	1065	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0025	Diphtheria IgG antibodies by EIA - Serum	Madhya Pradesh	D	2675	Fri : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0026	DMD/BMD (18 Exons) for Males Only	Madhya	D	6315	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
D0027	DMD/BMD-79 exons Deletion/Duplication Analysis MLPA - Blood	Madhya Pradesh	D	15900	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
D0028	DMD-Test Duchenne Muscular Dystrophy Amniotic Fluid*	Madhya Pradesh	D	18190	Daily : 9:00am	7th day : 6:00pm	Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
D0028_c	DMD/ BMD-79 exons Deletion/ Duplication Analysis (Dystrophin gene - Xp21)	Madhya Pradesh	D	18190	Daily : 9:00am	7th day : 6:00pm	Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition 3-5 ml EDTA Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							sample of Mother is mandatory. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
D0030	DNA (Double Strand) Antibody By IFA - Reflex end point titre Serum	Madhya Pradesh	C	1850	Tue,Fri: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0032	DNA (Double Strand) Antibody NcX DNA Serum	Madhya Pradesh	C	1850	Daily: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0033	DNA IsolationBlood	Madhya Pradesh	C	1230	Daily : 9:00am	Next day : 5:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
D0034	DNA Paternity 16 markers *	Madhya Pradesh	E	10000	Daily : 9:00am	10th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Duly signed Paternity Consent form and documents mandatory. Photo Identity proof of father and child 2 Passport size photograph attached to consent form, Birth certificate of child In case of alleged father mother declaration required.
D0036	DNPH (Dinitro Phenyl hydralazine) test, urine	Madhya Pradesh	B	380	Daily : 11:30am	After 6hrs	10 ml of Spot urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
D0037	Dopamine ELISA Plasma	Madhya Pradesh	D	2675	Wed : 11:30am	Next day : 7:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
D0038	Dopamine, HPLC, 24 hrs urine	Madhya Pradesh	D	3315	Wed : 11:30am	Next day : 7:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
D0039	Maternal screen-1st Trimester Dual Marker test serum	Madhya Pradesh	C	2500	Daily 07:30 am	same day 09:00 pm	3 ml of Serum
D0039_reflex	Pregascreen Dual (First trimester test-VII, serum Dual Marker test reflex Confirmation Either NIPT or Karyotyping/FISH of CVS)	Madhya Pradesh	D	3500	Daily 07:30 am	same day 09:00 pm	3 ml of serum , If screen is positive reflex confirmation Either NIPT or Karyotyping and FISH of amniotic fluid
D0040	Maternal screen-2nd Trimester- Quadruple Marker test, Serum	Madhya Pradesh	C	3000	Daily 07:30 am	same day 09:00 pm	3 ml of Serum
D0040_reflex	Pregascreen Quadruple (Second trimester test-IV, serum Quadruple marker test reflex confirmation Either NIPT or Karyotyping/ FISH of amniotic fluid)	Madhya Pradesh	C	4000	Daily 07:30 am	same day 09:00 pm	3 ml of serum , If screen is positive reflex confirmation Either NIPT or Karyotyping and FISH of amniotic fluid
D0041	Maternal screen (Triple test) 2nd trimester	Madhya Pradesh	C	3135	Daily: 9am to 9pm	Next day at 6pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D0042	DPD-Deoxyypyridinoline (Pyrilinks D)Urine	Madhya Pradesh	C	1960	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot Urine (Protect From Light)
D0043	Dihydropyrimidine Dehydrogenase (DPD)	Madhya	D	7490	Tue & Fri : 9:00am	5th day : 6:00pm	3-5 ml of EDTA whole blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Gene Mutation Analysis	Pradesh					required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
D0044	DPT Antibodies by EIA - Serum	Madhya Pradesh	D	6955	Fri : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D0045	Drugs Of Abuse Panel-5 Drug Panel, Urine	Madhya Pradesh	C	2520	Daily: 9am to 9pm	After 8 hrs	10 ml of Spot urine
D0046	Drugs Of Abuse Panel-6 Drug Panel, Urine	Madhya Pradesh	C	3305	Daily: 9am to 9pm	After 8 hrs	10 ml of Spot urine
D0047	Drugs Of Abuse Panel - 9 Drug, Urine	Madhya Pradesh	C	4590	Daily: 9am to 9pm	After 8 hrs	10 ml of Spot urine
D0048	DVT-Deep Vein Thrombosis Panel	Madhya Pradesh	D	12840	Mon and Thurs: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
D0049	Dysmorphic RBC detection - Spot Urine	Madhya Pradesh	B	165	Daily : 9:00am to 9:00pm	After 4 hrs	15ml of Spot Urine in Sterile leak proof container
D0050	D3 Hydroxybutarate (Ranbut), Serum	Madhya Pradesh	C	1010	Daily: 9am to 9pm	After 8 hrs	3ml of serum
D0053	Drugs Of Abuse Panel-12 Drug Panel, Urine	Madhya Pradesh	C	6830	Daily: 9am to 9pm	After 10 hrs	10 ml of Spot urine
D0054	DRPLA Gene Analysis	Madhya Pradesh	E	6220	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
D0059	Diabetes Complications assessment	Madhya Pradesh	E	2850	Mon and Thu: 9am	Same day 5pm	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma (F, PP), 10 ml of spot urine
D0060	Diabetes Risk Assessment Profile	Madhya Pradesh	E	2900	Wed: 9am	Next day 5pm	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of flouride plasma
D0061	Diabetes type I profile	Madhya Pradesh	E	11000	Monday & Thursday 9am	Next day 5pm	5 ml of serum
D0125	Dystonia gene panel	Madhya Pradesh	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
D0127	DFS-70 Antibody	Madhya Pradesh	C	2350	Wed : 11:30am	Same day : 9:00pm	2 ml serum required
D0128	Deafness panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
D0129	DLBCL (Diffuse large B-cell lymphoma) FISH PANEL Tissue	Madhya Pradesh	E	14000	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
D0130	11 Deoxycorticosterone**	Madhya Pradesh	R	15300	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
D0138	DERMATOPHAGOIDES FARINAE HOUSE DUST MITE (SPECIFIC IGG4)	Madhya Pradesh	C	1680	Mon,Thu: 9am	Next day: 5pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
D0139	DERMATOPHAGOIDES PTERONYSSINUS HOUSE DUST MITE (SPECIFIC IGG4)	Madhya Pradesh	C	1680	Mon,Thu: 9am	Next day: 5pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
D0140	Deletion/ Duplication Analysis	Madhya Pradesh	R	10810	Daily : 11:30am	After 4 weeks : 6:00pm	for 7 days and at Ambient temperature for a day. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
D0141	Dengue RNA Detection by Real Time PCR - Serum	Madhya Pradesh	C	4145	Daily : 10:00am	24 Hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
D1802	Digital X- Ray- Ankle - AP / Lateral view	Madhya Pradesh	O	800			
D1803	Digital X- Ray- Cervical Spine - Ap / Lateral	Madhya Pradesh	O	800			
D1804	Digital X- Ray- Chest - Decubits	Madhya Pradesh	O	400			
D1805	Digital X- Ray- Chest - Lateral	Madhya Pradesh	O	400			
D1809	Digital X- Ray- Elbow joint - AP view	Madhya Pradesh	O	400			
D1810	Digital X- Ray- Femur bone - AP / Lateral	Madhya Pradesh	O	800			
D1812	Digital X- Ray- Foot - AP / Oblique	Madhya Pradesh	O	800			
D1814	Digital X- Ray- Knee - Both Skyline View	Madhya Pradesh	O	400			
D1816	Digital X- Ray- Knee - joint AP / Lat	Madhya Pradesh	O	800			
D1817	Digital X- Ray- Knee Joints (both) - AP/ Lateral view	Madhya Pradesh	O	1200			
D1818	Digital X- Ray- KUB	Madhya Pradesh	O	800			
D1819	Digital X- Ray- L S Spine - AP / Lat	Madhya	O	800			

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
D1820	Digital X- Ray- L S spine - AP / Lat / Coned down	Madhya Pradesh	O	1200			
D1821	Digital X- Ray- L S Spine - Lateral	Madhya Pradesh	O	400			
D1822	Digital X- Ray- L S spine both - oblique	Madhya Pradesh	O	800			
D1823	Digital X- Ray- Mastoids Schuller s View	Madhya Pradesh	O	400			
D1824	Digital X- Ray- Nasal Bone	Madhya Pradesh	O	800			
D1825	Digital X- Ray- Pelvis with both hip	Madhya Pradesh	O	400			
D1827	Digital X- Ray- PNS Caldwell	Madhya Pradesh	O	400			
D1829	Digital X- Ray- Sacro coccyx Lat	Madhya Pradesh	O	400			
D1830	Digital X- Ray- Sacro iliac joint AP	Madhya Pradesh	O	400			
D1831	Digital X- Ray- Shoulder Joint AP /Lat	Madhya Pradesh	O	800			
D1832	Digital X- Ray- Skull Lat	Madhya Pradesh	O	400			
D1833	Digital X- Ray- Skull AP view	Madhya Pradesh	O	400			
D1834	Digital X- Ray- Tibia / Fibula AP / Lat	Madhya Pradesh	O	800			
D1835	Digital X- Ray- Wrist / Hand AP	Madhya Pradesh	O	400			
D1836	Digital X- Ray- Wrist joint AP / Lat	Madhya Pradesh	O	800			
D1854	Digital X- Ray - Right Mastoid	Madhya Pradesh	O	400			
D1901	DIC (Coagulation) Profile	Madhya Pradesh	C	2465	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
D5612	PKG D5612 (Diabetes Platinum)	Madhya Pradesh	E	3400	Daily	Next day 5.00 pm	4 ml Serum ,3 ml EDTA,3 ml Flouride.1Urine container
D5613	D5613 (Diabetes Gold)	Madhya Pradesh	E	1400	Daily	Next day 5.00 pm	4 ml Serum ,3 ml EDTA,3 ml Flouride.1Urine container
D8252	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR and Karyotyping, CVS - Prenatal	Madhya Pradesh	E	15500	Daily : 9:00am	15th day : 6:00pm	<p>â€¢ Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition â€¢ 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius .</p>
D8254	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR and Karyotyping, Amniotic Fluid-Prenatal	Madhya Pradesh	E	15500	Daily : 9:00am	15th day : 6:00pm	<p>â€¢ Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition â€¢ 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius</p>
D8259	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR and Chromosomal Array CGH, Amniotic Fluid - Prenatal	Madhya Pradesh	E	18000	Daily : 9:00am	15th day : 6:00pm	<p>â€¢ Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated</p>

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
D8263	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR and Chromosomal Array CGH, POC - Prenatal	Madhya Pradesh	E	18000	Daily : 9:00am	15th day : 6:00pm	<p>condition – 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius</p> <p>– POC - 1)Upto 12 weeks of gestation: 50-100 mg placental villi (fetal side) in 20 ml normal saline with 2-3 drops of gentamycin in a sterile container 2) In 12-18 weeks of gestation: a) 50-100 mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container b) and or fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container 3)IUFD >20 weeks of gestation: fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container. if possible in 20ml normal saline with 2-3 drops of gentamycin – 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius</p>
D8264	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR, Karyotyping and Chromosomal Array CGH, Amniotic Fluid – Prenatal	Madhya Pradesh	E	19000	Daily : 9:00am	15th day : 6:00pm	<p>– Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition – 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius</p>
D8266	ALLERGY DRUG PANEL - (COMPREHENSIVE - 35 DRUG)	Madhya Pradesh	E	16500	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
D8267	ALLERGY DRUG PANEL -(ANTIBIOTICS - 19 DRUG)	Madhya Pradesh	E	10900	Tue & Fri - 11:30 am	next day 7 pm	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D8268	ALLERGY DRUG PANEL - (ANESTHETIC - 4 DRUG)	Madhya Pradesh	E	4600	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
D8269	ALLERGY DRUG PANEL-(PAIN KILLER - 5 DRUG)	Madhya Pradesh	E	5400	Tue & Fri - 11:30 am	next day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
DL001	Amphetamine, Qualitative	Madhya Pradesh	B	450	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL002	Barbiturates, Qualitative	Madhya Pradesh	B	750	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL003	Benzodiazepine Diazepam	Madhya Pradesh	B	630	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL005	Cannabis Marijuana, THC, Qualitative	Madhya Pradesh	B	630	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL006	Cocaine Benzoyllecgonine, Qualitative	Madhya Pradesh	B	630	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL013	Dengue profile-	Madhya Pradesh	C	1250	Daily: 9am to 9pm	After 6 hrs	
DL015	Methadone Dolophine, Qualitative	Madhya Pradesh	B	450	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL016	Opiates Morphine, Qualitative, Urine Spot	Madhya Pradesh	B	750	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL017	PCP Phencyclidine Phosphate, Qualitative, Urine Spot_DL017	Madhya Pradesh	B	450	Daily : 9:00am to 9:00pm	After 10 mins	Urine
DL019	Troponin I (Qualitative)	Madhya Pradesh	C	1120	Daily: 9am to 9pm	After 4 hrs	
E0001	E Coli K1 / N. Meningitidis B Antigen detection by Agglutination - CSF	Madhya Pradesh	C	1745	Daily : 9:00am to 3:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
E0002	E Coli K1 / N. Meningitidis B Antigen detection by Agglutination - Serum	Madhya Pradesh	C	1715	Daily : 9:00am to 3:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0003	E Coli K1 / N. Meningitidis B Antigen detection by Agglutination - Urine	Madhya Pradesh	C	1345	Daily : 9:00am to 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
E0004	E2 Estradiol Serum	Madhya Pradesh	B	690	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0005	E3, unconjugated Estriol Serum	Madhya Pradesh	C	1570	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0006	EBV DNA Detection by Real Time PCR - EDTA Blood	Madhya Pradesh	D	4280	Mon, Thursday: 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
E0007	EBV DNA Detection by Real Time PCR - CSF	Madhya Pradesh	D	4280	Mon, Thursday: 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
E0008	EBV DNA Detection by Real Time PCR - Tissue	Madhya Pradesh	D	4280	Mon, Thursday: 9:00am	Next day : 5:00pm	Tissue with Normal Saline in Sterile leak proof container
E0009	EBV EA IgG Antibody by CLIA - Serum	Madhya Pradesh	C	2130	Mon : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0010	EBV NA IgG Antibody by CLIA - Serum	Madhya Pradesh	C	2130	Mon : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0011	EBV VCA IgG Antibody by CLIA - Serum	Madhya Pradesh	C	2130	Mon : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0012	EBV NA IgM Antibody by EIA - Serum	Madhya Pradesh	C	2130	Tuesday/Friday 9 am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0013	EBV VCA IgM Antibody by CLIA - Serum	Madhya Pradesh	C	1880	Daily : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0014	EBV encoded RNA (EBER) by CISH - FFPE	Madhya Pradesh	D	3425	Daily : 7:30pm	Same day : 5pm	Tissue In 10% Formalin Or Block / FFPE .
E0016	Echinococcus IgG Antibody by EIA - Serum	Madhya Pradesh	C	1815	Tue & Fri: 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
E0021	Electrolytes Serum	Madhya Pradesh	B	450	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0022	Electrolytes Urine 24H	Madhya Pradesh	B	655	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
E0023	Electrolytes Urine Spot	Madhya Pradesh	B	655	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
E0024	Electrolytes Dialysis Fluid	Madhya Pradesh	B	720	Daily: 9am to 9pm	After 6 hrs	2 ml of Dialysis Fluid
E0027	Endomysial Antibody IgA antibody Serum	Madhya Pradesh	C	2350	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0028	Endomysial Antibody IgG antibody Serum	Madhya Pradesh	C	2800	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0030	Enteroviruses RNA Detection by PCR - CSF**	Madhya Pradesh	R	5915	Daily : 11:30am	5th day : 6:00pm	2 ml of CSF in Sterile leak proof tube or container
E0031	Enteroviruses RNA Detection by PCR - Plasma**	Madhya Pradesh	R	6030	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA Plasma in Sterile Plain Vacutainer
E0032	Eosinophil Detection, URINE	Madhya Pradesh	B	360	Daily: 9am to 9pm	After 8 hrs	5 ml urine sample
E0033	Eosinophil detection (SPUTUM)	Madhya Pradesh	B	360	Daily: 9am to 9pm	After 8 hrs	Sputum In Sterile Container
E0035	Erythropoietin Serum	Madhya Pradesh	C	2350	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
E0038	ESR Automated Blood	Madhya Pradesh	A	110	Daily: 9am to 9pm	After 6 hrs	3 ml EDTA Whole Blood
E0039	Everolimus, Blood by LCMS method	Madhya Pradesh	D	5670	Daily: 7.30am	Same Day 7.30pm	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
E0042	Estrone, Serum	Madhya Pradesh	R	7140	Daily : 11:30am	After 4 weeks : 6:00pm	Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample. 5 ml Serum
E0043	EGFR MUTATION DETECTION (EXONS 18, 19, 20, 21)	Madhya Pradesh	D	8560	Mon, Wed & Fri: 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>10% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
E0045	EWS gene rearrangement (:Translocation t(11;22) /t(21;22) translocation/ rearrangement **	Madhya Pradesh	R	7955	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
E0051	EBV DNA Viral load by Real Time PCR - EDTA Blood	Madhya Pradesh	D	6740	Mon, Thursday: 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
E0053	Extended RAS Test (KRAS, BRAF, NRAS, PIK3CA) - FFPE Tissue	Madhya Pradesh	D	16050	Mon, Wed & Fri: 9:00am	10th day : 6:00pm	Tumor containing FFPE block. Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
E0055	Epileptic encephalopathy gene panel**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
E0056	EGFR T790M Mutation (Acquired EGFR TKI Resistance)*	Madhya Pradesh	E	8000	Daily : 9:00am	7th day : 6:00pm	refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. Tumor containing FFPE block (>10% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
E0057	EXPANDED LUNG CANCER PANEL BY FISH	Madhya Pradesh	E	24000	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
E0064	Allergy panel -Common Eczema	Madhya Pradesh	D	7760	Daily : 11:30am	Next day : 5:00pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
E0065	Allergy panel - Eczema Comprehensive Veg	Madhya Pradesh	D	17655	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
E0066	Allergy panel - Eczema Comprehensive Veg + non Veg	Madhya Pradesh	D	21775	Daily : 11:30am	Next day : 5:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							for 7 days and at Ambient temperature for a day.
E0068	ERA-Endometrium Receptivity Assay, Endometrial biopsy**	Madhya Pradesh	R	42840	Daily : 11:30am	15th Day	CytogeneticsandFISH.pdf
E0069	EWSR1-FLI1 **	Madhya Pradesh	R	12240	Daily : 11:30am	15th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
E0070	EWSR1(22q12)/NTRK3 (15q25) translocation **	Madhya Pradesh	R	12240	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
E0071	EBV DNA Viral load by Real Time PCR - CSF	Madhya Pradesh	D	6740	Mon & Thurs : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
E0072	EBV DNA Viral load by Real Time PCR - Tissue	Madhya Pradesh	D	6740	Mon & Thurs : 9:00am	Next day : 5:00pm	Tissue with Normal Saline in Sterile leak proof container
E0075	Endometritis Panel	Madhya Pradesh	D	7705	Wed : 9:00am	Next day : 5:00pm	Genital discharge, Ulcer, Ulcer base, High Vaginal swab, Endocervical swab, Endometrial Pus, Endometrial Fluid, Endometrial Tissue.
E1019	E Histolytica Antigen Detection by Immunochromatography -Stool **	Madhya Pradesh	R	2905	Daily : 11:30am	3rd day : 6:00pm	1gm of Stool specimen in Sterile leak proof container
E1020	Entamoeba histolytica detection by Wet Mount - Stool	Madhya Pradesh	B	380	Daily : 9:00am to 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
E1021	Kidney Biopsy by Electron Microscopy	Madhya Pradesh	R	9690		15 day	Biopsy in 3% Gluteraldehyde along with histopathology report & relevent clinical details
E1902	E2 estradiol, Serum	Madhya Pradesh	B	710	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
E1903	ESR Automated, EDTA whole blood	Madhya Pradesh	A	110	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
E8141	Endometrial Cancer panel by NGS	Madhya Pradesh	R	45900	Daily 7.30 am	4 Weeks	Paraffin Block
F0002	Factor V Leiden-Mutant Detection (G1691A Mutation) by Real Time PCR - Blood	Madhya Pradesh	D	5670	Wed & Sat : 9:00am	Next day : 5:00pm	3 ml of EDTA Whole Blood in EDTA Vacutainer
F0003	Factor V Level	Madhya Pradesh	D	5885	Fourth Mon: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
F0004	Factor VII Activity Proconvertin Citrated plasma	Madhya Pradesh	R	3980	Daily : 11:30am	10th day : 6:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0005	Factor VIII Activity Anti Haemophilia factor A, Citrated plasma	Madhya Pradesh	D	2355	Mon: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
F0006	Factor IX Activity Anti Haemophilic factor B (Christmas Disease), Citrated plasma	Madhya Pradesh	D	2035	Mon : 9:00am	Same day : 8pm	1 ml of citrated platelet poor plasma in frozen condition.
F0007	Factor X Activity Stuart Power Citrated plasma	Madhya Pradesh	R	4795	Daily : 11:30am	10th day : 6:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0008	Anti Factor X A Assay**	Madhya Pradesh	R	5405	Daily : 11:30am	4th day : 6:00pm	Sample 3 ml Citrate Plasma Frozen with medication history.
F0009	Factor XI Activity Rosenthal Citrated plasma	Madhya Pradesh	R	3875	Daily : 11:30am	10th day : 6:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0010	Factor XII Activity Hageman Citrated plasma	Madhya Pradesh	R	3960	Daily : 11:30am	5th day : 6:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0011	Factor XIII Activity, Citrated plasma	Madhya Pradesh	R	2580	Daily : 11:30am	Next day : 5:00pm	1 ml of citrated platelet poor plasma in frozen condition.
F0012	Fatty Acids With Very Long Chains C22 To	Madhya	R	15300	Daily : 11:30am	After 4 weeks :	3 ml of EDTA Plasma Or Heparin

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	C26 Plasma**	Pradesh				6:00pm	plasma
F0014	FDP DR 70	Madhya Pradesh	R	3670	Daily : 11:30am	7th day : 6:00pm	3 ml of Serum and 3 ml of citrate plasma
F0015	FDP Fibrin degradation products Reflex to end point titre.	Madhya Pradesh	C	1345	Daily: 9am and 3pm and 9pm	After 6 hrs	1 ml of citrated platelet poor plasma in frozen condition.
F0016	Febrile by Agglutination - Serum	Madhya Pradesh	B	1495	Daily : 7:30pm	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
F0017	Ferric Chloride, Urine	Madhya Pradesh	C	390	Daily : 11:30am	Next day : 5:00pm	10 ml of Spot urine
F0018	Ferritin Serum	Madhya Pradesh	B	500	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0019	Fertility Endocrine Panel-Female	Madhya Pradesh	B	2440	Daily: 9am to 9pm	After 8 hrs	5 ml of serum
F0020	Fertility Endocrine Panel-Male	Madhya Pradesh	B	3325	Tue and Thu and Sat: 7:30am	After 8 hrs	5 ml of serum
F0021	Fever Profile (Maxi)	Madhya Pradesh	C	2910	Daily: 9am to 9pm	Same Day	1.4 ml of EDTA whole blood 3 ml of serum 10 ml of spot urine
F0022	Fever Profile (Mini)	Madhya Pradesh	C	860	Daily: 9am to 9pm	Report Same Day: 1. Complete Blood Count(CBC Haemogram) 2. ESR - Erythrocyte Sedimentation Rate 3. Malaria Parasite Detection 4. Routine Examination Urine 4 ml of EDTA whole blood 3 ml of serum 10 ml of spot urine F5604 Fever Profile (Super Maxi) Madhya Pradesh C 4705 Daily : 9:00am	4ml of serum 3 ml of EDTA whole blood 20 Ml of spot urine Collect blood 2-3 ml in Bactec FX aerobic bottle

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
						Report Same Day: 1. Complete Blood Count(CBC Haemogram) 2. Routine Examination Urine 3. CRP - C Reactive Protein 4. Malaria Parasite Detection 5. SGPT 6. Typhi Dot IgM 7. Dengue	
F0023	Fibrinogen Plasma (F0023)	Madhya Pradesh	C	1065	Daily: 9am to 9pm	After 6 hrs	1 ml of citrated platelet poor plasma in frozen condition.
F0028	Fibromax Fibrosis, Activity, STEATO, ASH and NASH of Liver.	Madhya Pradesh	E	13900	Daily: 9am to 9pm	5th day	5 ml of serum and 2 ml of fluoride plasma . Fasting sample
F0029	Filaria Antibody detection by Immunochromatography - Serum	Madhya Pradesh	C	1120	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
F0030	Filaria Antigen detection by Immunochromatography - EDTA blood	Madhya Pradesh	C	1120	Daily : 9:00am to 9:00pm	After 6 hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer
F0031	Filaria Profile	Madhya Pradesh	B	1310	Daily: 9am and 3pm	After 6 hrs	3 ml EDTA whole blood and serum sample
F0032	FLT-3 (AML) by PCR ITD, D835 - Blood	Madhya Pradesh	D	5885	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.
F0033	FLT-3 (AML) by PCR ITD, D835 - Bone Marrow	Madhya Pradesh	D	5885	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.
F0034	FMC 7 by Flow cytometry	Madhya Pradesh	D	2620	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							condition(preferably). Clinical history and clinician number required.
F0039	By Conventional method	Madhya Pradesh	C	1290		3rd day	Prefixed smears along with air dry smears, relevant clinical & radiological details mandatory
F0040	By LBC	Madhya Pradesh	C	1570		3rd day	More than 1 MI fluid in a leak proof plain container along with relvent clinical & radiological details
F0043	Procedure charge	Madhya Pradesh	C	560			NA
F0046	Folic Acid Serum	Madhya Pradesh	C	1230			3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0048	Food intolerance Test	Madhya Pradesh	D	12000	Daily : 11:30am	3rd day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0050	Fragile X (FMR1)CGG Repeat Analysis	Madhya Pradesh	D	7385	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
F0051	Friedreich Ataxia Mutation Analysis	Madhya Pradesh	D	5885	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
F0052	Free Androgen Index FAI	Madhya Pradesh	D	2835	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0053	Free Beta HCG Serum	Madhya Pradesh	C	1570	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0054	Free Kappa Freelite chains Serum	Madhya	D	2780	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					ml serum for profiles and reflex tests.
F0055	Free Lambda Freelite chains Serum	Madhya Pradesh	D	2780	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0056	Fructosamine Serum	Madhya Pradesh	C	1065	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0057	Fructose biochemical by Biochemical - Semen	Madhya Pradesh	B	220	Daily : 9:00am to 9:00pm	After 4 hrs	2ml of Semen in Sterile leak proof container
F0058	FSH-Follicle Stimulating Hormone ECLIA Serum	Madhya Pradesh	B	540	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0059	FSH-LH-Prolactin Serum	Madhya Pradesh	C	1600	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0060	FSH-LH-Prolactin-Testosterone	Madhya Pradesh	C	1800	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0061	FSH-LH-Prolactin-TSH	Madhya Pradesh	C	1680	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0062	FSH-LH-Testosterone	Madhya Pradesh	C	1455	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
F0065	Fungal Antibigram (Yeast) MIC Panel - Pure culture	Madhya Pradesh	C	2800	Daily : 9:00am to 9:00pm	5th day : 5pm	Pure culture Growth on Agar slant/ Plate
F0066	Fungal Identification - Pure culture	Madhya Pradesh	C	2240	Daily : 9:00am	10th day : 6pm	Pure culture on SDA plate/ slant
F0067	Fungal stain by KOH & Calcofluor White stain - BAL, ET Secretion Respiratory Sample	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of BAL / Other Respiratory samples in Sterile leak proof container
F0068	Fungal stain by KOH & Calcofluor White stain- Specify type Body fluid	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	3ml of Fluid in Sterile leak proof container
F0069	Fungal stain by KOH & Calcofluor White stain - CSF	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
F0070	Fungal stain by KOH & Calcofluor White	Madhya Pradesh	B	545	Daily : 9:00am &	After 6 hrs	Ear specimens should be collected

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	stain - Ear Sample	Pradesh			3:00pm		by a medical specialist, an experienced technologist or nurse. Aspirate discharge and collect in a sterile container or collect with a sterile swab. Transport at the earliest. Amies transport medium can be used for transportation.
F0071	Fungal stain by KOH & Calcofluor White stain- Eye Sample	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	Eye specimens in Sterile leak proof container
F0072	Fungal stain by KOH & Calcofluor White stain - Hair	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	Hair sample in Sterile leak proof container
F0073	Fungal stain by KOH & Calcofluor White stain- Nail	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	Nail sample in Sterile leak proof container
F0074	Fungal stain by KOH & Calcofluor White stain- Other samples	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Vaginal swab/ Vaginal discharge in Sterile leak proof container.
F0075	Fungal stain by KOH & Calcofluor White stain- Pus	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of PUS in Sterile leak proof container
F0076	Fungal stain by KOH & Calcofluor White stain- Skin	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	Skin sample in Sterile leak proof container
F0077	Fungal stain by KOH & Calcofluor White stain- Sputum	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	2ml of Sputum in Sterile leak proof container
F0078	Fungal stain by KOH & Calcofluor White stain- Tissue	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	Tissue with Normal Saline in Sterile leak proof container
F0079	Fungal stain by KOH & Calcofluor White stain- Urine	Madhya Pradesh	B	545	Daily : 9:00am & 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
F0080	Fungal culture and identification by MALDITOF (Vitek MS) - BAL, E.T Secretion Respiratory Sample	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	2ml of BAL / Other Respiratory samples in Sterile leak proof container
F0081	Fungal culture and identification by MALDITOF (Vitek MS) - Blood	Madhya Pradesh	C	2800	Daily : 9:00am to 9:00pm	10th day : 6pm	3ml of Blood in MycoF/lytic Bactec Bottle or Sodium Heparin
F0082	Fungal culture and identification by MALDITOF (Vitek MS) - Body Fluid	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	2ml of Body Fluid in Sterile leak proof container
F0083	Fungal culture and identification by MALDITOF (Vitek MS) - Bone Marrow	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	3ml of Bone Marrow in MycoF/lytic Bactec Bottle or Sodium Heparin
F0084	Fungal culture and identification by MALDITOF (Vitek MS) - CSF	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	2ml of CSF in Sterile leak proof tube or container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
F0085	Fungal culture and identification by MALDITOF (Vitek MS) - Ear Sample	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	Ear specimens in Sterile leak proof container
F0086	Fungal culture and identification by MALDITOF (Vitek MS) - Eye Sample	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	Eye specimens in Sterile leak proof container
F0087	Fungal culture and identification by MALDITOF (Vitek MS) - Hair	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	Hair sample in Sterile leak proof container
F0088	Fungal culture and identification by MALDITOF (Vitek MS) - Nail	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	Nail sample in Sterile leak proof container
F0089	Fungal culture and identification by MALDITOF (Vitek MS) - Other Sample	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	2ml of Vaginal swab/ Vaginal discharge in Sterile leak proof container.
F0090	Fungal culture and identification by MALDITOF (Vitek MS) - Pus	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	2ml of PUS in Sterile leak proof container
F0091	Fungal culture and identification by MALDITOF (Vitek MS) - Skin	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	Skin sample in Sterile leak proof container
F0092	Fungal culture and identification by MALDITOF (Vitek MS) - Sputum	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	2ml of Sputum in Sterile leak proof container
F0093	Fungal culture and identification by MALDITOF (Vitek MS) - Stool	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	5gm of Stool specimen in Sterile leak proof container
F0094	Fungal culture and identification by MALDITOF (Vitek MS) - Tissue	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	Tissue with Normal Saline in Sterile leak proof container
F0095	Fungal culture and identification by MALDITOF (Vitek MS) - Urine	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	21st day : 5pm	15ml of Spot Urine in Sterile leak proof container
F0099	Fever Profile Detection by Real Time Multiplex PCR - EDTA Blood	Madhya Pradesh	D	5500	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
F0103	FIP1L1-PDGFR A GENE REARRANGEMENT **	Madhya Pradesh	R	11220	Tue : 11:30am	12th day : 6:00pm	4 ML EDTA whole blood / Bone marrow
F0105	Flu Profile-1 Detection by Real Time PCR - Nasopharangeal Swab	Madhya Pradesh	D	6955	Daily : 7:00am & 3:00pm	Same Day : 3:00pm & 11:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
F0106	Fibroblast Growth Factor 23, Plasma	Madhya Pradesh	D	11115	Wed 12:00 pm	same day 08:00 pm	3 ml of EDTA plasma to be sent in aliquot tube at proper temperature.
F0107	Flow for Lymphoma/ CLPD / lymphoma MRD**	Madhya Pradesh	R	14790	Monday to Thursday 09.00 am	3rd day	3 ml of EDTA Bone Marrow
F0109	Flow for B ALL MRD**	Madhya Pradesh	R	14790	Monday to Thursday 09.00 am	5th day	3 ml of EDTA Bone Marrow
F0110	Flow for T ALL MRD**	Madhya Pradesh	R	18870	Monday to	6th day	3 ml of EDTA Bone Marrow

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
F0112	Food Intolerance (Vegetarian panel)	Pradesh Madhya Pradesh	D	11500	Thursday 09.00 am Daily : 11:30am	3rd day : 6:00pm	Require serum (5ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
F0115	First trimester test-V, Serum Quad Test Free beta HCG, PAPPa, AFP and PLGF	Madhya Pradesh	D	5245	Daily cut off 7:30 AM	same day 9:00 PM	3 ml of serum is required. Mandatory Requirements: Latest Ultrasound (USG) report with a completely filled Maternal screening TRF form including CRL/BPD, nasal bone status, ductus venosus status, and NT value in millimeters. Please note that this test is not valid for triplets and above Preeclampsia: USG findings of mean arterial pressure (MAP) and uterine artery pulsatility index (UAPI), the patients blood pressure (BP), and height should be included as per maternal TRF requirements .
F0116	First trimester test-VI, serum Penta Screen Free beta HCG, PAPPa, AFP, Inhibin A, PLGF and NT 10 to 13.6 TBA	Madhya Pradesh	R	6120	Daily : 11:30am	10th day : 6:00pm	3 ml of serum is required. Mandatory Requirements: Latest Ultrasound (USG) report with a completely filled Maternal screening TRF form including CRL/BPD, nasal bone status, ductus venosus status, and NT value in millimeters. Please note that this test is not valid for triplets and above Preeclampsia: USG findings of mean arterial pressure (MAP) and uterine artery pulsatility index (UAPI), the patients blood pressure (BP), and height should be included as per maternal TRF requirements

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
F0117	Fragile X mutation analysis+Chromosomal Array profile	Madhya Pradesh	E	21000	Tue & Fri : 9:00am	7th day : 6:00pm	3 ml of Blood In Na-Heparin +3 ml of EDTA Vacutainer
F0118	FOXO1 (13q14) translocation **	Madhya Pradesh	R	16320	Daily : 11:30am	12th Day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
F0119	COVI-INFLUENZA VIRUS PANEL BY REAL TIME PCR (Inf A, Inf B, H1N1(Swine Flu), H3N2, RSV, Covid-19)	Madhya Pradesh	D	5885	Daily : 7:00am & 3:00pm	Same Day : 3:00pm & 11:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
F0120	Mould/Filamentous fungi by microbroth dilution - Agar plate or slant	Madhya Pradesh	C	9520	Daily : 9:00am to 9:00pm	15th day : 5pm	Pure culture Growth of Filamentous fungi on Agar slant/ Plate.
F0122	Factor VIII inhibitor assay**	Madhya Pradesh	D	6955	Daily : 11:30am	3rd day : 6:00pm	3 ml Citrate Plasma (Fasting & Frozen mandatory)
F0126	Foundation one CDX	Madhya Pradesh	R	300900	Daily 7.30 am	5 weeks	Paraffin Block (Carcinoma)
F1003	Free Fatty Acid (Non-Esterified)	Madhya Pradesh	D	13375	Daily 7.30 am	10th Days	3 ml Serum
F1901	Second opinion services, Slides	Madhya Pradesh	C	1120	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1902	Lymphnode	Madhya Pradesh	C	1230	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1903	CANINE-FNAC-Others, conventional	Madhya	C	1065	Daily: 9am	Report after 48 hrs	Please contact lab before sending

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1904	CANINE-FNAC-Thyroid, conventional	Madhya Pradesh	C	1065	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1905	Procedure Charges by appointment	Madhya Pradesh	C	1120	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F1906	Fungal stain	Madhya Pradesh	C	560	Daily: 9am	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
F5605	Fatty Liver	Madhya Pradesh	D	6170	Daily: 9am to 9pm	After 6hrs	3 ml of EDTA whole blood, 2 ml of Flouride plasma, 3 ml of serum
F5606	Fever Profile (Super Mini)	Madhya Pradesh	C	2090	Daily : 9:00am	Same Day	4 ml of EDTA whole blood 3 ml of serum 10 ml of spot urine
F8194	Foundation one Heme**	Madhya Pradesh	R	362100	Daily 7.30 am	5 weeks	Paraffin Block (Sarcoma)
F8215	FIT-Fecal Immunochemical Test, quantitative	Madhya Pradesh	D	4200	Daily : 9:00AM	Same day : 7:00 PM	100 gms of stool is required

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0001	G6PD Neonatal Screen Dried blood spot	Madhya Pradesh	C	505	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
G0002	G6PD Qualitative Blood	Madhya Pradesh	B	655	Daily: 9am to 3pm	After 8 hrs	3 ml EDTA Whole Blood
G0003	G6PD Quantitative Reflex Qualitative	Madhya Pradesh	C	1200	Daily: 9am to 9pm	After 8 hrs	3 ml EDTA Whole Blood
G0004	GAD Total Antibody Type I Diabetes Serum	Madhya Pradesh	D	5500	Tue,Thu & Sat: 1 pm	Same day 7pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0005	Galactose 6 Sulphate Sulphatase Blood	Madhya Pradesh	R	6630	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
G0006	Galactose Classical Transferase Blood**	Madhya Pradesh	R	2415	Daily : 11:30am	8th day : 6:00pm	6 ml Sodium heparin whole blood required as soon as possible after collection
G0008	Galactose screening,Blood **	Madhya	R	980	Daily : 11:30am	7th day : 6:00pm	4 ml of Blood In Sodium Heparin

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0011	Galactose-Total Neonatal Screen Dried blood spot	Pradesh Madhya Pradesh	C	505	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
G0012	Ganglioside IgG Antibody,serum	Madhya Pradesh	D	4065	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0013	Ganglioside IgM Antibody,serum	Madhya Pradesh	D	4065	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0016	Gastrin Serum	Madhya Pradesh	C	1570	Daily: 9am to 9pm	After 6 hrs	2 ml of Serum.Fasting Required
G0017	GBM-Antibody Glomerular Basement Membrane serum	Madhya Pradesh	C	1790	Wed,Sat : 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0018	Genetic Hair Loss Test Alopecia Blood	Madhya Pradesh	E	4500	Daily : 7:30am	15th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0019_A	GFR With Creatinine (Adult - Above 18 Years) Serum	Madhya Pradesh	B	340	Daily: 9am to 9pm	After 6 hrs	condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0019_C	GFR With Creatinine (Upto 18 Years) Serum	Madhya Pradesh	B	350	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0020	GGTP Gamma GT Serum	Madhya Pradesh	A	295	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0021	Giardia Antigen detection by Immunochromatography - Stool	Madhya Pradesh	B	1255	Daily : 9:00am & 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
G0022	Gliadin-IgA Antibodies De-amidated (DGP)Serum	Madhya Pradesh	C	2130	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0023	Gliadin-IgG Antibodies De-amidated (DGP)Serum	Madhya Pradesh	C	2130	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0024	Glucagon Stimulation Test - For C-peptide, Serum	Madhya Pradesh	D	3745	Daily: 9am to 9pm	After 8 hrs	3 ml of serum (0min, 5, 10, 15 min)
G0025	Glucosamine N Acetyl Transferase Blood	Madhya Pradesh	R	7650	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
G0026	Glucose Urine spot	Madhya Pradesh	A	135	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
G0027	Fasting Glucose, plasma	Madhya Pradesh	A	60	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0027_100	Glucose Post glucose (100 GMS) Plasma	Madhya Pradesh	A	65	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0027_120	GLUCOSE 120 MIN	Madhya Pradesh	A	65	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0027_180	GLUCOSE 180 MIN	Madhya Pradesh	A	65	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0027_240	GLUCOSE 240 MIN	Madhya Pradesh	A	65	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0027_300	GLUCOSE 300 MIN	Madhya Pradesh	A	65	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0027_50	Glucose Post glucose (50 GMS) Plasma	Madhya Pradesh	A	65	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0027_60	GLUCOSE 60 MIN	Madhya Pradesh	A	65	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0027_75	Glucose plasma(75 GMS)	Madhya Pradesh	A	65	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0027_P	Glucose post prandial plasma	Madhya Pradesh	A	60	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0027_R	Random, Glucose	Madhya Pradesh	A	60	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma
G0028	Glucose CSF	Madhya Pradesh	A	135	Daily: 9am to 9pm	After 6 hrs	2 ml of Body Fluid
G0028_B	Glucose Body Fluid	Madhya Pradesh	A	135	Daily: 9am to 9pm	Glucose Body Fluid	2 ml of Body Fluid
G0028_S	Glucose Synovial Fluid	Madhya Pradesh	A	135	Daily: 9am to 9pm	After 6 hrs	2 ml of Body Fluid
G0029	Glucose Tolerance Test 2 hours 3 sample	Madhya Pradesh	B	220	Daily: 9am to 9pm	After 6 hrs	0 min, 60 , 120 min -2ml Fluoride plasma
G0030	Glucose Tolerance Test 2.5 hours 5 samples	Madhya Pradesh	B	350	Daily: 9am to 9pm	After 6 hrs	0 min, 30, 60, 120, 150 min - 2ml Fluoride plasma
G0031	Glucose Tolerance Test 5 hours 8 samples	Madhya Pradesh	B	600	Daily: 9am to 9pm	After 6 hrs	0 min, 30, 60, 90, 120, 180, 240, 300 min - 2ml Fluoride plasma
G0032	Glucose Tolerance Test 3 samples for glucose and insulin	Madhya Pradesh	C	2240	Daily: 9am to 9pm	After 6 hrs	0 min, 60 , 120 min -2ml Fluoride plasma and 3ml serum
G0033	Glutamic Acid Quantitative Urine 24H**	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
G0034	Glutamic Acid Quantitative Plasma*	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
G0035	Glutamine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					plasma to be sent in aliquot tube at proper temperature.
G0036	Glutamine Quantitative Urine 24H**	Madhya Pradesh	R	6325	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
G0037	Glycine Quantitative CSF	Madhya Pradesh	D	6635	Daily : 11:30am	7th day : 6:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
G0038	Glycine Quantitative Urine 24H**	Madhya Pradesh	R	6325	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
G0039	Glycine Quantitative Plasma*	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
G0041	Glycophorin-A	Madhya Pradesh	D	2460	Daily: 9am	Next day 5pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
G0042	GNRH (LHRH) Stimulation Test For FSH and LH	Madhya Pradesh	B	3925	Daily: 9am to 9pm	Next day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0044	Goodpasture`s Syndrome Panel	Madhya Pradesh	D	3480	Fri 7:30 am	Next day 5 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0045	Gram stain For Neisseria sp by Microscopy - Urethral discharge	Madhya Pradesh	B	325	Daily : 9:00am to 3:00pm	After 6 hrs	2 ml of Urethral Discharge in Sterile leak proof container
G0046	Gram stain Nugent s score by Microscopy – Vaginal discharge	Madhya Pradesh	B	435	Daily : 9:00am to 3:00pm	After 6 hrs	Discharge on sterile swab in Amies transport medium
G0047	Gram stain by Microscopy - Other samples	Madhya Pradesh	B	325	Daily : 9:00am to 3:00pm	After 6 hrs	2 ml of Specimen in Sterile leak proof container
G0048	Growth Disorder Panel	Madhya Pradesh	D	5780	Daily: 9am to 9pm	After 8 hrs	5 ml of serum
G0049	Growth Hormone HGH Serum	Madhya Pradesh	C	895	Daily: 9am to 9pm	After 6 hrs	2 ml of Serum In Fasting And Resting
G0050	Growth Hormone Stimulation Test- Arginine, Serum	Madhya Pradesh	D	3530	Daily: 9am to 9pm	After 6 hrs	3 ml of serum (0min, 30, 60, 90, 120 min)

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0051	Growth Hormone Stimulation Test - Clonidine, Serum	Madhya Pradesh	D	3530	Daily: 9am to 9pm	2nd day	3 ml of serum (0min, 30, 60, 90, 120 min)
G0052	Growth Hormone Stimulation Test-Exercise, Serum	Madhya Pradesh	D	3530	Daily: 9am to 9pm	After 8 hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0053	Growth Hormone Stimulation Test - GHRH, Serum	Madhya Pradesh	D	3530	Daily: 9am to 9pm	After 8 hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0054	Growth Hormone Stimulation Test-Insulin	Madhya Pradesh	D	3530	Daily: 9am to 9pm	After 8 hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0055	Growth hormone stimulation test for L Dopa, Serum	Madhya Pradesh	D	3530	Daily: 9am to 9pm	After 8 hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0056	Growth Hormone Suppression Test-Glucose, Serum	Madhya Pradesh	D	3530	Daily: 9am to 9pm	After 8 hrs	3 ml of serum (0min, 30, 60, 90, 120 min)
G0058	Anti GAD-65 Paraneoplastic Syndrome Serum	Madhya Pradesh	D	6500	Wed: 1:00pm	Next day 5:00 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0059	GFR With Cystatin C Serum	Madhya Pradesh	C	1345	Daily : 11:30am	After 6hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0059_A	GFR With Cystatin C Serum (Adult - Above 18 Years)	Madhya Pradesh	C	1345	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0059_C	GFR With Cystatin C Serum (Upto 18 Years)	Madhya Pradesh	C	1345	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0061	Glycosaminoglycan,Urine	Madhya Pradesh	R	3160	Daily : 11:30am	15th day : 6:00pm	20 ml Spot Urine Sample In Leak Proof Container
G0065	GTT-Gestational ,Two Hour(75 grams glucose)	Madhya Pradesh	B	275	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma Each
G0066	GTT-Gestational ,Three Hour(100 grams glucose)	Madhya Pradesh	B	390	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma Each
G0073	G6pd Gene Mutation**	Madhya Pradesh	R	4490	Tue : 11:30am	12th day : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							at Ambient or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0074	GALACTOSEMIA GENE MUTATION**	Madhya Pradesh	R	8670	Tue : 11:30am	12th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0075	Gardenella Vaginalis - Vaginal Discharge/ cervical secretions**	Madhya Pradesh	R	1325	Daily : 9:00am to 9:00pm	7th day : 6pm	Collect vaginal discharge on a sterile cotton swab and place in into Amies transport medium with aseptic precautions. 2. Make vaginal discharge smears on clean glass slide for Gram staining.
G0078	Councelling Genetics	Madhya Pradesh	E	1200	Daily : 9:00am	Next day : 5:00pm	Prior Appointment
G0079	GIST PDGFRA gene mutation**	Madhya Pradesh	R	10200	Thu & Sat: 11:30am	15th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
G0080	Glucose Tolerance Test 2 hours 2 samples	Madhya Pradesh	B	155	Daily: 9am to 9pm	After 6 hrs	0 min, 120 min -2ml Fluoride plasma
G0082	GLIOMA PANEL (1P/19Q, IDH1/2,MGMT)**	Madhya Pradesh	D	26750		7th Day	(FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							hours and paraffin embedded. Clinical Diagnosis with Histopathology report, HistologyandIHC.pdf
G0083	Gene analysis Alpha-1 antitrypsin deficiency - SERPINA1**	Madhya Pradesh	R	21420	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0087	Gene analysis Achondroplasia FGFR3 G380R common mutation screening**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0088	Agammaglobulinemia (BTK)-Gene analysis, EDTA Blood	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0089	Gene analysis Alexander disease (GFAP)**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	condition. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0091	Alstrom syndrome (ALMS1) gene analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0092	Androgen receptor (AR)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0093	Gene analysis Ataxia-telangiectasia (ATM)**	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0094	Gene analysis ATRX**	Madhya Pradesh	E	19000	Daily: 9am	4 weeks	Genetci history is mandatory. 3 ml EDTA Whole Blood
G0095	Brown Vialetto-Van Laere syndrome-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0096	C8orf37 gene-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0097	Canavan disease (ASPA)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0098	Cantu syndrome (ABCC9) gene-Gene Analysis	Madhya Pradesh	R	22440	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0099	Carney complex (PRKAR1A)-Gene Analysis	Madhya Pradesh	R	17340	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0100	Chediak-Higashi syndrome (LYST)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0101	Cohen s syndrome (VPS13B)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0102	Congenital adrenal hyperplasia CYP21A2(21-0H)**	Madhya Pradesh	R	12240	Daily : 11:30am	After 6 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0103	Gene analysis Congenital adrenal hyperplasia CYP21A2 (21-OH) by NGS**	Madhya Pradesh	R	21420	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0104	Gene analysis Cystic fibrosis (CFTR)-Whole gene**	Madhya Pradesh	R	19260	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0105	Gene analysis Cystic megalencephaly (MLC1)**	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0106	Gene analysis DOCK8 (Hyper IgE Syndrome) gene**	Madhya Pradesh	R	17850	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0108	Dravet syndrome (SCN1A) gene panel,NGS**	Madhya Pradesh	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0109	Gene analysis duchenne muscular dystrophy (DMD)**	Madhya Pradesh	R	17000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0110	Fabry disease (GLA)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0111	Factor VII deficiency-Gene analysis, EDTA Blood	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0112	Fanconi bickel syndrome (SLC2A2)-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	at Amibent or referigerated condition based on sample stability. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0113	Gene analysis(Farber lipogranulomatosis - ASAH1 -**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0114	Fraser syndrome genes-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0115	Fucosidosis - Lysosomal storage disorder- Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0116	Galactosialidosis-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0117	GLB1 (GM1 gangliosidosis; mucopolysaccharidosis type IVB) gene analysis**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0118	Gene analysis GLUT1 deficiency (SLC2A1)**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0119	Haemophilia A-Gene analysis, EDTA Blood	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0120	Gene analysis Haemophilia B (F9)**	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	at Amibent or referigerated condition based on sample stability. Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.
G0121	Hereditary fructose intolerance-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0122	Hyperargininemia (ARG1)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0123	Krabbe disease (GALC)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0124	Lysosomal acid lipase deficiency (LIPA gene)-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0126	McArdle disease (PYGM)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0127	Menkes disease (ATP7A)-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0128	Mondini defect (SLC26A4)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0129	Gene analysis MYD88 gene L265P mutation**	Madhya Pradesh	R	18360	Daily : 11:30am	12th Day : 6:00pm	3 ml EDTA Whole Blood
G0130	Neurofibromatosis type 1 (NF1) Gene Analysis **	Madhya Pradesh	R	17340	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0131	Gene analysis Neurofibromatosis type 2 (NF2)**	Madhya Pradesh	R	17340	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0132	NOTCH3 (CADASIL) gene analysis -NGS**	Madhya Pradesh	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0133	Pallister Hall syndrome (GLI3)-Gene Analysis	Madhya Pradesh	R	17340	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0134	Phenylketonuria (PAH)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0135	PMP22 Gene analysis**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0136	Polycystic Liver Disease-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0137	Gene analysis Pompe Disease (GAA)**	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0138	PTEN-Gene Analysis	Madhya Pradesh	R	18360	Daily : 11:30am	After 4 weeks : 6:00pm	Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0139	Pulmonary alveolar microlithiasis (SLC34A2)-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0142	Gene analysis RB1*	Madhya Pradesh	R	17340	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0143	Gene analysis RBC membrane disorders and Enzymopathies**	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability. Amniotic Fluid or Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition.
G0144	Sialidosis-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0145	Sjogren-Larsson syndrome (ALDH3A2)-Gene Analysis	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0147	Gene analysis Smith-Lemli-Opitz syndrome (DHCR7)**	Madhya Pradesh	R	17340	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0149	Tay-Sachs disease (HEXA)-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0150	Thrombotic Thrombocytopenic Purpura (TTP) - ADAMTS13 gene analysis**	Madhya Pradesh	R	18900	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.
G0151	TSC1-Gene Analysis	Madhya Pradesh	R	17340	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0152	Gene analysis TSC1 & TSC2**	Madhya Pradesh	R	17340	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0153	TSC2-Gene Analysis	Madhya Pradesh	R	17340	Daily : 11:30am	After 4 weeks : 6:00pm	ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0154	Tyrosinemia-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0155	Gene analysis Von Hippel-Lindau syndrome (VHL)**	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.
G0156	Gene analysis Von Willebrand disease** (VWF)	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0157	WAS (Wiskott Aldrich Syndrome)-Gene analysis, EDTA Blood	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	condition based on sample stability. Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability.
G0158	Wilson disease (ATP7B) gene analysis,NGS**	Madhya Pradesh	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0159	Xanthinuria (XDH)-Gene Analysis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0162	4H syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0163	Aicardi-Goutieres syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0164	Alagille syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0165	Alport syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0166	Gene Panel Aplastic anaemia**	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0167	Arthrogryposis & congenital myasthenic syndrome Gene Panel**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0168	Bardet-Biedl syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0169	Gene Panel Bartter syndrome**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0170	Benign infantile epilepsy-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0171	Cardiac channelopathy-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0173	Gene Panel Charcot-Marie-Tooth and sensory neuropathies**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0174	Ciliopathy-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0175	Citrullinemia-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0176	Cockayne syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0177	Congenital afibrinogenemia-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0178	Congenital dyserythropoietic anemia gene panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0180	Congenital hypopituitarism-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0181	Gene Panel Congenital neutropenia	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0182	Congenital stationary night blindness-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0183	Gene Panel Connective tissue disorder**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0184	Corneal dystrophy-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0185	Cornelia de Lange Syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0186	Craniosynostosis-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0187	Cutis-laxa-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0188	Diamond blackfan anemia-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0189	Gene Panel Dyskeratosis congenita**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0190	Gene Panel Early-onset juvenile parkinsonism	Madhya Pradesh	R	22030	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0191	Ectodermal dysplasia-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0192	Gene Panel Ehler Danlos syndrome**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0193	Epidermolysis bullosa-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0194	Gene Panel Episodic ataxia**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0195	Familial female mental retardation/ epilepsy-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0196	Gene panel Familial hemiplegic migraine**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0197	Gene Panel Fanconi anemia**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0198	Fatty acid oxidation disorders-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	condition. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0199	Glycine encephalopathy-Gene Panel	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0200	Gene Panel Glycogen storage disorder**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0201	Glycosylation (CDG) disorders-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0202	Gene Panel Haemophilia (F8 & F9)Gene Panel, EDTA Blood	Madhya Pradesh	R	19950	Daily : 11:30am	After 4 weeks : 6:00pm	Genetci history is mandatory. Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability.
G0203	Gene Panel Hemophagocytic lymphohistiocytosis (HLH)**	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0204	Hereditary elliptocytosis-Gene Panel, EDTA Blood	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0205	Hereditary multiple exostoses Panel-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0206	Hereditary pancreatitis gene panel**	Madhya Pradesh	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0207	Hereditary Spherocytosis-Gene Panel, EDTA Blood	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0208	Gene Panel Homocystinuria**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	degree celsius. Sample to be preferably shipped at Referigerated condition. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0209	Gene Panel Hypercholesterolemia**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0210	Hyperekplexia-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0211	Hypomyelination syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0212	Icthyosis-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0213	Gene panel Joubert syndrome**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0214	Kallmann syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0215	Leber congenital amaurosis-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0216	Gene Panel Leigh syndrome & mitochondrial encephalopathy**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0217	Gene Panel Leukodystrophy**	Madhya Pradesh	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0218	Lissencephaly-Gene Panel, EDTA Blood	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0219	Maple syrup urine disease-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0220	Meckel Gruber syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0221	Gene Panel Methylmalonic aciduria (MMA)**	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0222	Microcephaly-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0223	Microphthalmia and anophthamia-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0224	Gene Panel Mucopolysaccharidosis**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0225	Myotonia congenita-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0226	Gene Panel Neuronal ceroid lipofuscinosis	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0227	Neuronal migration disorder-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0228	Niemann-Pick disease gene panel**	Madhya Pradesh	R	19380	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0229	Non-immune hydrops-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0230	Gene Panel Noonan syndrome**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0231	Oculocutaneous albinism-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0232	Optic atrophy-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0233	Organic acidemia-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0234	Gene Panel Osteogenesis imperfecta**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0235	Osteopetrosis-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0236	Pachydermoperiostosis & primary hypertrophic osteoarthropathy-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0237	Pantothenate kinase-associated neurodegeneration-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0238	Peroxisomal disorder-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0240	Gene Panel Polycystic kidney disease**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0241	Pontocerebellar hypoplasia-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0242	Gene Panel Primary hyperoxaluria**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0243	Gene Panel (Primary immunodeficiency)**	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0244	Progressive familial intrahepatic cholestasis-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0245	Progressive myoclonic epilepsy-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0246	Rasopathy-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0247	Retinal Degeneration-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0249	Seckel syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0250	Gene Panel Severe combined immunodeficiency (SCID)**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
G0251	Sideroblastic anaemia-Gene Panel, EDTA Blood	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Amibent or referigerated condition based on sample stability. Amniotic Fluidor Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0252	Gene Panel Skeletal dysplasia**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0253	Sotos syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0254	Gene Panel Urea cycle defects**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0255	Usher syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0256	Waardenburg syndrome-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0257	Xeroderma pigmentosum-Gene Panel	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
G0323	BTP-Beta Trace Protein GFR*	Madhya Pradesh	D	3530	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
G0324	Gastrointestinal panel by FilmArray (Biofire) - Stool	Madhya Pradesh	E	17500	Daily : 9:00am to 9:00pm	After 10 hrs	5gm of Stool specimen in Sterile leak proof container
G0333	Klippel-Feil syndrome-Gene panel	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0334	Nephrotic syndrome-Gene Panel	Madhya Pradesh	R	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0335	Genotypic phenotypic correlation / case review by International experts	Madhya Pradesh	R	5100	Daily : 11:30am	Next day : 6:00pm	Prior Appointment
G0336	Glioma Panel-1	Madhya Pradesh	D	20865		7 days	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
G0337	Gene panel ALL risk stratification - T-ALL**	Madhya Pradesh	E	17000	Daily : 11:30am	20th day : 6:00pm	4 ml EDTA Whole blood/bone marrow aspirate
G0338	Neurodegeneration with brain iron accumulation-Gene Panel	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0339	FKTN-Gene Sequencing	Madhya Pradesh	E	19000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0340	Gaucher disease (GBA)-Gene Analysis	Madhya Pradesh	E	19000	Daily: 7:30am	4 weeks	ice packs. Do not freeze. TRF, history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0341	FISH-MDM2 (12q15) GENE AMPLIFICATION**	Madhya Pradesh	E	9700	Daily : 11:30am	15th day : 6:00pm	Formalin fixed paraffin embedded tissue block. Ship at room temperature. Tissue block must contain at least 10% tumor tissue. Histopath report with History mandatory.
G0344	Clinical exome & whole mitochondrial genome sequencing	Madhya Pradesh	E	26000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0345	Giant axonal neuropathy-1 (GAN)-Gene analysis	Madhya Pradesh	E	19000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
G0348	Bone marrow failure syndrome gene panel**	Madhya Pradesh	R	21000	Daily : 11:30am	After 4 weeks : 6:00pm	4 ml EDTA Whole Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
G0349	GI mapping biopsies (Upto 5 containers)	Madhya Pradesh	C	5040	Daily : 9:00am	5th day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
G1901	CANINE-Gamma GT (GGTP)	Madhya Pradesh	A	295	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
G1903	Canine - Glucose fasting	Madhya Pradesh	A	100	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
G8175	Germline HRR gene panel	Madhya Pradesh	R	23460	Daily 7.30 am	25TH Days	3 ml of EDTA whole blood
H0001	H1N1 (Swine Flu) RNA Detection by Real Time PCR - Swab	Madhya Pradesh	C	5040	Daily : 7:00am & 3:00pm	Same Day : 3:00pm & 11:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
H0002	Hemochromatosis (HFE) Mutations Analysis (C282Y, H63D)	Madhya Pradesh	D	7490	Tue & Thu: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical history and TRF mandatory.
H0003	Haemoglobin-Free biochemical Plasma	Madhya Pradesh	R	1120	Daily : 11:30am	5th day : 6:00pm	2 ml of Citrated Plasma
H0004	Haemoglobin, free, urine	Madhya Pradesh	C	450	Daily : 11:30am	Same day : 9:00pm	10 ml of Spot urine
H0005	Haemophilia Profile	Madhya Pradesh	D	3745	Mon: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
H0007	Haemophilus Influenza B Antigen Antigen detection by Agglutination - CSF	Madhya Pradesh	D	1390	Daily : 9:00am to 3:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0008	Haemophilus Influenza B Antigen Antigen detection by Agglutination - Serum	Madhya Pradesh	D	1390	Daily : 9:00am to 3:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0009	Haemophilus Influenza B Antigen Antigen detection by Agglutination - Urine	Madhya Pradesh	D	1390	Daily : 9:00am to 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
H0011	Hanging Drop preparation for V Cholerae by Microscopy - Stool	Madhya Pradesh	B	435	Daily : 9:00am to 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
H0013	Haptoglobin genotype by PCR**	Madhya Pradesh	R	3895	Wed & Sat : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
H0014	Haptoglobin Serum	Madhya Pradesh	C	2350	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0015	HAV-IgG Antibody by CLIA - Serum	Madhya Pradesh	B	1200	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0016	HAV-IgM Antibody by CLIA - Serum	Madhya Pradesh	B	1200	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0017	HAV-RNA Detection by PCR - EDTA Blood	Madhya Pradesh	D	5245	Daily : 7:30am	7th day : 6pm	3 ml of Whole Blood in EDTA Vacutainer
H0018	HbA1C- Glycated Haemoglobin, by HPLC method	Madhya Pradesh	B	500	Daily: 9am and 3pm and 9pm	After 6 hrs	3 ml EDTA whole blood
H0026	Anti HBcAg total by CLIA - Serum	Madhya Pradesh	B	1090	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0027	Anti HBcAg-IgM by CLIA - Serum	Madhya Pradesh	B	1090	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0028	Anti HBeAg by CLIA - Serum	Madhya Pradesh	B	1090	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0029	Anti HbsAg Antibody by CLIA - Serum	Madhya Pradesh	B	1090	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0030	HBV DNA Detection by Taqman RT PCR - Serum or Plasma	Madhya Pradesh	C	5265	Daily : 07:30am	Same day : 11:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0031	HBV DNA Detection by PCR - Amniotic fluid**	Madhya Pradesh	D	4335	Daily : 11:30am	7th day : 6:00pm	EDTA Plasma in Plain Sterile Vacutainer 10 ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
H0032	HBV DNA Detection by PCR - Ascitic fluid**	Madhya Pradesh	D	4335	Daily : 11:30am	7th day : 6:00pm	3 ml of Fluid in Sterile leak proof container
H0033	HBV DNA Detection by PCR - CSF**	Madhya Pradesh	D	4335	Daily : 11:30am	7th day : 6:00pm	3 ml of Body Fluid in Sterile leak proof container
H0034	HBV DNA Detection by PCR - Liver tissue**	Madhya Pradesh	D	4335	Daily : 11:30am	7th day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
H0035	HBV DNA Genotype & Drug Resistance by PCR & Sanger Sequencing - Serum or Plasma	Madhya Pradesh	D	8560	Mon : 9:00am	Wed : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0037	HBeAg by CLIA - Serum	Madhya Pradesh	B	1090	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0038	HBsAg-Confirmation by CLIA - Serum	Madhya Pradesh	C	1090	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0039	HBsAg-Quantitative by CLIA - Serum	Madhya Pradesh	C	1230	Daily : 9:00am to 9:00pm	After 10 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0040	HBsAg-Screening by CLIA - Serum	Madhya Pradesh	B	400	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0041	HBV DNA Viral load by Taqman RT PCR - Serum or Plasma	Madhya Pradesh	C	7000	Daily : 7:30am	Same day : 11:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0042	HBV YMDD mutation by PCR & Sanger Sequencing - Serum	Madhya Pradesh	D	4495	Mon : 9:00am	Wed : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0043	HCG Beta Subunit, Serum	Madhya Pradesh	B	750	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0044	HCG Beta Subunit Urine 24H	Madhya Pradesh	B	870	Daily: 9am to 9pm	After 8 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
H0045	HCG Beta Subunit Urine Spot	Madhya Pradesh	B	260	Daily: 9am to 9pm	After 8 hrs	10 ml of Spot urine
H0045_U	PREGNANCY test by Immunochromatography - Urine	Madhya Pradesh	B	165	Daily : 9:00am to 9:00pm	After 4 hrs	15ml of Spot Urine in Sterile leak proof container
H0046	HCG Beta Subunit, CSF	Madhya Pradesh	B	870	Daily: 9am to 9pm	After 8 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
H0047	HCG Stimulation Test for DHT	Madhya Pradesh	D	3850	Daily: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0048	HCG stimulation test for FSH and LH	Madhya Pradesh	B	1745	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0049	HCG Stimulation test for Testosterone	Madhya Pradesh	B	1090	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0050	HCV Genofibro (Fibromax by CLIA , Genotyping and IL28B by PCR & Sequencing - Serum	Madhya Pradesh	C	26880	Daily : 9:00am	10th day : 6pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0051	HCV RNA Combo 1 (SGPT, Anti HCV by EIA, HCV-RNA Detection & Quantification by RT-PCR, HCV-Genotyping by Sanger Sequencing) - Serum	Madhya Pradesh	C	11200	Daily : 9:00am	4th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0052	HCV Genotyping by PCR & Sanger Sequencing - Serum or Plasma	Madhya Pradesh	D	7385	Tue & Fri : 9:00am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0053	HCV-IgM Antibody by EIA - Serum	Madhya Pradesh	C	1680	Mon & Thu : 9:00am	Same day : 5pm	Vacutainer 3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0054	HCV RNA Detection by Taqman RT PCR - Serum or Plasma	Madhya Pradesh	C	4370	Daily : 07:30am	Same day : 11:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0055	HCV-total Antibody by CLIA - Serum	Madhya Pradesh	B	800	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0056	HCV RNA Viral load by Taqman RT PCR - Serum or Plasma	Madhya Pradesh	C	7990	Daily : 07:30am	Same day : 11:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0057	HDL Cholesterol Direct serum	Madhya Pradesh	A	195	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0058	HDV-IgM Antibody by EIA - Serum	Madhya Pradesh	D	5135	Daily : 9:00am	7th day : 6pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0059	HDV-RNA Detection by PCR - EDTA Plasma**	Madhya Pradesh	D	6205	Daily : 7:30am	10th day : 6pm	3 ml of EDTA Plasma in Sterile Plain Vacutainer
H0060	HDV-total Antibody by EIA - Serum	Madhya Pradesh	C	2130	Mon & Thu : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0061	HE4 with ROMA index	Madhya Pradesh	C	2800	Daily:7:30 am	Same day 8:00 PM	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0062	Heinz BodiesBlood	Madhya Pradesh	B	240	Daily : 11:30am	After 8hrs	3 ml EDTA Whole Blood
H0063	Helicobacter Pylori Antigen detection by Immunochromatography - Stool	Madhya Pradesh	C	1680	Daily : 9:00am to 3:00pm	After 8 hrs	5gm of Stool specimen in Sterile leak proof container
H0064	Helicobacter Pylori Culture by Vitek/ MALDITOF - Tissue	Madhya Pradesh	C	1345	Daily : 9:00am to 9:00pm	15th day : 5pm	Tissue with Normal Saline in Sterile leak proof container
H0065	Helicobacter Pylori IgA antibodies by EIA - Serum	Madhya Pradesh	C	2130	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0066	Helicobacter Pylori IgG antibodies by CLIA - Serum	Madhya Pradesh	C	2130	Daily : 9:00am	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0067	Helicobacter Pylori IgM antibodies by EIA - Serum	Madhya Pradesh	C	2130	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0069	Hemoglobin	Madhya Pradesh	A	110	Daily: 9am to 9pm	After 6 hrs	3 ml EDTA Whole Blood
H0070	Hemosiderin, LBC	Madhya Pradesh	C	1735		2nd day	More than 1 MI fluid in a leak proof plain container along with relvent clinical & radiological details
H0071	Heparan Sulphate Sulphamidase Sanfillipo / MPS III A Blood**	Madhya Pradesh	R	10200	Daily : 11:30am	After 4 weeks : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
H0072	HBV Immunity Profile by CLIA - Serum	Madhya Pradesh	C	2465	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0073	HBV Profile by CLIA - Serum	Madhya Pradesh	C	5040	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0074	Hepatitis Marker Profile - Acute by CLIA - Serum	Madhya Pradesh	C	3695	Mon, Wed & Fri : 7:00am	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0075	Hepatitis Profile, Comprehensive by CLIA - Serum	Madhya Pradesh	C	7450	Mon & Thu : 9:00am	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0076	Hepatitis - Atypical Virus Profile by CLIA - Serum	Madhya Pradesh	D	4550	Daily : 7:30am	6th day : 7pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0077	Hepcidin Serum	Madhya Pradesh	D	1925	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0078	Her-2/Neu(C-Er-B2)(ERBB2)by FISH	Madhya Pradesh	D	11235	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
H0079	HEV IgG Antibody by EIA - Serum	Madhya Pradesh	B	1415	Mon & Thu : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0080	HEV IgM Antibody by EIA - Serum	Madhya Pradesh	B	1415	Daily : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0081	HGV-RNA Detection by PCR - EDTA plasma**	Madhya Pradesh	R	5305	Daily : 11:30am	15th day : 6:00pm	3 ml of EDTA Plasma in Sterile Plain Vacutainer
H0082	Hirsutism Profile-1	Madhya Pradesh	C	4145	Tue and Fri: 7am	Same day 7pm	5 ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
H0083	Hirsutism Profile-2	Madhya Pradesh	C	6285	Tue and Fri: 7am	Same day 7pm	5 ml of serum
H0084	Hirsutism Profile-3	Madhya Pradesh	C	8400	Tue and Fri: 7am	Same day 7pm	5 ml of serum
H0085	Histamine EIA Plasma	Madhya Pradesh	R	6885	Daily : 11:30am	20th day : 6:00pm	3 ml of EDTA plasma
H0086	Histidine Quantitative Serum	Madhya Pradesh	R	6325	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
H0087	Histidine Quantitative Urine 24H**	Madhya Pradesh	R	6325	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
H0088	Histochemistry (special stain with reporting) Alcian Blue Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0089	Histochemistry (special stain with reporting) Bleach Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0090	Histochemistry (special stain with reporting) Congo Red Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0091	Histochemistry (special stain with reporting) Fite Farraco Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0092	Histochemistry (special stain with reporting) Fontana Masson Silver Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0093	Histochemistry (special stain with	Madhya	C	1455		3rd Day	Paraffin blocks (Formalin fixed

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	reporting) Giemsa Tissue/FFPE	Pradesh					paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0094	Histochemistry (special stain with reporting) Gomori Silver Methanamine (GMS) Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0095	Histochemistry (special stain with reporting) Hale S Colloidal Iron Tissue/ FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0096	Histochemistry (special stain with reporting) Histochemistry Unspecified Special stain Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0097	Histochemistry (special stain with reporting) Masson Trichrome Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0098	Histochemistry (special stain with reporting) Mucicarmine Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0099	Histochemistry (special stain with reporting) Orcein Tissue/FFPE	Madhya Pradesh	C	1455			Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0100	Histochemistry (special stain with reporting) PAS with and/or Without Diastase Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0101	Histochemistry (special stain with reporting) PASM Stain Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0102	Histochemistry (special stain with reporting) Prussian Blue Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0103	Histochemistry (special stain with reporting) PTAH Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0104	Histochemistry (special stain with reporting) Reticulin Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0105	Histochemistry (special stain with reporting) Van Gieson Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0106	Histochemistry (special stain with reporting) Verhoeff S Vangieson Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0107	Histochemistry (special stain with reporting) Von Kossa Stain Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
H0108	Histochemistry (special stain with reporting) ZN Tissue/FFPE	Madhya Pradesh	C	1455		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							detail
H0109	Histochemistry (Special Stain Only) Alcian Blue Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0110	Histochemistry (Special Stain Only) Bleach Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0111	Histochemistry (Special Stain Only) Congo Red Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0112	Histochemistry (Special Stain Only) Fite Farraco Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0113	Histochemistry (Special Stain Only) Fontana Masson Silver Method Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0114	Histochemistry (Special Stain Only) Giemsa Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0115	Histochemistry (Special Stain Only) Gomori Silver Methanamine Stain (GMS) Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0116	Histochemistry (Special Stain Only) Hale S Colloidal Iron Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0117	Histochemistry (Special Stain Only) Masson Trichrome Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0118	Histochemistry (Special Stain Only) Mucicarmine Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0120	Histochemistry (Special Stain Only) PAS With and/or Without Diastase Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0121	Histochemistry (Special Stain Only) PASM Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0122	Histochemistry (Special Stain Only) Prussian Blue Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0123	Histochemistry (Special Stain Only) PTAH Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0124	Histochemistry (Special Stain Only) Reticulin Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0125	Histochemistry (Special Stain Only) Unspecified Special Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0126	Histochemistry (Special Stain Only) Van Gieson Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0127	Histochemistry (Special Stain Only) Verhoeff S Vangieson Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0128	Histochemistry (Special Stain Only) Von Kossa Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0129	Histochemistry (Special Stain Only) ZN Stain Tissue/FFPE	Madhya Pradesh	C	785		3rd Day	Single Paraffin block / section on charged slide
H0130	Histology Frozen Section (Intra Operative) day frozen 8am to 8pm, working days	Madhya Pradesh	D	10165	Daily: 9am to 5pm	After 6 hours	Fresh tissue to be sent along with relevant clinical & radiological details
H0131	Second opinion	Madhya Pradesh	D	2675		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded)/ Stained unstained slides Please attach original Histopathology report and all relevant clinical and radiological detail
H0132	Histone Antibody Serum	Madhya Pradesh	C	1850	W,S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0133	Histopathological examination of-Amputation of Limb	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0135	Histopathological examination of-Appendix	Madhya Pradesh	C	1680		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0136	Histopathological examination of-Bone Marrow Trepine biospy	Madhya Pradesh	C	5040		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0137	Histopathological examination of-Bone Specimen	Madhya Pradesh	C	5040		7th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0138	Histopathological examination of-Brain Tissue	Madhya Pradesh	C	5040		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							proof container along with relevant clinical & radiological details
H0139	Histopathological examination of-Breast Lump	Madhya Pradesh	C	2240		4th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0141	Histopathological examination of- Breast MRM / Breast conservative surgery reflex FISH	Madhya Pradesh	D	8025		7th day for histo	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0142	Histopathological examination of-Buccal Mucosal Biopsy	Madhya Pradesh	C	1680		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0143	Histopathological examination of- Cervix Biopsy	Madhya Pradesh	C	1680		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0144	Histopathological examination of- Cervix Cone Biopsy	Madhya Pradesh	C	4145		6th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0145	Histopathological examination of- Colectomy / APR	Madhya Pradesh	C	8400		7th day for histo	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0147	Histopathological examination of- Commando With RND/MND	Madhya Pradesh	C	6720		6th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0148	Histopathological examination of - Endobronchial biopsy	Madhya Pradesh	C	1680		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0149	Histopathological examination of- Endometrium (D&C Material)	Madhya Pradesh	C	1120		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0150	Histopathological examination of- Eye ball Specimen	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0151	Histopathological examination of - Fistulectomy / Sinus tract includes special stains	Madhya Pradesh	C	1680		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0152	Histopathological examination of- Gall Bladder	Madhya Pradesh	C	1680		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0153	Histopathological examination of- Gastrectomy	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0154	Histopathological examination of- GI Endoscopic Biopsy / Colonic Biopsy	Madhya Pradesh	C	1680		4th day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0157	Histopathological examination of-Hepatic Resection	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0158	Histopathological examination of- Hirschsprungs Disease	Madhya Pradesh	C	7840		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0159	Histopathological examination of -Kidney Biopsy	Madhya Pradesh	C	2350		4th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0160	Histopathological examination of- Kidney Biopsy with Immunofluorescence IgG, IgM, IgA. C3, C1q etc	Madhya Pradesh	D	4815		7th Day	HistologyandIHC.pdf
H0161	Histopathological examination of-Liver Biopsy	Madhya Pradesh	C	5040		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0162	Histopathological examination of-Lung (Lobectomy)	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0163	Histopathological examination of-Lymph Node	Madhya Pradesh	C	1680		4th day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0164	Histopathological examination of-Mediastinal Mass	Madhya Pradesh	C	5040		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0165	Histopathological examination of-Muscle Biopsy	Madhya Pradesh	C	2015		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0166	Histopathological examination of-Nasal Polyp	Madhya Pradesh	C	1680		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0167	Histopathological examination of-Nephrectomy (Non Tumorous)	Madhya Pradesh	C	2800		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0168	Histopathological examination of-Oesophagectomy	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0169	Histopathological examination of-Orchidectomy	Madhya Pradesh	C	2015		4th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0170	Histopathological examination of-Ovarian Mass (Large)	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0171	Histopathological examination of-Ovarian Mass (Small To Medium)	Madhya Pradesh	C	2015		4th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0172	Histopathological examination of-Penis	Madhya Pradesh	C	3360		4th Day	clinical & radiological details Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0173	Histopathological examination of-Prostate Chips (Tur-P)	Madhya Pradesh	C	2240		4th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0174	Histopathological examination of-Prostate Trucut Needle Biopsy	Madhya Pradesh	C	5040		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0175	Histopathological examination of-Radical Cystectomy	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0177	Histopathological examination of-Radical Nephrectomy	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0178	Histopathological examination of-Radical Prostatectomy	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0179	Histopathological examination of-Sebaceous cyst	Madhya Pradesh	C	1680		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0180	Histopathological examination of-Skin Biopsy	Madhya Pradesh	C	1825		4th day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0181	Histopathological examination of- Skin Biopsy With Immunofluorescence IgG, IgM, IgA. C3, C1q	Madhya Pradesh	D	4815		7th Day	HistologyandIHC.pdf
H0182	Histopathological examination of-Small Intestinal Resection	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0183	Histopathological examination of-Soft Tissue Tumor	Madhya Pradesh	C	8400		5th Day	proof container along with relevant clinical & radiological details Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0184	Histopathological examination of-Spleen	Madhya Pradesh	C	2240		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0185	Histopathological examination of-Synovial Biopsy	Madhya Pradesh	C	1680		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0186	Histopathological examination of-Thyroidectomy	Madhya Pradesh	C	3640		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0187	Histopathological examination of- Trucut Biopsy (CT/USG Guided) Special stains & Reflex IHC	Madhya Pradesh	C	5040		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0188	Histopathological examination of-TUR Urinary Bladder	Madhya Pradesh	C	2240		4th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0189	Histopathological examination of-Unspecified specimen-Large	Madhya Pradesh	C	2700		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0190	Histopathological examination of-Unspecified specimen-Large & Complex	Madhya Pradesh	C	3920		6th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0191	Histopathological examination of-Unspecified specimen-Medium	Madhya Pradesh	C	1500		4th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0192	Histopathological examination of-Unspecified specimen-Small	Madhya Pradesh	C	1000		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0193	Histopathological examination of-Uterus With/Without Adnexae	Madhya Pradesh	C	2350		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0194	Histopathological examination of-Werthiems Hysterectomy	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0195	Histopathological examination of-Whipples Resection	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0198	HIT Antibody Heparin Induced Thrombocytopenia / PF4 Serum**	Madhya Pradesh	R	7500	Daily : 11:30am	4th day : 6:00pm	2ml serum required
H0201	HIV-1 Proviral DNA Detection by RT PCR - EDTA Blood	Madhya Pradesh	C	4200	Daily : 07:30am	Same day : 11:00pm	3 ml of Whole Blood in EDTA Vacutainer
H0202	HIV-1 RNA Detection by Taqman RT PCR - EDTA Plasma	Madhya Pradesh	C	3470	Daily : 07:30am	Same day : 11:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0205	HIV-1&2 Western blot by EIA - Serum	Madhya Pradesh	C	3025	Daily : 7:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0206	HIV-1 RNA Combo 1 - EDTA Blood & Plasma	Madhya Pradesh	C	7800	Daily : 7:30am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer & 3ml of EDTA Plasma in Sterile Plain Vacutainer
H0207	HIV1/2 - Human Immunodeficiency Virus - confirmatory serology by immuno-transfer technique - serum	Madhya Pradesh	R	10710	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0208	HIV-2 RNA Detectionby Real Time PCR - EDTA Plasma	Madhya Pradesh	C	4255	Daily : 7:30am	3rd day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0209	HIV-2 RNA Quantitative Viral load by Real Time PCR - EDTA Plasma	Madhya Pradesh	C	6720	Daily : 7:30am	3rd day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0210	HIV-2 RNA Combo by Flow Cytometry & Real Time PCR - EDTA Blood & Plasma	Madhya Pradesh	C	8400	Daily : 7:30am	3rd day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer & 3ml of EDTA Plasma in Sterile Plain Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0211	HIV-DUO by CLIA - Serum	Madhya Pradesh	B	430	Daily : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0212	HIV Early Detection Profile-1 - EDTA Plasma & Serum	Madhya Pradesh	C	5040	Daily : 7:30am	3rd day : 5pm	3ml of EDTA Plasma in Sterile Plain Vacutainer & 3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0213	HIV Early Detection Profile-2 - EDTA Blood & Serum	Madhya Pradesh	C	5040	Daily : 7:30am	3rd day : 5pm	3 ml of Serum in Gel Top Vacutainer or Plain Sterile Vacutainer & 3 ml of Whole Blood in EDTA Vacutainer
H0214	Donor Specific Antibody (DSA)	Madhya Pradesh	D	9525	Daily: 9am	Same day 6pm	Recipient: 5 ml serum + 6 ml Heparin/ACD/EDTA whole blood, Donor: 20 ml Heparin/ACD/EDTA Blood
H0215	HLA Typing A Locus	Madhya Pradesh	D	9525	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0216	HLA typing A,B,C,DR,DQ loci (ClassI & Class II)	Madhya Pradesh	D	20865	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0217	HLA Typing A,B,C Locus (ClassI)	Madhya Pradesh	D	13910	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0218	HLA Typing B Locus	Madhya Pradesh	D	9525	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0221	HLA Typing B5 Locus	Madhya Pradesh	D	9525	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0222	HLA Typing B51 Locus	Madhya Pradesh	D	9525	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0224	HLA Typing C Locus	Madhya Pradesh	D	9525	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0225	HLA Typing DQ Locus	Madhya Pradesh	D	7260	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0226	HLA Typing DR Locus	Madhya Pradesh	D	6420	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0227	HLA Typing DR & DQ Loci (ClassII)	Madhya Pradesh	D	11235	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0233	Single Antigen Bead Assay (Class I & Class II)*	Madhya Pradesh	E	40000	Daily: 9am	Same day 6pm	Recipient: 5 ml serum
H0234	Auto Antibody Detection Test*	Madhya Pradesh	D	1390	Daily: 9am	Same day 6pm	Recipient: 5 ml serum + 6 ml

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0235	HLA (Lymphocyte) Cross Match Total (T+B) cells*	Pradesh Madhya Pradesh	D	2460	Daily: 9am	Same day 6pm	Heparin/ACD/EDTA Blood Recipient: 5 ml serum + 6 ml Heparin/ACD/EDTA Blood, Donor: 6 ml Heparin/ACD/EDTA Blood
H0249	HLA-B*5701 Associated HCP5 rs2395029 Genotyping for Abacavir Hypersensitivity	Madhya Pradesh	D	3425	Tue & Thu: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
H0250	HLA B27 (Human Leucocyte Antigen B27)	Madhya Pradesh	C	2450	Daily: 9am	Next day 5pm	3 ml of EDTA whole blood
H0252	HLA-B*27 PCR	Madhya Pradesh	C	3470	Mon, Wed & Fri: 9:00am	Next day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
H0253	HOMA Index Insulin Resistance Test	Madhya Pradesh	C	1230	Daily: 9am to 9pm	Next day 5pm	3 ml of serum and 2 ml of fluoride plasma. Fasting sample.
H0254	Homocysteine Serum	Madhya Pradesh	C	1290	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0256	Homogentisic Acid (ALKAPTANURIA)	Madhya Pradesh	C	895	Daily : 11:30am	Same day : 9:00pm	10 ml of Spot urine
H0258	HPV-DNA Detectionby PCR and Typing by PCR & Sanger Sequencing - Vaginal swab	Madhya Pradesh	C	2910	Mon & Thurs : 9:00am	3rd day : 6:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
H0259	HPV-DNA Detectionby PCR and Typing PCR & Sanger Sequencing - Cervical scrapping	Madhya Pradesh	C	2910	Mon & Thurs : 9:00am	3rd day : 6:00pm	Thin Prep Container
H0260	HPV-DNA Detectionby PCR and Typing by PCR & Sanger Sequencing - Tissue	Madhya Pradesh	C	3135	Mon & Thurs : 9:00am	3rd day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
H0261	HPV-DNA Detectionby PCR and Typing by PCR & Sanger Sequencing - LBC	Madhya Pradesh	C	2910	Mon & Thurs : 9:00am	3rd day : 6:00pm	Thin Prep Container
H0264	HRT Profile	Madhya Pradesh	C	1370	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0265	HsCRP -High Sensitivity CRPSerum	Madhya Pradesh	B	830	Mon: 7.30am	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0266	HSV-1 IgG and IgM Antibody by EIA -	Madhya Pradesh	B	1310	Daily : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Serum	Pradesh					in Plain Sterile Vacutainer
H0268	HSV-1 IgG Antibody by EIA - Serum	Madhya Pradesh	B	720	Daily : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0269	HSV-1 IgM Antibody by EIA - Serum	Madhya Pradesh	B	720	Daily : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0270	HSV 1&2 DNA Detection by Real Time PCR - Blister fluid	Madhya Pradesh	C	4705	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	3ml of Fluid in Sterile leak proof container
H0271	HSV 1&2 DNA Detection by Real Time PCR - CSF	Madhya Pradesh	C	4705	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	2ml of CSF in Sterile leak proof tube or container
H0272	HSV 1&2 DNA Detection by Real Time PCR - Serum	Madhya Pradesh	C	4705	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0273	HSV 1&2 DNA Detection by Real Time PCR - Swab	Madhya Pradesh	C	4705	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
H0274	HSV 1&2 IgG and IgM Antibody by EIA - Serum	Madhya Pradesh	B	1310	Daily : 7:30am	5th day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0276	HSV 1+2 IgG Antibody by CLIA - Serum	Madhya Pradesh	B	720	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0277	HSV 1&2 IgG Antibody by western blot - Serum	Madhya Pradesh	D	2355	Daily: 11:30am	Next day 8pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0278	HSV 1+2 IgM Antibody by CLIA - Serum	Madhya Pradesh	B	720	Daily : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0279	HSV 1&2 IgM Antibody by western blot - Serum	Madhya Pradesh	D	2355	Daily: 11:30am	Next day 8pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0280	HSV-2 IgG and IgM Antibody by EIA - Serum	Madhya Pradesh	B	1310	Daily : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0282	HSV-2 IgG Antibody by EIA - Serum	Madhya Pradesh	B	720	Daily : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0283	HSV-2 IgM Antibody by EIA - Serum	Madhya Pradesh	B	720	Daily : 7:30am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0284	HTLV 1 And 2 Virus Total Antibody by CLIA - Serum	Madhya Pradesh	D	5245	Mon : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0285	Hu AntibodyNeuronal ag in Paraneoplastic syndromesSerum	Madhya Pradesh	D	4815	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0286	Huntington Disease Mutation Analysis	Madhya Pradesh	D	5135	Mon & Wed:	7th day : 6:00pm	3-5 ml of EDTA whole blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh			9:00am		required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
H0287	Hydatid Cyst Detection by Microscopy - Fluid	Madhya Pradesh	C	860	Daily : 9:00am to 3:00pm	After 6 hrs	3ml of Fluid in Sterile leak proof container
H0288	Hypertension Profile-Routine	Madhya Pradesh	C	2015	Daily: 9am to 9pm	After 8 hrs	3 ml of EDTA whole blood, 3 ml of serum, 2 ml of fluoride plasma, 10 ml of spot urine
H0289	Heart failure profile- maxi	Madhya Pradesh	D	5245	Mon: 7am	4th day	5 ml of serum
H0295	Extra slide set Histopathology cases	Madhya Pradesh	C	560		4th day	Stained slide set is issued
H0300	HBV mutation Profile by PCR & Sanger Sequencing - Serum**	Madhya Pradesh	R	11220	Daily : 11:30am	10th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0301	HBV precore mutations by PCR & Sanger Sequencing - Serum**	Madhya Pradesh	R	5610	Wed : 11:30am	10th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0302	HCV Reflex Profile - Serum	Madhya Pradesh	B	1525	Tue & Fri : 9:00am	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0304	Histopathological examination of - Fallopian tube	Madhya Pradesh	C	1680		3rd Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0306	Histopathological examination of-Placenta	Madhya Pradesh	C	2240		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0307	Histopathological examination of- Testicular biopsy	Madhya Pradesh	C	1680		4th day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0308	Single Antigen Bead Assay (Class I)*	Madhya Pradesh	E	23200	Daily: 9am	Same day 6pm	Recipient: 5 ml serum
H0309	Single Antigen Bead Assay (Class II)*	Madhya Pradesh	E	20000	Daily: 9am	Same day 6pm	Recipient: 5 ml serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0310	Homocysteine reflex B12-folate Serum	Madhya Pradesh	C	1905	Daily: 9am to 9pm	Next day 1pm	4ml of Serum
H0311	HIV-1 RNA Drug Resistance for integrase inhibitor by PCR & Sanger Sequencing - EDTA Plasma	Madhya Pradesh	D	5350	1st & 3rd Mon : 9:00am	15th day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0312	HIV-1 RNA Coreceptor Tropism (Genotype) by PCR & Sanger Sequencing - EDTA Plasma	Madhya Pradesh	D	7060	1st & 3rd Mon : 9:00am	15th day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0313	HBV Basal Core Promoter mutation by PCR & Sanger Sequencing - Serum**	Madhya Pradesh	R	5610	Wed : 11:30am	8th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0314	Homovanillic acid urine	Madhya Pradesh	R	6120	Daily : 11:30am	12th Day : 6:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 30 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
H0316	HIV-1 RNA Drug Resistance-Extended by PCR & Sanger Sequencing - EDTA Plasma	Madhya Pradesh	D	19795	1st & 3rd Mon : 9:00am	15th day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0317	HPFH (Hereditary persistence of fetal haemoglobin) 5 deletion**	Madhya Pradesh	R	13260	Daily : 11:30am	After 4 weeks : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability. Amniotic Fluidor

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Chorionic villus sample: 10 ml of Amniotic Fluid or CVS sample in sterile container or Folin wu tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition.
H0318	HNPCC/ Lynch Syndrome Panel	Madhya Pradesh	R	18000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
H0320	Breast Carcinoma Profile I	Madhya Pradesh	D	9630		10th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0321	Breast carcinoma Profile II	Madhya Pradesh	D	11770		10th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0322	Lung Carcinoma (NSCLC) Profile	Madhya Pradesh	D	12305	Daily : 9:00am	After 2 weeks : 6:00pm	(FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical Diagnosis with Histopathology report, HistologyandIHC.pdf
H0323	Colorectal Carcinoma Profile	Madhya Pradesh	D	6955		10th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0324	Head and neck cancer (HNSCC profile)	Madhya Pradesh	D	6955		5th Day	HistologyandIHC.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0326	Heavy Metals profile-1, Hair**	Madhya Pradesh	R	5610	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
H0327	Heavy Metals profile-2, Hair**	Madhya Pradesh	R	11730	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
H0328	Heavy Metals profile-3, Hair**	Madhya Pradesh	R	17850	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
H0329	Heavy metals for ayurvedic tablets test	Madhya Pradesh	R	5610	Daily : 11:30am	10th day : 6:00pm	Medicine sample (10 grm), Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clincial symptoms if any. Reported in 7 days
H0357	Histopathological examination of- GI Endoscopic Biopsy with reflex IHC	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0793	H3 Viral marker Profile by CLIA - Serum	Madhya Pradesh	C	2140	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
H0794	Heavy and Trace Metals profile-1	Madhya Pradesh	D	11500	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
H0795	Heavy and Trace Metals profile-2	Madhya Pradesh	D	5885	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
H0796	CELL BLOCK	Madhya Pradesh	C	4480		5th Day	HistologyandIHC.pdf
H0797	HIV-1 RNA Viral load by Taqman RT PCR - EDTA Plasma	Madhya Pradesh	C	5800	Daily : 07:30am	Same day : 11:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0798	HIV-1 RNA Quantitative Viral load by Taqman RT PCR - CSF	Madhya Pradesh	C	6700	Daily : 07:30am	Same day : 11:00pm	2ml of CSF in Sterile leak proof tube or container
H0799	Panel Reactive Antibody (PRA) Identification by Luminex* (ClassI)	Madhya Pradesh	D	11480	Daily: 9am	Same day 6pm	Recipient: 5 ml serum
H0800	Panel Reactive Antibody	Madhya Pradesh	D	11480	Daily: 9am	Same day 6pm	Recipient: 5 ml serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0801	(PRA)Identification by Luminex* (ClassII) Panel Reactive Antibody (PRA) Screening by Luminex (ClassI & Class II)	Pradesh Madhya Pradesh	D	6315	Daily: 9am	Same day 6pm	Recipient: 5 ml serum
H0808	HEV RNA Detection by Real Time PCR - EDTA Blood	Madhya Pradesh	D	5350	Daily : 9:00am	5th day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer
H0810	Hereditary Cancer panel **	Madhya Pradesh	R	18360	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
H0813	HLA typing A,B (ClassI) and DR loci (Class II)	Madhya Pradesh	D	13500	Daily: 9am	Same day 6pm	6 ml EDTA/ACD Blood
H0824	Any specimen with reflex IHC	Madhya Pradesh	C	8400		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0825	Histopathological examination of- Nasopharyngeal Biopsy	Madhya Pradesh	D	4815		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0826	Histopathological examination of-Salivary Gland	Madhya Pradesh	C	4030		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0827	Histopathological examination of- Laryngeal Biopsy	Madhya Pradesh	C	2465		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0828	Histopathological examination of-Nasal Polyp with Reflex IHC	Madhya Pradesh	D	2355		5th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0829	Rush Biopsy service	Madhya	D	3210		5th Day	NA

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
H0841	HPV-DNA Detection by Hybrid Capture reflex Typing by PCR & Sanger Sequencing - LBC or Cervical sampler	Madhya Pradesh	C	2910	Mon & Thurs : 9:00am	3rd Day : 6:00pm	Thin Prep Container/ Digene Sampler
H0849	High Resolution HLA Typing (A, B, C, DR, DQ, DP) - Solid Organ Transplant	Madhya Pradesh	E	14980	Mon & Thur : 9.00 am	7 Days	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
H0850	High Resolution HLA Typing (A, B, C, DR, DQ, DP) - BMT/ Stem cell Transplant	Madhya Pradesh	E	14980	Mon & Thur : 9.00 am	7 Days	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
H0856	Sub speciality expert opinion	Madhya Pradesh	D	10485		7th day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container Or paraffin blocks along with relevant clinical & radiological details
H0860	Histology - Frozen Section (Intra Operative) Night frozen 8 pm to 8 am, week off & public holidays	Madhya Pradesh	D	17120			Fresh tissue to be sent along with relevant clinical & radiological details
H0861	HLA-B* 1502 for carbamazepine toxicity**	Madhya Pradesh	D	9095	Daily: 7:30am	15 days	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
H0862	HLA - T and B Cell Flow Cytometry Cross Match (FCXM)	Madhya Pradesh	D	7100	Daily : 11:30am	3rd day : 6:00pm	Recipient: 5 ml serum, Donor: 15 ml Heparin
H0863	High Resolution HLA Typing B**	Madhya Pradesh	R	6020	Tue & Fri: 11:30am	10th day : 6:00pm	6 ml EDTA Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H0864	HCG as Tumour Marker, Serum	Pradesh Madhya Pradesh	C	895	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
H0865	High Resolution HLA Typing DQB1**	Madhya Pradesh	R	6020	Tue & Fri: 11:30am	10th day : 6:00pm	6 ml EDTA Blood
H0866	Pre Transplant Workup panel	Madhya Pradesh	D	37450	Daily: 9am	Same day 6pm	Recipient : 3 tubes EDTA/ ACD blood & 1 Gel tube, Donor-6 tubes EDTA/ACD Blood
H0867	PRA Identification Class I and Class II*	Madhya Pradesh	D	19260	Daily: 9am	Same day 6pm	Recipient: 5 ml serum
H0868	Histopathological examination of- Eyelid Biopsy	Madhya Pradesh	C	1680		4th Day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0869	HBV DNA Combo RT-PCR & Sanger Sequencing - Serum or Plasma	Madhya Pradesh	C	15120	Mon : 9:00am	Wed : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0870	HCV RNA Combo 2 by RT-PCR & Sanger Sequencing - Serum or Plasma	Madhya Pradesh	C	11200	Tue & Fri : 9:00am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer/ 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H0871	HIV-1 RNA Combo 2 by PCR & Sanger Sequencing - EDTA Plasma	Madhya Pradesh	D	22470	1st & 3rd Mon : 9:00am	15th day : 6:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
H0872	Anti HMGCR antibodies - serum	Madhya Pradesh	R	10200	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml Serum
H0874	NextGen HRR Gene Panel	Madhya Pradesh	E	25805	Daily : 11:30am	After 4 weeks : 6:00pm	FFPE block with a minimum of 10% (as assessed by the pathologist) of tumor cells should be present
H0876	Hydroxychloroquine (HCQ)	Madhya Pradesh	D	1390	Tue & Fri: 11:30am	Next day : 5:00pm	3 ml EDTA Whole Blood
H0877	Histopathology for fungal infection	Madhya Pradesh	C	1680	Daily : 9:00am	4th day	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
H0879	Hemophagocytic lymphohistiocytosis (T-	Madhya	R	9690	Daily : 11:30am	7th day : 6:00pm	Sample required: EDTA Blood 3ml

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H1030	cell Activation HLH Panel)** HEREDITARY SPASTIC PARAPLEGIA GENE PANEL	Pradesh Madhya Pradesh	R	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
H1901	CANINE-Histo-Large specimen (unspecified) More than 5 cm	Madhya Pradesh	C	2800	Daily: 9am	Report after 8 days	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
H1902	Second opinion service, Sldies and blocks	Madhya Pradesh	C	2800	Daily: 9am	Report after 8 days	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
H1903	CANINE-Histo-Medium specimen (unspecified) 2 to 5 cm	Madhya Pradesh	C	1230	Daily: 9am	Report after 8 days	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
H5671	Hepatitis B Monitoring profile	Madhya Pradesh	D	10410	Daily : 7:30am	Next day 5pm	3ml of serum, 3 ml of EDTA Plasma in Plain Sterile Vacutainer
H5682	Tru Heart Basic H5682	Madhya Pradesh	W	1399	Daily	Next day 5.00 pm	4 ml Serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
H5692	Hrd track stand-alone test**	Madhya Pradesh	R	132600	Daily : 11:30am	After 4 weeks : 6:00pm	FFPE blocks with 20% tumor content and HP report.
H8036	HIV-Screening by Immunochromatography	Madhya Pradesh	B	325	Daily: 9am to 9pm	After 6 hrs	
H8037	HBsAg-Screening by Immunochromatography	Madhya Pradesh	B	325	Daily: 9am to 9pm	After 6 hrs	
H8038	HCV-Screening by Immunochromatography	Madhya Pradesh	B	870	Daily: 9am to 9pm	After 6 hrs	
H8335	[Combo] Haemophilia A (F8) gene analysis (along with analysis of the F8-intron 22 inversion**	Madhya Pradesh	R	24480	Daily : 11:30am	After 4 weeks : 6:00pm	4 ml EDTA Whole Blood
H8395	Hypersensitive pneumonitis panel-Extended	Madhya Pradesh	E	10500	Mon, Wed & Fri : 11:30am	Next day : 5:00pm	3ml serum sample required
I0001	IA2 Insulin Autoantibodies serum	Madhya Pradesh	D	2890	Mon & Thu 9 am	Next day 5 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0002	IBD Profile	Madhya Pradesh	C	6050	Fri: 7am	Next day 5pm	3 ml of serum, 100 gms of stool in sterile container
I0004	Iduronate 2 Sulphate Sulphatase Hunter/ MPS II Blood	Madhya Pradesh	R	10100	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
I0005	IgA CSF	Madhya Pradesh	C	1570	Daily: 9am to 9pm	6th day	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
I0006	IgA Serum	Madhya Pradesh	C	785	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0007	IgD by Flow cytometry Blood*	Madhya Pradesh	D	3425	Daily : 11:30am	Next day 5pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
I0008	IGF BP3 Serum	Madhya Pradesh	C	3810	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0009	IGF-1 Somatomedin C Serum	Madhya Pradesh	C	3650	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0010	IgG by Flow cytometry Blood**	Madhya Pradesh	R	3010	Daily : 11:30am	4th day : 6:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
I0011	IgG Total antibody CSF	Madhya Pradesh	C	1570	Daily: 9am to 9pm	After 8 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
I0012	IgG Total antibody Serum	Madhya Pradesh	C	785	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0013	IgG4 Sub Class Serum	Madhya Pradesh	D	5500	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0014	IgM by Flow cytometry Blood**	Madhya Pradesh	R	3060	Daily : 11:30am	4th day : 6:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
I0015	IgM total CSF	Madhya Pradesh	C	1570	Daily: 7.30am	4th day	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
I0016	IgM total Serum	Madhya Pradesh	C	785	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0017	IHC Panels with reporting- Bone Marrow	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0018	IHC Panels with reporting-Brain Tumor	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							original Histopathology report and all relevant clinical and radiological detail
I0019	IHC Panels with reporting Breast IV Comprehensive diagnosis	Madhya Pradesh	D	9095		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0020	IHC panel-Breast-I, tissue	Madhya Pradesh	C	4255		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0021	IHC Panels with reporting Breast-II ER, PR, c-erb-B2	Madhya Pradesh	C	5375		3rd Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0022	IHC Panels with reporting Breast-III ER, PR, c-erb-B2, Ki 67	Madhya Pradesh	C	6945		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0023	IHC Panels with reporting - Cervical Cancer	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0025	IHC Panels with reporting - GastroIntestinal neoplasm	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0026	IHC Panels with reporting - Liver Neoplasm	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0027	IHC Panels with reporting - Lung Neoplasm	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0028	IHC Panels with reporting - Lymph node	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0032	IHC Panels with reporting Ovarian Neoplasm	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0033	IHC Panels with reporting Pancreatic Neoplasm	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0034	IHC Panels with reporting Peritoneal Adenocarcinoma Vs Mesothelioma	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0035	IHC Panels with reporting Pituitary Marker Profile - ACTH, LH, FSH, TSH, PROLACTIN	Madhya Pradesh	C	15120		7th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0037	IHC Panels with reporting - Round Cell Tumor(Adult/Pediatric)	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0038	IHC Panels with reporting - Carcinoma of unknown primary (CUP)	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0039	IHC Panels with reporting - Salivary Gland Neoplasm	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0040	IHC Panels with reporting - Soft Tissue Tumor	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0041	IHC Panels with reporting - Testicular Neoplasm	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0042	IHC Panels with reporting - Thyroid Neoplasm	Madhya Pradesh	C	9520		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0043	IHC Panels with reporting - 34 Beta E12 (HMWCK)	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0046	IHC single marker with reporting Alpha Feto Protein	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0048	IHC single marker with reporting Bcl-2 Oncoprotein	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0049	IHC single marker with reporting Bcl-6 Oncoprotein	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0050	IHC single marker with reporting Beta Catenin	Madhya Pradesh	D	2355		4th Day	detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0051	IHC single marker with reporting Beta HCG	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0053	IHC single marker with reporting C4D	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0054	IHC single marker with reporting CA 125	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0055	IHC single marker with reporting CA-19.9 **	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0056	IHC single marker with reporting Calcitonin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0057	IHC single marker with reporting Calponin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0058	IHC single marker with reporting Carletinin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0061	IHC multiplex marker with reporting CD 20 + CD 3 Cocktail	Madhya Pradesh	D	3210		4th Day	all relevant clinical and radiological detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0062	IHC single marker with reporting CD10	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0063	IHC single marker with reporting CD117	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0064	IHC single marker with reporting CD138	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0065	IHC single marker with reporting CD15	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0066	IHC single marker with reporting CD1a	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0067	IHC single marker with reporting CD20	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0068	IHC single marker with reporting CD21	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0070	IHC single marker with reporting CD23	Madhya Pradesh	D	2355		4th Day	original Histopathology report and all relevant clinical and radiological detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0071	IHC single marker with reporting CD3	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0072	IHC single marker with reporting CD30	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0073	IHC single marker with reporting CD31	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0074	IHC single marker with reporting CD34	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0075	IHC multiplex marker with reporting CD4 + CD8 Cocktail	Madhya Pradesh	D	3210		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0076	IHC single marker with reporting CD43	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0077	IHC single marker with reporting CD45 LCA	Madhya	D	2355		4th Day	Paraffin blocks (Formalin fixed

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0078	IHC single marker with reporting CD5	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0079	IHC single marker with reporting CD56	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0080	IHC single marker with reporting CD57	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0081	IHC single marker with reporting CD68	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0082	IHC single marker with reporting CD79a	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0083	IHC single marker with reporting CD99 Mic2	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0084	IHC single marker with reporting CDX2	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0085	IHC single marker with reporting CEA	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0086	IHC-CERB2 (Her 2/neu)-IHC single marker, tissue	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0087	IHC single marker with reporting C-erb-B2 for Gastric carcinoma Her2 Neu	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0088	IHC single marker with reporting Chromogranin A	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0089	IHC single marker with reporting CK10	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0090	IHC single marker with reporting CK19	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0091	IHC single marker with reporting CK20	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0092	IHC single marker with reporting CK5/6	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0093	IHC single marker with reporting CK7	Madhya Pradesh	D	2355		4th Day	detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0094	IHC single marker with reporting CK8	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0095	IHC single marker with reporting p16INK4a	Madhya Pradesh	D	4815		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0096	IHC single marker with reporting Cyclin D1	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0097	IHC single marker with reporting Desmin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0098	IHC multiplex marker with reporting Desmoglein 3 + Napsin A Cocktail	Madhya Pradesh	D	3210		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0099	IHC single marker with reporting DOG-1	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0100	IHC single marker with reporting E Cadherin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0101	IHC single marker with reporting EGFR	Madhya Pradesh	D	2355		4th Day	all relevant clinical and radiological detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0102	IHC single marker with reporting EMA Epithelial Membrane Antigen	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0103	IHC-Estrogen Receptor (ER)-single marker, tissue	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0104	IHC single marker with reporting Factor VIII	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0106	IHC single marker with reporting GCDFP15	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0107	IHC single marker with reporting GFAP	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0108	IHC single marker with reporting Glypican-3	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0110	IHC single marker with reporting H.Pylori	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							original Histopathology report and all relevant clinical and radiological detail
I0113	IHC single marker with reporting Hep-Par 1	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0114	IHC single marker with reporting HMB45	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0115	IHC single marker with reporting HPV	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0119	IHC single marker with reporting Inhibin Alpha	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0120	IHC single marker with reporting Insulin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0121	IHC multiplex marker with reporting Kappa + Lambda Cocktail	Madhya Pradesh	D	3210		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0122	IHC single marker with reporting Ki 67	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0126	IHC single marker with reporting Lysoyme	Madhya	D	2355		4th Day	Paraffin blocks (Formalin fixed

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0127	IHC single marker with reporting Melan A	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0128	IHC single marker with reporting Mesothelin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0129	IHC single marker with reporting MOC 31	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0130	IHC single marker with reporting MPO Myeloperoxidase	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0131	IHC single marker with reporting Mum1 Protein	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0132	IHC single marker with reporting My F4 Myogenin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0133	IHC single marker with reporting Myoglobin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0134	IHC single marker with reporting NSE Neuron Specific Enolase	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0135	IHC single marker with reporting Oct.2	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0137	IHC single marker with reporting p53	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0138	IHC single marker with reporting p63	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0139	IHC single marker with reporting Pan CK AE1/AE3	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0140	IHC single marker with reporting Pax 5	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0141	IHC single marker with reporting Pax 8	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0142	IHC multiplex marker with reporting PIN 4 Cocktail Ck5 + Ck14 +P63 + P504S	Madhya Pradesh	D	3210		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0143	IHC single marker with reporting PLAP Placental Alkaline Phosphatase	Madhya Pradesh	D	2355		4th Day	detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0144	IHC-Progesterone Receptor (PR)single marker,tissue	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0146	IHC single marker with reporting Prostatic Acid Phosphatase	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0147	IHC single marker with reporting PSA Prostate Specific Antigen	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0148	IHC single marker with reporting S100 Protein	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0149	IHC single marker with reporting Smooth Muscle Actin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0150	IHC single marker with reporting Synaptophysin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0151	IHC single marker with reporting TdT	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0152	IHC single marker with reporting Thyroglobulin	Madhya Pradesh	D	2355		4th Day	all relevant clinical and radiological detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0154	IHC single marker with reporting TTF1	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0155	IHC single marker with reporting Vimentin	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0156	IHC single marker with reporting WT1 Wilm S Tumor Protein	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0157	IHC stain only 34 Beta E12	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0158	IHC stain only ALK1	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0159	IHC stain only Alpha Feto Protein	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0160	IHC stain only Bcl-2 Oncoprotein	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0161	IHC stain only Bcl-6 Oncoprotein	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0162	IHC stain only Beta Catenin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0163	IHC stain only Beta HCG	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0164	IHC stain only CA125	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0165	IHC stain only Calcitonin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0166	IHC stain only Calponin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0167	IHC stain only Calretinin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0170	IHC stain only CD10	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0171	IHC stain only CD117	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0172	IHC stain only CD138	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0173	IHC stain only CD15	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0174	IHC stain only CD1a	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0175	IHC stain only CD20	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0176	IHC stain only CD21	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0178	IHC stain only CD23	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0179	IHC stain only CD3	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0180	IHC stain only CD30	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0181	IHC stain only CD31	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0182	IHC stain only CD34	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0184	IHC stain only CD43	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0185	IHC stain only CD45 LCA	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0186	IHC stain only CD5	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0187	IHC stain only CD79a	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0188	IHC stain only CD99 Mic2	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0189	IHC stain only CDX2	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0190	IHC stain only CEA	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0191	IHC stain only c-erb-B2 Her 2 Neu	Madhya Pradesh	D	1605		3rd Day	Single Paraffin block / section on charged slide
I0192	IHC stain only Chromogranin A	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0194	IHC stain only CK 19	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0195	IHC stain only CK 20	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0196	IHC stain only CK 5/6	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0197	IHC stain only CK 7	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0198	IHC stain only CK8	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0200	IHC stain only Cyclin D1	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0201	IHC stain only Desmin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0202	IHC stain only DOG-1	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0203	IHC stain only E Cadherin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0204	IHC stain only EGFR	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0205	IHC stain only EMA Epithelial Membrane Antigen	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0206	IHC stain only ER Estrogen Receptor	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0207	IHC stain only Factor VIII	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0208	IHC stain only GCDFP15	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0209	IHC stain only GFAP	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0210	IHC stain only Glypican-3	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0211	IHC stain only H.Pylori	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0212	IHC stain only Hep-Par 1	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0213	IHC stain only HMB45	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0215	IHC stain only Inhibin Alpha	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0216	IHC stain only Insulin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0217	IHC stain only Ki 67	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0218	IHC stain only Melan A	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0219	IHC stain only Mesothelin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0220	IHC stain only MOC 31	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0221	IHC stain only MPO Myeloperoxidase	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0222	IHC stain only Multiplex CD 20 + CD3 Cocktail	Madhya Pradesh	D	2140		3rd Day	Single Paraffin block / section on charged slide
I0223	IHC stain only Multiplex CD4 + CD8	Madhya Pradesh	D	2140		3rd Day	Single Paraffin block / section on charged slide
I0224	IHC stain only Multiplex Desmoglein 3 + Napsin A	Madhya Pradesh	D	2140		3rd Day	Single Paraffin block / section on charged slide
I0225	IHC stain only Multiplex Kappa + Lambda	Madhya Pradesh	D	2140		3rd Day	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0226	IHC stain only Multiplex Pin-4 Cocktail Ck5 + Ck14 + P63 + P504S	Madhya Pradesh	D	2140		3rd Day	Single Paraffin block / section on charged slide
I0227	IHC stain only Mum1 Protein	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0228	IHC stain only My F4 (Myogenin)	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0230	IHC stain only NSE Neuron Specific Enolase	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0231	IHC stain only OCT.2	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0233	IHC stain only p53	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0234	IHC stain only p63	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0235	IHC stain only Pan CK AE1/AE3	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0236	IHC stain only Pax 5	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0237	IHC stain only Pax 8	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0238	IHC stain only PLAP	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0239	IHC stain only PR Progesterone Receptor	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0240	IHC stain only Prostatic Acid Phosphatase	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0241	IHC stain only PSA Prostate Specific Antigen	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0242	IHC stain only S100 Protein	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0243	IHC stain only Single Marker (Any other)	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0244	IHC stain only SMA Smooth Muscle Actin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0245	IHC stain only Synaptophysin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0246	IHC stain only TdT	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0247	IHC stain only Thyroglobulin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0248	IHC stain only TTF1	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0249	IHC stain only Vimentin	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0250	IHC stain only WT1 Wilm S Tumor Protein	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0251	IHC Comprehensive panel for diagnosis	Madhya Pradesh	C	8900		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0258	IL28 B SNP Genotyping by PCR & Sanger Sequencing - EDTA Blood	Madhya Pradesh	D	5350	Daily : 9:00am	10th day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer
I0259	Immune Status Panel	Madhya Pradesh	C	6830	Daily: 9am	After 8 hrs	5 ml of serum
I0260	Immunofixation Qualitative Serum	Madhya Pradesh	C	6270	Tue, Thu & Sat: 11:30am	Next day : 5:00pm	4 ml of Serum
I0261	Immunofixation Qualitative Urine Profile	Madhya Pradesh	C	7170	Thu : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0262	Immunofixation- Quantitative Serum	Madhya Pradesh	C	7980	Tue, Thu & Sat: 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0263	Immunofixation-Quantitative Urine	Madhya Pradesh	C	8960	Mon : 11:30am	Next day : 5:00pm	10 ml 24 Hrs Urine (Preferred) OR Spot urine
I0264	Immunoglobulin Profile IgG, IgM and IgA,CSF	Madhya Pradesh	C	4030	Daily: 7.30am	Next day 5pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
I0265	Immunoglobulin Profile IgG, IgM and IgA, Serum	Madhya Pradesh	C	1680	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0266	Impotence Profile	Madhya Pradesh	C	3695	Daily: 9am to 9pm	After 8 hrs	5 ml of serum
I0267	India Ink Preparation by Microscopy - CSF	Madhya Pradesh	B	380	Daily : 9:00am to 3:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
I0268	Infectious Mononucleosis Profile by IFA - Serum	Madhya Pradesh	C	5600	Mon & Thu : 7:30am	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
I0269	Infertility Profile , Female	Madhya Pradesh	C	7055	Daily: 9am	5th day	3 ml of EDTA whole blood, 5 ml of serum, 2 ml of flouride plasma, 10 ml of spot urine
I0270	Infertility Profile, Male	Madhya Pradesh	C	5880	Daily: 9am	5th day	3 ml of EDTA whole blood, 5 ml of serum, 2 ml of flouride plasma, 10 ml of spot urine, Semen sample
I0273	Inhibin A Serum	Madhya Pradesh	C	1345	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0274	Inhibin B Serum	Madhya Pradesh	C	2240	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0275	Insulin (Fasting) Serum	Madhya Pradesh	B	900	Daily: 9am to 9pm	After 6 hrs	2 ml of Serum Collected In Fasting.
I0275_PG	Insulin, Post Glucose Serum	Madhya Pradesh	B	895	Daily: 9am to 9pm	After 6 hrs	2 ml of serum collected 2 hours post Glucose Monohydrate 82.5 gm
I0275_PP	Insulin, Post Prandial Serum	Madhya Pradesh	B	895	Daily: 9am to 9pm	After 6 hrs	2 ml of serum collected 2 hours post lunch
I0275_R	Insulin (Random) Serum	Madhya Pradesh	B	895	Daily: 9am to 9pm	After 6 hrs	2 ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
I0276	Insulin Antibody Type I & II diabetes serum	Madhya Pradesh	D	2140	Monday & Thursday 9am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0277	Insulin Suppression Test - C-peptide Serum	Madhya Pradesh	C	7730	Daily: 9am to 9pm	After 8 hrs	3 ml of Serum (basal, 10,20,30,40,50,60 Minutes)
I0279	Interleukin -6 CLIA Plasma	Madhya Pradesh	C	1120	Daily 9 am - 9 pm	After 6 hrs	3 ml of EDTA plasma to be send in aliquot tube at proper temperature.
I0281	Intrinsic Factor Antibody by IFA Serum	Madhya Pradesh	C	2690	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0282	Iron biochemical Serum	Madhya Pradesh	B	350	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0284	Iron Urine 24H**	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
I0285	Iron Serum	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
I0286	Iron Studies, Serum	Madhya Pradesh	B	700	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0287	Islet Cell Antibody ICA-512 serum	Madhya Pradesh	D	2245	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
I0288	Isoleucine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
I0289	Isoleucine Quantitative Urine 24H**	Madhya Pradesh	R	6325	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
I0290	Isospora by smear examination - Stool	Madhya Pradesh	B	600	Daily : 9:00am to 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
I0291	IHC single marker with reporting Unspecified stain	Madhya Pradesh	D	2355		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0292	Reticulocyte count profile	Madhya Pradesh	B	435	Daily: 9am to 9pm	After 6 hrs	3 ml EDTA Whole Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					
I0293	IHC stain only CD 56_I0293	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0294	IHC stain only CD 57_I0294	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0295	IHC stain only CD68_I0295	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0296	IHC stain only Lysozyme_I0296	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0297	IHC p16INK4a stain only	Madhya Pradesh	D	3105		3rd Day	Single Paraffin block / section on charged slide
I0300	IgVH GENE MUTATION**_I0300	Madhya Pradesh	R	9325	Wed : 11:30am	10th day : 6:00pm	3 ml of EDTA whole blood
I0301	BCR ABL1 Mutation Analysis	Madhya Pradesh	E	10200	Daily : 9:00am	8th day : 6:00pm	5-10 ml of EDTA whole blood required. Sample stable for 3 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. Sample should reach to the lab within 48 to 72 hours.
I0305	IGH/CCND1:t(11;14) translocation Mantle Cell Lymphoma FISH, TISSUE	Madhya Pradesh	D	8560	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
I0312	IGH gene rearrangement Translocations involving 14q32 locus	Madhya Pradesh	D	8560	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
I0315	IDH1 and IDH2 gene mutation - FFPE	Madhya	E	7000	Mon, Wed & Fri:	5th day : 6:00pm	Tumor containing FFPE block

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Tissue	Pradesh			9:00am		(>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
I0316	IGH/BCL2 :t(14;18) translocation Follicular Lymphoma FISH, BONE MARROW	Madhya Pradesh	D	4815	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
I0317	IGH/BCL2 :t(14;18) translocation Follicular Lymphoma FISH, BLOOD	Madhya Pradesh	D	4815	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
I0318	IGH/BCL2 fusion	Madhya Pradesh	D	8560	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
I0319	IHC single marker with reporting p40	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0320	IHC single marker with reporting Uroplakin II	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0321	IHC stain only p40	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0322	IHC stain only Uroplakin II	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0323	Immature Platelet fraction (IPF) Count	Madhya Pradesh	B	490	Daily: 9am to 9pm	After 6 hrs	3 ml EDTA whole Blood with direct smear
I0324	Influenza A virus IgG Antibody by EIA - Serum	Madhya Pradesh	C	2015	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0325	Influenza A virus IgM Antibody by EIA - Serum	Madhya Pradesh	C	2015	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
I0326	Influenza B virus IgG Antibody by EIA - Serum	Madhya Pradesh	C	2015	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
I0327	Influenza B virus IgM Antibody by EIA - Serum	Madhya Pradesh	C	2015	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
I0331	ICG Clearance Test Serum	Madhya Pradesh	D	5885	Daily	Same Day after 6 hrs	5 ml Serum Sample(Sample should be collected in plain red vaccutainer) 2 tubes (1 tube for ICG clearance test & 1 tube for creatinine) 0 mins(4 red top vaccutainers-5 ml of blood in each tube) 5 mins(1 red top vaccutainer-2 ml of blood) 10 mins (1 red top vaccutainer-2 ml of blood) 15 mins (1 red top vaccutainer-2 ml of blood) 20 mins (1 red top vaccutainer-2 ml of blood) Sample should reach within 2 hrs of collection at 2-8 °C. Sample tubes to be wrapped in aluminium foil to avoid exposure to light. Left over dye solution to be sent across to the lab along with the samples.
I0332	IHC Androgen Receptor with interpretation	Madhya Pradesh	D	2355		6th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							all relevant clinical and radiological detail
I0337	PDL1 (SP263) by IHC - Ventana	Madhya Pradesh	D	8560		8th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0338	IMMUNOHISTOCHEMISTRY C-MYC ANTIBODY	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0339	IMMUNOHISTOCHEMISTRY EBV	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0342	IHC single marker with reporting TFE3	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0347	IHC stain only C-myc antibody	Madhya Pradesh	D	1285		4th Day	Single Paraffin block / section on charged slide
I0348	IHC stain only EBV	Madhya Pradesh	D	1285		3rd Day	Single Paraffin block / section on charged slide
I0350	IHC stain only CD7	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0352	IHC stain only IgG4	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0353	IHC stain only GATA3	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0354	IHC stain only H caldesmon	Madhya Pradesh	D	2355		4th Day	all relevant clinical and radiological detail Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0355	IHC stain only ERG	Madhya Pradesh	D	2355		4th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0363	IHC stain only Androgen Receptor	Madhya Pradesh	D	1605		3rd Day	Single Paraffin block / section on charged slide
I0366	Hypersensitive pneumonitis panel	Madhya Pradesh	D	7200	Mon, Wed & Fri : 11:30am	Next day : 5:00pm	3ml serum sample required
I0366_1	Penicillium Chrysogenum(Specific IgG)	Madhya Pradesh	C	1905	Mon & Thu : 11:30am	Next day : 6:00pm	Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
I0366_2	Cladosporium herbarum (Specific IgG)	Madhya Pradesh	C	1905	Mon & Thu : 11:30am	Next day : 6:00pm	Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
I0366_3	Aspergillus fumigatus(Specific IgG)	Madhya Pradesh	C	1905	Mon & Thu : 11:30am	Next day : 5:00pm	Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
I0366_4	Mucor Racemosus(Specific IgG)	Madhya Pradesh	C	1905	Mon & Thu : 11:30am	Next day : 6:00pm	Require serum (3ml) sample collected in fasting status

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0366_5	Alternaria alternata(Specific IgG)	Madhya Pradesh	C	1905	Mon & Thu : 11:30am	Next day : 6:00pm	preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day. Require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
I0366_6	Pigeon serum proteins, feathers and droppings(Specific IgG)	Madhya Pradesh	C	1905	Mon,Thu: 9am	Next day: 5pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
I0368	IGRAS (QFT Plus) by ELISA - Lithium Heparin Blood	Madhya Pradesh	D	3500	Mon, Wed, Fri : 7:30am	Next day : 5pm	5 ml of Lithium Heparin Blood
I0371	Iron, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
I0372	IDH1 and IDH2 gene mutation, AML, Blood	Madhya Pradesh	E	6500	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.
I0372_BM	IDH1 and IDH2 gene mutation, AML, Bone Marrow	Madhya Pradesh	E	6500	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0373	IGRAS (Quantiferon TB Gold) - 3tube by ELISA - Lithium Heparin Blood	Madhya Pradesh	C	3200	Mon, Wed, Fri : 7:30am	Next day : 5pm	5 ml of Lithium Heparin Blood
I0378	Immunofluorescence Analysis in Michels Medium (Kidney / Skin Biopsy)	Madhya Pradesh	D	3800	Daily : 9:00am	7 th day	Biopsy to be sent in Michelâ€™s Medium in a Leak proof container along with relevant clinical & radiological details
I0380	Isavuconazole **	Madhya Pradesh	R	6120	Daily : 7:30am	7th day : 6:00pm	3 ml of EDTA Plasma.Please indicate history of the drug dose taken, time of sample collection and clinical History along with the sample.
I0384	IHC stain only ATRX	Madhya Pradesh	D	3745	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0385	IHC single marker with reporting ATRX	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0386	IHC stain only IDH 1	Madhya Pradesh	D	3745	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0387	IHC single marker with reporting IDH 1	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0388	IHC stain only SOX 10	Madhya Pradesh	D	3745	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0389	IHC single marker with reporting SOX 10	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0391	IHC single marker with reporting HHV 8	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0394	IHC stain only SV 40	Madhya	D	3745	Daily : 9:00am	3rd day	Single Paraffin block / section on

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					charged slide
I0395	IHC single marker with reporting SV 40	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0397	IGLON5 ANTIBODY Serum	Madhya Pradesh	D	12840	Tue & Fri: 11:30am	Next day : 6:00pm	3 ml Serum
I0399	IHC stain only PLA 2 R1	Madhya Pradesh	D	3745	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0400	IHC single marker with reporting PLA 2 R1	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0401	IHC stain only STAT 6	Madhya Pradesh	D	2675	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0402	IHC single marker with reporting STAT 6	Madhya Pradesh	D	3210	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0403	IHC stain only INI 1	Madhya Pradesh	D	3745	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0404	IHC single marker with reporting INI 1	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0405	IHC stain only Sall 4	Madhya Pradesh	D	3745	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0406	IHC single marker with reporting Sall 4	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0407	IHC stain only NKX 3.1	Madhya Pradesh	D	3745	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I0408	IHC single marker with reporting NKX 3.1	Madhya Pradesh	D	4815	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0409	IHC stain only CD 4	Madhya Pradesh	D	2140	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0410	IHC single marker with reporting CD 4	Madhya Pradesh	D	3210	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0411	IHC stain only CD 8	Madhya Pradesh	D	2140	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
I0412	IHC single marker with reporting CD 8	Madhya Pradesh	D	3210	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I0413	IHC single marker with reporting CMV	Madhya Pradesh	D	2355	Daily : 9:00am	7th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
I1005	IgG Subclass 1(IgG1) serum by Nephelometry**	Madhya Pradesh	R	5100	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum (Age & gender manadatory)
I1006	IgG Subclass 2(IgG2) serum by Nephelometry**	Madhya Pradesh	R	5100	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum (Age & gender manadatory)
I1007	IgG Subclass 3(IgG3) serum by Nephelometry**	Madhya Pradesh	R	5100	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum (Age & gender manadatory)
I1025	IODINE, Serum by ICP-MS method**	Madhya Pradesh	R	9425	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml Serum
I1026	IODINE, Urine by ICP-MS method**	Madhya Pradesh	R	9180	Daily : 11:30am	After 4 weeks : 6:00pm	Individual Metals and Minerals.pdf
I1901	Impression Smears	Madhya Pradesh	D	1100	Daily: 9am	Report after 48 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
I1902	Insulin, Serum by CMIA	Madhya Pradesh	B	925	Daily: 9am to 9pm	Report after 24 hrs	available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
I5073	Isoniazid**	Madhya Pradesh	R	7140	Daily 7.30 am	10 days	2 ml Na-Heparin Plasma
I5649	FLU-XPRT VIRAL PANEL	Madhya Pradesh	D	4815	Daily : 7:00am & 3:00pm	Same Day : 3:00pm & 11:00pm	Swab with Sufficient Material in Viral Transport Medium (VTM)
I8151	Immunofixation image	Madhya Pradesh	C	560	Tue, Thu & Sat: 11:30am	Next day : 5:00pm	Pre-requisite: Kindly mention the VID for gel electrophoresis test asked.
I8155	PD-L1 Dako 22C3 IHC	Madhya Pradesh	D	15160	Daily 7.30 am	10th days	Paraffin Block
I8335	ILD Complete	Madhya Pradesh	D	18725	Tue , Fri: 9am	Next day 5:00 pm	5 ML OF SERUM
I8336	ILD Screen	Madhya Pradesh	D	7490	Mon,Thu: 9am	Next day 5:00 pm	3 ml EDTA WHOLE BLOOD,5 ML OF SERUM,
I8337	ILD Autoimmune Maxi	Madhya Pradesh	D	14445	Tue , Fri: 9am	Next day 5:00 pm	5 ML OF SERUM
I8338	ILD Autoimmune Profile Mini	Madhya Pradesh	D	8025	Tue , Fri: 9am	Next day 5:00 pm	5 ML OF SERUM
J0001	JAK-2 Mutation (CMPD) Exon 12 Mutation Bone Marrow	Madhya Pradesh	D	6955	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
J0002	JAK-2 Mutation (CMPD) (V617F) by Real Time PCR - Bone Marrow	Madhya Pradesh	D	6955	Tue & Fri : 9:00am	Next day : 7:00pm	3 ml of Bone Marrow in EDTA Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
J0003	JAK-2 Mutation (CMPD) Exon 12 Mutation Blood	Madhya Pradesh	D	6955	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
J0004	JAK-2 Mutation (CMPD) (V617F) by Real Time PCR - Blood	Madhya Pradesh	D	6800	Tue & Fri : 9:00am	Next day : 7:00pm	3 ml of Whole Blood in EDTA Vacutainer
J0005	JC-BK Virus DNA Detection by Real Time PCR - CSF	Madhya Pradesh	D	6100	Mon & Thurs : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
J0006	JC-BK Virus DNA Detection by Real Time PCR - EDTA Plasma	Madhya Pradesh	D	6100	Mon & Thurs : 9:00am	Next day : 5:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
J0007	JC-BK Virus DNA Detection by Real Time PCR - Urine	Madhya Pradesh	D	6100	Mon & Thurs : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
J0008	JEV RNA Detection by Real Time PCR - EDTA Blood	Madhya Pradesh	D	4300	Fri : 7:30am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
J0009	JEV RNA Detection by Real Time PCR - CSF	Madhya Pradesh	D	4300	Fri : 7:30am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
J0010	JEV RNA Detection by Real Time PCR - EDTA Plasma	Madhya Pradesh	D	4300	Fri : 7:30am	Next day : 5:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
J0011	JEV RNA Detection by Real Time PCR - Serum	Madhya Pradesh	D	4300	Fri : 7:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
J0012	Jo-1 Antibody Serum	Madhya Pradesh	C	1735	Wed and Sat 9:00 am	Same day 5:00 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
J0013	JC DNA Detection by Real Time PCR - CSF	Madhya Pradesh	D	4920	Mon & Thurs : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
J0014	JC DNA Detection by Real Time PCR - EDTA Plasma	Madhya Pradesh	D	4920	Mon & Thurs : 9:00am	Next day : 5:00pm	3ml of EDTA Plasma in Sterile Plain Vacutainer
J0015	JC DNA Detection by Real Time PCR - Urine	Madhya Pradesh	D	4920	Mon & Thurs : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container
J0016	JAK2 V617 Mutation, with Reflex to Jak2 Ex- 12, CALR Ex-9 Mutation and MPL W515, S505 Mutation	Madhya Pradesh	D	12305	Wed & Sat : 9:00am	10th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
J0017	JAK2 V617 Mutation, with Reflex to Jak2 Ex- 12, CALR Ex-9 Mutation and MPL W515, S505 Mutation,Bone Marrow	Madhya Pradesh	D	12305	Wed & Sat : 9:00am	10th day : 6:00pm	referring clinicians contact number and TRF mandatory. 3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.
J0020	JC DNA Detection by Real Time PCR - EDTA Blood	Madhya Pradesh	D	4920	Mon & Thurs : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
J0021	JC-BK Virus DNA Detection by Real Time PCR - EDTA Blood	Madhya Pradesh	D	6100	Mon & Thurs : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
K0001	Kappa And Lambda-Freelite, Urine	Madhya Pradesh	D	5350	Daily: 9am to 9pm		10 ml of 24 hrs urine (Preferred) / Spot urine
K0002	Kappa And Lambda-Freelite, Serum	Madhya Pradesh	D	5000	Daily: 9am to 9pm		3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
K0003	Kappa Light Chains by Flow cytometry	Madhya Pradesh	D	2890	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
K0004	Kappa Light Chains-Total Urine 24H	Madhya Pradesh	D	3745	Daily: 9am to 9pm		24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
K0005	FISH for Aneuploidy in POC - Chromosomes - 13,16,18, 21, 22, X and Y	Madhya Pradesh	D	9630	Daily : 9:00am	6th day : 6:00pm	avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab. POC in sterile container with normal saline with few drops of antibiotic. In case of expected delay store at 4 degrees. Never add formalin or freeze. Note: Whole Fetus is not accepted.
K0006	Karyotyping reflex FISH-13,18,21,X,Y AmnioticFluid	Madhya Pradesh	D	11575	Daily : 9:00am	15th day : 6:00pm	20 ml of Amniotic Fluid In sterile centrifuge Falcon 15 ml tubes.
K0007	Karyotyping reflex FISH(13,18,21,X,Y)Chorion Villi	Madhya Pradesh	D	10595	Daily : 9:00am	15th day : 6:00pm	20-30 mg Chorionic Villus in Sterile nutrient medium in 15 ml Falcon tubes(Medium will be provided by Metropolis).
K0010	FISH for Aneuploidy in POC - Chromosomes 13,18,21, X,Y	Madhya Pradesh	D	8560	Daily : 9:00am	6th day : 6:00pm	POC in sterile container with normal saline with few drops of antibiotic. In case of expected delay store at 4 degrees. Never add formalin or freeze.Note: Whole Fetus is not accepted.
K0012	Karyotyping by G-Banding Peripheral Blood	Madhya Pradesh	D	4300	Daily : 9:00am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
K0014	Karyotyping by G-Banding Bone marrow	Madhya Pradesh	D	4280	Daily : 9:00am	10th day : 6:00pm	2-3 ml Bone Marrow/3 ml Blood in Na-heparin vacutainer
K0017	Karyotyping by high resolution banding (HRB) Blood	Madhya Pradesh	D	5350	Daily : 9:00am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
K0018	Kaolin Clotting Time	Madhya Pradesh	B	1365	Daily: 9am	After 8 hrs	1 ml of citrated platelet poor plasma in frozen condition.
K0019	KRAS-Exon 2 and 3 Mutation Analysis - FFPE Tissue	Madhya Pradesh	D	8025	Mon, Wed & Fri: 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>10% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
K0021	Kennedy disease (CAG repeats) by PCR	Madhya Pradesh	R	7240	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
K0022	Karyotyping from cell pellet /slides	Madhya Pradesh	D	1605	Daily : 9:00am	10th day : 6:00pm	Cell pellet / Slides
K0023	Kala azar	Madhya Pradesh	R	1895	Daily : 11:30am	5th day : 6:00pm	3ml of serum
K0026	Karyotyping G Banding and FISH (Chr 13,18, 21, X and Y),Amniotic Fluid	Madhya Pradesh	D	16585	Daily : 9:00am	15th day : 6:00pm	20 ml of Amniotic Fluid In Sterile Falcon centrifuge tube
K0030	Karyotyping by G-Banding for neonates	Madhya Pradesh	D	4815	Daily : 9:00am	7th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
K0031	Karyotyping by G-Banding for 50 metaphase cells	Madhya Pradesh	D	4815	Daily : 9:00am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
K0033	Karyotyping by G-Banding reflex FISH	Madhya Pradesh	D	10595	Daily : 9:00am	15th day : 6:00pm	3 ml of Prenatal Cord Blood In Na-Heparin Vacutainer
K1901	Canine - Kidney panel	Madhya Pradesh	D	1125	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
K8112	Kidney biopsy for Medical Renal Disease	Madhya Pradesh	D	5350	Daily : 9:00am	7th day	1) Tissue in 10% NB formalin for histopathology 2) Tissue in Michel Medium for Immunofluorescence. Please use this code when you are sending 2 sample vials ; One in Formaldehyde & other in Michels medium.
K8113	Kidney biopsy for Transplant Pathology	Madhya Pradesh	D	10595	Daily : 9:00am	14th day	1) Tissue in 10% NB formalin for histopathology 2) Tissue in Michel

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Medium for Immunofluorescence. Please use this code when you are sending 2 sample vials ; One in Formaldehyde & other in Michels medium.
L0001	Lactate CSF	Madhya Pradesh	C	1175	Daily: 9am to 9pm	After 6 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
L0002	Lactate Biochemical Plasma	Madhya Pradesh	C	1175	Daily: 9am to 9pm	After 6 hrs	2 ml of Fluoride Plasma Labile Analyte
L0003	Lactose Tolerance Test 5 Samples for Glucose, Plasma	Madhya Pradesh	C	670	Daily: 9am to 9pm		2 ml of flouride plasma (5Samples)
L0004	Lambda Light Chains by Flow cytometry	Madhya Pradesh	D	2890	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
L0005	Lambda Light Chains-Total Urine	Madhya Pradesh	D	3745	Daily: 7.30am	Same day 5pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
L0006	Lamotrigine Level Plasma**	Madhya Pradesh	D	4385	Tue : 11:30am	Next day : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0007	LAP Score	Madhya Pradesh	B	1200	Daily : 11:30am	3rd day : 6:00pm	5 Peripheral Blood Smears and EDTA Blood. Clinical history and clinician no is mandatory
L0008	Liver cytosolic antigen type-1(LC-1)Antibody, Serum	Madhya Pradesh	C	3135	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0009	LD (Leigh`s Disease),Blood **	Madhya Pradesh	R	14790	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
L0010	LD bodies Leishmania Dodovani bodies by Microscopy - Blood	Madhya Pradesh	B	1035	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA Bone Marrow in EDTA Vacutainer
L0011	LD bodies Leishmania Dodovani bodies by Microscopy - Bone Marrow	Madhya Pradesh	C	1065	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA Whole Blood in EDTA Vacutainer
L0012	LDH Lactate Dehydrogenase Serum	Madhya Pradesh	A	380	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0013	LDH Lactate Dehydrogenase Ascitic Fluid	Madhya Pradesh	A	490	Daily: 9am to 9pm	After 6 hrs	2 ml of Ascitic Fluid
L0014	LDH Lactate Dehydrogenase Peritoneal Fluid	Madhya Pradesh	A	490	Daily: 9am to 9pm	After 6 hrs	2 ml of Peritoneal Fluid
L0015	LDH Lactate Dehydrogenase Pleural Fluid	Madhya Pradesh	A	490	Daily: 9am to 9pm	After 6 hrs	2 ml of Pleural Fluid
L0016	LDH- Isoenzyme by Elctrophoresis Serum	Madhya Pradesh	D	4920	Wed : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0017	LDL Cholesterol -Direct serum	Madhya Pradesh	A	350	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0019	Lead, Blood	Madhya Pradesh	C	2130	Daily:7:30 am	Same day 8:00 PM	HeavyMetalSamplingguidelines.pdf
L0020	Lead, Urine 24 Hrs	Madhya Pradesh	C	2130	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0021	Lead,Urine Spot	Pradesh Madhya Pradesh	C	2410	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
L0022	Legionella Culture by Vitek/ MALDITOF - Sputum / BAL**	Madhya Pradesh	R	1225	Daily : 9:00am to 9:00pm	15th day : 5pm	2 ml of Sputum in Sterile leak proof container
L0026	Legionella Pneumophila Antigen detection Urine	Madhya Pradesh	C	3305	Daily : 9:00am to 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
L0027	Legionella Pneumophila IgG antibodies by EIA - Serum	Madhya Pradesh	C	3080	Mon & Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0028	Legionella Pneumophila IgM antibodies by EIA - Serum	Madhya Pradesh	C	3080	Mon & Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0029	Leishmania Antibody-IgG Antibody for Kala Azar by EIA - Serum	Madhya Pradesh	C	1990	Tuesday/Friday 9 am	Same day 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0034	Leptin Human Serum	Madhya Pradesh	D	4495	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0035	Leptospira Detection by smear examination - Blood	Madhya Pradesh	B	655	Daily : 9:00am to 3:00pm	After 6 hrs	3 ml of Whole Blood in EDTA Vacutainer
L0036	Leptospira Detection by smear examination - CSF	Madhya Pradesh	B	765	Daily : 9:00am to 3:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
L0037	Leptospira Detection by smear examination - Urine	Madhya Pradesh	B	765	Daily : 9:00am to 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
L0038	Leptospira DNA Detection by Real Time PCR - Blood	Madhya Pradesh	D	2355	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
L0039	LeptospiraDNA Detection by Real Time PCR - Urine	Madhya Pradesh	D	2355	Daily : 10:00am	After 24hrs	15ml of Spot Urine in Sterile leak proof container
L0040	Leptospira IgG antibodies by EIA - Serum	Madhya Pradesh	C	1290	Daily : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0041	Leptospira IgM antibodies by EIA - Serum	Madhya Pradesh	C	1290	Daily : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0042	Leptospira Profile	Madhya Pradesh	C	2715	Daily : 11:30am	same day 11:00 pm	3 ml of Whole Blood in EDTA Vacutainer, 3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0043	Leucine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0044	Leucine Quantitative Urine 24H**	Madhya Pradesh	R	6325	Daily : 11:30am	7th day : 6:00pm	proper temperature. 10 ml of Spot urine
L0046	Leukemia panel by FISH - ALL	Madhya Pradesh	D	14980	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0047	Leukemia panel by FISH - AML	Madhya Pradesh	D	16050	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0051	Leukemia panel by FISH - CLL	Madhya Pradesh	D	14980	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0052	Leukemia panels by FISH - Eosinophilia CBF-Beta, PDGFR-Alpha & beta, FGFR-1 Bone Marrow/Leukemic Blood	Madhya Pradesh	D	12840	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0053	Myelodysplastic Syndrome (MDS) Panel by FISH, Bone marrow	Madhya Pradesh	D	9630	Daily : 9:00am	5th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0054	Comprehensive Multiple Myeloma Panel by FISH	Madhya Pradesh	D	17120	Daily : 9:00am	7th day : 6:00pm	3 ml of Heparin Bone Marrow
L0055	Leukemia -Acute Leukemia Panel, Blood	Madhya Pradesh	D	11235	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA or heparin blood and Direct smears required
L0056	Leukemia - Acute Leukemia Panel, Bone Marrow	Madhya Pradesh	D	11235	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA and heparin bone marrow. Direct smears required
L0057	Leukemia - Chronic Lymphoproliferative Disorder Panel, Blood	Madhya Pradesh	D	12840	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA or heparin blood and Direct smears required
L0058	Leukemia - Chronic Lymphoproliferative Disorder Panel, Bone Marrow	Madhya Pradesh	D	14765	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA or heparin bone marrow and Direct smears required
L0065	Myeloma panel by flowcytometry, Bone marrow**	Madhya Pradesh	D	9845	Daily : 11:30am	5th day : 6:00pm	3 ml of EDTA or heparin bone marrow and Direct smears required
L0066	LH-Luteinizing Hormone ECLIA Serum	Madhya Pradesh	B	540	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0067	LHON (Leber's hereditary optic neuropathy) by PCR**	Madhya Pradesh	R	16320	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0068	Lipase Serum	Madhya Pradesh	B	640	Daily: 9am to 9pm	After 6 hrs	ice packs. Do not freeze. TRF, history is mandatory. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0068_Fluid	Fluid Lipase	Madhya Pradesh	B	720	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0069	Lipase Urine	Madhya Pradesh	B	720	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
L0070	Lipid Profile - Maxi	Madhya Pradesh	C	1500	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0071	Lipid Profile-Mini	Madhya Pradesh	B	550	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0071_NF	Lipid profile 2 (Non- fasting)	Madhya Pradesh	B	545	Daily: 9am to 9pm	Lipid profile 2 (Non- fasting)	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0072	Lipid ScreenCholesterol and TriglyceridesSerum	Madhya Pradesh	B	490	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0072_NF	Lipid profile 3 (Non- fasting)	Madhya Pradesh	B	490	Daily: 9am to 9pm	Lipid profile 3	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh				(Non- fasting)	ml serum for profiles and reflex tests.
L0073	Total Lipids, serum	Madhya Pradesh	A	415	Daily: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0076	Lipoprotein Electrophoresis, Serum	Madhya Pradesh	C	3025	Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0077	Lithium Serum	Madhya Pradesh	C	560	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0079	Liver Function Test-1(Maxi)	Madhya Pradesh	B	1390	Daily: 9am to 9pm	After 6 hrs	5 ml of serum and 3 ml of citrate plasma
L0080	Liver Function Test-2(Mini)	Madhya Pradesh	B	600	Daily: 9am to 9pm	After 6 hrs	5 ml of serum
L0081	Liver Onco Marker Profile	Madhya Pradesh	D	1605	Daily: 9am to 9pm	After 8 hrs	5 ml of serum
L0082	Liver Function Test-3 (LIVER SCREEN)	Madhya Pradesh	B	435	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0083	LKM1 by IFA Liver Kidney Microsomes Reflex to end point titre Serum	Madhya Pradesh	C	2130	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0084	Lp(a)Lipoprotein(a) serum	Madhya Pradesh	C	1175	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0085	Lung Marker Profile	Madhya Pradesh	D	7650	Wed/Sat- 7:30 am	Next day 5:00 pm	3 ml of Serum and 3 ml of EDTA plasma
L0086	Lupus Anticoagulants (LAC) profile	Madhya Pradesh	D	2350	Daily: 9am and 5pm	After 6 hrs	1 ml of citrated platelet poor plasma in frozen condition.
L0087	LymeBorrelia Burgdorferi IgG antibodies by EIA - Serum	Madhya Pradesh	C	1905	Sat : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0089	LymeBorrelia Burgdorferi IgM antibodies by EIA - Serum	Madhya Pradesh	C	1905	Sat : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
L0091	Lysine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0092	Lysine Quantitative Urine	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
L0097	Leukemia Characterisation Panel for CLPD, Blood	Madhya Pradesh	D	21185	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA and heparin blood. Direct smears required
L0098	Leukemia Characterisation Panel for CLPD , Bone Marrow	Madhya Pradesh	D	21185	Daily : 11:30am	7th day : 6:00pm	3 ml of EDTA and heparin bone marrow and Direct smears required
L0099	Leukemia Comprehensive workup panel- Acute Leukemia, Blood	Madhya Pradesh	D	26750	Daily : 11:30am	12th day : 6:00pm	3 ml of EDTA and heparin blood. Direct smears required
L0100	Leukemia Comprehensive workup panel- Acute Leukemia, Bone Marrow	Madhya Pradesh	D	26750	Daily : 11:30am	12th day : 6:00pm	3 ml of EDTA and heparin bone marrow. Direct smears required
L0101	Leukemia CMPD Panel: JAK2 by PCR and BCR-ABL by FISH,Blood	Madhya Pradesh	D	9415	Daily : 9:00am	5th day : 6:00pm	3 ml of EDTA and Na-heparin blood and Direct smears required
L0102	Leukemia CMPD Panel: JAK2 by PCR and BCR-ABL by FISH, Bone Marrow	Madhya Pradesh	D	9415	Daily : 9:00am	5th day : 6:00pm	3 ml of EDTA Bone Marrow
L0103	Lp- PLA2 (PLAC)Cardiac risk marker serum	Madhya Pradesh	D	2675	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
L0105	Leukemia panel by FISH - ALL, Leukemic Blood	Madhya Pradesh	D	14980	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0106	Leukemia AML panel by FISH, Leukemic Blood	Madhya Pradesh	D	16050	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0107	Leukemia panel by FISH - CLL, Leukemic Blood	Madhya Pradesh	D	14980	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0108	Leukemia panels by FISH - Eosinophilia CBF-Beta, PDGFR-Alpha & beta, FGFR-1 Leukemic Blood	Madhya Pradesh	D	12840	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0109	Myelodysplastic Syndrome (MDS) Panel by FISH	Madhya Pradesh	D	9630	Daily : 9:00am	5th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0111	Lymphocyte subset panel with NK cells	Madhya Pradesh	D	7385	Daily : 11:30am	Next day : 5:00pm	6 ml EDTA whole Blood
L0112	LDH Lactate Dehydrogenase CSF	Madhya Pradesh	A	490	Daily: 9am to 9pm	After 6 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
L0113	Lyme-Borrelia Burgdorferi IgM by EIA - CSF	Madhya Pradesh	C	5375	Sat : 11:30am	Next day : 5:00pm	2 ml of CSF in Sterile leak proof tube or container
L0114	Lyme-Borrelia Burgdorferi IgG by EIA - CSF	Madhya Pradesh	C	5375	Sat : 11:30am	Next day : 5:00pm	2 ml of CSF in Sterile leak proof

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0117	Leukemia translocation Panel (ALL)* - t(1;19), t(4;11), t(9;22), t(12;21)	Pradesh Madhya Pradesh	D	8025	Daily : 9:00am	7th day : 6:00pm	tube or container 5-10 ml of EDTA whole blood required. Sample stable for 3 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. Sample should reach to the lab within 48 to 72 hours.
L0118	Leukemia translocation Panel (AML)* - t(8;21), t(15;17), Inv(16)	Madhya Pradesh	D	8025	Daily : 9:00am	7th day : 6:00pm	5-10 ml of EDTA whole blood required. Sample stable for 3 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. Sample should reach to the lab within 48 to 72 hours.
L0120	LDH	Madhya Pradesh	A	390	Daily	After 6 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
L0123	LDH Body fluid	Madhya Pradesh	A	490	Daily: 9am to 9pm	After 6 hrs	2 ml body fluid
L0133	LUNG CANCER TARGETED THERAPY PANEL (EGFR+ALK+ROS1)	Madhya Pradesh	E	19000	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical Diagnosis with Histopathology report, HistologyandIHC.pdf
L0134	Lung Panel by NGS [EGFR, BRAF, MET (includes Ex-14 skipping), RET, ERBB2, KRAS (SNVs & small INDELs); ALK, ROS1, RET and NTRK1, NTRK2, NTRK3 (Fusions)]	Madhya Pradesh	E	32000	Daily 9:00 am	25th Days	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							FFPE block to be shipped at Ambient temperature.
L0135	Lead, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
L0136	Lead, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
L0138	Lithium, Spot Urine**	Madhya Pradesh	C	3080	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
L0139	Leukemia Panel-Ph like B-ALL, Bone Marrow **	Madhya Pradesh	R	19380	Daily : 11:30am	10th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0139_BLOOD	Leukemia Panel-Ph like B-ALL, Blood**	Madhya Pradesh	R	19380	Daily : 9:00am	8th day : 6:00pm	3 ml BLOOD in Na-heparin vacutainer
L0140	Leukemia Panel-T-ALL, Blood**	Madhya Pradesh	R	16320	Daily : 11:30am	8th day : 6:00pm	3 ml Blood in Na-heparin vacutainer
L0140_BM	Leukemia Panel-T-ALL, Bone marrow**	Madhya Pradesh	R	16320	Daily : 11:30am	10th day : 6:00pm	3 ml Bone Marrow in Na-heparin vacutainer
L0141	Leukemia Panel-AML Mutation, Blood	Madhya Pradesh	D	12520	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
L0141_BM	Leukemia Panel-AML Mutation, Bone marrow	Madhya Pradesh	D	9630	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
L0142	NextGen Oncomine Myeloid 69 Gene Panel	Madhya Pradesh	E	35000	DAILY : 9 : 00 AM	20 Days	6 ml EDTA whole Blood
L0143	Leptin deficiency (LEP) -Gene Analysis	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L0144	Myeloma MRD Flow cytometry Bone marrow **	Madhya Pradesh	R	18870	Daily : 11:30am	5th day : 6:00pm	condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
L0145	Flow for AML MRD**	Madhya Pradesh	R	18870	Daily : 11:30am	5th day : 6:00pm	sample : 3 ml of EDTA Bone Marrow & direct smear
L0146	Lip biopsy for Sjogrens syndrome (Focus score) with WSI	Madhya Pradesh	C	2465	Daily : 9:00am	4th day	3 ml of EDTA Bone Marrow & General Instruction Local Samples reaching in 24 Hours can be send in Ambient condition and outstation samples to be shipped at Refrigerated condition. Do not keep EDTA or Heparin blood samples in contact with ice packs. Do not freeze. Medical history / Heamato-pathological reports with referring doctors number.
L0152	Lung Cancer NGS panel (SNVs, InDels and Fusions)- Liquid Biopsy**	Madhya Pradesh	E	46200	Daily : 11:30am	25th day : 6:00pm	Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
L1008	Levetiracetam, serum	Madhya Pradesh	D	5350	Wed : 11:30am	Next day : 5:00pm	Peripheral Blood in Streck tube (2 tubes of 10 ml each) (Tube supplied upon request)
L1015	LISTERIA MONOCYTOGENES - SEROLOGY	Madhya Pradesh	R	7140	Daily: 7.30am	1 month	3 ml of serum(sample collect in, plain vacutainer, not in gel tube)
							3 ml of serum for single test and 5 ml serum for profiles and reflex tests. International outsource : Category R for referred and I for international. Plz esquire before sending sample for arrangement of proper courier to international lab. Additional sample tubes may be needed for sample adequacy. TAT may vary based on information by concerned lab of delay. Courier

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
L1901	Canine - Lipase	Madhya Pradesh	B	655	Daily: 9am to 9pm	Report after 24 hrs	charges have to be paid extra based on sample temperature condition and shipment. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
L1902	Canine - Liver profile	Madhya Pradesh	D	910	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
L8015	LH/FSH Ratio Serum	Madhya Pradesh	B	925	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0001	Anti MAG (Myelin associated Glycoprotein) Ab (M0001)	Madhya Pradesh	D	4600	Tue & Fri : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0002	Magnesium Serum	Madhya Pradesh	B	520	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0003	Magnesium, Spot Urine	Madhya Pradesh	B	620	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
M0004	Magnesium Urine 24H	Madhya Pradesh	B	620	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
M0006	Malaria Parasite detection by Direct Smear examination by Microscopy - Blood	Madhya Pradesh	B	130	Daily : 9:00am, 1:00pm, 5:00pm	After 6 hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer
M0007	Malaria DNA Detection by Real Time PCR - Blood	Madhya Pradesh	C	3705	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
M0008	Malaria Detection by QBC - Blood	Madhya Pradesh	B	380	Daily : 9:00am, 1:00pm, 5:00pm	After 6 hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer
M0009	Malaria Profile reflex Malaria antigen by Immunochromatography - Blood	Madhya Pradesh	C	425	Daily : 9:00am, 1:00pm, 5:00pm	After 4 hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer
M0010	Malarial Antigen Vivax & Falciparum by Immunochromatography - Blood	Madhya Pradesh	B	540	Daily : 9:00am, 1:00pm, 5:00pm	After 6 hrs	3 ml of EDTA Whole Blood in EDTA Vacutainer
M0011	Manganese Blood	Madhya Pradesh	D	3380	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0012	Manganese Serum	Madhya Pradesh	D	3380	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0013	Manganese Urine	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0014	Manganese Urine 24H**	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0015	Mantoux test Tuberculin test (10TU)	Madhya Pradesh	O	225	Daily : 9:00am to 9:00pm	3rd day : 5pm	Mycobacteriology.pdf
M0015_2TU	Mantoux test Tuberculin test (2TU)	Madhya Pradesh	O	225	Daily : 9:00am to 9:00pm	3rd day : 5pm	Mycobacteriology.pdf
M0015_5TU	Mantoux test Tuberculin test (5TU)	Madhya Pradesh	O	225	Daily : 9:00am to 9:00pm	3rd day : 5pm	Mycobacteriology.pdf
M0019	MDMA (Qualitative) Ecstasy, XTC Urine Spot	Madhya Pradesh	D	855	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
M0020	Measles Virus IgG Antibody by CLIA - Serum	Madhya Pradesh	C	1905	Daily : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0021	Measles Virus IgG Antibody by EIA - CSF	Madhya Pradesh	C	3585	Tue & Fri : 7:00am	Next day : 5pm	2 ml of CSF in Sterile leak proof tube or container & 2 mL Serum
M0022	Measles Virus IgM Antibody by CLIA - Serum	Madhya Pradesh	C	1905	Daily : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0023	Medullary Thyroid Carinoma Mutation Detection Sequencing method Blood**	Madhya Pradesh	R	15690	Daily : 11:30am	10th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0024	MELAS-Mutation Detection Blood**	Madhya Pradesh	R	6120	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0026	Meningitis (viral) Profile IgG by CLIA - Serum	Madhya Pradesh	D	5135	Wed : 9:00am	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0027	Meningitis (viral) Profile IgG & IgM by CLIA - Serum	Madhya Pradesh	D	8455	Wed : 9:00am	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0028	Meningitis (viral) Profile IgM by CLIA - Serum	Madhya Pradesh	D	5190	Wed : 9:00am	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0029	Menopause	Madhya Pradesh	C	3250	Daily: 9am to 9pm	After 8 hrs	5 ml of serum
M0030	Menopause diagnostic profile	Madhya Pradesh	C	2015	Daily: 9am to 9pm	After 6 hrs	5 ml of serum
M0031	Menopause monitoring profile	Madhya Pradesh	C	2240	Daily: 9am to 9pm	After 6 hrs	5 ml of serum
M0032	Mercury Blood	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0033	Mercury Serum	Madhya Pradesh	D	3380	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0034	Mercury Urine	Madhya Pradesh	D	3640	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0035	Mercury Urine 24H**	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0036	MERRF Mutation Detection Blood**	Madhya Pradesh	R	10200	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0037	Neoxpert by GCMS**	Madhya Pradesh	R	7040	Daily : 11:30am	4th day : 6:00pm	Spot Urine:10 ml of Spot urine required with clinical history and clinician number.
M0038	Metabolic Screen	Madhya Pradesh	R	7140	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
M0039	Metanephrine- Free Plasma	Madhya Pradesh	D	5500	Tue & Fri : 11:30am	Next day : 5:00pm	Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
M0040	Metanephrine level, HPLC, 24 hrs urine	Madhya Pradesh	D	3500	Mon & Thu : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
M0041	Methadone Dolophine Urine Spot	Madhya Pradesh	B	740	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
M0042	Methamphetamine Qualitative Urine Spot	Madhya Pradesh	B	740	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
M0042	Methamphetamine Qualitative Urine Spot	Madhya Pradesh	B	680	Daily : 9:00am to 9:00pm	After 10 mins	Urine
M0043	Methaqualone Qualitative Urine Spot	Madhya Pradesh	B	820	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
M0044	Meth-Haemoglobin Blood	Madhya Pradesh	C	1120	Daily : 11:30am	Next day : 5:00pm	3 ml of Heparin Blood. Sample required within 12 to 24 hrs after collection & avoid holiday & weekends.
M0045	Methionine Quantitative Plasma*	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
M0046	Methionine Quantitative Urine SPOT**	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
M0047	Methotrexate Serum	Madhya Pradesh	C	3585	W, S: 7.30am	Same day at 6pm	3 ml of serum(sample collect in, plain vacutainer, not in gel tube)
M0048	Methyl Malonic Acid Qualitative Urine Spot	Madhya Pradesh	R	715	Daily : 11:30am	5th day : 6:00pm	15 ml of Spot Urine
M0049	Methyl Malonic Acid Quantitative Urine SPOT**	Madhya Pradesh	R	13260	Daily : 11:30am	5th day : 6:00pm	15 ml of Spot Urine
M0050	MGMT PCR Temozolomide Resistance Tissue**	Madhya Pradesh	R	16625	Daily : 11:30am	25th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
M0051	MICM By PCR Maternally Inherited Cardiomyopathy Blood**	Madhya Pradesh	R	7140	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0052	Microalbumin Urine 24H	Madhya Pradesh	B	600	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
M0053	Microalbumin / Creatinine Ratio Urine Spot	Madhya Pradesh	B	590	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
M0054	Microsporidium Detection by modified ZN Stain - Stool	Madhya Pradesh	B	600	Daily : 9:00am to 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
M0055	Mitochondrial (M2) Antibody by IFA Reflex to end point titre Serum	Madhya Pradesh	C	2130	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0059	Mixing Studies (APTT)	Madhya Pradesh	B	2725	Daily: 9am and 4pm	After 6 hrs	1 ml of citrated platelet poor plasma in frozen condition.
M0060	MMC-Mutation Detection Maternal Myopathy With Cardiomyopathy Blood**	Madhya Pradesh	R	7140	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0061	MM-Mutation Detection Mitochondrial Myopathy	Madhya Pradesh	R	8670	Daily : 11:30am	8th day : 6:00pm	Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0063	MPO by EIA, p-ANCASerum	Madhya Pradesh	C	1570	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0064	(cMPO)Myeloid cell marker	Madhya Pradesh	D	2460	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition (preferably). Clinical history and clinician details required.
M0065	MPO (Myeloperoxidase)Stain, Blood	Madhya Pradesh	C	895	Daily : 11:30am	3rd day : 6:00pm	3 Direct Peripheral smear slides, EDTA blood. Clinical history and clinician no is mandatory
M0066	MPO (Myeloperoxidase) Stain, Bone Marrow	Madhya Pradesh	C	1120	Daily : 11:30am	3rd day : 6:00pm	3 Direct Bone Marrow Aspirate slides. Clinical history and clinician no is mandatory
M0068	MRSA Screen - Nasal swab	Madhya Pradesh	C	1010	Daily : 9:00am to 9:00pm	2nd day : 5pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
M0069	MSUD-Maple Syrup Urine Disorder Neonatal Screen Blood	Madhya Pradesh	D	535	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
M0070	MTHFR C677T/A1298C Mutations	Madhya Pradesh	D	7275	Wed & Fri: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
M0071	Mucopolysaccharides (MPS) Type-I Blood**	Madhya Pradesh	R	5610	Daily : 11:30am	7th day : 6:00pm	10 ml sodium Heparin OR EDTA Blood
M0072	Mucopolysaccharides (MPS) Type-VI Blood**	Madhya Pradesh	R	2245	Daily : 11:30am	7th day : 6:00pm	10 ml sodium Heparin OR EDTA Blood
M0074	Mucopolysaccharides, urine	Madhya Pradesh	D	590	Daily : 11:30am	Next day : 5:00pm	10 ml of Spot urine
M0075	Multiple Sclerosis Profile, Maxi	Madhya Pradesh	D	10165	Tue and Fri: 9am	Next day 7pm	2 ml of CSF in leakproof container (foolin vu conical tube is recommended to prevent leakage). 3 ml of Serum.
M0076	Multiple Sclerosis Profile, Mini	Madhya Pradesh	D	6410	Tue and Fri: 9am	Next day 7pm	2 ml of CSF in leakproof container (foolin vu conical tube is recommended to prevent leakage). 3 ml of Serum.
M0078	Mumps virus IgG Antibody by CLIA - Serum	Madhya	C	1905	Daily : 9:00am to	After 6 hrs	3ml of Serum in Gel Vacutainer or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0079	Mumps virus IgM Antibody by CLIA - Serum	Pradesh Madhya Pradesh	C	1905	Daily : 11:30am	Next day : 5:00pm	in Plain Sterile Vacutainer 3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0080	MUSK Antibody Myasthenia gravis Serum	Madhya Pradesh	C	6945	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0081	Myasthenia Gravis Profile	Madhya Pradesh	D	11450	Mon: 9am	5th day	5 ml of serum
M0086	Mycoplasma Pneumoniae IgG Antibody by EIA - Serum	Madhya Pradesh	C	2410	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0087	Mycoplasma Pneumoniae IgM Antibody by EIA - Serum	Madhya Pradesh	C	2410	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
M0091	Myoglobin Serum	Madhya Pradesh	C	2630	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0092	Myoglobin Urine Spot	Madhya Pradesh	C	1345	Daily: 9am to 9pm	After 8 hrs	10 ml of Spot urine
M0095	Myotonic Dystrophy Type 1 (DM1) Mutation Analysis	Madhya Pradesh	E	7350	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
M0099	MRSA screen Methicillin Resistant Staphylococcus aureus - PCR	Madhya Pradesh	C	2240	Daily : 9:00am	Same day : 5:00pm	One swab from Nasal, Groin and Axilla.
M0100	Multiple myeloma workup panel	Madhya Pradesh	D	21345	Daily : 9:00am	7th day : 6:00pm	3 ml of EDTA Bone Marrow, 4ml of heparin Bone marrow, Serum, Direct smears
M0101	MPL W515/S505 Mutation Analysis, Blood	Madhya Pradesh	E	6150	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.
M0102	MPL W515/S505 Mutation Analysis, Bone marrow	Madhya Pradesh	E	6150	Wed & Sat : 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0103	Microsatellite Instability (MSI), HNPCC/ Lynch Synd	Madhya Pradesh	D	14980	Mon, Wed & Fri: 9:00am	10th day : 6:00pm	at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
M0104	MYCN Neuroblastoma (n MYC gene amplification)**	Madhya Pradesh	R	13260	Daily : 11:30am	12th Day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
M0108	Metachromatic Leukodystrophy gene**	Madhya Pradesh	R	20400	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
M0132	MISMATCH REPAIR (MMR) PROTEIN IMMUNOHISTOCHEMISTRY ONLY	Madhya Pradesh	D	9095		5th Day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological detail
M0136	Marfan s syndrome By NGS**	Madhya Pradesh	R	21420	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0140	Muscular Dystrophy & Congenital Myopathy gene panel**	Madhya Pradesh	R	17000	Daily : 11:30am	After 4 weeks : 6:00pm	Genetci history is mandatory. 3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
M0141	Methylmalonic acids Quantitative Serum**	Madhya Pradesh	R	3430	Thu : 11:30am	4th day : 6:00pm	3ml of serum
M0142	MET Amplification FISH	Madhya Pradesh	E	10000	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
M0145	Maturity-Onset Diabetes of the Young (MODY) & neonatal diabetes gene panel**	Madhya Pradesh	R	17000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
M0147	MLH1 Promoter Hypermethylation	Madhya Pradesh	E	7100	Mon, Wed & Fri: 9:00am	5th day : 6:00pm	â€¢ Tumor containing FFPE block (>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable.
M0149	MSH2 Sequencing and deletion duplication analysis	Madhya Pradesh	E	18000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0150	MMR IHC Reflex to BRAF CD600 and MHL1 promoter Hypermethylation	Madhya Pradesh	D	11770	Daily : 9:00am	After 2 weeks : 6:00pm	for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory. (FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical Diagnosis with Histopathology report, HistologyandIHC.pdf
M0151	Metanephrine- Fractionated Plasma	Madhya Pradesh	D	7490	Tue & Fri : 11:30am	Next day : 5:00pm	Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
M0153	Combine screening of Pre-eclampsia & Trisomies	Madhya Pradesh	D	3745	Daily cut off 7:30 AM	same day 9:00 PM	3 ml of serum is required. Mandatory Requirements: Latest Ultrasound (USG) report with a completely filled Maternal screening TRF form including CRL/BPD, nasal bone status, ductus venosus status, and NT value in millimeters. Please note that this test is not valid for triplets and above Preeclampsia: USG findings of mean arterial pressure (MAP) and uterine artery pulsatility index (UAPI), the patients blood pressure (BP), and height should be included as per maternal TRF requirements
M0154	Meningitis/Encephalitis panel by FilmArray (Biofire) - CSF	Madhya Pradesh	E	18500	Daily : 9:00am to 9:00pm	After 10 hrs	2ml of CSF in Sterile leak proof tube or container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0160	MMR IHC Reflex to MSI PCR	Madhya Pradesh	D	8455	Daily : 9:00am	After 2 weeks : 6:00pm	(FFPE block) Tumour containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical Diagnosis with Histopathology report, HistologyandIHC.pdf
M0161	Myositis Profile-IgG(16 Antigen)	Madhya Pradesh	D	10500	Daily: 11:30am	Next day 8pm	4 ml of serum
M0163	Magnesium, Hair tests by ICPMS**	Madhya Pradesh	R	5610	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
M0164	Manganese, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
M0165	Manganese, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
M0166	Marijuana confirmation by LCMS, Spot Urine**	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
M0167	Mercury, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
M0169	Metals profile-1, EDTA whole blood **	Madhya Pradesh	D	6420	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0170	Metals profile-1, Spot urine	Madhya Pradesh	D	6420	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0172	Metals profile-2, Spot urine**	Madhya Pradesh	R	11730	Daily : 11:30am	10th day : 6:00pm	HeavyMetalSamplingguidelines.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0174	Metals profile-3, Spot urine**	Madhya Pradesh	R	15810	Daily : 11:30am	10th day : 6:00pm	HeavyMetalSamplingguidelines.pdf
M0175	Methadone confirmation by LCMS, Spot Urine**	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
M0176	Methaqualone confirmation by LCMS, Spot Urine**	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
M0178	Minerals profile, Hair**	Madhya Pradesh	R	11730	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
M0179	Minerals profile, Spot Urine**	Madhya Pradesh	R	11730	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
M0181	Molybdenum, Hair tests by ICPMS**	Madhya Pradesh	R	5610	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
M0182	Molybdenum, Spot Urine**	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
M0183	Mucorales (Mucormycosis) DNA Detection by Real Time PCR - Body fluids/CSF/EDTA Blood/Pure culture/Sputum/Tissue	Madhya Pradesh	C	7280	Daily : 7:30am	3rd Day : 6:00pm	3ml of Body fluids/ CSF/ EDTA Blood/ Pure culture/ Sputum/ Tissue in Sterile leak proof container
M0193	MetroHealth Instant Glucose Fasting	Madhya	L0	100	Daily : 9:00am to	After 5 mins	Whole Blood

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M0194	MetroHealth Instant Glucose Post Prandial	Pradesh Madhya Pradesh	L0	100	Daily : 9:00am to 9:00pm	After 5 mins	Whole Blood
M0195	MetroHealth Instant Glucose Random	Madhya Pradesh	L0	100	Daily : 9:00am to 9:00pm	After 5 mins	Whole Blood
M0196	MetroHealth Instant Haemoglobin	Madhya Pradesh	L0	100	Daily : 9:00am to 9:00pm	After 5 mins	Whole Blood
M0197	MetroHealth Instant Urine	Madhya Pradesh	L0	150	Daily : 9:00am to 9:00pm	After 5 mins	Urine
M0823	Macro TSH (TSH, Free TSH,% Recovery) Serum	Madhya Pradesh	C	560	Daily: 9am to 9pm	Next day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
M0824	Anti diabetic drug contents	Madhya Pradesh	R	5610	Daily : 11:30am	10th day : 6:00pm	Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clinical symptoms if any. Reported in 7 days
M0825	Medicine test profile (steroids, heavy metals, anti diabetic, NSAIDS)	Madhya Pradesh	R	20400	Daily : 11:30am	10th day : 6:00pm	Medicine sample (10 grm), Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clinical symptoms if any. Reported in 7 days
M0826	NSAID contents	Madhya Pradesh	R	3060	Daily : 11:30am	10th day : 6:00pm	Medicine sample (10 grm), Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clinical symptoms if any. Reported in 7 days
M0829	Monogenic and syndromic obesity-Gene Panel	Madhya Pradesh	E	20000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
M0830	MTB Antibiogram (TB NGS) by Next Generation Sequencing - Pure Culture	Madhya Pradesh	R	8670	Daily : 11:30am	15th day : 6:00pm	Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
M0831	MPA (Mycophenolic Acid)	Madhya Pradesh	C	2520	Daily: 9am to 9pm	After 6 hrs	require serum (3ml) sample collected in fasting status preferably but not mandatory.
M5085	MPN Comprehensive Panel (Blood)	Madhya Pradesh	D	16000	Tue & Fri : 9:00am	7th day : 6:00pm	3ml whole blood in EDTA Vacutainer & 3ml Heparin blood
M5086	MPN Comprehensive Panel (Bone Marrow)	Madhya Pradesh	D	14445	Tue & Fri : 9:00am	7th day : 6:00pm	3ml Bone Marrow in EDTA Vacutainer & 3ml Bone Marrow
M5097	MTB triple combo (For pulmonary samples)	Madhya Pradesh	D	7170	Daily : 7:30am	After 6 weeks : 6:00pm	2ml of Pulmonary specimen in Sterile leak proof container
M5098	MTB triple combo (For extra pulmonary samples)	Madhya Pradesh	D	7170	Daily : 7:30am	After 6 weeks : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
M5768	MetroHealth Instant 5 min Screen Pkg	Madhya Pradesh	L0	350	Daily : 9:00am to 9:00pm	After 5 mins	Whole Blood-Finger Prick , 30 to 60 ml of Urine in a leak proof container.
M8027	Myotonic Dystrophy Type 2 (3q 21)	Madhya Pradesh	R	9180	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
M8228	MIC OF COLISTIN BY MICRO BROTH DILUTION	Madhya Pradesh	C	560	Daily : 9:00am to 9:00pm	3rd day : 6:00pm	Pure culture isolate is required with identification upto species level.
M8268	MTB Xpert XDR Rapid (For Pulmonary Samples)	Madhya Pradesh	D	3940	Daily : 7:30am	Same day : 6:00pm	2ml of Pulmonary specimen in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
M8269	MTB Xpert XDR Rapid (For Extra pulmonary Samples)	Madhya Pradesh	D	3940	Daily : 7:30am	Same day : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
M8346	MELD package	Madhya Pradesh	R	1225	Daily 9.00am to 9.00pm	Same day at 5.00pm	
N0001	N- Acetyl Glucosamine 6 Sulphate Sulphatase Blood**	Madhya Pradesh	R	8670	Daily : 11:30am	After 4 weeks : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
N0002	NARP Neurogenic Ataxia Retinitis Pigmentosa Blood**	Madhya Pradesh	R	6630	Daily : 11:30am	8th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
N0003	Natural Killer cells	Madhya Pradesh	D	4815	Daily : 11:30am	Next day : 5:00pm	6 ml EDTA whole Blood
N0006	Neisseria Gonorrhoeae Culture & Identification by Vitek/ MALDITOF - Urethral discharge	Madhya Pradesh	C	1290	Daily : 9:00am to 9:00pm	3rd day : 5pm	2 ml of Urethral Discharge in Sterile leak proof container
N0007	Neisseria Gonorrhoeae Culture & Identification Vitek/ MALDITOF - Vaginal discharge	Madhya Pradesh	C	1230	Daily : 9:00am to 9:00pm	3rd day : 5pm	2 ml of Vaginal Discharge in Sterile leak proof container
N0008	Neisseria Gonorrhoea DNA Detection by Real Time PCR - Genital swab	Madhya Pradesh	C	3810	Wed & Sat : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
N0009	Neisseria Gonorrhoea DNA Detection by Real Time PCR - Urethral discharge	Madhya Pradesh	C	3810	Wed & Sat : 9:00am	Next day : 5:00pm	2ml of Urethral Discharge in Sterile leak proof container
N0010	Neisseria Gonorrhoea DNA Detection by Real Time PCR - Urine	Madhya Pradesh	C	3810	Wed & Sat : 9:00am	Next day : 5:00pm	15ml of Spot Urine in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
N0011	Neisseria Meningitidis A,C,Y W135 Antigen detection by Agglutination - CSF	Madhya Pradesh	C	2855	Daily : 9:00am to 3:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
N0012	Neoxpert IAP (Newborn Screening 3 Analytes)	Madhya Pradesh	C	950	Daily 11:30am	Next Day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
N0013	Neoxpert 5 (Newborn Screening 5 Analytes)	Madhya Pradesh	C	1735	Daily 9:00am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
N0014	Neoxpert 8 (Newborn Screening 8 Analytes)	Madhya Pradesh	C	2800	Daily 9:00am	Next day : 5:00pm	<p>the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.</p> <p>Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.</p>
N0015	Neoxpert by TMS (Newborn Screening 40+ Analytes)	Madhya Pradesh	C	5490	Daily : 11:30am	5th day : 6:00pm	<p>Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between</p>

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
N0017	Neuronal Antibody Profile	Madhya Pradesh	E	9000	Daily: 11:30am	Next day 8pm	3 ml of serum
N0018	Nickel EDTA Whole Blood	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
N0019	Nickel Serum	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
N0020	Nickel, Urine Spot	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
N0021	Nicotine Metabolite Quantitative Serum	Madhya Pradesh	C	1570	Daily: 7.30am	Same Day 7.30pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0022	Nicotine Metabolite Quantitative Urine Spot	Madhya Pradesh	C	1570	Daily: 7.30am	Same Day 7.30pm	2 ml of Spot Urine
N0023_NIPS	NEXTGEN NON-INVASIVE PRENATAL SCREENING TEST	Madhya Pradesh	E	16000	Daily : 9:00am	8th day : 6:00pm	Sample collection require in special tube - Streak Tube - 10 ml. We Required filled TRF which include clinician stamp, Ultra Sound report (USG), Dual and Quatraple Marker Report, Photo

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Identification Proof.
N0027	NMDA Receptor Antibody (NR1), CSF	Madhya Pradesh	D	6635	Mon & Thu : 11:30am	Same day : 9:00pm	2 ml of CSF in sterile screw capped container.
N0028	NMDA Receptor Antibody (NR1) Autoimmune encephalitis Serum	Madhya Pradesh	D	6635	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0029	NMO (Aquaporin 4)Neuromyelitis Optica Antibodies Serum	Madhya Pradesh	C	5040	Tue & Fri : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0030	NMO (Aquaporin 4)Neuromyelitis Optica Antibodies CSF	Madhya Pradesh	C	5040	Tue & Fri : 11:30am	Same day : 9:00pm	2 ml of CSF in sterile screw capped container.
N0032	Nocardia Detection by smear examination - Sputum	Madhya Pradesh	B	740	Daily : 9:00am to 3:00pm	After 6 hrs	2ml of Sputum in Sterile leak proof container
N0033	Nocardia Detection by smear examination - Tissue	Madhya Pradesh	B	600	Daily : 9:00am to 3:00pm	After 6 hrs	Tissue with Normal Saline in Sterile leak proof container
N0034	Nor-Adrenaline (Nor-Epinephrine) Plasma	Madhya Pradesh	D	3315	Wed : 11:30am	Next day : 7:00pm	Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
N0035	Nor-Adrenaline (Nor-Epinephrine) Urine-HPLC	Madhya Pradesh	D	3850	Wed : 11:30am	Next day : 7:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
N0036	Nor-Metanephrine-Free Plasma	Madhya Pradesh	D	4800	Tue & Fri : 11:30am	Next day : 7:00pm	to be conveyed to lab. Plasma : 3 ml of EDTA plasma , Do not consume bananas, pineapple, chocolates, coffee, icecream, high diet in cereals and vanilla, potatoes and B complex 48 hours before collection. Provide CT scan , USG report and medication history.
N0037	Nor-Metanephrine, HPLC, urine 24 hrs	Madhya Pradesh	D	2460	Mon & Thu : 11:30am	Next day : 7:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
N0039	NPM1 mutations AML - Blood	Madhya Pradesh	D	4280	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.
N0040	NPM1 mutations AML - Bone Marrow	Madhya Pradesh	D	4280	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
N0041	NRAS mutations Exon-2 & 3 FFPE Tissue	Madhya Pradesh	D	7490	Mon, Wed & Fri: 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>10% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
N0042	NSE Neuron Specific Enolase CSF	Madhya Pradesh	D	3960	Tue & Fri: 11:30am	After 8hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
N0043	NSE Neuron Specific Enolase Serum	Madhya Pradesh	D	3425	Tue & Fri: 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0044	NT-Pro BNP Serum	Madhya Pradesh	C	2900	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
N0045	NTxN-Tellopeptide (Collagen Crosslinked)Urine 24H	Madhya Pradesh	C	3585	Wed : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
N0047	NMDA- CSF & Serum	Madhya Pradesh	D	11770	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests. 2 ml of CSF in sterile screw capped container.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
N0048	NMO- CSF & Serum	Madhya Pradesh	D	8025	Tue & Fri : 11:30am	Same day : 9:00pm	2 ml of CSF in sterile screw capped container.
N0056	Neuronal Antibody Profile-3	Madhya Pradesh	E	13900	Daily: 11:30am	Next day 8pm	3 ml of serum
N0057	Neoxpert Expanded Screening Panel (Newborn Screening 48+ Analytes)	Madhya Pradesh	D	6420	Daily : 11:30am	5th day : 6:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion history and sample to be collected post 72 hrs after transfusion.
N0058	Comprehensive Neurology panel**	Madhya Pradesh	R	22440	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
N0059_Serum	NMO with MOG Antibody Profile	Madhya	D	6900	Daily : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					ml serum for profiles and reflex tests.
N0060_CSF	NMO with MOG Antibody Profile	Madhya Pradesh	D	7060	Daily : 11:30am	Same day : 9:00pm	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
N0061	Nocardia Susceptibility by Microbroth dilution - Pure Culture	Madhya Pradesh	C	3640	Daily : 9:00am	15th day : 5pm	Pure growth of Nocardia spp. on Blood/Chocolate Agar plate. Kindly provide Previous identification report of Nocardia.
N0062	Nitrosonaphthol Test Qualitative, Spot Urine	Madhya Pradesh	R	20400	Daily : 11:30am	5th day : 6:00pm	10 ml of Spot sample in Leakproof container. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
N0064	Nickel, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
N0065	Nickel, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
N0066	NR4A3 (TEC) 9p22.3 translocation**	Madhya Pradesh	R	12240	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
N0067	NTRK3 and ETV6 translocation**	Madhya Pradesh	R	14280	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
N0068	Neurology comprehensive panel & whole mitochondrial genome sequencing	Madhya Pradesh	E	25000	Daily: 7:30am	5 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
N0069	Neurotransmitter disorders-Gene Panel	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
N0070	New born screening panel by NGS	Madhya Pradesh	E	16000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
N0096_C	Integrated test, standard	Madhya Pradesh	E	5500	Daily : 11:30am	7th day : 6:00pm	3 ml of Serum for all test. Serum to be shipped at 2 to 8 degree Refrigerated. Serum sample to be stored at 2 to 8 degree for 2 days and post 2 days to be send in frozen condition. Latest USG is recommended mentioning CRL or BPD, Single or twin fetus (test not valid for triplets or above), NT value in mm. Maternal screening TRF is mandatory with history of date of

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
N0096_S	Integrated test, sequential	Madhya Pradesh	E	5500	Daily : 11:30am	7th day : 6:00pm	<p>birth, weight, Race, Smoking, type 1 diabetes. In IVF pregnancy Egg donor DOB, egg extraction date and embryo transfer date. in Preeclampsia cases history of BP or family histroy additionally required with USG findings of MAP, UAPI values etc as per TRF requirements. Reflex test will be done for positive screen cases and will be conveyed by lab for sample requirements further based on consent with Gynecologist or pateint and TAT of report will be based on the test added and registered. For NIPS 10 ml of Whole blood in strek tube to be send in ambinet condition upto 1 day or refrigerated conditio upto 3 days.</p> <p>3 ml of Serum for all test. Serum to be shipped at 2 to 8 degree Refrigerated. Serum sample to be stored at 2 to 8 degree for 2 days and post 2 days to be send in frozen condition. Latest USG is recommended mentioning CRL or BPD, Single or twin fetus (test not valid for triplets or above), NT value in mm. Maternal screening TRF is mandatory with history of date of birth, weight, Race, Smoking, type 1 diabetes. In IVF pregnancy Egg donor DOB, egg extraction date and embryo transfer date. in Preeclampsia cases history of BP or family histroy additionally required with USG findings of MAP, UAPI values etc as per TRF requirements. Reflex test will be done for positive screen cases and</p>

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							will be conveyed by lab for sample requirements further based on consent with Gynecologist or pateint and TAT of report will be based on the test added and registered. For NIPS 10 ml of Whole blood in strek tube to be send in ambinet condition upto 1 day or refrigerated conditio upto 3 days.
N5543	NextGen NIPS Advanced	Madhya Pradesh	E	18500	Daily : 9:00am	8th day : 6:00pm	Sample collection require in special tube - Streak Tube - 10 ml. We Required filled TRF which include clinician stamp, Ultra Sound report (USG), Dual and Quatraple Marker Report, Photo Identification Proof.
N8209	NextGen Precision Assay (500 genes +TMB+MSI)**	Madhya Pradesh	R	65000	11.30 am	20 day 6:00pm	EDTA SAMPLE + Paraffin Block
O0001	Obesity Profile	Madhya Pradesh	C	6890	1st and 3rd Wed	Next day 5 pm	5 ml of serum, 2 ml of flouride plasma , 20 ml of spot urine
O0002	Occult blood by chemical test - Stool	Madhya Pradesh	A	165	Daily : 9:00am to 9:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
O0003	Occult blood by chemical test - Urine	Madhya Pradesh	A	100	Daily : 9:00am to 9:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
O0004	Oligoclonal band,CSF	Madhya Pradesh	C	5150	Tue & Fri: 11:30am	Next day : 5:00pm	3 ml of serum & 2 ml of CSF in sterile screw capped container.
O0005	Opiates Morphine Urine Spot	Madhya Pradesh	B	740	Daily: 9am to 9pm	After 6 hrs	5 ml of Spot Urine
O0006	Oppportunistic infection panel by Modified ZN smear - Stool	Madhya Pradesh	C	1790	Daily : 9:00am to 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
O0008	Ornithine Quantitative Plasma*	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
O0009	Ornithine Quantitative Urine 24H**	Madhya Pradesh	R	6325	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
O0010	Osmolality freezing point Serum	Madhya Pradesh	C	1010	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tests.
00011	Osmolality freezing point Urine	Madhya Pradesh	C	1010	Daily: 9am to 9pm	After 6 hrs	5 ml of Spot Urine
00012	Osmotic Fragility Test	Madhya Pradesh	C	950	Daily: 9am	Same day 5pm	3 ml Sodium heparin sample
00013	Osteocalcin BGP-Bone Gla Protein Serum	Madhya Pradesh	C	3695	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
00014	Osteomon Profile	Madhya Pradesh	D	3200	Mon and Thu: 8am	Same day 5pm	3 ml of serum and 3 ml of EDTA plasma
00015	Osteoporosis Profile-I	Madhya Pradesh	C	2185	Daily: 9am to 9pm	After 8 hrs	2 ml of EDTA Plasma
00016	Osteoporosis Profile-II	Madhya Pradesh	C	6550	Mon and Thu: 8am	Same day 5pm	3 ml of serum, 2 ml of EDTA plasma and 10 ml of spot urine
00017	Osteoporosis Profile-III	Madhya Pradesh	C	8850	Mon and Thu: 8am	Same day 5pm	3 ml of serum, 2 ml of EDTA plasma and 10 ml of spot urine
00018	OT Microbiological Surveillance Panel Bacterial & Fungal by Vitek/ MALDITOF - Per swab	Madhya Pradesh	C	560	Daily : 9:00am to 9:00pm	3rd day : 5pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
00020	Ovarian Antibody AOA serum	Madhya Pradesh	C	3975	Mon : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
00021	Ovarian Marker Profile	Madhya Pradesh	D	4335	Tue,Fri: 7:30am	After 8 hrs	5 ml of serum
00022	Ovarian Reserve Profile-Maxi	Madhya Pradesh	C	6890	Mon: 9am	Next day 5pm	5 ml of serum
00023	Ovarian Reserve Profile-Mini	Madhya Pradesh	C	3695	Mon and Thu: 9am	Same day 7pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
00024	Oxalate, Urine	Madhya Pradesh	C	2185	Daily : 11:30am	Same day : 9:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
00025	OT Microbiological Surveillance Panel Bacterial & Fungal - 5 swabs	Madhya Pradesh	C	2520	Daily : 9:00am to 9:00pm	3rd day : 5pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
00026	Oxcarbazepine, Serum**	Madhya Pradesh	D	5085	Fri : 11:30am	Next day : 5:00pm	3 ml of serum(sample collect in, plain vacutainer, not in gel tube)
00027	Ova and Parasite Examination,Other Samples	Madhya Pradesh	B	325	Daily : 9:00am to 3:00pm	Next day : 5pm	5gm of Stool specimen in Sterile leak proof container
00039	Open Neural tube defect screening, AFP MOM, Second trimester - Serum	Madhya Pradesh	C	1065	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
00040	Opiates confirmation by LCMS, Spot Urine**	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
00041	Oncotrack - Ultima [Liquid biopsy for 56 theranostic genes]**	Madhya Pradesh	R	49400	Daily : 11:30am	After 4 weeks : 6:00pm	Peripheral blood in 2 Streck tubes, (tube supplied upon request)
00042	NextGen Oncomine Focus Assay - 52 Gene Panel	Madhya Pradesh	E	40290	Daily 9:00 am	4 Weeks	Tumor containing FFPE block (>10% of tumor). Clinical history with Histopathology report.
01801	Other diagnostic services- 2 D Echo	Madhya Pradesh	O	2500			
01806	Other diagnostic services- ECG	Madhya Pradesh	O	275			

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
O1810	Other diagnostic services- TMT (Stress test)	Madhya Pradesh	O	1200			
O1818	OPG - Orthopantomogram	Madhya Pradesh	O	500			
O1901	CANINE-Occult blood Stool	Madhya Pradesh	A	130	Daily: 9am, 2pm, 8pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
O1902	Occult blood	Madhya Pradesh	A	110	Daily: 9am, 2pm, 8pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P0001	P1NP Intact N-Terminal Propeptide of Type I Procollagen Serum	Madhya Pradesh	C	1905	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0002	Pacer-26	Madhya Pradesh	C	2575	Daily: 9am to 9pm	After 8 hrs	5 ml of serum and 2 ml of F flouride plasma
P0003	Palmitoyl Protein Thioesterase1	Madhya Pradesh	R	5710	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
P0004	Pancreatic (Acute) Profile	Madhya Pradesh	D	1710	Daily: 9am to 9pm	After 8 hrs	3 ml of serum and 2 ml of EDTA whole blood
P0005	Pancreatic Marker Profile	Madhya Pradesh	D	3905	Daily: 9am to 9pm	After 8 hrs	5 ml of serum
P0006	Panfungal-DNA detection by PCR & Sanger Sequencing - Blood	Madhya Pradesh	D	10860	Daily : 7:30am	3rd day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer
P0007	Panfungal-DNA detection by PCR & Sanger Sequencing - Body fluids	Madhya Pradesh	D	10305	Daily : 7:30am	3rd day : 6:00pm	2ml of Body Fluid in Sterile leak proof container
P0008	Panfungal-DNA detection by PCR & Sanger	Madhya Pradesh	D	9630	Daily : 7:30am	3rd day : 6:00pm	2ml of CSF in Sterile leak proof

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Sequencing - CSF	Pradesh					tube or container
P0009	Panfungal-DNA detection by PCR & Sanger Sequencing - Sputum	Madhya Pradesh	D	10730	Daily : 7:30am	3rd day : 6:00pm	2ml of Sputum in Sterile leak proof container
P0010	Panfungal-DNA detection by PCR & Sanger Sequencing - Tissue	Madhya Pradesh	D	10730	Daily : 7:30am	3rd day : 6:00pm	Tissue with Normal Saline in Sterile leak proof container
P0011	PAP smear by Conventional method	Madhya Pradesh	C	560		2nd day	Pre fixed smears along with P/S findings
P0012	PAP smear by Liquid based cytology (LBC)	Madhya Pradesh	C	1200		3rd day	Cervical scraping in a LBC vial
P0013	PAPPa-Pregnancy Associated Plasma Protein Serum	Madhya Pradesh	D	1500	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0014	Paracetamol Acetaminophen Serum	Madhya Pradesh	D	1175	W, S: 7.30am	Same day at 6pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0015	Parainfluenza 1,2 & 3 IgG Antibody by EIA - Serum	Madhya Pradesh	C	2015	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0016	Parainfluenza 1,2 & 3 IgM Antibody by EIA - Serum	Madhya Pradesh	C	2015	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0017	Parathyroid Panel	Madhya Pradesh	C	3190	Daily: 9am to 9pm	After 8 hrs	2 ml of EDTA plasma, 3 ml of serum and spot sample urine
P0018	Parietal Cell Antibody GPCA by IFA Serum	Madhya Pradesh	C	2350	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0019	Parvovirus B19 DNA Detection by Real Time PCR - Amniotic fluid	Madhya Pradesh	D	3305	Mon, Thursday: 9:00am	Next day : 5:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
P0020	Parvovirus B19 DNA Detection by Real Time PCR - EDTA Blood	Madhya Pradesh	D	3305	Mon, Thursday: 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
P0021	Parvovirus B19 DNA Detection by Real Time PCR - CSF	Madhya Pradesh	D	3305	Mon, Thursday: 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
P0022	Parvovirus B19 IgG Antibody by EIA - Serum	Madhya Pradesh	C	2910	Wed : 9:00am	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0023	Parvovirus B19 IgM Antibody by EIA - Serum	Madhya Pradesh	C	2910	Wed : 9:00am	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0024	PAS Stain, Blood	Madhya	C	785	Mon & Thu :	3rd day : 6:00pm	3 Direct Peripheral smear slides,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh			11:30am		EDTA blood. Clinical history and clinician no is mandatory
P0025	PAS Stain, Bone Marrow	Madhya Pradesh	C	785	Mon & Thu : 11:30am	3rd day : 6:00pm	3 Direct Bone Marrow Aspirate slides. Clinical history and clinician no is mandatory
P0026	Paul bunnel test by Agglutination - Serum	Madhya Pradesh	B	1090	Daily : 9:00am & 3:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0028	PCP Phencyclidine Phosphate Urine Spot	Madhya Pradesh	B	740	Daily: 9am to 9pm	After 6 hrs	5 ml of Spot Urine
P0032	Peripheral smear examination (PBS)	Madhya Pradesh	B	200	Daily: 9am and 3pm	After 6 hrs	3 ml EDTA Whole Blood
P0033	pH by Indicator Strip - Stool	Madhya Pradesh	A	110	Daily : 9:00am to 9:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
P0034	pH by Indicator Strip - Urine	Madhya Pradesh	A	80	Daily : 9:00am to 9:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
P0035	Phenobarbitone Gardinal Serum**	Madhya Pradesh	C	1110	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0036	Phenyl Alanine Neonatal Screen Dried blood spot	Madhya Pradesh	B	490	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. DBS card to be send in Ambient temperature only . Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age, or as close to discharge as possible if discharged prior to 24 hours of age. Baby DOB, sample collection date and relevant clincial history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper.The blood should fill all printed circles and is applied evenly on one side of the filter

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Plz provide blood transfusion histroy and sample to be collected post 72 hrs after transfusion.Spot Urine:10 ml of Spot urine required with clinical history and clinician number.Amino acid Plasma: 3 ml of separated EDTA plasma required . Sample to be shipped at Frozen condition . Sample stable at refrigerated condition 2 to 8 degree for 1 day. Clinical history, metabolic disorders or inborn errors of metabolsim suspicion or provisional diagnosis of clincian required fasting recommended not mandatory. Amino acid supplementation or medications needs to stoped based on physician recommendations
P0037	Phenylalanine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
P0038	Phenylalanine Quantitative Spot Urine*	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
P0039	Phenytoin Eptoin, Dilantin Serum	Madhya Pradesh	C	1110	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0040	Pheochromocytoma Marker Profile-I	Madhya Pradesh	D	7385	Daily: 9am	7th day	20 ml of 24 hrs urine in HCL preservative. Quote Volume.
P0041	Pheochromocytoma Marker Profile-II	Madhya Pradesh	D	17655	Daily: 9am	7th day	4 ml of EDTA plasma and 20 ml of urine (Urine 24H)
P0050	Phospholipid profile IgG	Madhya Pradesh	C	5825	Thu: 1pm	Next day 5:00 pm	5 ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P0051	Phospholipid profile IgM	Madhya Pradesh	C	5825	Thu: 1pm	Next day 5:00 pm	5 ml of serum
P0052	Phospholipid Syndrome Profile	Madhya Pradesh	D	4500	Tuesday, Thursday and Saturday 7:30 am	same day 7pm	3 ml of Citrate plasma, 3 ml of serum
P0053	Phosphorus-Inorganic Serum	Madhya Pradesh	A	175	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0054	Phosphorus-Inorganic Urine Spot	Madhya Pradesh	A	325	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
P0055	Phosphorus-Inorganic Urine 24H	Madhya Pradesh	A	325	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
P0056	Phosphotidylserine IgG antibodySerum	Madhya Pradesh	C	1600	Thu: 1pm	Next day 5:00 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0057	Phosphotidylserine IgM antibodySerum	Madhya Pradesh	C	1600	Thu: 1pm	Next day 5:00 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0058	PIK3CA Gene Mutations - FFPE Tissue	Madhya Pradesh	E	6800	Mon, Wed & Fri: 9:00am	5th day : 6:00pm	Tumor containing FFPE block (>40% of tumor). Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours. If tissue provided then

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							tissue processing and FFPE block making charges will be applicable. Clinical history with Histopathology report.
P0059	Pituitary Marker Profile	Madhya Pradesh	C	4200	Daily: 9am to 9pm	After 8 hrs	3 ml of EDTA plasma, 3 ml of serum
P0060	PLA2 receptor antibody Primary membranous Glomerulopathy Serum	Madhya Pradesh	D	4920	Tue : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0061	Plate/Slant For Identification & Susceptibility - Automated Plate/Slant	Madhya Pradesh	C	1400	Daily : 9:00am to 9:00pm	3rd day : 5pm	Pure culture isolate is required
P0062	Platelet Antibodies by IFA Serum	Madhya Pradesh	D	9460	Mon : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0063	Platelet Count ,EDTA Blood	Madhya Pradesh	A	175	Daily: 9am to 9pm	After 6 hrs	3 ml EDTA Whole Blood
P0066	Pneumocystis Carinii Detection by IFA - BAL / Sputum	Madhya Pradesh	C	3305	Daily : 9:00am to 3:00pm	Next day : 5pm	2 ml of BAL / Sputum in Sterile leak proof container
P0070	Pneumonia Viral-III Panel by EIA - Serum	Madhya Pradesh	D	6955	Thurs: 7:30 am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
P0071	PNH confirmation test, Blood	Madhya Pradesh	D	7115	Daily: 9am	Next day 5pm	3 ml EDTA whole Blood
P0072	PNMA2 (MA2/TA) Antibody Serum	Madhya Pradesh	D	5135	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0073	Synovial Fluid by Polarizing microscopy	Madhya Pradesh	C	1010		4th day	0.5 ml Synovial fluid in a plain container
P0074	Pompe Disease Enzyme Assay Blood**	Madhya Pradesh	R	7855	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
P0075	Porphobilinogen (Quantitative), urine 24 Hrs	Madhya Pradesh	C	4255	Thu : 11:30am	Next day : 5:00pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
P0076	Porphyria Profile, 24 hrs Urine	Madhya Pradesh	D	6065	Thu : 11:30am	Next day : 5:00pm	Add 10 ml of 6 M HCL preservative in 24 hrs urine container. Collect 24 hrs urine sample and send 10 ml of 24 Hrs Urine In Dark Container.
P0077	Porphyria Profile, Urine Spot	Madhya Pradesh	D	6155	Thu: 9am	Next day 5pm	10 ml of 24 Hours Urine OR spot urine In Dark Container
P0078	Potassium ISE Serum	Madhya Pradesh	A	170	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0079	Potassium Urine 24H	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
P0080	Potassium, urine by ISE_P0080	Madhya Pradesh	A	295	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P0081	PR3 - c-ANCA Serum	Madhya Pradesh	C	1625	Tue,Thu, S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0082	PRA- Plasma Renin Activity	Madhya Pradesh	C	5040	Tue & Fri: 11:30am	Next day : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
P0083	Pre Albumin Serum	Madhya Pradesh	D	3530	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0086	Pre-Operative Profile (Maxi)	Madhya Pradesh	C	2015	Daily: 7pm	After 8 hrs	3 ml of EDTA whole blood, 2 ml of citrate plasma, 3 ml of serum, 20 ml of spot urine
P0087	Pre-Operative Profile (Mini)	Madhya Pradesh	C	1175	Daily: 8pm	After 8 hrs	3 ml of EDTA whole blood, 3 ml of serum, 20 ml of spot urine
P0088	Procalcitonin Serum	Madhya Pradesh	C	1500	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0089	Progesterone (P4) Serum	Madhya Pradesh	B	650	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0090	Prolactin Serum	Madhya Pradesh	B	600	Daily: 9am to 9pm	After 8 hrs	4 ml of Serum
P0091	Propoxyphene Qualitative Urine Spot	Madhya Pradesh	C	1010	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
P0092	Prostate Profile	Madhya Pradesh	D	2140	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0094	Protein C Activity Functional ,Citrate plasma	Madhya Pradesh	C	4650	Mon and Thurs: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
P0095	Protein C Antigen Citrated plasma (P0095)	Madhya Pradesh	C	5320	Daily: 7:30am	Same day : 9:00pm	1 ml of citrated platelet poor plasma in frozen condition.
P0096	Serum, Protein Electrophoresis	Madhya Pradesh	C	830	Daily:9.00 am	Same day 8:00 PM	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0097	Protein Electrophoresis, Urine 24 Hrs	Madhya Pradesh	C	950	Mon, Thu: 9am	Next day 5:00 pm	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
P0098	Protein Electrophoresis reflex immunofixation Serum	Madhya Pradesh	C	3350	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0099	Protein S activity, Plasma	Madhya Pradesh	C	4030	Mon and Thurs: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
P0100	Protein S Antigen (Free)Citrated plasma	Madhya Pradesh	C	6160	Mon and Thurs: 9am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
P0101	Proteins Serum	Madhya Pradesh	B	165	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0102	Proteins Urine 24H	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea ,

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
P0103	Proteins Ascitic Fluid	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	2 ml of Ascitic Fluid
P0104	Proteins CSF	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
P0105	Proteins Pedicardial Fluid	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	2 ml of Pericardial Fluid
P0106	Proteins Pleural Fluid	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	2 ml of Pleural Fluid
P0107	Proteins Synovial Fluid	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	2 ml of Synovial Fluid
P0108	Proteins/Creatinine Ratio	Madhya Pradesh	B	590	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
P0109	Prothrombin G20210A (Factor II) Mutation	Madhya Pradesh	D	5885	Wed & Fri: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
P0110	Prothrombin Time (PT)	Madhya Pradesh	B	290	Daily: 9am to 9pm	After 6 hrs	Local Labs can send whole Blood within 24 hrs of collection and for Outstations Frozen Plasma required.
P0111	PSA-Free Prostate Specific Antigen Serum	Madhya Pradesh	C	1010	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0112	PSA-Total Prostate Specific Antigen Serum	Madhya Pradesh	B	840	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0113	Pseudo Cholinesterase OPC poisoning Serum	Madhya Pradesh	C	1120	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0114	PTH-Intact Molecule Parathyroid hormone Serum/Plasma	Madhya Pradesh	B	1660	Daily: 9am to 9pm	After 6 hrs	2 ml of Serum / EDTA Plasma
P0115	PUO Profile	Madhya Pradesh	C	5600	Daily: 9am to 9pm	5th day	2 Blood culture bottle (Aerobic and anaerobic), 3 ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P0117	Pyruvate by Biochemical Method	Madhya Pradesh	C	2350	Daily: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0119	PSA with Free PSA ratio	Madhya Pradesh	C	1680	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0122	Pro Gastrin Releasing Peptide	Madhya Pradesh	R	12035	Daily: 7.30am	1 month	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0123	Prolactin, Active & macroprolactin serum	Madhya Pradesh	C	1010	Daily : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0124	PSA reflex Free PSA, Serum	Madhya Pradesh	C	1175	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0125	Porphobilinogen (Quantitative), Urine Spot	Madhya Pradesh	C	4255	Thu : 11:30am	Next day : 5:00pm	10 ml of Spot urine
P0126	Protein Body Fluid	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	2 ml body fluid
P0127	Potassium Body Fluid	Madhya Pradesh	A	295	Daily: 9am to 9pm	After 6 hrs	2 ml body fluid
P0128	Proteins, Peritoneal Fluid	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	2 ml body fluid
P0135	Neurotropic Virus panel by real time PCR	Madhya Pradesh	D	6420	Daily : 9:00am	3rd day : 6:00pm	2ml of CSF in Sterile leak proof tube or container
P0139	Pancreatic Elastase Test	Madhya Pradesh	D	6955	Wednesday-9.00AM	same day 5.00 PM	100 gms of stool in leak proof container. History of medicines, morning sample preferred.
P0149	PV panel (Jak-2 V617F Reflex to Jak-2 Exon 12 Mutation)-Blood	Madhya Pradesh	D	8560	Tue & Fri : 9:00am	8th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.
P0152	Pre-eclampsia screening Serum	Madhya Pradesh	D	3425	Mon & Thu : 11:30am	Next day : 6:00pm	3 ml of Serum for all test. Serum to be shipped at 2 to 8 degree Refrigerated. Serum sample to be

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							stored at 2 to 8 degree for 2 days and post 2 days to be send in frozen condition. Latest USG is recommended mentioning CRL or BPD, Single or twin fetus (test not valid for triplets or above), NT value in mm. Maternal screening TRF is mandatory with history of date of birth, weight, Race, Smoking, type 1 diabetes. In IVF pregnancy Egg donor DOB, egg extraction date and embryo transfer date. in Preeclampsia cases history of BP or family histroy additionally required with USG findings of MAP, UAPI values etc as per TRF requirements. Reflex test will be done for positive screen cases and will be conveyed by lab for sample requirements further based on consent with Gynecologist or pateint and TAT of report will be based on the test added and registered. For NIPS 10 ml of Whole blood in strek tube to be send in ambinet condition upto 1 day or refrigerated conditio upto 3 days.
P0154	Posaconazole Serum	Madhya Pradesh	R	6120	Mon & Thu : 11:30am	Next day : 7:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
P0155	Plasma Acylcarnitine profile**	Madhya Pradesh	R	5610	Daily : 11:30am	7th day : 6:00pm	3 ml of separated EDTA plasma required . Sample to be shipped at Frozen condition . Sample stable at refrigerated condition 2 to 8 degree for 1 day. Clinical history, metabolic disorders or inborn errors of metabolsim suspicion or provisional diagnosis of clinician required fasting recommended not mandatory. Amino acid

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							supplementation or medications needs to stoped based on physician recommendations
P0156	PV panel (Jak-2 V617F Reflex to Jak-2 Exon 12 Mutation)-Bone Marrow	Madhya Pradesh	D	8560	Tue & Fri : 9:00am	8th day : 6:00pm	3-5 ml of EDTA bone marrow required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription referring clinicians contact number and TRF mandatory.
P0158	PML/RARA - t(15;17) Qualitative by PCR*	Madhya Pradesh	D	5245	Daily : 9:00am	4th day : 6:00pm	5-10 ml of EDTA whole blood required. Sample stable for 3 Days at 2 to 8 degree celsius. Clinical diagnosis, test prescription, referring clinicians contact number and TRF mandatory. Sample should reach to the lab within 48 to 72 hours.
P0159	Phencyclidine confirmation by LCMS, Spot Urine**	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
P0160	Platinum, EDTA Whole Blood**	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
P0161	Platinum, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
P0162	Platinum, Spot Urine**	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
P0163	Propoxyphine confirmation by LCMS, Spot Urine**	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	10ml of spot urine sample required to be shipped at 2-8C. The sample is stable for 1day, 1 day & 30days

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P0165	POLG-Gene Analysis	Madhya Pradesh	E	17000	Daily: 7:30am	4 weeks	<p>at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.</p> <p>3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.</p>
P0166	Pulmonary surfactant metabolism dysfunction-Gene Panel	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	<p>3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.</p>
P0167	PGS (Preimplantation genetic screening), Embryo	Madhya Pradesh	R	22440	Daily : 11:30am	15th day : 6:00pm	CytogeneticsandFISH.pdf
P0168	Pre screening for Bleeding Disorders	Madhya Pradesh	C	4705	Daily : 9 am	Next day 5pm	1 ml of citrated platelet poor plasma in frozen condition.
P0169	Panfungal-DNA detection by PCR & Sanger Sequencing - Pure Culture	Madhya Pradesh	D	9630	Daily : 7:30am	3rd Day : 6:00pm	Pure culture on SDA plate/ slant
P1002	Plasminogen Activator Inhibitor-1 (PAI-1) / SERPINE-1, 4G/5G Genotyping	Madhya Pradesh	R	14585	Wed : 11:30am	10th day : 6:00pm	<p>Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped</p>

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							at Amibent or referigerated condition based on sample stability.
P1003	PCOS Profile - 1	Madhya Pradesh	D	3600	Daily : 7:30am	Same day : 5:00pm	5 ml of serum and 2 ml of F, PP flouride plasma
P1004	PCOS Profile - 2	Madhya Pradesh	D	7060	Tue and Fri: 9am	Same day 7pm	2 ml of Fluoride plasma (F, PP), 5 ml of fasting serum sample
P1005	PCOS profile - 3	Madhya Pradesh	D	10595	Tue and Fri: 9am	Same day 7pm	2 ml of Fluoride plasma (F, PP), 6 ml of fasting serum sample
P1008	Prader-Willi / Angelman Syndrome by MS-PCR	Madhya Pradesh	D	7490	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
P1009	Prolactin (Tripooled Sample) serum	Madhya Pradesh	D	640	Daily: 9am to 9pm	After 8 hrs	3 ml of serum(0, 20, 40 min)
P1017	Phenol, Urine**	Madhya Pradesh	R	2925	Daily : 11:30am	20th day : 6:00pm	31 ml Urine (collect Urine after end of work shift on last working day of shift)
P1901	CANINE-Peripheral smear examination	Madhya Pradesh	B	380	Daily: 9am, 3pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1902	Canine - Phosphorous	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1903	Platelet Count, EDTA whole blood	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1904	Canine - Potassium	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1905	Pre ? Operative Panel	Madhya Pradesh	C	1010	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1906	Progesterone (P4), Serum by CMIA	Madhya Pradesh	B	710	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1907	Protein Electrophoresis, Serum	Madhya Pradesh	C	860	Daily: 7am	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1908	Canine - Proteins, Serum	Madhya Pradesh	C	245	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
P1909	Prothrombin Time (PT), Citrate plasma	Madhya Pradesh	B	415	Daily: 9am to 9pm	Report after 24 hrs	for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P1910	Canine-PTH	Madhya Pradesh	B	1745	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
P5629	Prega-Plus	Madhya Pradesh	E	8599	Daily 9 am	Upto 72 hours	2 ml of EDTA whole blood in 2 vacutainers, , 5 ml of serum, 2 ml of flouride plasma , 20 ml of spot urine, 10 ml of urine in Sterile Container .(Mid-stream or early morning) . Upon 1st visit sample for all other tests is collected Patient has to be then given 82.5 grams of commercial glucose (75 gram anhydrous Glucose) and Sample has to be collected exactly after two hours.Sample for Pregascreen Dual is collected between 11-13 weeks gestation. Tests for partner to be registered separately under P0630.)
P5631	Prega-Pro	Madhya Pradesh	E	5599	Daily 9 am	Next day 5pm	2 ml of EDTA whole blood in 2 vacutainers, , 5 ml of serum, 2 ml

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							of flouride plasma , 20 ml of spot Urine 1st visit sample for all other tests is collected Patient has to be then given 82.5 grams of commercial glucose (75 gram anhydrous Glucose) and Sample has to be collected exactly after two hours.Sample for Pregascreen Dual is collected between 11-13 weeks gestation...
P8220	Pneumocystis jirovecii (PCP) Detection	Madhya Pradesh	D	7500	Every Tuesday, Thursday and Saturday. 9AM	Same day 11pm	3ml Sputum or BAL
P8229	PLA2 receptor antibody Quantitative , Serum	Madhya Pradesh	D	5000	Tue : 11:30am	Next day : 5:00pm	2 ml Serum
P8277	Liquid Biopsy PIC3CA NGS Mutation	Madhya Pradesh	E	45000	Daily 7.00 am	25 th days	Peripheral Blood in streck tube
P8530	PDL-1 (22C3) Head & Neck squamous cell carcinoma	Madhya Pradesh	D	15020	Daily 9 am	7th day	Paraffin Block
P8531	PDL-1 (22C3) Breast	Madhya Pradesh	D	15020	Daily 9 am	7th day	Paraffin Block
P8532	PDL-1 (22C3) Cervix	Madhya Pradesh	D	15020	Daily 9 am	7th Day	Paraffin Block
P8533	PDL-1 (22C3) Lung	Madhya Pradesh	D	15020	Daily 9 am	7th Day	Paraffin Block
P8534	PDL-1 (22C3) Urothelial	Madhya Pradesh	D	15020	Daily 9 am	7th Day	Paraffin Block
Q8101	Detection of Fetal Aneuploidies (13, 18, 21 and XY) by QF PCR	Madhya Pradesh	E	5100	Daily : 9:00am	Next day : 6:00pm	â€¢ Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition or â€¢ Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Refrigerated condition â€ 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . â€ çPOC - 1)Upto 12 weeks of gestation: 50-100 mg placental villi (fetal side) in 20 ml normal saline with 2-3 drops of gentamycin in a sterile container 2) In 12-18 weeks of gestation: a) 50-100 mg placental villi (fetal side) in 20ml normal saline with 2-3 drops of gentamycin in a sterile container b) and or fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container 3)IUFD >20 weeks of gestation: fetal Muscle tissue preferably medial aspect of thigh muscle or a fetal toe in sterile container. if possible in 20ml normal saline with 2-3 drops of gentamycin
R0001	RA test Rheumatoid Arthritis Serum	Madhya Pradesh	B	450	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0002	RA test Rheumatoid Arthritis Synovial Fluid	Madhya Pradesh	B	765	Daily: 9am to 9pm	After 6 hrs	1ML of Synovial Fluid
R0004	Rabies virus antibodies by CLIA - Serum	Madhya Pradesh	D	3850	Mon: 7:30 am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0005	RA-IgA Rheumatoid Arthritis Serum	Madhya Pradesh	R	865	Daily : 11:30am	7th day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0007	RBC Folate EDTA Whole Blood	Madhya Pradesh	C	2295	Daily : 9:00am	Same day after 8 hrs	EDTA Blood (2 ml) And Serum (2 ml)
R0009	Reducing Substances Paediatric - Urine	Madhya	A	110	Daily : 9:00am to	After 6 hrs	15ml of Spot Urine in Sterile leak

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh			9:00pm		proof container
R0010	Reducing Substance by chemical test - Stool	Madhya Pradesh	A	220	Daily : 9:00am to 9:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
R0011	Renal (Kidney) Function Tests-1 RFT Maxi	Madhya Pradesh	C	2000	Daily: 9am to 9pm	After 6 hrs	3 ml of EDTA whole blood, 2 ml of Flouride plasma, 3 ml of serum, 20 ml of spot urine
R0012	Renal (Kidney) Function Tests	Madhya Pradesh	C	600	Daily: 9am to 9pm	After 6 hrs	3 ml of serum, 20 ml of spot urine
R0012_WithouteGFR	Renal (Kidney) Function Tests-2 RFT Mini	Madhya Pradesh	C	1175	Daily: 9am to 9pm	After 4 hrs	3 ml of serum, 20 ml of spot urine
R0013	Renin Direct Plasma	Madhya Pradesh	C	3810	Daily 9:00 am	Same day 5:00 pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
R0014	Pneumonia (Atypical) panel-IgG Respiratory Tract Infection Profile-IgG by EIA - Serum	Madhya Pradesh	D	6955	Thurs: 7:30 am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0015	Pneumonia (Atypical) panel-IgM Respiratory Tract Infection Profile-IgM by EIA - Serum	Madhya Pradesh	D	6955	Thurs: 7:30 am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0018	Rett Syndrome MECP2 Sequencing Blood	Madhya Pradesh	E	8150	Daily : 9:00am	10th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
R0020	Rh Antibody, Serum	Madhya Pradesh	B	980	Daily: 9am and 5pm	After 6 hrs	3 ml Serum sample
R0021	Rheumatic Fever Panel	Madhya Pradesh	C	1850	Daily: 9am to 9pm	3rd day	3 ml of EDTA blood and serum. And Throat swab
R0022	Rheumatoid Arthritis Panel-1	Madhya Pradesh	C	4030	Wed and Sat: 7:30am	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0023	Rheumatoid Arthritis Panel-2	Madhya Pradesh	C	5710	Wed and Sat: 7:30am	Same day 7 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0024	Rheumatoid Arthritis Panel-3	Madhya	C	8680	Wed and Sat:	Same day 7 pm	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh			7:30am		ml serum for profiles and reflex tests.
R0025	Rheumatology Profile	Madhya Pradesh	C	6610	Wed and Sat: 7:30am	Same day 7 pm	5 ml of serum
R0026	Ri Antibody, Serum	Madhya Pradesh	D	4815	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0027	RNP-Sm Antibody Serum	Madhya Pradesh	C	1960	W,S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0028	Rota virus Antigen detection - Stool	Madhya Pradesh	C	1455	Daily : 9:00am & 3:00pm	After 6 hrs	5gm of Stool specimen in Sterile leak proof container
R0030	Routine examination Comprehensive examination by Microscopy - Stool	Madhya Pradesh	A	220	Daily : 9:00am & 3:00pm	Next day : 5pm	5gm of Stool specimen in Sterile leak proof container
R0031	Routine examination - Body fluid	Madhya Pradesh	C	450	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Body Fluid in Sterile leak proof container
R0032	Routine examination - Ascitic fluid	Madhya Pradesh	C	450	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Ascitic Fluid in Sterile leak proof container
R0032_2	Routine examination - Ascitic fluid - 2	Madhya Pradesh	C	450	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Ascitic Fluid in Sterile leak proof container
R0033	Routine examination - CSF	Madhya Pradesh	C	450	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of CSF in Sterile leak proof tube or container
R0033_2	Routine examination - CSF	Madhya Pradesh	C	1400	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of CSF in Sterile leak proof tube or container
R0034	Routine examination (Pericardial Fluid)	Madhya Pradesh	C	450	Daily : 9:00am & 3:00pm	Next day : 10am	2 ml of Pericardial Fluid in Sterile leak proof container
R0035	Routine examination (Peritoneal Fluid)_R0035_N	Madhya Pradesh	C	450	Daily : 9:00am & 3:00pm	Next day : 10am	2 ml of Peritoneal Fluid in Sterile leak proof container
R0036	Routine examination (Pleural Fluid)_R0036_N	Madhya Pradesh	C	450	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Pleural Fluid in Sterile leak proof container
R0036_2	Routine examination (Pleural Fluid)-2	Madhya Pradesh	C	1555	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Pleural Fluid in Sterile leak proof container
R0037	Routine examination (Prostatic Fluid)	Madhya Pradesh	C	450	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Prostatic Fluid in Sterile leak proof container
R0038	Semen Analysis by Microscopy - Semen	Madhya Pradesh	B	600	Daily : 9:00am to 9:00pm	After 6 hrs	2ml of Semen in Sterile leak proof container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
R0039	Routine examination Sputum	Madhya Pradesh	B	435	Daily : 9:00am & 3:00pm	Next day : 5pm	2ml of Sputum in Sterile leak proof container
R0040	Routine Examination of Synovial Fluid	Madhya Pradesh	C	450	Daily : 9:00am & 3:00pm	Next day : 10am	2ml of Synovial Fluid in Sterile leak proof container
R0041	Routine Examination Urine	Madhya Pradesh	A	110	Daily : 9:00am & 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
R0042	RSV-IgG Antibody by EIA - Serum	Madhya Pradesh	C	2350	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0043	RSV-IgM Antibody by EIA - Serum	Madhya Pradesh	C	2350	Thurs : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0044	Rubella Virus Avidity test by EIA - Serum	Madhya Pradesh	C	1455	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0045	Rubella Virus IgG and IgM Antibody by CLIA - Serum	Madhya Pradesh	C	1345	Daily : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0047	Rubella Virus IgG Antibody by CLIA - Serum	Madhya Pradesh	B	720	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0048	Rubella Virus IgM Antibody by CLIA - Serum	Madhya Pradesh	B	720	Daily : 9:00am to 9:00pm	After 8 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
R0049	Rubella Virus RNA Detection by PCR - Amniotic Fluid	Madhya Pradesh	C	3695	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
R0050	Rubella Virus RNA Detection by PCR - EDTA Blood	Madhya Pradesh	C	3695	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	3 ml of Whole Blood in EDTA Vacutainer
R0051	Rubella Virus RNA Detection by PCR - CSF	Madhya Pradesh	C	3695	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	2ml of CSF in Sterile leak proof tube or container
R0052	Rubella Virus RNA Detection by PCR - CVS	Madhya Pradesh	C	3695	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	2ml of CVS in Sterile leak proof tube or container
R0053	Rubella Virus RNA Detection by PCR - Urine	Madhya Pradesh	C	3695	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	15ml of Spot Urine in Sterile leak proof container
R0055	RVVT-Russell Viper Venom Time	Madhya Pradesh	C	1120	Daily: 9am and 3pm	After 6 hrs	1 ml of citrated platelet poor plasma in frozen condition.
R0056	Retinitis Profile by Real Time PCR - Aqueous/Vitreous Fluid	Madhya Pradesh	D	10700	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	3ml of Aqueous/ Vitreous Fluid in Sterile leak proof container
R0059	Retinol binding protein**	Madhya Pradesh	R	2855	Daily : 11:30am	15th day : 6:00pm	3ml of serum
R0060	Ricketssia DNA Detection by Real Time	Madhya Pradesh	D	5670	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
R0063	PCR - Blood ROS1 Gene Rearrangement by FISH	Pradesh Madhya Pradesh	D	12305	Daily : 9:00am	7th day : 6:00pm	Vacutainer (FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
R0066	Respiratory panel by FilmArray (Biofire) - Nasopharyngeal swab	Madhya Pradesh	E	16800	Daily : 9:00am to 9:00pm	After 10 hrs	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
R0067	RET Gene Rearrangement by FISH, Tissue	Madhya Pradesh	E	9000	Daily : 9:00am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
R0069	RNA Polymerase III Antibodies, IgG**	Madhya Pradesh	C	3695	Tue : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
R0077	Allergy panel - Common Rhinitis/ Wheeze	Madhya Pradesh	D	7490	Daily : 11:30am	Next day : 5:00pm	All tests require serum (3ml) sample collected in fasting status preferably but not mandatory. The samples should be shipped at R (2-8C) and can be stored at 2-8C for 7 days and at Ambient temperature for a day.
R0080	Renal biopsy profile (LM, IF, EM)	Madhya Pradesh	R	11500	Daily : 9:00am	20 day	1) Tissue in 10% NB formalin 2) Tissue in Michel Medium. Please use this code when you are sending 2 sample vials ; One in Formaldehyde & other in Michels medium. For Electron microscopy please send 3% glutaraldehyde
R0086	Rifampicin**	Madhya Pradesh	D	6420	Daily : 11:30am	10th day : 6:00pm	2 ml K2-EDTA plasma
R1004	Reverse T3**	Madhya	R	22950	Daily: 7.30am	1 month	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					ml serum for profiles and reflex tests. International outsource : Category R for referred and I for international. Plz esquire before sending sample for arrangement of proper courier to international lab. Additional sample tubes may be needed for sample adequacy. TAT may vary based on information by concerned lab of delay. Courier charges have to be paid extra based on sample temperature condition and shipment.
R1008	Ribosome antibodies, serum**	Madhya Pradesh	R	7140	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
R1901	Reticulocyte panel - Automated, EDTA whole blood	Madhya Pradesh	B	435	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
R1902	CANINE-Routine Examination Stool	Madhya Pradesh	A	165	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
R1903	Routine examination reflex to cytology, Body fluids	Madhya Pradesh	C	895	Daily: 9am, 3pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
R1904	Canine - Routine Examination Urine	Madhya Pradesh	A	130	Daily : 9:00am to	Report after 24 hrs	10 ml of Spot urine

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
R5588	Reassure BOH Mini	Pradesh Madhya Pradesh	D	6315	Daily 7.00am	Reporting 3rd day	5 ml serum, 3 ml citrate plasma
R5589	Reassure BOH Maxi	Madhya Pradesh	D	11770	Daily 7 pm,	Reported 15th day	5 ml serum, 4 ml Heparin whole Blood (Husband and wife), 3 ml citrate plasma
S0001	Spinocerebellar Ataxia PANEL (SCA1, 2, 3, 6, 7, 12 and 17)	Madhya Pradesh	E	15000	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0002	Spinocerebellar Ataxia (SCA1) DNA TEST	Madhya Pradesh	D	3265	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0003	Spinocerebellar Ataxia 2 (SCA2) DNA TEST	Madhya Pradesh	D	3265	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0004	SCA3 (Machado-Joseph Disease) DNA Test	Madhya Pradesh	D	3265	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0005	Spinocerebellar Ataxia 6 (SCA6) DNA TEST	Madhya Pradesh	D	3265	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0006	Spinocerebellar Ataxia TYPE 7 (SCA7) DNA TEST	Madhya Pradesh	D	3265	Mon & Wed: 9:00am	7th day : 6:00pm	clinicians contact number and TRF mandatory. 3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0007	Spinocerebellar Ataxia Type 12 (SCA12) DNA Test	Madhya Pradesh	D	3305	Mon & Wed: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0008	Scl-70 Antibody Serum	Madhya Pradesh	C	1850	W,S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0009	Scleroderma Panel	Madhya Pradesh	C	6610	Wed/Sat - 7:30 am	Next day 5pm	5 ml of serum
S0010	Selenium, Blood	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
S0011	Selenium, Serum	Madhya Pradesh	D	2945	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
S0012	Selenium, Urine	Madhya Pradesh	D	3640	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
S0013	Selenium Urine 24H**	Madhya Pradesh	D	3640	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
S0014	Serine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
S0015	Serine Quantitative Urine 24H**	Madhya Pradesh	D	6635	Daily : 11:30am	7th day : 6:00pm	10 ml of Spot urine
S0016	Serotonin5-Hydroxy Tryptamine Serum	Madhya Pradesh	D	4710	Wed : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0017	Serotonin5-Hydroxy Tryptamine,Urine	Madhya	D	5415	Wed : 11:30am	Next day : 5:00pm	10 ml of 24 hrs urine (Preferred) /

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0018	SGOT AST Serum	Pradesh Madhya Pradesh	A	160	Daily: 9am to 9pm	After 6 hrs	Spot urine 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0019	SGPT ALT Serum	Madhya Pradesh	A	170	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0020	SHBG Sex hormone binding globulin Serum	Madhya Pradesh	C	2575	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0021	Sickle Cell Anemia Mutation Analysis	Madhya Pradesh	D	3745	Wed & Fri: 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Clinical history, referring clinicians contact number, HPLC Report, CBC findings, clinical/family history and TRF mandatory.
S0022	Sickling Test	Madhya Pradesh	B	325	Daily : 9:00am to 9:00pm	after 6 hrs	3 ml EDTA whole blood
S0024	Silver ICPMS Urine	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
S0026	SIROLIMUS, BLOOD BY LC MSMS	Madhya Pradesh	C	6160	Tue & Fri: 11:30am	Next day : 5:00pm	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
S0027	Sjogrens Syndrome Profile	Madhya Pradesh	D	2675	Wed and Sat 7:30am	Next day 5 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0028	Anti-SLA/LP, IgG	Madhya Pradesh	C	3135	Tues/fri 9am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0029	SLE profile	Madhya Pradesh	C	5375	Sat - 7:30 am	Next day 5 pm	3 ml of EDTA whole blood , 3 ml of serum
S0030	Sm Antibody Smith Antibody Serum	Madhya Pradesh	C	1735	W,S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0031	SMA 12+2	Madhya Pradesh	D	1660	Daily: 9am to 9pm	After 8 hrs	2 ml of flouride plasma , 3 ml of serum
S0032	Sodium ISE Serum	Madhya Pradesh	A	200	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0033	Sodium Urine	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	2 ml of Spot Urine
S0034	Sodium Urine 24H	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
S0035	Sodium Body Fluid	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	2 ml of Body Fluid
S0036	Soluble Transferrin Receptor Serum	Madhya Pradesh	C	1175	Daily: 9am to 9pm	3rd day	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0037	Specific gravity by Fully Automated Urine chemistry & Microscopy - Urine	Madhya Pradesh	A	80	Daily : 9:00am to 9:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
S0038	Spectral Karyotyping (SKY) (m-FISH) 24 Color Painting Blood**	Madhya Pradesh	D	18725	Daily: 9am	15th day	3 ml of Blood In Na-Heparin Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0039	Spectral Karyotyping (SKY) (m-FISH) 24 Color Painting Bone Marrow**	Madhya Pradesh	R	16320	Daily : 11:30am	15th day : 6:00pm	2-3ml Bone marrow/3ml Blood in Sodium Heparin Vacutainer
S0040	Sperm Antibody-Total ASAB serum	Madhya Pradesh	C	1345	Tue,Fri: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0041	Sphingomyelinase Niemann Pick Disease Blood**	Madhya Pradesh	R	7650	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
S0044	Spinal Muscular Atrophy Microdeletion By PCR SMN1 Gene Deletion	Madhya Pradesh	D	3425	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0045	Sry Deletion FISH	Madhya Pradesh	D	5670	Daily : 9:00am	12th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
S0046	SSA-Antibody Ro Serum	Madhya Pradesh	C	1850	W,S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0047	SSB-Antibody La Serum	Madhya Pradesh	C	1850	W,S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0048	STD Profile -I	Madhya Pradesh	D	3200	Tue,Thurs & Sat : 7:30am	5th day : 5pm	3 ml of serum
S0049	STD Profile -II	Madhya Pradesh	D	5350	Daily : 8:00pm	5th day : 5pm	4 ml Serum
S0050	STD Profile -III	Madhya Pradesh	D	8880	Daily : 9:00am	5th day : 6:00pm	6 ml of EDTA whole blood, 15ml Urine in sterile container. For HPV : Swab, Tissue for male and LBC, Swab, Tissue for female
S0051	Steroid analysis From medicine **	Madhya Pradesh	R	11220	Daily : 11:30am	10th day : 6:00pm	Medicine sample (10 grm)
S0052	Stomach Marker Profile	Madhya Pradesh	D	3800	Tue and Fri: 9am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0053	Calculus (Stone) Analysis by automated	Madhya	D	1390	Daily : 11:30am	Next day : 5:00pm	Collect sample in dry plain

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	FTIR (Gall bladder calculus)	Pradesh					container, mention stone type and history
S0054	Calculus (Stone) Analysis by automated FTIR (Kidney (Renal) calculus)	Madhya Pradesh	D	1390	Daily : 11:30am	Next day : 5:00pm	Collect sample in dry plain container, mention stone type and history
S0055	Calculus (Stone) Analysis by automated FTIR (Urine calculus)	Madhya Pradesh	D	1390	Daily : 11:30am	Next day : 5:00pm	Collect sample in dry plain container, mention stone type and history
S0056	Stone screening profile	Madhya Pradesh	D	2300	Daily: 9:00am	5th day : 6:00pm	24 hours urine sample with appropriate preservative.
S0058	Streptococcus group B Antigen detection by Agglutination - CSF	Madhya Pradesh	C	1455	Daily : 9:00am to 3:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
S0059	Streptococcus group B Antigen detection by Agglutination - Serum	Madhya Pradesh	C	1455	Daily : 9:00am to 3:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0060	Streptococcus pneumoniae Antigen detection by Agglutination - CSF	Madhya Pradesh	C	1455	Daily : 9:00am to 3:00pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
S0061	Streptococcus pneumoniae Antigen detection by Agglutination - Serum	Madhya Pradesh	C	1455	Daily : 9:00am to 3:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0062	Streptococcus pneumoniae Antigen detection by Agglutination - Urine	Madhya Pradesh	C	1455	Daily : 9:00am to 3:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
S0063	Sucrose Lysis Test (PNH Screening tests) (S0063)	Madhya Pradesh	B	870	Daily: 9am	Next day 5pm	3 ml of EDTA Whole Blood And 3 ml Serum
S0066	Glucose Sugar - Urine	Madhya Pradesh	A	80	Daily : 9:00am to 9:00pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
S0067	Syphilis Antibodies Treponema by EIA - Serum	Madhya Pradesh	C	1345	Tue & Fri: 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0070	Scrub typhus profile Weil felix and scrub typhus-IgM by EIA - Serum	Madhya Pradesh	C	1230	Tue & Fri: 11:30am	Next day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0071	Scrub typhus IgM by EIA - Serum	Madhya Pradesh	C	1065	Tue, Thu & Sat: 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S0074	SMA Carrier Detection Genetic Disorders**	Madhya Pradesh	R	9180	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0075	SCA17 PCR-Fragment analysis	Madhya Pradesh	D	3425	Mon & Wed: 9:00am	7th day : 6:00pm	EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory. 3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical/family history, referring clinicians contact number and TRF mandatory.
S0077_Sperm-DFI	SPERM DNA FRAGMENTATION (DFI)	Madhya Pradesh	D	5300	Daily : 11.30 am	3rd day 1 pm	2-3 ml of freshly collected Semen sample in a sterile leak proof container. If reaching the lab within 4hrs of collection, transported at ambient temperature.If more than 4 hrs - then transported in frozen condition.3 days abstinence is mandatory.
S0078	Calculus (Stone) Analysis by automated FTIR (other)	Madhya Pradesh	D	1390	Daily : 11:30am	Next day : 5:00pm	Collect sample in dry plain container, mention stone type and history
S0079	SYT gene rearrangement (Translocation: t(X;18) translocation) **	Madhya Pradesh	R	5305	Daily : 11:30am	7th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
S0083	Campylobacter Culture and Identification by Vitek/ MALDITOF - Stool	Madhya Pradesh	C	740	Daily : 9:00am to 9:00pm	5th day : 5pm	5gm of Stool specimen in Sterile leak proof container
S0084	SUCCINYLACETONE, BLOOD**	Madhya Pradesh	R	4795	Daily : 11:30am	8th day : 6:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. DBS card to be send in Ambient temperature only . Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age, or as close to discharge as possible if discharged

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0085	Succinylacetone Urine**	Madhya Pradesh	R	4795	Daily : 11:30am	7th day : 6:00pm	<p>prior to 24 hours of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Plz provide blood transfusion history and sample to be collected post 72 hrs after transfusion. Spot Urine: 10 ml of Spot urine required with clinical history and clinician number. Amino acid Plasma: 3 ml of separated EDTA plasma required . Sample to be shipped at Frozen condition . Sample stable at refrigerated condition 2 to 8 degree for 1 day. Clinical history, metabolic disorders or inborn errors of metabolism suspicion or provisional diagnosis of clinician required fasting recommended not mandatory. Amino acid supplementation or medications needs to be stopped based on physician recommendations</p> <p>Spot Urine: 10 ml of Spot urine required with clinical history and clinician number. Baby DOB, sample collection date and</p>

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
S0096	STRONTIUM ICPMS Blood	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	relevant clinical history mentioned . 3 ml EDTA Blood / Metal free EDTA Blood
S0097	Sexually Transmitted Infections (STI) Urethritis Extended Panel	Madhya Pradesh	D	6955	Wed : 9:00am	Next day : 5:00pm	Urethral discharge, Genital discharge, Blister fluid, Ulcer, Ulcer base, High Vaginal swab, Urethral Swab, Endocervical swab
S0103	Soft Tissue Sarcoma Panel**	Madhya Pradesh	R	40800	Daily 7.30 am	10th days	Paraffin Block
S0547	sFlt-1 (Soluble fms-like tyrosine kinase-1) Serum	Madhya Pradesh	D	4065	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0548	PIGF (Placental Growth Factor) Serum	Madhya Pradesh	D	4065	Mon & Thu : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0549	sFlt-1/PIGF ratio serum	Madhya Pradesh	D	6205	Daily cut off 7:30 AM	same day 9:00 PM	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S0551	Sanger Sequencing / Additional family member screening**	Madhya Pradesh	R	9180	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
S0552	Selenium, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
S0553	Selenium, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
S0554	Silver, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							before sending
S0555	Silver, Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
S0556	Steroids test by MS method	Madhya Pradesh	R	11220	Daily : 11:30am	10th day : 6:00pm	Require clinician provisional diagnosis, homeopathic, allopathic, ayurvedic or unani medicine history, duration of intake, dose and clincial symptoms if any. Reported in 7 days
S0558	SARS- CoV-2 (COVID 19) Detection by Real Time PCR - Nasopharangeal & Oropharangeal Swab	Madhya Pradesh	O	700	24hr registration	After 24hrs	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
S0559	SARS CoV2 Antigen Detection by Immunchromatography - Nasopharangeal Swab	Madhya Pradesh	O	300	Daily : 9:00am	Same day : 5:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
S0560	SARS CoV2 Antigen Reflex to SARS CoV2 real time RT PCR (COVID 19 reflex test)* - Nasopharangeal & Oropharangeal Swab	Madhya Pradesh	O	3900	Daily: 9:00 am Cut-off	24 hrs	For antigen test :at ambient temperature within an hour or if transported at 2 to 8 degree Celsius for 4 to 6 hours.
S0562	SARS- CoV-2 (COVID 19) viral load by Real Time PCR*	Madhya Pradesh	E	2000	Daily : 9:00am	Same day : 5:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
S0565	Soluble CD25 (IL 2 Receptor)	Madhya Pradesh	D	12485	Daily : 11:30am	5th day : 6:00pm	3 ml Serum
S1005	SCC-TA4 ANTIGEN - SERUM **	Madhya Pradesh	R	8975	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
S1006	Schistosomiasis - confirmation by EIA - Serum**	Madhya Pradesh	R	9690	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S1007	Schistosomiasis Ab by Hemagglutination - Serum**	Madhya Pradesh	R	5865	Daily : 9:00am	After 30 days : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S1020	Strongyloidiasis serology, Serum by EIA** - Serum	Madhya Pradesh	R	5815	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
S1038	Sickle Cell Anemia Mutation Prenatal Diagnosis	Madhya Pradesh	D	7060	Wed & Fri: 9:00am	7th day : 6:00pm	Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							condition or Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS) specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
S1039	Spinal Muscular Atrophy (SMA) Microdeletion	Madhya Pradesh	D	6420	Daily : 9:00am	7th day : 6:00pm	Amniotic Fluid (15-20 ml) in a sterile screw cap tube. Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition 3-5 ml EDTA Blood sample of Mother is mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius . Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G - FORM OF CONSENT (For invasive techniques). Clinical indication and family history mandatory.
S1039_c	Spinal Muscular Atrophy (SMA)-SMN1 Gene deletion	Madhya Pradesh	D	6420	Daily : 9:00am	7th day : 6:00pm	Chorionic Villi Sampling (CVS) 20 mg of chorionic villus (CVS)

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							specimen in a sterile, screw-top container filled with tissue culture transport medium or in sterile saline Sample stable for 1 day at Ambient temperature and 3 Days at 2 to 8 degree celsius. Sample to be preferably shipped at Referigerated condition Duly signed consent form as per PCPNDT must which include SAMPLE FORWARDING DOCTORS COVERING LETTER and FORM G -FORM OF CONSENT (For invasive techniques) Mother 3-5 ml EDTA Blood sample mandatory. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius .
S1040	Stone screening profile-1	Madhya Pradesh	D	3960	Daily: 9:00am	5th day : 6:00pm	5 ml of serum and 20 ml of 24 Hr Urine sample. Appropriate preservative needed.
S1041	Systemic sclerosis Profile	Madhya Pradesh	D	4280	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
S1901	Canine - SGOT (AST)	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
S1902	Canine - SGPT (ALT)	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
S1903	Skin Scraping Examination	Madhya	C	1230	Daily: 9am	Report after 24 hrs	Please contact lab before sending

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
S1904	Canine - Sodium	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
S1905	CANINE-Calculus (Stone) Analysis by automated FTIR	Madhya Pradesh	C	1455	Daily: 7am	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
S8322	SARS-CoV-2 (COVID 19) Detection (Qualitative) by Rapid	Madhya Pradesh	E	1975	Daily : 9:00am	Same day : 5:00pm	1 Swab with Sufficient Material in Viral Transport Medium (VTM)
S8399	Streptococcus Group A Antigen Detection	Madhya Pradesh	L4	1200	Daily: 9am to 9pm	After 6 hrs	Throat Swab
S8514	Salivary Cortisol	Madhya Pradesh	E	2800	Tuesday 9 am	next day 5 pm	4-5 ml saliva in sterile lolli sponge tube
T0004	t(1 19)by FISH E2A(ALL)	Madhya Pradesh	D	4175	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0005	t(1 19)By FISH E2A (ALL) Bone Marrow	Madhya Pradesh	D	4175	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0006	t(11;14) by FISH IGH/CCND1 (ALL, NHL) Bone Marrow**	Madhya Pradesh	D	5255	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0007	t(11;14) by FISH IGH/CCND1 (ALL, NHL) Blood**	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0008	t(11q23) by FISH MLL gene rearrangement	Madhya Pradesh	D	4385	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	(ALL, AML) Blood	Pradesh					vacutainer
T0009	t(11q23) by FISH MLL gene rearrangement (ALL, AML) Bone Marrow	Madhya Pradesh	D	4385	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0010	TEL/AML1(ALL) t(12 21) by FISH	Madhya Pradesh	D	4175	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0011	t(12 21) by FISH TEL/AML1 (ALL) Bone Marrow	Madhya Pradesh	D	4175	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Blood in Na-heparin vacutainer
T0012	PML Ra Ra t(15 17) (AML-M3) by FISH	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0013	t(15 17) by FISH PML Ra Ra (AML-M3) Bone Marrow	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0014	PML Ra Ra (Quanti) by PCR BLOOD**	Madhya Pradesh	D	5885	Daily: 7:30 am	7th day	3 ml of EDTA Whole Blood
T0015	PML Ra Ra (Quanti) by PCR BONE MARROW**	Madhya Pradesh	D	5885	Daily: 7:30 am	7th day	2-3 ml Bone Marrow in EDTA vacutainer
T0016	PDGFR Alpha rearrangement t(4q12) (HES) by FISH	Madhya Pradesh	D	5455	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0017	t(4q12) by FISH PDGFR Alpha Rearrangement (HES) Bone Marrow**	Madhya Pradesh	D	5455	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0018	PDGFR Beta rearrangement t(5q33) (HES) Blood by FISH	Madhya Pradesh	D	5455	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0019	PDGFR Beta rearrangement t(5q33) (HES) Bone marrow by FISH	Madhya Pradesh	D	5455	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0020	ETO/AML1 (AML-M2) t(8;21) by FISH	Madhya Pradesh	D	4600	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0021	t(8 21) by FISH ETO/AML1 (AML-M2) Bone Marrow	Madhya Pradesh	D	4600	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0024	Bcr-Abl (PHILADELPHIA) t(9 22) by FISH, BLOOD	Madhya Pradesh	D	4710	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0025	Bcr-Abl (PHILADELPHIA) t(9 22) by FISH, Bone Marrow	Madhya Pradesh	D	4710	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0028	T3-Free Free Tri-Iodothyronine Serum	Madhya Pradesh	B	300	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0029	T3-Total Tri Iodothyronine Serum	Madhya Pradesh	B	220	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0030	T4-Free Free Thyroxine Serum	Madhya Pradesh	B	300	Daily: 9am to 9pm	After 6 hrs	tests. 3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0031	T4-Total Thyroxine Serum	Madhya Pradesh	B	220	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0032	Tacrolimus By ECLIA Blood	Madhya Pradesh	C	4500	Daily: 1pm	Same Day 7.30pm	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
T0033	Tacrolimus By LCMS, Blood	Madhya Pradesh	C	4480	Tue , Thur and Sat: 9: 00 am	Same day 9:00 pm	Blood to be collected in 3ml EDTA vacutainer tube and shipped at Ambient temperature. The sample is stable for 6 hours and 7days at Ambient & 2-8C respectively. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
T0035	IGRAS Interferon gamma release assay by ELISA - Lithium Heparin Blood	Madhya Pradesh	C	2990	Mon, Wed, Fri : 7:30am	Next day : 5pm	5 ml of Lithium Heparin Blood
T0036	TCR-Alpha Beta	Madhya Pradesh	D	2460	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
T0037	TCR-Gamma Delta	Madhya Pradesh	D	2460	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							to be transported in Ambient condition(preferably). Clinical history and clinician number required.
T0038	TDT (Nuclear)	Madhya Pradesh	D	2460	Daily : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
T0039	Testicular Marker Profile	Madhya Pradesh	D	1390	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0040	Testosterone Bioavailable	Madhya Pradesh	C	2690	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0041	TestosteroneFreeSerum	Madhya Pradesh	C	1770	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0042	Testosterone Total Serum	Madhya Pradesh	B	700	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0043	Testosterone Profile	Madhya Pradesh	C	3920	Tue, Thu & Sat: 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0044	Tetanus Toxoid-IgG Quantitative test by EIA - Serum	Madhya Pradesh	C	4705	Fri : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0045	Thalassemia Studies	Madhya Pradesh	C	1680	Daily: 7.30am	After 9 Hrs	3 ml EDTA whole blood and serum sample
T0046	Thallium Blood	Madhya Pradesh	D	3265	Daily : 11:30am	5th day : 6:00pm	HeavyMetalSamplingguidelines.pdf
T0047	Thallium Serum	Madhya Pradesh	D	3265	Daily : 11:30am	5th day : 6:00pm	HeavyMetalSamplingguidelines.pdf
T0048	Thallium Urine Spot	Madhya Pradesh	D	3265	Daily : 11:30am	5th day : 6:00pm	HeavyMetalSamplingguidelines.pdf
T0049	Theophylline Serum	Madhya Pradesh	C	1155	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0050	Threonine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	tests. 3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
T0051	Threonine Quantitative Urine**	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
T0052	Thrombin Time (TT)	Madhya Pradesh	B	1035	Daily: 9am and 3pm and 9pm	After 6 hrs	1 ml of citrated platelet poor plasma in frozen condition.
T0053	Thrombotic Risk DNA Panel (FV, FII, MTHFR)	Madhya Pradesh	D	12840	Wed & Fri: 9:00am	10th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
T0054	Thrombophilia Profile-Maxi	Madhya Pradesh	D	16500	Mon and Thurs: 9am	Next day 5pm	Serum and 1 ml of citrated platelet poor plasma in frozen condition.
T0055	Thrombophilia Profile-Mini	Madhya Pradesh	D	9500	Mon and Thurs: 9am	Next day 5pm	Serum and 1 ml of citrated platelet poor plasma in frozen condition.
T0056	Thrombotic Profile	Madhya Pradesh	D	16585	Mon and Thurs: 9am	Next day 5pm	Serum and 1 ml of citrated platelet poor plasma in frozen condition.
T0057	Thyroglobulin Serum	Madhya Pradesh	C	1790	Daily: 9am to 9pm	After 8 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0058	Thyroglobulin Antibody ATA CSF	Madhya Pradesh	C	1510	Daily: 9am to 9pm	After 6 hrs	2 ml of CSF in leakproof container (foilin vu conical tube is recommended to prevent leakage).
T0059	Thyroglobulin Antibody ATA Serum	Madhya Pradesh	B	1310	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0060	Thyroid Antibodies-AtAb TPO and ATA,CSF	Madhya Pradesh	C	2465	Daily: 9am to 9pm	After 6 hrs	2 ml of CSF in leakproof container (foilin vu conical tube is recommended to prevent leakage).
T0061	Thyroid Antibodies-TPO and ATG,Serum	Madhya Pradesh	C	2130	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0062	Thyroid Comprehensive Profile-1	Madhya Pradesh	B	2945	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0063	Thyroid Comprehensive Profile-2	Madhya Pradesh	B	1745	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					ml serum for profiles and reflex tests.
T0064	Thyroid Marker Profile	Madhya Pradesh	D	4270	Wed and Sat: 9am	Next day 5pm	5 ml of serum
T0065	Thyroid panel-1 (T3/T4/TSH)	Madhya Pradesh	B	550	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0066	Thyroid panel - 2(FT3, FT4, TSH)	Madhya Pradesh	B	800	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0067	Thyroid Panel-3 (Reflex)	Madhya Pradesh	B	435	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0068	Thyroid Panel-4	Madhya Pradesh	B	470	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0069	Thyroxine Binding Globulin (TBG) Biochemical Serum	Madhya Pradesh	C	5935	Daily : 11:30am	5th day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0071	TNF AlphaTumour Necrosis FactorSerum	Madhya Pradesh	C	4370	Wed & Sat : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0072	TORCH PCR (Toxoplasma, CMV, Rubella and HSV 1&2) by PCR - Amniotic Fluid	Madhya Pradesh	C	8850	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
T0073	TORCH PCR (Toxoplasma, CMV, Rubella and HSV 1&2) by PCR - Blood	Madhya Pradesh	C	8850	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	3 ml of Whole Blood in EDTA Vacutainer
T0074	TORCH by PCR **	Madhya Pradesh	C	8850	Mon, Wed & Fri : 9:00am	Same day : 9:00pm	Chorionic Villus Sample In Sterile saline container
T0075	TORCH PCR (Toxoplasma, CMV, Rubella and HSV 1&2) by PCR - CSF	Madhya Pradesh	C	8850	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	2ml of CSF in Sterile leak proof tube or container
T0076	TORCH PCR (Toxoplasma, CMV, Rubella and HSV 1&2) by PCR - Urine	Madhya Pradesh	C	8850	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	15ml of Spot Urine in Sterile leak proof container
T0077	TORCH-10	Madhya Pradesh	C	3360	Daily: 9am	Same day 5pm	5 ml of serum
T0078	TORCH-13	Madhya Pradesh	C	3920	Mon: 9am	Next day 5pm	5 ml of serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0079	TORCH-4 IgG	Madhya Pradesh	C	1625	Daily: 9am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0080	Tacrolimus CYP3A5 genotyping**	Madhya Pradesh	R	8160	Tue & Thu: 11:30am	7th day : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
T0081	TORCH-4 IgM	Madhya Pradesh	C	1625	Daily: 9am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0083	TORCH-5 IgG	Madhya Pradesh	C	1790	Daily: 9am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0084	TORCH-5 IgM	Madhya Pradesh	C	1790	Daily: 9am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0085	TORCH-8	Madhya Pradesh	C	2910	Daily: 9am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0086	Toxoplasma Avidity test by EIA- Serum	Madhya Pradesh	C	1455	Mon & Thu : 11:30am	Next day : 5:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0087	Toxoplasma DNA Detection by PCR - Amniotic Fluid**	Madhya Pradesh	C	3695	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
T0088	Toxoplasma DNA Detection by PCR - Blood	Madhya Pradesh	C	3695	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	3 ml of Whole Blood in EDTA Vacutainer
T0089	Toxoplasma DNA Detection by PCR - Chorionic Villus	Madhya Pradesh	C	3695	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	2ml of CVS in Sterile leak proof tube or container
T0090	Toxoplasma DNA Detection by PCR - Cord Blood**	Madhya Pradesh	C	3695	Mon & Thu : 7:30am	Next day : 5pm	3ml of Cord Blood in EDTA Vacutainer
T0091	Toxoplasma DNA Detection by PCR - CSF	Madhya Pradesh	C	3695	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	2ml of CSF in Sterile leak proof tube or container

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0092	Toxoplasma DNA Detection by PCR - Urine	Madhya Pradesh	C	4200	Mon, Wed & Fri : 9:00am	Same day: 9:00pm	15ml of Spot Urine in Sterile leak proof container
T0094	Toxoplasma-IgG IgG Antibody by CLIA - Serum	Madhya Pradesh	B	720	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0095	Toxoplasma-IgG & IgM by CLIA - Serum	Madhya Pradesh	B	1415	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0096	Toxoplasma-IgM IgM Antibody by CLIA - Serum	Madhya Pradesh	B	720	Daily : 9:00am to 9:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0097	Treponema Pallidium Haemagglutination Assay by Agglutination - Serum	Madhya Pradesh	B	840	Daily : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0098	TPMT Enzyme ActivityThiopurine Methyl TransferasePlasma	Madhya Pradesh	D	4750	W,S: 7:30am	Same day 5pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
T0099	Thiopurine S-Methyltransferase (TPMT) Genotype	Madhya Pradesh	D	5670	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius .
T0100	TPO (Thyroid Peroxidase) Antibody(AMA) Microsomal Antibody CSF	Madhya Pradesh	C	1400	Daily: 9am to 9pm	After 6 hrs	2 ml of CSF in leakproof container (folin vu conical tube is recommended to prevent leakage).
T0101	TPO (Thyroid Peroxidase) Antibody (AMA) Microsomal Antibody Serum	Madhya Pradesh	B	1230	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0102	TransferrinSerum	Madhya Pradesh	B	1200	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0103	TRH STIMULATION TEST FOR PROLACTIN, SERUM	Madhya Pradesh	C	1510	Daily: 9am to 9pm	After 6 hrs	3 ml of serum (0, 30, 60 min)
T0104	TRH stimulation test For TSH, Serum	Madhya Pradesh	C	1230	Daily: 9am to 9pm	After 6 hrs	3 ml of serum (0, 30, 60 min)
T0105	Triglycerides Urine	Madhya Pradesh	A	410	Daily: 9am to 9pm	After 6 hrs	2 ml of Spot Urine
T0106	Triglycerides Serum	Madhya Pradesh	A	220	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0107	Triglycerides Body Fluid	Madhya Pradesh	A	410	Daily: 9am to 9pm	After 6 hrs	2 ml of Body Fluid

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0108	Tripeptidyl Peptidase 1 Blood	Madhya Pradesh	R	5710	Daily : 11:30am	10th day : 6:00pm	sample type - 5 ml Heparin Blood (sample required within 12 - 24 hrs of collection, Avoid collection on weekend & Holidays)
T0109	Trisomy 12 (CLL) by FISH**	Madhya Pradesh	D	4175	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0110	Trisomy 12 by FISH CLL Bone Marrow	Madhya Pradesh	D	4175	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0111	Trisomy 8 by FISH AML, MDS Blood	Madhya Pradesh	D	4175	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
T0112	Trisomy 8 by FISH AML, MDS Bone Marrow	Madhya Pradesh	D	4175	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
T0113	Trisomy (Aneuploidy) for Chromosome 13 & 21 By FISH Amniotic Fluid	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	20 ml of Amniotic Fluid In Sterile Falcon centrifuge tube
T0114	Trisomy (Aneuploidy) Chromosome 13 & 21 by FISH	Madhya Pradesh	D	5245	Daily : 9:00am	5th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
T0115	Trisomy (Aneuploidy) for Chromosome 13 & 21 By FISH Chorionic Villus	Madhya Pradesh	D	5245	Daily : 9:00am	3rd day : 6:00pm	20-30 mg Chorionic Villus in Sterile nutrient medium in 15 ml Falcon tubes(Medium will be provided by Metropolis).
T0116	Trisomy for Chromosome 13 & 21 By FISH Cord blood	Madhya Pradesh	D	5245	Daily : 9:00am	5th day : 6:00pm	3 ml of Prenatal Cord Blood In Na-Heparin Vacutainer
T0118	Trisomy(Aneuploidy)Chromosome13,18,21, X,Y by FISH	Madhya Pradesh	D	10970	Daily : 9:00am	3rd day : 6:00pm	20 ml of Amniotic Fluid In Sterile Falcon centrifuge tube
T0119	Trisomy for Chromosome 13, 18, 21, X, Y By FISH Blood	Madhya Pradesh	D	9415	Daily : 9:00am	5th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
T0120	Trisomy (Aneuploidy) for Chromosome 13, 18, 21, X, Y By FISH Chorionic Villus	Madhya Pradesh	D	9415	Daily : 9:00am	3rd day : 6:00pm	20-30 mg Chorionic Villus in Sterile nutrient medium in 15 ml Falcon tubes(Medium will be provided by Metropolis).
T0121	Trisomy for Chromosome 13, 18, 21, X, Y By FISH Cord blood	Madhya Pradesh	D	9415	Daily : 9:00am	5th day : 6:00pm	3 ml of Prenatal Cord Blood In Na-Heparin Vacutainer
T0123	Trisomy (Aneuploidy) Chromosome 18, X,Y by FISH	Madhya Pradesh	D	5670	Daily : 9:00am	3rd day : 6:00pm	20 ml of Amniotic Fluid In Sterile Falcon centrifuge tube
T0124	Trisomy for Chromosome 18, X & Y By FISH Blood	Madhya Pradesh	D	5670	Daily : 9:00am	5th day : 6:00pm	3 ml of Blood In Na-Heparin Vacutainer
T0125	Trisomy (Aneuploidy) for Chromosome 18,	Madhya	D	5670	Daily : 9:00am	3rd day : 6:00pm	20-30 mg Chorionic Villus in Sterile

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	X & Y By FISH Chorionic Villus	Pradesh					nutrient medium in 15 ml Falcon tubes(Medium will be provided by Metropolis).
T0126	Trisomy for Chromosome 18, X & Y By FISH Cord blood	Madhya Pradesh	D	5670	Daily : 9:00am	5th day : 6:00pm	3 ml of Prenatal Cord Blood In Na Heparin Vacutainer
T0128	Troponin-I, Serum	Madhya Pradesh	C	1770	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0129	Troponin-T, Serum	Madhya Pradesh	C	1850	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0130	TSH (Ultrasensitive) ECLIA Serum	Madhya Pradesh	B	380	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0131	TSH Receptor Antibody LATS, TSI - Serum	Madhya Pradesh	C	4600	W,S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0132	TSH-Thyroid Stimulating Hormone Neonatal Screen Dried blood spot	Madhya Pradesh	B	490	Daily : 11:30am	Next day : 5:00pm	Dried blood spot: 1 drop of heel prick or finger prick blood on special filter paper DBS card. Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age. Baby DOB, sample collection date and relevant clinical history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper. The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Provide blood transfusion

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							history and sample to be collected post 72 hrs after transfusion.
T0133	tTG Antibody-IgA Tissue Transglutaminase,Celiac disease Serum	Madhya Pradesh	C	1600	Daily : 11:30am	After 8hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0134	tTG Antibody-IgG Tissue Transglutaminase,Celiac disease Serum	Madhya Pradesh	C	1455	Mon and Thu: 7.30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
T0135	Typhi Dot IgM IgM antibody for Typhoid by Immunochromatography - Serum	Madhya Pradesh	B	640	Daily : 9:00am	Same day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0136	Typhoid DNA by Real Time PCR - Blood	Madhya Pradesh	C	6160	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
T0137	Typhoid DNA by Real Time PCR - CSF	Madhya Pradesh	C	6160	Daily : 10:00am	After 24hrs	2ml of CSF in Sterile leak proof tube or container
T0138	Typhoid Profile	Madhya Pradesh	D	2055	Daily: 9am to 9pm	5th day	3 ml of EDTA whole blood, 3 ml of serum and 1 blood culture bottle
T0139	Tyrosine Quantitative Plasma	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
T0140	Tyrosine Quantitative Urine**	Madhya Pradesh	R	6325	Thu : 11:30am	Tue : 5:00pm	10 ml of Spot urine
T0141	Tzanck smear (Microscopic examination)	Madhya Pradesh	B	435	Daily: 9am	After 8 hrs	Air Dried Smears From Vesicular fluid
T0143	TBE-TICKBORNE ENCEPHALITIS Antibodies by EIA - Serum	Madhya Pradesh	D	7810	Daily : 11:30am	After 4 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0145	Trypanosoma Cruzi IgG by EIA - Serum**	Madhya Pradesh	R	2295	Daily : 11:30am	20th day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
T0158	t(12;16):FUS rearrangement **	Madhya Pradesh	R	16320	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature.
T0160	t(X;18):SSX18 translocation **	Madhya Pradesh	R	16320	Daily : 11:30am	12th day : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0162	Tumor marker panel (270 genes), comprehensive, FFPE	Madhya Pradesh	D	47000		4 weeks	formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath/IHC report along with doctors number. FFPE block to be shipped at Ambient temperature. (FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report. 4 ml of EDTA whole blood or collection in special tubes as per test. Sample to be shipped at Ambient temperature for local samples and refrigerated condition for outstation samples. do not freeze or keep EDTA blood samples in contact with ice packs. EDTA samples stable for 2 days at Ambient temperature and 5 days at Refrigerated condition
T0168	NextGen Tumor Mutation Burden	Madhya Pradesh	E	43000	Daily 9:00 am	25th Days	(FFPE block) Tumor containing tissue fixed in 20% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
T0169	Thallium, hair by ICPMS	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
T0170	Tin, EDTA Whole Blood**	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
T0171	Tin, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
T0172	Tin, Spot Urine**	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0173	TUNGSTEN, EDTA WHOLE BLOOD	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
T0174	Tungsten, Hair tests by ICPMS	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
T0175	TUNGSTEN, SPOT URINE	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	Individual Metals and Minerals.pdf
T0182	TB CHECK - Pulmonary	Madhya Pradesh	D	10700	Daily : 7:30am	After 6 weeks : 6:00pm	Pulmonary Smear positive samples or Culture positive specimen
T0184	TB CHECK - Extra Pulmonary	Madhya Pradesh	D	10700	Daily : 7:30am	After 6 weeks : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
T0192	TB Spectrum - Pulmonary	Madhya Pradesh	D	18725	Daily : 7:30am	After 6 weeks : 6:00pm	Pulmonary Smear positive samples or Culture positive specimen
T0200	TRIO Clinical Exome sequencing - 26MB (80-100x) **	Madhya Pradesh	E	56000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
T0201	TB Spectrum - Extra Pulmonary	Madhya Pradesh	D	18725	Daily : 7:30am	After 6 weeks : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
T0502	TBC1D24-Gene Analysis	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
T0504	TRIO-Whole exome sequencing Sure SelectV5-45-50MB (80-100x)	Madhya Pradesh	E	56000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
T0505	TTN-Gene sequencing	Madhya Pradesh	E	20000	Daily: 7:30am	4 weeks	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambient and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, history is mandatory.
T0506	Trisomy for Chromosome 18, X & Y by FISH, Ovarian follicular fluid **	Madhya Pradesh	D	5670	Daily : 11:30am	5th day : 6:00pm	CytogeneticsandFISH.pdf
T0507	Thrombophilia Super Maxi	Madhya Pradesh	D	22000	Monday : 9am	Wednesday : 7 pm	2 ml of citrated platelet poor plasma in frozen condition ,3 ml serum and 3 ml of EDTA whole blood is required.
T0938	TPMT and NUDT15 gene analysis**	Madhya Pradesh	R	14995	Daily : 11:30am	20th day : 6:00pm	4 ml EDTA Whole blood
T0987	Tricyclic antidepressants Drug(TCA) (Qualitative), Urine	Madhya Pradesh	D	695	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
T0988	Taurine DBS**	Madhya Pradesh	D	6635	Thu : 11:30am	Tue : 5:00pm	ried blood spot: 1 drop of heel prick or finger prick blood on special

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T0989	2 step procedure for Screening and diagnosis of Gestational Diabetes	Madhya Pradesh	E	300	Daily: 9am to 9pm	After 6 hrs	<p>filter paper DBS card. DBS card to be send in Ambient temperature only . Newborn Screening sample is obtained in the hospital between 48 hours and 5 days of age, or as close to discharge as possible if discharged prior to 24 hours of age. Baby DOB, sample collection date and relevant clincial history mentioned on DBS card to be filled mandatory. Positive screen cases may require repeat testing or further workup in terms of confirmation by plasma or LCMS or GCMS method. No foreign substances should contaminate the filter paper.The blood should fill all printed circles and is applied evenly on one side of the filter paper only, free of layering and clots. The specimen should be dried in a horizontal position for at least three hours, away from direct contact with surfaces, heat, and sunlight. Plz provide blood transfusion histroy and sample to be collected post 72 hrs after transfusion. Sample to be shipped at Frozen condition . Sample stable at refrigerated condition 2 to 8 degree for 1 day. Clinical history, metabolic disorders or inborn errors of metabolsim suspicion or provisional diagnosis of clincian required fasting recommended not mandatory. Amino acid supplementation or medications needs to stoped based on physician recommendations</p> <p>Glucose Screen in pregnancy. ADA recommends 2 step procedure for</p>

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	Mellitus(ADA)						Screening and diagnosis of Gestational Diabetes Mellitus (GDM). Step 1: A 50 g glucose challenge test (GCT) is used for screening without regard to the time of last meal or time of the day.Step 2: If 1 hour GCT value is more than 140 mg/dl, 100g Oral Glucose Tolerance Test (OGTT) is recommended and plasma glucose is estimated at 0, 1, 2 and 3 hours.G0027_50(post glucose 1 hr 50 gms) and G0066 (GTT-Gestational 3 hrs 100 gms glucose). If Post Glucose 1hr 50 gms value Above 140 then reflex G0066 GTT-Gestational ,Three Hour(100 grams glucose) taken.
T1032	Tobramycin**	Madhya Pradesh	R	6020	Daily : 11:30am	5th day : 6:00pm	3ml of Serum sample required to be shipped at 2-8C. Do not use gel tubes. Sample is stable for 1day, 7 days & 30days at A, R and F temperatures respectively. General Instructions: 1. Please indicate history of the drug dose taken, time of sample collection (peak or trough levels, e.g C0 or C2) and clinical History along with the sample.
T1053	Tropheryma whipplei - Direct Diagnosis by PCR - Blood**	Madhya Pradesh	R	24890	Daily : 9:00am	After 30 days : 5pm	3 ml of EDTA Whole Blood in EDTA Vacutainer
T1170	BCR-ABL IS (International Scale) Transcript fusion Quantification by Real Time PCR - Blood	Madhya Pradesh	D	7250	Tue & Fri : 9:00am	3rd day : 6:00pm	3 ml of Whole Blood in EDTA Vacutainer
T1171	BCR-ABL IS (International Scale) Transcript fusion Quantification by Real Time PCR - Bone Marrow	Madhya Pradesh	D	6420	Tue & Fri : 9:00am	3rd day : 6:00pm	3 ml of Bone Marrow in EDTA Vacutainer
T1901	Canine - T3 (Total)	Madhya Pradesh	B	185	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T1902	Canine-T4-Free, Serum	Madhya Pradesh	B	275	Daily: 9am to 9pm	Report after 24 hrs	range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician. Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1903	Canine-T4 (Total)	Madhya Pradesh	B	220	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1904	Testosterone Total, Serum	Madhya Pradesh	B	765	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1905	CANINE-Thyroid panel-1	Madhya Pradesh	C	615	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1906	CANINE-Thyroid panel - 2	Madhya Pradesh	C	840	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1908	Canine-TSH (Ultrasensitive), Serum	Madhya Pradesh	B	380	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1910	Canine - Triglyceride	Madhya Pradesh	A	240	Daily: 9am to 9pm	Report after 24 hrs	Please contact lab before sending samples. Kits used are validated for Human samples. Reference range will be provided where available for canine or any other animal. Please correlate clinically and to be interpreted by the clinician.
T1911	Thyroid cancer mutation panel, FFPE Tissue	Madhya Pradesh	E	35000	Daily : 11:30am	After 4 weeks : 6:00pm	(FFPE block) Tumor containing tissue fixed in 10% neutral buffered formalin with duration of 6 to 72 hours and paraffin embedded. Clinical history with Histopath report.
T5041	TruHealth Screen	Madhya Pradesh	W	499			See Individual Test Details
T5585	TruHealth Youth	Madhya Pradesh	W	4000			See Individual Test Details
T5586	TruHealth Master	Madhya Pradesh	W	4500			See Individual Test Details
T5587	TruHealth Senior	Madhya Pradesh	W	5000			See Individual Test Details
T5588	TruHealth Women	Madhya Pradesh	W	5000			See Individual Test Details
T5589	TruHealth Senior Women	Madhya Pradesh	W	5000			See Individual Test Details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
T5610	Tru Health Heart Advance	Madhya Pradesh	W	6000	Mon and Thu: 7.30am	Same day 5pm	5 ml of serum
T5611	Tru Health Heart Basic	Madhya Pradesh	W	1500	Daily: 9am to 9pm	After 6 hrs	5 ml of serum,2 ml of flouride plasma
T5612	TruHealth Advance	Madhya Pradesh	W	2800			See Individual Test Details
T5628	Tru Health Heart Pro	Madhya Pradesh	W	3500	Daily: 9am to 9pm	After 6 hrs	5 ml of serum
T5705	T5705 : TruHeart Gold	Madhya Pradesh	W	2499	Daily	Next day 5.00 pm	4 ml Serum
T5709	TB Viram screening	Madhya Pradesh	D	6540	Daily : 9:00am	After 6 weeks : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer. 2ml of Fluoride Plasma. 2ml of Pulmonary or Extra Pulmonary specimen in Sterile leak proof container. Also note, kindly refer to individual test for the sample guideline.
T5710	TB Viram susceptibility	Madhya Pradesh	D	20490	Daily : 9:00am	After 4 weeks : 6:00pm	2ml of Sputum/BAL/other respiratory samples in Sterile leak proof container Pure Culture (Mtb Complex) Bottle Or MGIT Tube or LJ slant
T8006	Tissue processing	Madhya Pradesh	C	560			Tissue to be sent in adequate 10% Neutral Buffered Formalin in a Leak proof container along with relevant clinical & radiological details
T8500	Transplant Panel for Viral Infection	Madhya Pradesh	D	15000	Every Tuesday, Thursday and Saturday 9:00 AM	Same day 11:00 PM	2ml EDTA Plasma or Serum
U0001	U1-SnRNP (68 KDa)Serum	Madhya Pradesh	C	1850	W,S: 7:30am	Same day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
U0002	UGT1A1 Gene Polymorphism (TA Repeat)	Madhya Pradesh	D	6420	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Amibent Temperature and 7 Days at 2 to 8 degree celsius .
U0003	Urea, Serum	Madhya Pradesh	A	200	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
U0004	Urea Urine 24H	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	tests. 24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample. Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea , coffee , cold drink avoidance for accurate results and such history to be conveyed to lab.
U0005	Urea Breath Test For Helicobacter Pylori	Madhya Pradesh	C	2240	Daily : 9:00am	Same day : 5pm	NA
U0006	Urea Clearance Test Serum	Madhya Pradesh	A	600	Daily: 9am to 9pm	After 6 hrs	10 ml of 24 Hrs Urine & 2 ml of serum
U0007	Uric Acid Serum	Madhya Pradesh	A	170	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
U0008	Uric Acid Spot urine	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	10 ml of Spot urine
U0009	Uric Acid Urine 24H	Madhya Pradesh	A	300	Daily: 9am to 9pm	After 6 hrs	24 H urine: 24 hour Urine should be collected in mug and transferred to relevant Urine container with or without preservative. 20 ml of aliquot to be send to lab in refrigerated condition along with clinical history and mention Urine volume of 24 Hr collected sample for calculation purpose. Time of collection, start time to end time and exclude first urine sample and include last voided urine sample.

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							Many analytes like VMA, 5 HIAA, adrenaline may require drug avoidance as well as bananas, tea, coffee, cold drink avoidance for accurate results and such history to be conveyed to lab.
U0010	Urobilinogen Urine	Madhya Pradesh	A	110	Daily: 9am to 9pm	After 6 hrs	15ml of Spot Urine in Sterile leak proof container
U0014	Urea Urine Spot	Madhya Pradesh	A	305	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
U1801	Ultrasound (USG) by Appointment- Biophysical profile	Madhya Pradesh	O	2500			
U1803	Ultrasound (USG) by Appointment- Follicular Study	Madhya Pradesh	O	1500			
U1804	Ultrasound (USG) by Appointment- Full Abdomen with pelvis	Madhya Pradesh	O	2950			
U1805	Ultrasound (USG) by Appointment- K.U.B	Madhya Pradesh	O	1500			
U1806	Ultrasound (USG) by Appointment- Musculoskeletal (Shoulder, Elbow, Ankle, Knee, etc)	Madhya Pradesh	O	1500			
U1807	Ultrasound (USG) by Appointment- Pelvis	Madhya Pradesh	O	1500			
U1808	Ultrasound (USG) by Appointment- Pregnancy - 2nd trimester	Madhya Pradesh	O	2000			
U1809	Ultrasound (USG) by Appointment- Pregnancy - 3rd trimester	Madhya Pradesh	O	1750			
U1810	Ultrasound (USG) by Appointment- Pregnancy - early	Madhya Pradesh	O	1500			
U1811	Ultrasound (USG) by Appointment- Pregnancy - level 3 Scan	Madhya Pradesh	O	2000			
U1812	Ultrasound (USG) by Appointment- Small parts (Thyroid, Scrotum, etc)	Madhya Pradesh	O	1500			
U1813	Ultrasound (USG) by Appointment- Trans rectal with K.U.B	Madhya Pradesh	O	1600			
U1814	Ultrasound (USG) by Appointment- Upper	Madhya Pradesh	O	1500			

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	abdomen	Pradesh					
U8003	Urine Sugar Fasting	Madhya Pradesh	A	25	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
U8004	Urine sugar, 2 hrs post prandial	Madhya Pradesh	A	25	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
U8004_PP	Urine sugar, post glucose	Madhya Pradesh	A	25	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
U8005	Urine sugar, random	Madhya Pradesh	A	25	Daily : 9:00am to 9:00pm	After 4hrs	15ml of Spot Urine in Sterile leak proof container
V0001	Valproic Acid Sodium Valporate Serum	Madhya Pradesh	C	980	Daily: 9am to 9pm	After 6 hrs	3 ml of Serum (Mention Time of Drug)
V0002	Vancomycin Serum	Madhya Pradesh	R	6630	Daily : 11:30am	4th day : 6:00pm	3 ml of Serum (Mention Time of Drug)
V0003	VDRL RPR by Flocculation - Serum	Madhya Pradesh	B	200	Daily : 9:00am to 3:00pm	After 6 hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
V0004	VDRL by Flocculation - CSF	Madhya Pradesh	C	480	Daily : 9:00am, 12:00pm, 4:00pm, 8:00 pm	After 6 hrs	2ml of CSF in Sterile leak proof tube or container
V0005	Voltage Gated Calcium Channel (VGCC) Antibody**	Madhya Pradesh	R	10200	Daily : 11:30am	After 4 weeks : 6:00pm	5 ml Serum
V0006	VGKC Antibody, CSF	Madhya Pradesh	D	9845	Mon & Thu : 11:30am	Same day : 9:00pm	2 ml of CSF in sterile screw capped container.
V0007	VGKC Antibody, Serum	Madhya Pradesh	D	9500	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0008	Vitamin A (Retinol), Serum	Madhya Pradesh	D	4280	Mon & Thu : 11:30am	Next day : 5:00pm	2 ml Serum Or EDTA Plasma In aluminium Foil
V0009	Vitamin B1 (Thiamine), EDTA Whole blood by HPLC (V0009)	Madhya Pradesh	D	3530	Wed & Sat : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood covered In Aluminium Foil
V0010	Vitamin B12 Cyanocobalamin Serum	Madhya Pradesh	B	1100	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0011	Vitamin B12 Active Holo Transcobalamin Serum	Madhya Pradesh	C	1455	Daily: 9am to 9pm	After 6 hrs	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0012	Vitamin B2 Riboflavin EDTA Whole Blood	Madhya	D	3425	Wed : 11:30am	Next day : 5:00pm	3 ml of EDTA whole blood covered

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
		Pradesh					In Aluminium Foil
V0013	Vitamin B6 Pyridoxine EDTA Whole Blood	Madhya Pradesh	D	3425	Sat : 11:30am	Mon : 7:00pm	3 ml of EDTA whole blood covered In Aluminium Foil
V0014	Vitamin C, Serum/Heparin plasma by HPLC	Madhya Pradesh	D	3750	Tue & Fri: 11:30am	Next day : 5:00pm	2 ml Lithium heparin Plasma or serum In Fasting condition covered In Aluminium Foil
V0015	Vitamin D Total-25 Hydroxy Plasma/Serum	Madhya Pradesh	B	1600	Daily: 9am to 9pm	After 6 hrs	2 ml of serum / EDTA Plasma
V0016	Vitamin D plus profile	Madhya Pradesh	A	1800	Daily: 9am to 9pm	After 6 hrs	2 ml Serum / Edta Plasma
V0017	Vitamin D-1,25-Dihydroxy Serum	Madhya Pradesh	C	3810	Daily: 7.30am	Next day 5pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0018	Vitamin E Tocopherol Serum	Madhya Pradesh	D	4280	Mon & Thu : 11:30am	Next day : 5:00pm	2 ml Serum Or EDTA Plasma In aluminium Foil
V0019	Vitamin profile	Madhya Pradesh	C	3135	Daily: 9am to 9pm	After 8 hrs	3 ml Serum and 3 ml of EDTA plasma
V0020_LC	Vanillyl Mandelic Acid (VMA), LC MSMS, 24hrs urine	Madhya Pradesh	C	3975	Mon and Thu: 1pm	Next day 5:00 pm	10 ml 24 Hrs Urine
V0021	VRE Screen Vancomycin Resistant Enterococci by Vitek/ MALDITOF - Perianal swab	Madhya Pradesh	C	1120	Daily : 9:00am to 9:00pm	3rd day : 5pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube
V0022	Vw-Von Willebrand Antigen, Citrated plasma	Madhya Pradesh	D	9095	Mon : 9:00am	Same day : 8pm	1 ml of citrated platelet poor plasma in frozen condition.
V0023	VZV DNA Detection by Real Time PCR - Amniotic fluid	Madhya Pradesh	C	3695	Mon & Thurs : 9:00am	Next day : 5:00pm	10ml of Amniotic Fluid in Sterile Leak proof Falcon Tube/ Sterile Container
V0024	VZV DNA Detection by Real Time PCR - Blood	Madhya Pradesh	C	3695	Mon & Thurs : 9:00am	Next day : 5:00pm	3 ml of Whole Blood in EDTA Vacutainer
V0025	VZV DNA Detection by Real Time PCR - CSF	Madhya Pradesh	C	3695	Mon & Thurs : 9:00am	Next day : 5:00pm	2ml of CSF in Sterile leak proof tube or container
V0026	VZV DNA Detection by Real Time PCR - CVS	Madhya Pradesh	C	3695	Mon & Thurs : 9:00am	Next day : 5:00pm	2ml of CVS in Sterile leak proof tube or container
V0027	VZV DNA Detection by Real Time PCR - Swab	Madhya Pradesh	C	3695	Mon & Thurs : 9:00am	Next day : 5:00pm	1 Swab in Normal Saline with Sufficient Material in Sterile Falcon Tube

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
V0029	VZV IgG Antibody by CLIA - Serum	Madhya Pradesh	C	1905	Daily : 9:00am	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
V0030	VZV IgM Antibody by CLIA - Serum	Madhya Pradesh	C	1905	Daily : 9:00am	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
V0031	VDRL reflex Syphilis Total ABS - Serum	Madhya Pradesh	B	710	Mon & Thu : 11:30am	Next day : 6:00pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
V0032	Voriconazole EDTA Plasma	Madhya Pradesh	D	5885	Tue & Fri: 11:30am	Next day : 5:00pm	3 ml of EDTA plasma or heparin plasma to be sent in aliquot tube at proper temperature.
V0034	VGKC - CSF & Serum	Madhya Pradesh	E	13500	Mon & Thu : 11:30am	Same day : 9:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.2 ml of CSF in sterile screw capped container.
V0036	Vitamin B Complex profile	Madhya Pradesh	E	5000	Daily: 9am	7th day	3 ml of serum and 3 ml of EDTA whole bloodIn aluminium Foil
V0037	Vitamin K1**	Madhya Pradesh	R	4490	Daily : 11:30am	15th day : 6:00pm	5 ml serum without gel tube. Wrap in aluminium foil. Frozen.
V0038	VANADIUM ICPMS BLOOD	Madhya Pradesh	R	2805	Daily : 11:30am	10th day : 6:00pm	3 ml EDTA Blood / Metal free EDTA Blood
V0144	VITAMIN D GOLD(covers Vitamin D2 , vitamin D3 & D total)	Madhya Pradesh	D	2675	Daily : 11:30am	3rd day : 6:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0148	Vasculitis profile, Serum	Madhya Pradesh	D	4710	Tue, Thur, Sat 7:30 am	Next day 5 pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V0149	Vesicobullous Disorder profile	Madhya Pradesh	E	9000	Mon : 11:30am	Next day : 5:00pm	3 ml of serum for single test and 5 ml serum for profiles and reflex tests.
V5023	Vascular Endothelial Growth Factor (VEGF) â€” serum**	Madhya Pradesh	D	6420	Daily 7.30 am	4 Weeks	3 ml Serum
V7005	Vascular Endothelial Growth Factor (VEGF) serum**	Madhya Pradesh	R	6120	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml Serum
V8123	Vitamin B3 (Niacin)	Madhya Pradesh	D	3850	Daily 7.30 am	7th Days	4 ml Serum
V8124	Vitamin B7 (Biotin)	Madhya Pradesh	D	3850	Daily 7.30 am	7th Days	4 ml Serum
V8125	Vitamin B5	Madhya	D	3850	Daily 7.30 am	7th Days	4 ml Serum

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
W0001	Warfarin Sensitivity CYP2C9*2,CYP2C9*3 &VKORC1 Mutation Blood**	Pradesh Madhya Pradesh	R	14690	Daily : 11:30am	10th day : 6:00pm	Clinical history and TRF mandatory. Blood : 3 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius . Do not send sample in frozen condition. Sample can be shipped at Ambient or refrigerated condition based on sample stability. below details are mandatory along with sample for the test "Genotyping for warfarin dose determination" 1.Clinical details 2.Medication details (Dose) 3.Height and Weight
W0002	Water analysis Chemical**	Madhya Pradesh	R	1410	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
W0003	Water Culture for coliforms Water	Madhya Pradesh	C	670			50ml of Water in Sterile leak proof container
W0006	WBC:Total & Differential Counts, EDTA blood.	Madhya Pradesh	A	220	Daily: 9am to 9pm	After 6 hrs	3 ml EDTA whole Blood
W0007	WBC:Total & Differential Counts, Body Fluid	Madhya Pradesh	A	275	Daily: 9am to 9pm	After 6 hrs	Body fluid to be send in sterile EDTA container in refrigerated condition. Clinical history and clinician number required.Body fluids(CSF, Pleural, Peritoneal, Pericardial, synovial fluid)
W0009	WBC-Total Counts, Body Fluid.	Madhya Pradesh	A	220	Daily: 9am to 9pm	After 6 hrs	Body fluid to be send in sterile EDTA container in refrigerated condition. Clinical history and clinician number required. Body fluids (CSF, Pleural, Peritoneal, Pericardial, synovial fluid)
W0010	Weil Felix Test For Rickettsia by Agglutination - Serum	Madhya Pradesh	B	1070	Daily : 7:30pm	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
W0011	Widal Test for Typhoid by Agglutination - Serum	Madhya Pradesh	B	200	Daily : 7:30pm	Next day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
W0013	WilsonS Disease Panel	Madhya Pradesh	D	4495	Mon and Thu: 9am	Same day 5pm	3 ml of EDTA whole blood, 2 ml of flouride plasma random, 3 ml of serum, 20 ml aliquot of 24 hrs urine (Volume reqd)
W0014	West nile virus IgG by EIA - Serum	Madhya Pradesh	D	3745	1st & 3rd Fri : 9:00am	4th day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
W0015	West nile virus IgM by EIA - Serum	Madhya Pradesh	D	3745	1st & 3rd Fri : 9:00am	4th day : 5pm	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
W0017	West nile virus DNA Detection by Real Time PCR - EDTA Blood	Madhya Pradesh	D	6100	Daily : 10:00am	After 24hrs	3 ml of Whole Blood in EDTA Vacutainer
W0018	Whole mitochondrial genome sequencing**	Madhya Pradesh	R	14000	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
W0019	WT1 MUTATION DETECTION**	Madhya Pradesh	R	7700	Daily : 11:30am	12th Day : 6:00pm	4 ML EDTA whole blood
W0020	Whole Exome panel**	Madhya Pradesh	R	25500	Daily : 11:30am	After 4 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory.
W0021	Whole genome sequencing Mean 30x	Madhya Pradesh	R	95550	Daily : 11:30am	After 6 weeks : 6:00pm	3 ml of EDTA whole blood required to be send at Ambient temperature for local samples and refrigerated condition for outstation samples. Samples stable for 2 days at Ambeint and 2 to 3 days for refrigerated condition. Do not keep

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
W0022	Whole Slide Imaging per slide	Madhya Pradesh	C	560		2nd Day	EDTA whole blood in contact with ice packs. Do not freeze. TRF, Genetci history is mandatory. Well mounted H &E stained Or IHC stained slide, Doctors TRF & E mail ID
W0023	Water chemistry panel - 3**	Madhya Pradesh	R	4590	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
W0024	Water test profile, Routine	Madhya Pradesh	R	2550	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
W0025	Water test profile, Comprehensive	Madhya Pradesh	R	9180	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
W0026	Water chemistry panel - 1**	Madhya Pradesh	R	610	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
W1802_F	Wellness packages (Body watch) by Appointment- Fitness check-Female Below 45 years	Madhya Pradesh	W	4500			See Individual Test Details
W1802_M	Wellness packages (Body watch) by Appointment- Fitness check-Male Below 45 years	Madhya Pradesh	W	4500			See Individual Test Details
W1803_F	Wellness packages (Body watch) by Appointment- Life style check -Female Below 45 years	Madhya Pradesh	W	8500			See Individual Test Details
W1803_M	Wellness packages (Body watch) by Appointment- Life style check -Male Below 45 years	Madhya Pradesh	W	8500			See Individual Test Details
W1804_F	Wellness packages (Body watch) by Appointment- Vitality check-Female Below 45 years	Madhya Pradesh	W	16500			See Individual Test Details
W1804_M	Wellness packages (Body watch) by Appointment- Vitality check-Male Below 45 years	Madhya Pradesh	W	16500			See Individual Test Details
W1805_F	Wellness packages (Body watch) by Appointment- Fitness check-Female Above	Madhya Pradesh	W	6500			See Individual Test Details

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
	45 years						
W1805_M	Wellness packages (Body watch) by Appointment- Fitness check-Male Above 45 years	Madhya Pradesh	W	6500			See Individual Test Details
W1806_F	Wellness packages (Body watch) by Appointment- Life style check-Female Above 45 years	Madhya Pradesh	W	12500			See Individual Test Details
W1806_M	Wellness packages (Body watch) by Appointment- Life style check-Male Above 45 years	Madhya Pradesh	W	12500			See Individual Test Details
W1807_F	Wellness packages (Body watch) by Appointment- Vitality check- Female Above 45 years	Madhya Pradesh	W	20500			See Individual Test Details
W1807_M	Wellness packages (Body watch) by Appointment- Vitality check-Male Above 45 years	Madhya Pradesh	W	20500			See Individual Test Details
X0001	X And Y (Sex Mismatch) by FISH	Madhya Pradesh	D	4710	Daily : 9:00am	3rd day : 6:00pm	3 ml Blood in Na-heparin vacutainer
X0002	X And Y (Sex Mismatch) by FISH Bone Marrow Transplant Bone Marrow	Madhya Pradesh	D	4710	Daily : 9:00am	3rd day : 6:00pm	2-3 ml Bone Marrow in Na-heparin vacutainer
X8011	XDR Rapid plus (For extra pulmonary samples)	Madhya Pradesh	D	6525	Daily : 7:30am	Same day : 6:00pm	Tissue, Lymph Node, Endometrium, Biopsy, Cervical Lymph Node, FNAC, Biopsy 2ml sample in sterile container (Ascitic fluid, Pericardial fluid, Urine, Gastric lavage, Pleural fluid, Pus, Peritoneal Fluid, CSF, Synovial fluid)
X8012	XDR Rapid plus (For pulmonary samples)	Madhya Pradesh	D	6525	Daily : 7:30am	Same day : 6:00pm	2ml of Pulmonary specimen in Sterile leak proof container
Y0001	Y-Chromosome Microdeletion Analysis	Madhya Pradesh	D	7275	Tue & Fri : 9:00am	7th day : 6:00pm	3-5 ml of EDTA whole blood required. Sample stable for 3 Days at Ambient Temperature and 7 Days at 2 to 8 degree celsius. Clinical history and TRF mandatory.
Y0002	Yersinia Enterocolitica Culture by Vitek/ MALDITOF - Stool	Madhya Pradesh	C	2410	Daily : 9:00am to 9:00pm	3rd day : 5pm	5gm of Stool specimen in Sterile leak proof container
Y0005	Yo Antibody Neuronal antibody for Paraneoplastic syndromes Serum	Madhya Pradesh	D	4815	Daily: 11:30am	Next day 8pm	3 ml of serum for single test and 5 ml serum for profiles and reflex

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
Z0001	ZAP-70	Madhya Pradesh	D	5885	Daily : 11:30am	4th day : 6:00pm	tests. 3 ml of EDTA whole blood required. Sample stable for 2 days at Ambient temperature and sample to be transported in Ambient condition(preferably). Clinical history and clinician number required.
Z0005	Zinc Urine	Madhya Pradesh	D	2355	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
Z0006	Zinc, Blood	Madhya Pradesh	D	2245	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
Z0007	Zinc, Serum	Madhya Pradesh	D	2300	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
Z0009	Zinc Urine 24H**	Madhya Pradesh	D	2245	Daily : 11:30am	3rd day : 6:00pm	HeavyMetalSamplingguidelines.pdf
Z0011	Zika Virus RNA Detection by Real Time PCR - Serum	Madhya Pradesh	D	6955	Daily : 10:00am	After 24hrs	3ml of Serum in Gel Vacutainer or in Plain Sterile Vacutainer
Z0012	Zinc , Water tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	2 litter of water in leakproof sterile container. Storage and shipping at 2 to 8 degree celsius.
Z0013	Zinc, Hair tests by ICPMS**	Madhya Pradesh	R	1835	Daily : 11:30am	10th day : 6:00pm	1 gram of hair in leak proof container. Ambient temperature storage and shipping. Plz enquire before sending
Z0014	Zika Virus RNA Detection by Real Time PCR - Urine	Madhya Pradesh	D	6955	Daily : 10:00am	24 Hrs	15ml of Spot Urine in Sterile leak proof container
Z0015	Zinc Transporter 8 (ZnT8) Antibody	Madhya Pradesh	D	3210	Mon & Thu : 9:00 am	Next day : 5:00pm	3 ml Serum
Â I0390	IHC stain only HHV 8	Madhya Pradesh	D	3500	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
Â I0392	IHC stain only MDM 2	Madhya Pradesh	D	3500	Daily : 9:00am	3rd day	Single Paraffin block / section on charged slide
Â I0393	IHC single marker with reporting MDM2	Madhya Pradesh	D	4500	Daily : 9:00am	4th day	Paraffin blocks (Formalin fixed paraffin embedded) Please attach original Histopathology report and all relevant clinical and radiological

Test Code	Test Name	location	Category	Price(Rs.)	Testschedule	Reportedon	Sample Guideline
							detail